DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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
ANNEX I: WIRING DIAGRAM

(This diagram is just for you reference. Manufacturer reserves the right to modify or improve it).

provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the dry bath (burns, hurts...) or damages in nearby equipments.

In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the dry bath and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

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INSTRUCTIONS ON ENVIRONMENT PROTECTION

Do not dispose of this product in the usual household garbage at the end of it's life cycle; hand it over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal.

The materials are recyclable as mentioned in it's marking. By recycling, material recycling or other forms or re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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
If you are not using the equipment for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Security

- ◆ The dry bath must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the dry bath in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the equipment near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the dry bath for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the dry bath by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the dry bath (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The dry bath is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the equipment by yourself, since you will lose the warranty and may

The system will save the data and increase to 100 °C automatically. The cursor will disappear as showed in chart.

◆ When the system temperature reaches 100 °C, the temperature of PV is actual temperature and the temperature of "AdjTemp" is 100 °C with the cursor. The "Adj" and "*" is flickering.

◆ After 20 min, the actual temperature of thermometer is 98.0 °C. Press buttons remove to the left or right to make the cursor rest on the digital you want to change. Press increasing or decreasing to adjust.

◆ After the finish of adjustment, the screen displays as chart.

◆ After the finish of adjustment, press button removes to the right and ADJ at the same time and keep off at once. After 8s the cursor display and the system accepts the adjustment. The screen returns to the display as the chart.

After the adjustment, the temperature displayed is the same the actual temperature

NOTE: During the adjustment, press the button removes to the right and ADJ contemporary to cancel the adjustment. The system keeps the last adjustment, so don't press casually.

5. Exchange of the metal block

Open the transparent lid and pull out the four screws which fix the block to the heating board with the screwdriver, close the lid, and put out the block from the main engine. Take another model of block, steadily lays aside on the main engine. Put in the instalment holes the screws, and fix with the spanner the metal block on the instrument.

The bath doesn't include any block. The following blocks are available:

Code	Capacity
90603021	96 wells of 0.2 mL
90603022	54 x 0.5 mL
90603023	35 x 1.5 mL
90603024	35 x 2.0 mL

PV: 40.1 C	Adj *
AdjTemp=	40.1 C

PV: 100.1 C	Adj *
AdjTemp=	100.0 C

PV: 100.1 C	Adj *
AdjTemp=	98.0 C

PV: 100.1 C	Adj *
AdjTemp=	100.0 C

PV: 98.9c	T:1
SV: 40.0	1:00M

PV: 40.0C	T:1
SV: 40.0	1:00M

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System-Testing
.....

◆ Press the button STOP to stop the operation. The flicker of the "*" and ":" will be stopped. The operation will re-run after pressing the START button and the symbol of "*" and ":" will flicker again.

◆ Pressing the STOP button and holding 3 seconds, the program operation will be stopped. The display will return to the beginning picture. The symbol of "*" will disappear.

3. Temperature error adjustment

The temperature of the Instrument has been adjusted before it is sold out. But if there is deviation between the actual temperature and the displayed temperature due to some reasons, you can do as follows to correct the error.

Notes: The Instrument uses double temperatures adjustment to ensure its veracity. This means it is linearly adjusted on 40 °C and 100 °C two points. The temperature veracity will be within ± 0.5 °C after the double temperature adjustment.

Both circumstances and the block temperature should be lower than 35 °C.

4. Adjustment methods as follows:

◆ After the start up of the Instrument, it enters waiting interface. Make sure the temperature in display is below 35 °C. If the temperature is higher than 35 °C, you should wait until the temperature is below of 35 °C.

◆ Inject olefin oil into one of the cone-shaped wells, and then put a thermometer into this well (Make sure the precision of the thermometer should be within 0.1 °C and the temperature ball should be absolutely immersed into the cone-shaped well). Heat insulation material is needed on the block to separate it from the circumstance.

NOTE: Please read the actual value after 20minutes' constant temperature to ensure the adjustment precision.

◆ Press button removes to the right and ADJ button contemporary to enter the temperature adjustment interface as in the drawing. The temperature begins to rise to 40.0 °C automatically and "*" is flickering in the right corner.

PV: 20.5 C	Adj *
AdjTemp=	20.5 C

The temperature of "PV:" and "AdjTemp" is current temperature 20.5 °C. When the temperature reaches 40 °C, the temperature of "PV" is current temperature of block while the temperature of "AdjTemp" is still 40 °C with cursor. The "Adj" and "*" is flickering.

PV: 40.1 C	Adj *
AdjTemp=	40.0 C

◆ After 20min, the actual temperature of thermometer is 38.8 °C. Press buttons removes to the left or right to make the cursor rest on the digital you want to change. Press increasing or decreasing button to adjust and press START to confirm.

PV: 40.1 C	Adj *
AdjTemp=	38.8 C

◆ Press the button of SET and keep off at once. Now, after the temperature set (SV) there is a cursor "_". Press button removes to the right or to the left to move the cursor, and press increasing or decreasing button to change the digital on the cursor into the one required. T:1 means controlling temperature for the first segment. The time also can be set using the former way. (For example: as show in the chart, the current temperature is 30 °C, the set temperature is 37 °C, the first temperature control segment and the temperature control time is 20 h).

PV: 30.0C	T:1
SV: _37.0C	20:00M

◆ Press the button of SET again, and can be set using the former way. T:2 means controlling temperature for the second segment. (for example: as show in the chart, the current temperature is 30 °C, the set temperature is 99 °C, the second temperature control segment and the temperature control time is 1h 30 min).

PV: 30.0C	T:2
SV: _99.0C	1:30M

NOTE: After the first temperature control segment is completed, it steps into the second temperature control segment automatically.

◆ After the completion of the set, press the SET again. About 8 seconds later, the cursor will disappear. The system accepts the set and operates under the set. The picture returns to T:1 temperature-control segment.

PV: 30.0C	T:1
SV: 37.0C	20:00M

NOTE: After pressing the button of START, the operation starts according to the set finally. If the time set is 0 h 0 min, the operation starts after pressing START button. When the temperature reaches its set temperature, the instrument will keep this temperature forever.

2. Run and stop functions

◆ After pressing START and keeping off at once, the operation runs at the set. The symbol of "*" will be displayed and flickering in the circle of 2 s. After the first segment is completed, it steps into the second segment automatically.

PV: 32.5C	T:1 *
SV: 37.0 C	20:00M

PV: 37.0 C	T:2 *
SV: 99.0 C	01:30M

NOTES:

1. After the start of operation, the temperature PV displays the current temperature of the module. When the current temperature reaches the set temperature, the time is counted backwards and the symbol of ":" will flicker in the circle of 2 s.

2. When the second temperature-control is finished, the buzzer alarms and displays the signal of "PROGRAM END".

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equipment to the net, make sure that it does not present any damage because of the shipment. In case the dry bath presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Net wire
- Instruction manual
- Warranty certificate

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT

Please put the dry bath on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the dry bath near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

Avoid inflammable or toxic substances in the working area.

The dry bath is supplied with a Schuko standard wire

Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz ± 10% provided with earth wire and to the other end to the dry bath connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz ± 10%.

Setting up

1. Temperature and time set


◆ The LCD will display the picture as the chart when the instrument powers on and the instrument goes into the initial state with the sound of "du...";

System-Testing

◆ About 6s later, the figure 30.0 (P) is the block's current temperature; the figure 37.0 (S) is the former set temperature. T: 1 means controlling the temperature for the first segment. Time display 10:00 is the former set time.

PV: 30.0c	T: 1
SV: 37.0c	10: 00M

- 3.1 Display of the current temperature
- 3.2 Display of the set temperature
- 3.3 The first segment
- 3.4 Operation remind sigma
- 3.5 Display of the set time

3. TECHNICAL SPECIFICATIONS

1. Normal working condition:

- The using environmental temperature: 5 C - 30 C
- The relative humidity: ≤ 70%
- The using power: AC 220V, 50Hz

2. The basic parameters and performance:

Code	50603010
Model	603/10
Temperature range	RT +5 °C - 100 °C
Timer set	1min - 99h 59 min
Accuracy of the temperature	≤ 0.3 °C
Display accuracy	0.1 °C
Heating time (from 20 °C to 100 °C)	12 min
Power supply	120 W
Fuse	250V 2.5A ?5×20
Dimensions (mm)	300 x 212 x 180
Weight	2.5 kg

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4. INSTALLATION / SETTING UP

Preliminary inspection

Unwrap the dry bath, take off the involving plastic and take off the polispan protection in which it comes fitted. Take off all the protective items and, without connecting the

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
- 1.1 Standard block
- 1.2 Operation panel

- 1.3 Power switch
- 1.4 Power socket
- 1.5 Fuse socket

Operation panel and display diagram:

Button instruction:

2.1 ADJ: Temperature adjustment button: Used together with button removes to the right to adjust the error between actual temperature and the displayed temperature.

2.2 SET: Set button: used together with increasing and decreasing button to set the operation temperature and time required

2.3 START: Run button: Start the operation.


2.4 STOP: Stop / pause button: when you press the button once, the operation is paused; if you press it again, the operation will stop.

2.5 ▲ : Increasing button of the temperature or time

2.6 ▼ : Decreasing button of the temperature or time

2.7 ► : Button removes to the right: Used together with ADJ button to adjust the error between actual temperature and the displayed temperature.

2.8 ◀ : Button removes to the left.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish	2-12
English	13-23
CE certificate	22

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	14
2. DESCRIPTION	15
3. TECHNICAL SPECIFICATIONS	16
4. INSTALLATION / SETTING UP.....	16
5. MAINTENANCE AND CLEANING	21
ANNEX I: WIRING DIAGRAM	23

1. USES OF THE INSTRUMENT

The dry bath (Incubator) 603/10 is a new product controlled by the miniature computer, which can be widely used for the preservation and reaction of the samples, the amplification of DNA, the pre-denaturation of the electrophoresis and blood serum coagulation.

The dry bath is supplied without metallic block allowing each user to make up the equipment according to its own necessities. For this purpose, there are available metallic blocks with different capacities.

Features of this product as follows:

1. Temperature and time digital displayed in LCD and digital controlled
2. Display of the set and current temperature contemporary
3. The time of temperature-control can be set up
4. Using of metal modules to avoid product pollution
5. Easy replacement, cleaning and disinfecting of metal modules
6. Built-in over-temperature protection, more reliable for using
7. Temperature adjustment
8. Buzzer alarm in the end of temperature control duration
9. Pause function

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
Thank you for choosing this equipment. We sincerely wish that you enjoy your dry bath Nahita 603/10. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the dry bath Nahita 603/10, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the dry bath avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the dry bath to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the dry bath.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the dry bath warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the equipment warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.