DESTILADOR DE AGUA GLASS WATER STILL

Nahita

REFERENCIA/ CODE 59021105

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Thanks for choosing a NAHITA stainless steel water still, its high mechanical quality will allow you to have this equipment available for many years. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

In order to use this equipment, it is advisable to follow these points:

- ◆ This manual should be available for all the users of these equipments.
- ◆ Any doubt can be solved getting into touch with the NAHITA Technical Service.
- ◆ As an accuracy instrument, it should be used carefully avoiding the sudden movements and knocks.
- ◆ Do not use the equipment with fluids that can give off steams or produce explosive or inflammable mixtures.

ATTENTION!!
IF AN EQUIPMENT IS NOT PROPERLY CLEANED AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

CONTENTS

1- USES OF THE INSTRUMENT.....	5
2- DESCRIPTION.....	6
3- TECHNICAL SPECIFICATIONS.....	6
4- OPERATING INSTRUCTIONS.....	6
ANNEX I: WARRANTY INFORMATION.....	7

1. USES OF THE INSTRUMENT

NAHITA stainless steel water still provides you with water of high purity, able for common processes in laboratories requiring the use of water as a disolvent. Thanks to distillation, we obtain water of the best quality.

We will obtain water suitable for reactivities preparation in qualitative and semiquantitative tests, bacteriological cultive environments and the lastest wash of glassware.

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
2. DESCRIPTION

1. Water exit
2. Water exit of the refrigerant
3. Refrigerant water entrance in the deposit or water exit to the drain
4. Draining tank
5. On/ off light
6. Tank's capacity display
7. Distillated water exit
8. Water entrance to the refrigerant (knob of different material from the other exits and entrances)

3. TECHNICAL SPECIFICATIONS

<i>Continuous production equipment</i>	5 litres per hour
<i>Feeding</i>	220V/50Hz
<i>Dimensions</i>	35 x 74 cm
<i>Weight</i>	4 Kg

4. OPERATING INSTRUCTIONS

- ◆ Place the water still in a stable surface, keeping it away from shocks and vibrations.
- ◆ Please, check the rubbers used in the setting are in good conditions.
- ◆ Put the water entrance rubbers to the refrigerant (8). The draining of refrigerant's water takes place through the stainless steel tube to the funnel which can be introduced again in the tank to be distilled or to the drain if the tank is full.
- ◆ Place an exit rubber for spare water (3) to the drain.
- ◆ Fill the deposit with water with the aid of the funnel or opening it and adding water into it.
- ◆ From (7) pick the distilled water.
- ◆ Once the settling is done, plug the equipment to electrical current. The red light will light as a probe the equipment is on (5) and the equipment will turn on.

WARRANTY

AUXILAB, S.L. warrants this equipment be free from defects in material and workmanship for a period of one year from date of purchase from AUXILAB S.L. or an authorised AUXILAB S.L. dealer. This warranty is limited to such purchaser and is not transferable.

During the year warranty period, AUXILAB S.L. will repair or replace, at its option, any defective equipments or pieces at no additional charge, provided that the product is returned, shipping prepaid, to AUXILAB S.L., an authorised dealer or an authorised service location. The purchaser is responsible for insuring any returned product and assumes the risk of loss during shipping. All replaced parts become the property of AUXILAB.

This limited warranty does not extend to any products which have been damaged as a result of accident, misuse, abuse or as a result of service or modification by anyone other than AUXILAB, an authorised dealer or an authorised service location.

Except as expressly set forth above, no other warranties are expressed or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. AUXILAB expressly disclaims all expressed or implied warranties not stated herein.

In order to obtain a free reparation or a free piece-replacement by our Technical Service over the above-mentioned period of time, it is absolutely necessary to submit the Warranty Certificate. This warranty doesn't apply for the maintenance of the equipment or for damaged caused by a natural, everyday use.

Some states do not allow the exclusion or limitation of incidental or consequential damages for consumer products, and some states do not allow limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This limited warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

Dealer:			
Equipment:		Model:	
Serial Number:			
Date:			
Brief failure description:			

AUXILAB, S.L.

Polígono Morea Norte, 8 - 31191 BERRAIN (Navarra)

Tel(948) 310 513 - Fax(948) 310 500 - 312 071

e-mail: asistencia@auxilab.com <http://www.auxilab.es>