DECLARACIÓN DE CONFORMIDAD CE  
CE DECLARATION OF CONFORMITY  
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi spectrophotometers 4101 and 4111RS. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual is inseparable from the Zuzi spectrophotometers 4101 and 4111RS, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the spectrophotometer avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the spectrophotometer to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the spectrophotometer.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by email to [info@auxilab.com](mailto:info@auxilab.com) or by telephone: +34 807 117 040 (1.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

#### INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this spectrophotometer by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

**6. TROUBLESHOOTING**

TROUBLE	CAUSE	SOLUTION
The equipment does not turn on	No current reaches the equipment	Check the current line
	The switch is off	Turn the switch on
Very unstable readings	Stabilization period too short	Wait at least 30 min for stabilization after turning the equipment on
	Voltage fluctuations	Keep the voltage within the range 85-220 V +/- 10%, without variations
	Vibrations on the working area	Put the equipment on a stable surface
Wrong readings	Damaged lamp	Replace the lamp
	Equipment out of calibration	Perform dark current
	Cell wrongly placed	Put the cell in the adapter with the optical faces in the light path. Insert completely the adapter with the cell in the sample compartment
	The optical faces of the cell are scratched or dirty	Change or clean the cell When cleaning the cell do not make too much pressure. Always hold the cells by the opaque faces
	Sample with bubbles or particles in suspension	Check the sample preparation and the analytical process
	Not enough volume of sample	Fill the cell with sample

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**INDEX OF LANGUAGES**

English ..... 12-21  
 Spanish ..... 2-11  
 French ..... 22-31

**INDEX OF CONTENTS**


1. USES OF THE INSTRUMENT ..... 13  
 2. DESCRIPTION ..... 13  
 3. TECHNICAL SPECIFICATIONS ..... 15  
 4. INSTALLATION / SETTING UP ..... 15  
 5. MAINTENANCE AND CLEANING ..... 19  
 6. TROUBLESHOOTING ..... 20  
 ANNEX I: CE CERTIFICATE ..... 32

**1. USES OF THE INSTRUMENT**

Zuzi spectrophotometers 4101 and 4111RS are very useful and easy-to-use instruments that allows performing qualitative and quantitative (only model 4111RS) analysis in the visible range of the spectrum. Thus, they are really useful in several fields such as petrochemical and biochemical industry, hospitals, environmental conservation as well as in research and scientific development laboratories.

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

- 2.1 On/off switch.  
2.2 Socket.

- 2.3 RS232 digital output (only model 4111RS).  
2.4 Analogical output.

- 3.1 Power on/off indicator (in use).  
3.2 Digital display.  
3.3 Indicators of the selected mode.  
3.4 Mode: selection of the measurement mode.  
3.5 0A: blank (0% Abs in mode A or 100% T in mode T) / Only model 4111RS: increase of the reference concentration (mode C) or factor (mode F) value.

- 3.6 0% T: auto zero (mode T) / Only mode 4111RS: decrease of the reference concentration (mode C) or factor (mode F) value.  
3.7 Print (only model 4111RS): print connection (mode A and T) / Confirmation of the reference concentration (mode C) or factor (mode F) value.

## 5. MAINTENANCE AND CLEANING

*To get the best results and a higher duration of this equipment it is essential to follow the processes of use.*

*NOTE: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The spectrophotometer is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the spectrophotometer by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the spectrophotometer (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.
- ◆ If you are not using the spectrophotometer for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).
- ◆ Always hold the cells by the opaque faces avoiding the direct contact of the optical faces with fingers. After using, the cells must be washed with distilled water to avoid contaminations.

### Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the spectrophotometer and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**Auto zero**

The equipment can be calibrated at 0% T to assure reliable and accurate measurements. For this:

- ◆ Put the black cell into the sample compartment.
- ◆ Once the displayed value stabilizes, press 0% T (3.6).

**NOTA:** *auto zero must be performed periodically and whenever the transmittance value deviates from 0.00 when putting the black cell into the sample compartment.*

**Data output**

- ◆ Model 4110RS through the digital output (3.1) can be connected, using the adequate cable, to a printer with serial port or to the serial port of a computer.
- ◆ Through the analogical output (3.2) the equipment can be connected to an analogical register.

**Security**

- ◆ The spectrophotometer must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the spectrophotometer on top of a horizontal, plane, stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the spectrophotometer near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the spectrophotometer for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the spectrophotometer by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the spectrophotometer (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

**3. TECHNICAL SPECIFICATIONS**

<b>Code</b>	54101010	54111010
<b>Model</b>	4101	4111RS
<b>Wavelength adjustment</b>	Manual	
<b>Wavelength range</b>	330-1000 nm	
<b>Bandwidth</b>	8 nm	
<b>Optical system</b>	Single beam, grating 1200 lines /mm	
<b>Wavelength accuracy</b>	+/- 2 nm	
<b>Wavelength resolution</b>	+/- 1 nm	
<b>Photometric range</b>	0-100% T 0-1.999 A	0-100% T 0-1.999 A 0-1999 C; 0-1999 F
<b>Photometric accuracy</b>	± 2% T	
<b>Stability</b>	±0.002 A/hr at 500 nm	
<b>Spray light</b>	<1% at 400 nm	<0.5% at 340
<b>Light source</b>	Halogen lamp 6 V, 10W	
<b>Detector</b>	Silicon photodiode	
<b>Power</b>	85 – 240 V, 50 Hz ±10%%	
<b>Dimensions (LxWxH)</b>	28x32x12 cm	
<b>Weight</b>	4.5 Kg	

**4. INSTALLATION / SETTING UP****1. Preliminary inspection**

Unwrap the spectrophotometer, take off the involving plastic and take off the protection in which it comes fitted. Take off all the protective items and, without connecting the spectrophotometer to the net, make sure that it does not present any damage because of the shipment. In case the spectrophotometer presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied. Please check that all the accessories are enclosed with the equipment:

- Cell adapter.
- 2 glass cells.
- Black cell for auto zero.
- Schuko cable.
- RS232 cable (only model 4111RS).
- Software (only model 4111RS).
- Cover.
- User's manual.
- Warranty certificate.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

## 2. Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.**

- ◆ Please put the spectrophotometer on top of a horizontal, plane and stable table making a free space of at least 30 cm per side. Do not put the spectrophotometer near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ During operation, dangerous material such as flammable liquids or pathological material, must be out of the working area.
- ◆ The spectrophotometer is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 85-240 V, 50 Hz  $\pm$  10% provided with earth wire and to the other end to the spectrophotometer connector (2.2).

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 85-240 V, 50 Hz  $\pm$  10%.

If you are not using the spectrophotometer for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

## 3. Operation mode

- ◆ Turn the spectrophotometer on by pressing the on/off switch (2.1) and wait for about 30 min until it stabilizes.
- ◆ By pressing Mode (3.4), select the desired measurement mode; the indicator (3.3) corresponding to the selected mode will light up:
  - Absorbance (A).
  - Transmittance (T).
  - Calculation of concentration (C), only model 4111RS.
  - Factor value adjustment (F), only model 4111RS.
- ◆ By turning the manual knob (1.3) select the wavelength value to perform measurements.

**NOTA:** *If the wavelength value is modified considerably during operation of the equipment, wait for a few minutes before performing new measurements to let the detector stabilize.*

## Absorbance/Transmittance measurements

Once the Absorbance or Transmittance mode is selected, proceed as follows:

- ◆ Make the blank. For this, put the cell with the blank sample in the sample compartment and press 0A (3.5); "BLA" will be displayed while the equipment is blanking. Finally "0.000" (Absorbance mode) or "100.0" (Transmittance mode) will be displayed indicating that blanking has finished successfully.

**NOTE:** *The equipment is supplied with a cell adapter. Put the cell with the sample in the adapter and the adapter into the sample compartment. Make sure that both the adapter and the cell are in the correct position so as the light beam passes through the sample.*

- ◆ Take the cell with the blank out and put the cell with the sample. Depending on the selected mode, the absorbance or transmittance reading will be automatically displayed.

## Concentration measurement (only model 4111RS)

**NOTE:** *Before performing concentration calculations it is necessary to make the blank as previously described.*

The equipment allows carrying out quantitative analysis to calculate the sample concentration by using 2 methods: taking a sample of known concentration as reference or inputting the value of the F factor ( $Abs = F \cdot Conc$ ).

- ◆ Using as reference a simple of known concentration:
  - By pressing Mode (3.4) select the Concentration mode (C).
  - Put the reference simple into the simple compartment and by pressing 0A/Inc (3.5) and 0% T/Dec (3.6) adjust the displayed value to the concentration value of the reference sample.
  - Press Print/Ent (3.7) to confirm the concentration value of the reference sample.
  - Take the reference simple out and put the sample to be measured; its concentration value will be automatically displayed.
- ◆ Inputting the F factor value ( $Abs = F \cdot Conc$ ):
  - By pressing Mode (3.4) select the mode Factor (F).
  - By pressing 0A/Inc (3.5) and 0% T/Dec (3.6) input the value of the F factor.
  - Press Print/Ent (3.7) to confirm this value.
  - Put the simple to be measured and its concentration value will be automatically displayed.