ESTEREOMICROSCOPIOS STEREOMICROSCOPES

Zuzi

Serie / Series 234

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

ANNEX I: CE CERTIFICATE

AUXILAB S.L.

CE DECLARATION OF CONFORMITY

URA TECHNIC STEREOMICROSCOPES SERIES 234 of AUXILAB, S.L. for the Directive of Machines (89/392/CEE modified) and the regulations adopted for their transposition

NAME OF THE MANUFACTURER / IMPORTER:
AUXILAB, S.L.

ADDRESS: Polígono Morea Norte, 8
31191 Beriáin (Navarra)

WE STATE THAT:

URA TECHNIC STEREOMICROSCOPES SERIES 234
CODES 50234000, 50234001, 50234200, 50234201

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ◆ Directive 73/23/CEE modified over the electric security.
- ◆ Directive 89/336/CEE modified over the electromagnetic compatibility.

And that the following harmonized rules have been applied (or part of them):

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BERIAIN 7th July 2004

Signed by: ALFONSO AINCIBURU SANZ
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

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- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The stereomicroscope is supplied with a Schuko standard wire and it should be connected to a current wire provided with an earth wire, it should be handy to be disconnected in case of emergency.
- ◆ Never try to repair the stereomicroscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the stereomicroscope (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Ura Technic Technical Assistance Department.
- ◆ If the are lamps blown you should replace them by another original Ura Technic, being careful not to touch them with the naked hands. Do not use any lamp with higher power, as this could provoke overheating or any other malfunction.
- ◆ Were it essential to replace the fuse, you should change it as follows: unscrew the protective lid (FUSE) that you will find at the rear part of the base, replace the fuse locating the new one at the same position and screw it again once you have finished.
- ◆ Do please always use the plastic cover whenever the stereomicroscope is not being used to avoid dust on the optical parts.
- ◆ Please keep the original packaging so as to use it in transport, whenever you have to return it for overhaul or just to keep it in a safe place when it is not going to be used for a long time.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the stereomicroscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ Lenses must never be dismantled by the user. In case there is any dirt on the external parts of the lenses, you should clean it with a smooth piece of cloth, fluff-free, dampened with a bit of xylol or toluene.
- ◆ Blow the dust laid down on the lenses with a plastic bulb or clean it either with a soft brush made of natural hair or any special gauze for lenses.
- ◆ For cleaning the mechanical parts you should use non-corrosive lubricants, being careful not to touch the optical parts.

ATTENTION!!
IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

6. CHOOSING EYEPIECES AND OBJECTIVES

- ◆ The observed image loses surface and sharpness as magnification increases.
- ◆ The mentioned increase should be done by changing objectives and putting each time a more powerful one and not by changing eyepieces to a higher magnification, as eyepieces only magnify the image obtained with the objective. Thus, the more magnification the eyepiece has, the higher will be the loss of sharpness, clarity and surface of the resulting image.
- ◆ For routine observations you should use the eyepiece with lesser magnification with more powerful objectives. The eyepiece with the higher magnification should be kept back for particular occasions, bearing in mind that it decreases definition and does not increase resolution.

7. SECURITY

- ◆ The stereomicroscope must be used by previously qualified staff that knows how it works thanks to the user's manual.
- ◆ You should put the stereomicroscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the stereomicroscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the stereomicroscope for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ *It is essential to have the equipment switched off and unplugged from the net* before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the stereomicroscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the stereomicroscope (see *Maintenance*).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

8. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

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	<i>Modelo 234</i>	<i>Modelo 234/1</i>	<i>Modelo 235</i>	<i>Modelo 235/1</i>
<i>Referencia</i>	50234000	50234001	50234200	50234201
<i>Cabezal</i>	Binocular	Binocular	Triocular	Triocular
<i>Zoom</i>	1x - 4x	1x - 4x	1x - 4x	1x - 4x

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- ◆ If you are not using the stereomicroscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Setting up

- ◆ Choose the most appropriate light source (incident and/or transmitted). To adjust incident illumination's intensity use the adjusting knob in the base. These models are provided with a ground glass stage of 95 mm diameter for transmitted lighting, and a black & white contrast stage for incident lighting.
 - ◆ You should observe with both eyes at the same time and graduate the interpupillary distance until the separate images you see become one stereoscopic vision.
 - ◆ Centre the sample into the stage and move the focusing control slowly until you manage to see a sharp image.
 - ◆ By alternately closing first one eye then the other you will appreciate a clear difference in focusing. To correct this defective vision you should move the dioptic compensation ring of each eyepiece.
- These equipments have a clear advantage among others: you can change magnifications gradually. Thus, it helps you to see the whole topographic structure of the sample so as you can concentrate on the more appealing parts to observe them at higher magnifications at the same time you are changing from a certain field's diameter to a higher one without changing objectives. Zoom range goes from 1x up to 4x.
- ◆ To change the stage you only have to loose the front screw, replace the current stage for the one desired and tighten the screw again.

Installing optional objectives:

- ◆ If you need more accurate magnifications you can attach optional objectives as: 0.5x, 1.5 and 2x.
- ◆ To do so, remove the zoom protective cover and screw the desired objective. Zoom usage is exactly the same. See chart on page 12 to find out the magnifications you can obtain.

5. TAKING PHOTOGRAPHS (Models 235, 235/1 optional)

- ◆ Place the photographic adapter on the head's tube and screw it tightly. The photographic adapter is composed by the photographic tube, the photographic ocular and the viewfinder on the focusing eyepiece. The viewfinder is provided with a T mount where you place the camera's corresponding T ring and then you screw the camera on it.
- ◆ Pull outwards the observation/photo rod located on the stereomicroscope's head in order to choose *photo*. Thus, light will be directed to the photographic tube, as the left eyepiece is void.
- ◆ Turn the dioptic compensation ring until you can see a sharp image. The area you are about to photograph is contained on the eyepiece's photographic focusing area.
- ◆ Now you can take the photographs you need by pressing the camera's shut button.

Working distance

Zoom	Objectives	Working distance
1 X	1.0 X	90 mm
↓	0.5 X	150 mm
	1.5 X	49 mm
4 X	2.0 X	34 mm

4. INSTALLATION / SETTING UPPreliminary inspection

- ◆ Unwrap the stereomicroscope, take off the involving plastic and take off the polispán protection in which it comes fitted. Take off all the protective items and, without connecting the stereomicroscope to the net, make sure that it does not present any damage because of the shipment. In case the equipment presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.
- ◆ Please check that all the accessories are enclosed with the equipment:
 - Black&white contrast stage
 - Pair of rubber eyecups
 - Lamps
 - Manual instructions
 - Fuse
 - Warranty certificate
 - Plastic cover

We will not accept any equipment in return period unless it comes in its original wrapping.

Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ◆ Please put the stereomicroscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the stereomicroscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The stereomicroscope is supplied with a Schuko standard wire
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth connection and to the other end to the stereomicroscope connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220V 50Hz $\pm 10\%$.

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- Base: metallic and very stable.
- Switch: off/ incident lighting/reflected lighting.
- Light: incident and transmitted.
 - ♦ - Incident: halogen lamp 6V 15W.
 - ♦ - Transmitted: halogen lamp 6V 10W.
- Stage: ground glass of 95 mm of diameter for transmitted lighting and black & white contrast for incident lighting.
- Focusing controls: bilateral by pinion and rack, with safety stop at the end of the route and tension adjustment.
- Zoom control: bilateral
- Binocular head: inclined 45°. Adjustable interpupillary distance, with dioptic compensation rings in both eyepieces.
- Pair of rubber eyecups.

Optical specifications

- Eyepieces: WF (large field) 10X: 20 mm diameter field
- Zoom: It enables a continuous magnification adjust without a loss of image's sharpness. Zoom control changes continuously the magnifier's factor from 1X to 4X. Working distance is exactly the same in all zoom range. It is operational at both sides, being always completely parfocal and centred.
- Options:
 - ♦ WF 15X eyepiece (Code 90234146)
 - ♦ WF 20X eyepiece (Code 90234148)
 - ♦ 2.5X photographic eyepiece (Code 90230191)
 - ♦ 5X photographic eyepiece (Code 90234193)
 - ♦ 0.5X objective (Code 90234200)
 - ♦ 1.5X objective (Code 90234201)
 - ♦ 2.0X objective (Code 90234202)

ZOOM	OBJECTIVES	WF10X		WF15X		WF20X	
		magnification	Field of view (mm)	magnification	Field of view (mm)	magnification	Field of view (mm)
1 X	1.0X	10x-40x	26 - 5.8	15x -60x	20 - 4.4	20x - 80x	16 - 3.5
↓	0.5X	5x-20x	52 - 11.6	7.5x -30x	40 - 8.82	10x - 40x	32 - 7
	1.5X	15x-60x	17.3 -3.19	22.5x -90x	13.3 - 2.94	30x - 120x	10.6 - 2.35
4 X	2.0X	20x-80x	13 - 2.9	30x -120x	10 - 2.2	40x - 160x	8 - 1.76

optional

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	12
2. DESCRIPTION	12
3. TECHNICAL SPECIFICATIONS.....	13
4. INSTALLATION / SETTING UP.....	14
5. TALKING PHOTOGRAPHS	15
6. CHOSE EYEPIECE AND OBJECTIVES.....	16
7. SECURITY.....	16
8. MAINTENANCE AND CLEANING.....	16
ANNEX I: CE CERTIFICATE.....	18


1. USES OF THE INSTRUMENT

These instruments are useful to watch enlarged images of small elements. The image is not inverted but stereoscopic, it is seen through the eyepiece in the same colour as the original sample. Ura Technic equipment will help you in most difficult tasks related to controlling, researching and preparing documents, in the best work conditions.

We stand out some field where stereomicroscopes can be used:

- Education: primary, secondary, vocational training and university.
- Industry: electronics, clock making, textile, jewellery, painting, food industry, tool setup, casting (metallurgy), graphics, galvanising, etc. In every process of inspection, assembly and repairing of instruments, measuring plates or any other kind of precision piece, as well as in quality control processes.
- Agriculture and stockbreeding, biology, prosthodontists, artistic restoration, philately, mycology, graphology, mineralogy, optics, research, etc.

These equipments provide a great comfort in their use and an excellent quality of the observed image.

2. DESCRIPTION

1. Eyepiece
2. Dioptic adjustment ring
3. Head
4. Zoom control
5. Focusing adjustment handle
6. Base
7. Holding clips
8. Stage
9. On/off switch
10. Incident light / transmitted light / incident and transmitted light switch
11. Incident lighting

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Thank you for choosing this equipment. We sincerely wish that you enjoy your stereomicroscope Ura Technic series 234. We highly recommend looking after this equipment according to what is stated in this manual.

Ura Technic develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of *Maintenance*, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the stereomicroscope URA Technic series 234, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the stereomicroscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the stereomicroscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the stereomicroscope.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Ura Technic Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0,30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the stereomicroscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the stereomicroscope warranty.
- ◆ Neither fuses (2.5A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Ura Technic Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!!

IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish.....	3-10
English.....	11-18