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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita high vacuum pump. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita high vacuum pump, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the oscillating stirrer.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

| | |
|---------------|-------|
| Spanish | 2-9 |
| English | 10-17 |
| French | 18-26 |

INDEX OF CONTENTS

| | |
|------------------------------------|----|
| 1. USES OF THE INSTRUMENT | 11 |
| 2. DESCRIPTION | 11 |
| 3. TECHNICAL SPECIFICATIONS | 12 |
| 4. INSTALLATION / SETTING UP | 12 |
| 5. MAINTENANCE AND CLEANING | 15 |
| 6. TROUBLESHOOTING | 15 |
| ANNEX I: CE CERTIFICATE | 27 |


1. USES OF THE INSTRUMENT

These high vacuum pumps are suitable for using in laboratories and chemical industry, packing, refrigeration, lyophilization equipments, lamps, vacuum impregnation, etc.

2. DESCRIPTION

- 1.1 Handle
- 1.2 Switch on/off
- 1.3 Oil inlet
- 1.4 Oil level sight glass
- 1.5 Oil outlet
- 1.6 Gas ballast valve
- 1.7 Base with anti-skid legs
- 1.8 Vacuum connection
- 1.9 Protecting caps

3. TECHNICAL SPECIFICATIONS

| | | |
|-----------------------|---|----------------------------------|
| Code | 59001101 | 55001102 |
| Type of pump | Two stage rotary vane | |
| Two stage rotary vane | 15 micrones / 0.015 mm Hg (Torr) / 2x10 ⁻⁶ MPa / 0.02 mbar | |
| Flow | 35 L/min – 2.1 m ³ /h | 65 L/min – 3.9 m ³ /h |
| Connections | 1/4"-SAE | 1/4"-SAE; 3/8"-SAE |
| Consumption | 160 W | 190 W |
| Dimensions | 3.8 kg | 5.7 kg |
| Weight | 271x166x195 mm | 301x152x204 mm |
| Power | 220/240 V, 50/60 Hz | |
| Oil capacity | 125 mL | 250 mL |
| Motor speed | 2400 prm | |

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the high vacuum pump, take off the involving plastic and take off the polystyrene protection in which it comes fitted. Take off all protections and without connecting the pump to the power supply, make sure that it does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all

5.) The pump presents is too hot when operating

| CAUSE | SOLUTION |
|--|---|
| Gases sucked at high temperatures | Verify the process. Study the possibility of using a cooling system |
| There is no enough oil or it is too contaminated | Correct the level of oil. If it is contaminated, drain it and add new oil |

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the stirrer and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

2) *The pump does not make a good vacuum or spends more time than expected.*

| CAUSE | SOLUTION |
|--|--|
| The oil is contaminated | Drain the contaminated oil and add new oil |
| There are leaks or losses in the vacuum system | Review the vacuum hose to detect possible leaks |
| The level of oil is not enough | Stop the pump and add oil until reaching the correct level if the oil is in perfect conditions of use. If not, change the oil. |
| Loss of oil due to a defect of the pump | Contact the Nahita Technical Assistance Department |
| Internal parts of the pump damaged or worn out by the use or a bad lubrication | Contact the Nahita Technical Assistance Department |

3) *The pump starts but presents a high electric consumption (amperage).*

| CAUSE | SOLUTION |
|---|--|
| The oil is too dense (high viscosity) or environmental temperature is lower than 5 °C | Use a lower viscosity oil |
| The pump is overloaded with oil | Remove the excess of oil |
| Foreign bodies inside the pump; broken vanes; damaged or blocked bearings | Contact the Nahita Technical Assistance Department |
| Low tension | Verify the tension used |

4) *The pump is very noisy when operating*

| CAUSE | SOLUTION |
|-------------------------------|--|
| Connection worn out or loosen | Contact the Nahita Technical Assistance Department |
| Bearings damaged | Contact the Nahita Technical Assistance Department |

the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Oil bottle, 250 mL
- User's manual
- Warranty certificate

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using the vacuum pump, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT

- ◆ Put the pump on top of a horizontal, plane and stable surface so as vibrations produced during operation do not disturb nearby equipments. The back side of the pump must be perfectly ventilated and free of obstructions.
- ◆ Do not put the pump near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. During operation, flammable liquids must be out of the working area.
- ◆ Open the oil inlet cap (1.5) and fill the pump until the oil level reaches the middle of the sight glass (1.4).
- ◆ Make sure that the switch (1.2) is off and insert the Schuko wire that feeds the AC electric current in the base of current 220 V, 50 Hz ± 10% provided with earth wire.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz ± 10%.

Setting up

- ◆ High vacuum pumps are designed to offer highly satisfactory results for a long time and don't need a special assembly.
- ◆ The aeration hose must be flexible to avoid a possible transmission of vibrations, water-tight and with the same or higher calibre than the pump connection (depending on the distance to cover) to obtain an optimum output.

Note: model 59001102 presents double vacuum connection. During operation, the connection that is not in use must be perfectly protected by its cap (1.8) to avoid leaks.

◆ It is advisable to use a thin filter to avoid filings, dust or any other object to go into and damage the pump.

◆ Lubrication: it is very important to control periodically the oil level. It must be checked with the pump on, and in vacuum. The oil level must be approximately at the middle of the sight glass (1.4).

ATTENTION: *Verify always that the oil level is correct since an excess of oil will produce a wrong operation of the pump.*

Oil change

◆ The oil change must be performed whenever the colour of the oil has changed and must be always done with the pump off and warmed. For this purpose proceed as follows:

◆ Open the oil outlet (1.5) and drain the oil.

◆ Close again the oil outlet (1.5) and switch the pump on for 1 minute approximately, to let the remaining oil inside the pump to be removed.

◆ Stop the pump, repeat the process to drain the remaining oil and close the oil outlet (5).

◆ Add new oil through the oil inlet (1.3) until the level reaches the middle of the glass sight (1.4).

◆ With the pump on and in vacuum, verify that the oil level is correct.

Note: *Always use a type of oil suitable for vacuum pumps.*

Gas ballast valve

◆ Nahita high vacuum pumps are provided with a gas ballast valve (1.6) to avoid the oil emulsion with agents such as moisture, gases, condensed vapours, etc.

◆ If you are not using the pump for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Security

◆ The pump must be used by previously qualified staff that knows how the equipment works thanks to the user manual.

◆ Put the pump on top of a horizontal, plane and stable surface so as vibrations produced during operation do not disturb nearby equipments. The back side of the pump must be perfectly ventilated and free of obstructions.

◆ Do not put the pump near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid vibrations, dust and dry environments.

◆ During operation, flammable liquids must be out of the working area.

◆ When you are not using the pump for a long period of time please make sure it is unplugged in order to avoid possible accidents.

◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece.

◆ Never try to repair the pump by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).

◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the pump (see Maintenance).

◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: *All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

◆ Please follow the processes of use of this manual.

◆ This manual should be available for all users of this equipment.

◆ The brushless motor is maintenance-free.

◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.

◆ Always use a type of oil suitable or vacuum pumps; the use of any other type of oil will provoke damages to the equipment.

◆ The pump is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.

◆ Never try to repair the pump by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the pump (burns, hurts...) or damages in nearby equipments.

◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

6. TROUBLESHOOTING

1) *The pump does not start*

| CAUSE | SOLUTION |
|--|------------------------|
| The electric motor: - Is not properly connected - Does not have the properly power supply - Low tension - The starting is not correct according to its power | Check the power supply |