**PIPETAS AUTOMÁTICAS.
VOLUMEN FIJO Y VARIABLE
AUTOMATIC PIPETTES
FIXED AND VARIABLE VOLUME
PIPETTES AUTOMATIQUES.
VOLUME FIXE ET VARIABLE**

DIGIPETTE

**Serie Profesional / Serie ECO
Professional Series / ECO Series
Série Professionnelle / Série ECO**

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Thanks for choosing a micropipette digipette NAHITA. It is advisable to read carefully this manual and follow the processes of use.

NAHITA develops its products putting more emphasis on Ergonomics and Security for users. Its high mechanical and optical quality will allow you to have this equipment available for many years.

BEFORE CONNECTING OR OPERATING WITH THIS INSTRUMENT, READ THESE INSTRUCTIONS THOROUGHLY.

In order to use this equipment, bear in mind the following points:

- ◆ Any doubts can be clarified by contacting the NAHITA Technical Service.
- ◆ As it is a precision instrument, it should be handled carefully avoiding sudden movements and knocks.
- ◆ In order to prevent fire or electrical unloadings, avoid atmospheres excessive humidity.
- ◆ To prevent fire or shock hazard, do not expose the balance to rain or moisture.
- ◆ Any doubts can be clarified by contacting the NAHITA Technical Service, by e-mail asistencia@auxilab.es.
- ◆ Warranty certificate only includes Professional pipettes of variable volume, with a warranty period of 2 years.

ATTENTION!!

IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT REPAIRED BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

| | |
|--------------|-------|
| SPANISH..... | 2-11 |
| ENGLISH..... | 12-21 |
| FRENCH..... | 22-31 |

CONTENTS

| | |
|-----------------------------------|----|
| 1. INTRODUCTION..... | 13 |
| 2. DESCRIPTION | 13 |
| 3. TECHNICAL SPECIFICATIONS | 14 |
| 4. OPERATING INSTRUCTIONS | 15 |
| 5. PIPETTING RECOMENDATIONS..... | 17 |
| 6. MAINTENANCE | 18 |
| 7. PIPETTE CALIBRATION | 19 |
| 8. ACCESSORIES..... | 21 |
| ANNEX I: CE CERTIFICATE | 32 |

1. INTRODUCTION

The micropipettes Digipette are suitable for medicine and health, chemical industry and the scientific research departments.

They have a large precision and accuracy across the whole volume range, from 0.2 μl to 5.000 μl .


Their ergonomic design minimises errors due to fatigue or irregular handling.

A special tip ejecting button for removing the suction tip is set in the pipette to make the tip displacement quickly and easily, except models 66632350 and 66612350.

They have two stop positions for the plunger: the first stops discharging the fluid, depressing the second stops blowing out the residual fluid.

2. DESCRIPTION

Picture 1

3. TECHNICAL SPECIFICATIONS

Autoclavable Professional Series

Variable volume

| <i>Model</i> | <i>Volume</i> | <i>Accuracy %</i> | <i>Precision %</i> |
|--------------|---------------|-------------------|--------------------|
| 66632020 | 0,1µl | ± 6,0 | ≤ 6 |
| | 0,5µl | ± 5,0 | ≤ 2.5 |
| | 2µl | ± 1,5 | ≤ 0.7 |
| 66632110 | 1µl | ± 2,5 | ≤ 1.25 |
| | 5µl | ± 1,5 | ≤ 0.6 |
| | 10µl | ± 1,0 | ≤ 0.4 |
| 66632120 | 2µl | ± 4.0 | ≤ 1.5 |
| | 5µl | ± 4.0 | < 1.5 |
| | 10µl | ± 3.0 | < 1.0 |
| | 20µl | ± 2.5 | < 0.8 |
| 66632210 | 10µl | ± 2.5 | < 0.8 |
| | 50µl | ± 2.0 | < 0.5 |
| | 100µl | ± 1.3 | < 0.3 |
| 66632220 | 20µl | ± 2.5 | < 1.0 |
| | 100µl | ± 1.3 | < 0.7 |
| | 200µl | ± 1.0 | < 0.5 |
| 66632310 | 100µl | ± 2.0 | ≤ 0.8 |
| | 200µl | ± 1.5 | ≤ 0.6 |
| | 500µl | ± 1.0 | ≤ 0.4 |
| | 1000µl | ± 0.8 | ≤ 0.28 |
| 66632350 | 1000µl | ± 1.3 | ≤ 0.4 |
| | 2000µl | ± 0.8 | ≤ 0.3 |
| | 5000µl | ± 0.6 | ≤ 0.2 |

Fixed volume

| <i>Model</i> | <i>Volume</i> | <i>Accuracy %</i> | <i>Precision %</i> |
|--------------|---------------|-------------------|--------------------|
| 66630110 | 10 µl | ± 1.00 | < 0.50 |
| 66630120 | 20 µl | ± 1.00 | < 0.30 |
| 66630150 | 50 µl | ± 0.80 | < 0.30 |
| 66630210 | 100 µl | ± 0.80 | < 0.25 |
| 66630220 | 200 µl | ± 0.80 | < 0.15 |
| 66630250 | 500 µl | ± 0.80 | < 0.20 |
| 66630310 | 1.000 µl | ± 0.80 | < 0.13 |

Tips

| <i>Model</i> | <i>Code</i> | <i>Model</i> | <i>Code</i> |
|--------------|-------------|--------------|-------------|
| 66632020 | 46622010 | 66630110 | 46622010 |
| 66632110 | 46622010 | 66630120 | 46622022 |
| 66632120 | 46622010 | 66630150 | 46622022 |
| 66632210 | 46622022 | 66630210 | 46622022 |
| 66632220 | 46622022 | 66630220 | 46622022 |
| 66632310 | 46622030 | 66630250 | 46622030 |
| 66632350 | 46622035 | 66630310 | 46622030 |

Accuracy and precision of pipettes are calculated by experienced people that use an appropriate technique.

ECO Series

Variable volume

| <i>Model</i> | <i>Volume</i> | <i>Accuracy %</i> | <i>Precision %</i> |
|--------------|---------------|-------------------|--------------------|
| 66612020 | 0,2µl | ±1.5 | <10 |
| | 2µl | ±2 | ≤ 1 |
| 66612110 | 1µl | ±4 | ≤ 2 |
| | 10µl | ±1,2 | ≤ 0.7 |
| 66612120 | 2µl | ±5 | ≤ 3 |
| | 20µl | ±2.5 | < 0.8 |
| 66612210 | 20µl | ±3.7 | < 1.2 |
| | 100µl | ±1.9 | < 0.4 |
| 66612220 | 50µl | ±2.0 | < 1.0 |
| | 200µl | ±1.0 | < 0.4 |
| 66612310 | 200µl | ±1.5 | ≤ 0.6 |
| | 1000µl | ±1.0 | ≤ 0.3 |
| 66612350 | 1000µl | ±1.5 | ≤ 0.6 |
| | 5000µl | ±0.8 | ≤ 0.3 |

Fixed volume

| <i>Model</i> | <i>Volume</i> | <i>Accuracy %</i> | <i>Precision %</i> |
|--------------|---------------|-------------------|--------------------|
| 66610110 | 10 µl | ± 2.5 | <1.00 |
| 66610120 | 20 µl | ± 2 | < 1.00 |
| 66610150 | 50 µl | ± 2 | < 1.00 |
| 66610210 | 100 µl | ± 1.2 | < 0.50 |
| 66610220 | 200 µl | ± 1.0 | < 0.40 |
| 66610250 | 500 µl | ± 1.0 | < 0.40 |
| 66610310 | 1.000 µl | ± 1.0 | < 0.30 |

Tips

| <i>Model</i> | <i>Code</i> | <i>Model</i> | <i>Code</i> |
|--------------|-------------|--------------|-------------|
| 66612020 | 46622010 | 66610110 | 46622010 |
| 66612110 | 46622010 | 66610120 | 46622022 |
| 66612120 | 46622010 | 66610150 | 46622022 |
| 66612210 | 46622022 | 66610210 | 46622022 |
| 66612220 | 46622022 | 66610220 | 46622022 |
| 66612310 | 46622030 | 66610250 | 46622030 |
| 66612350 | 46622035 | 66610310 | 46622030 |

Accuracy and precision of pipettes are calculated by experienced people that use an appropriate technique.

4. OPERATING INSTRUCTIONS

Preliminary examination

◆ Remove the pipette from the packing material and examine it carefully to make sure no damage has occurred during shipping. If there is any damage, notify to your dealer/forwarding agency.

◆ Keep all packing material until you are sure that the instrument works correctly. All defective items must be returned in the original packing together with any other supplied accessories.

◆ Verify you have received all the accessories that come with the balance:

- Tips
- Piece for calibration (Professional series)
- Instruction manual
- Warranty

Installation

Before starting this instrument, read carefully this manual.

BEFORE CONNECTING OR OPERATING WITH THIS INSTRUMENT, READ THESE INSTRUCTIONS THOROUGHLY.


Volume setting

◆ Both the *Autoclavable Professional* and the *ECO Series* are divided into fixed and variable volume micropipette. In the micropipettes of variable volume, the volume of the dosage is

◆ Extract the adjustment-knob (2.1) towards out. Turning the volumen the required. Put the adjustment-knob.

◆ Press again the adjustment-knob inward to fix the selected volume.


Picture 2

¡ATTENTION! The thread turns always in an easy and soft way, otherwise don't push it and verify if the selected volume is correct.

Principle of use

Picture 3

- ◆ Place the tip on the cone of the pipette.
- ◆ Press the tip on firmly to ensure an airtight seal. Make sure the sealing surfaces of the pipette and the tip are clean.
- ◆ With the pipette held vertically, press the pushbutton to the first stop (3 a).
- ◆ Immerse the tip in the liquid (3 b) under test to the depth of 2-4 mm.
- ◆ Release the pushbutton slowly, 2-3 s. for the pipettes to 200 ml, 3-5 s. for the pipette 200-5000. Pause, then remove the tip from the liquid (3 c).
- ◆ Holding the pipettes with tip in contact with the vessel wall at an angle 10-45° empty the tip by pressing pushbutton slowly to the first stop (3 d). The speed of dispensing should be the same as speed of aspirating.
- ◆ Wait approximately 1 s.
- ◆ Empty completely the tip by pressing the pushbutton to the second stop and subsequently remove the tip from the laboratory vessel sliding it along the side of the vessel (3 e).
- ◆ Eject the tip to waste container by pressing the tip eject button (3 f).
- ◆ The pipette is ready for next pipetting.

5. PIPETTING RECOMENDATIONS

Do not aspirate liquids without a tip, since if the inner part of the pipette gets wet then the next values will be wrong. Do not use corrosive liquids since they can damage the inside walls of the pipette and it can cause the loss of its warranty.

- ◆ When the pipette is new or you have not used it for a long time, it is advisable to carry out several tests (with water) so that the piston assembly, the teflon seal and the o-ring grease correctly.
- ◆ Do not immerse the tip more than 5 mm below the surface of the liquid.
- ◆ Do not release the push button during the aspiration of the liquid since the measurement would be wrong.
- ◆ Do not turn down the pipette when the tip is wet or filled with liquid.
- ◆ In case of using liquids which wet the tip walls (e.g. sera, proteins, organic solvents...), it is advisable to aspirate and dispense the liquids more slowly than with other kind of liquids as well as to dry the tip carefully before dispensing the liquid.


6. MAINTENANCE

- ◆ This micropipette is a high precision instrument, so the following instructions about its maintenance must be followed carefully.
- ◆ When the pipette is used to dose acids or strong bases, after a long period of use or when there are irregularities in the measurement (such as liquid drip...) it is advisable to check the instrument. In order to carry out a checking of the pipette, you have to be sure of your pipette model and then follow these instructions:

Autoclavable Professional Series

- ◆ Extract the tip ejector, pressing the tip eject button and pulling downwards.
- ◆ Unscrew the shaft coupling that connects the shaft and the handle, being careful not to lose any piece. Check that the piston assembly is clean, if it is rusty or corroded, it should be replaced.
- ◆ Avoid the liquid going into the shaft so that the measurements are correct. Check that there aren't any leakages since it would also give place to wrong measurements.
- ◆ Check that the o-ring and the teflon seal are not damaged (dripping of the pipette tip) and replace

Picture 4

them if it is necessary. Assemble the piston assembly together with the teflon seal and the o-ring (pic. 4). Add a very small quantity of silicone grease in such a way that connects the teflon seal with the top of the shaft (1.6). It is important to add the exact quantity of grease since it could have a remarkable influence in the measurements.

- ◆ Reassemble the pipette being careful that the shaft coupling is set correctly at the top of the shaft. If the shaft coupling and the shaft were not placed properly, it would affect both the volume and the reproducibility of the measurement.

ECO Series

- ◆ Pull upwards the push button, when you have taken it away, cover the hole of the push button with your thumb and pull downwards the tip ejector, you have to cover the hole because the spring and a metal screw could jump outside the inner part of the pipette and could be lost.
- ◆ Unscrew the shaft coupling that connects the shaft and the handle, being careful not to lose any

piece. Check that the piston assembly is clean, if it is rusty or corroded, it should be replaced.

- ◆ At the same time check the shaft, it is the piece that is in continuous contact with the samples, so it should be completely clean. Avoid that liquid goes into the shaft so that the measurements are correct. Check that there aren't any leakages since it would also give place to wrong measurements.

- ◆ Reassemble the pipette being careful that the shaft coupling is set correctly at the top of the shaft. If the shaft coupling and the shaft were not placed properly, it would affect both the volume and the reproducibility of the measurement.

7. PIPETTE CALIBRATION

- ◆ The estimate time for the calibration of the micropipettes will depend on the user's requirements. All the micropipettes are calibrated after being manufactured, but its use can modify the original measurement.

- ◆ Prepare an analytical precision balance, you will need higher or lower precision of the balance depending on the range or volume of the pipette you are going to make the calibration. Check that the environmental temperature is $\pm 20^{\circ}\text{C}$. Place a vessel to pour the distilled water and work out its tare, you will have to choose the suitable volume depending on the pipette you are going to make the calibration.

- ◆ For fixed volume micropipettes, place the tip of the pipette and fill it with distilled water, then weigh this volume on the balance. Repeat this operation five times and note down all the results.

- ◆ For variable volume micropipettes, it is advisable to check at least three points of the total range of the pipette, the three points which are usually checked are: the maximum, the minimum and another intermediate point. Repeat this operation five times and note down the results.

- ◆ With these results, work out the average measurement, the repeatability and check the good functioning of the micropipette. Compare the results you have obtained with the technical specifications of the charts (Chapter 3).

Setting the volume

◆ We will only adjust the micropipette, when we have checked that there are not any damaged pieces (see maintenance) and when the pipette accuracy does not exceed 10% (in such case the pipette should be sent to DIGIPETTE Technical Service).

Fixed and variable volume pipettes, professional series

Proceed as follows:

- ◆ Make sure that the adjustment knob (2.1) is in its fixed position (not turn) (variable volume pipettes)
- ◆ Insert the calibration piece or an instrument with thin tip in the holes of the knob wheel.
- ◆ By turning right or left the pipetted volume will increase or decrease respectively.
- ◆ Check result by weighing again the aspirated volume (make tests with distilled water) and repeat these steps until obtaining the desired value.
- ◆ Make sure that the volume scale and the selected volume (in variable volume pipettes) are correct.

Fixed and variable volume pipettes, eco series

Proceed as follows:

- ◆ Hold the adjustment knob manually avoiding it to turn (variable volume pipettes).
- ◆ Insert the calibration piece or an instrument with thin tip in the holes of the knob wheel.
- ◆ By turning right or left the pipetted volume will increase or decrease respectively.
- ◆ Check result by weighing again the aspirated volume (make tests with distilled water) and repeat these steps until obtaining the desired value.
- ◆ Make sure that the volume scale and the selected volume (in variable volume pipettes) are correct.

8. ACCESSORIES

Pipette tips, made of polypropylene and autoclavable.

| Code | Capacity | Color | Units/ bag |
|----------|-------------------|--------|------------|
| 46622010 | 0.1-10 μ l | White | 1000 |
| 46622022 | 2-200 μ l | Yellow | 1000 |
| 46622030 | 100-1000 μ l | Blue | 5 |
| 46622035 | 1000-5000 μ l | White | 250 |

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INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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