POLARÍMETRO
POLARIMETER
POLARIMÈTRE

Modelo / Model / Modèle 418

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Thank you for choosing this equipment. We sincerely wish that you enjoy your equipo-marca-modelo. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user. The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of *Maintenance*, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the polarimeter, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the polarimeter avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the polarimeter to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the polarimeter.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0,30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the polarimeter warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the polarimeter warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!!
IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

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◆ When the D.C power supplying system is out of order for the sodium lamp, the instrument can also be used to perform the measurement with the sodium lamp under the condition of being power-supplied with A.C. But, the performance of the instrument will be slightly decreased.

◆ When the sample of smaller rotation (less than 0.5°) is put in, the indicating value will change probable. In such a case, only by pressing the *Repetition - Measurement* button, new value will be indicated.

Measurement of the concentration and content

◆ First dilute the standard or reference sample of a known purity into several test samples of different concentrations in accordance with a certain ratio, and respectively measure their optical rotations. Then, plot the optical rotational curve, using lateral axis to represent the concentrations and longitudinal axis to represent the optical rotation.

◆ Usually, the optical rotational curve should be tabulated into a look up table in accordance with the arithmetic inter-polation.

5. MAINTENANCE

◆ The instrument should be placed at a dry place with creculation of air, so as to prevent the instrument from being sostened and mildwed.

◆ If the light source is accumulated with dust or damaged, you may open the cover to wipe the light source clean or replace it with a new one.

◆ Alter turning on the power source, if the sodium lamp does not light up, you may check the 3A fuse.

INDEX OF LANGUAGES

Spanish 1-6
 English..... 7-11
 French 12-16

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT 8
 2. DESCRIPTION 8
 3. TECHNICAL SPECIFICATIONS10
 4. SETTING UP.....10
 5. MAINTAINANCE11

1. USES OF THE INSTRUMENT


Polarimeter Zuzi 418 is used in chemical industry, hospitals, schools, universities and research laboratories as taking optical rotation measurement of organic substances with the aid of the polarimeter, you can get to know the specific weight, the hardness grade, the concentration and content of thousands of optically active materials, such as sugar, turpentine, industrial oil, camphor, etc.

This polarimeter is widely used in:

- ◆ Food industry, to know sugar and condiments contents.
- ◆ Sugar refining plants, to control sugar concentration in industrial process.
- ◆ Pharmaceutical industries, to measure medicine's optical rotation and perfume's oils.
- ◆ High schools and universities can also make a good use of this polarimeter in teaching and demonstration tasks.

2. DESCRIPTION

Picture 1 - Circuit


- 1.1. Light source
- 1.2. Opening
- 1.3. Lens
- 1.4. Coloured filter
- 1.5. Polarizer
- 1.6. Magnetic rotation coil
- 1.7. Sample
- 1.8. Polarizer

- 1.9. Photomultiplier
- 1.10. Preamplifier
- 1.11. Automatic high voltage
- 1.12. Frequency amplifier
- 1.13. Amplifier's switch
- 1.14. Servomotor
- 1.15. Gear system
- 1.16. Counter

Polarimeters 418 use a lamp of 20W as light source. Objective's lens and the opening, compose a bright beam parallel to the tube (Picture 1)

After passing through the polarizer (1.5) the parallel light starts polarizing. The vibration plane is 00 (Picture 2.a). When the polarized light pass through the magnetic rotation Faraday coil, the plane vibrates creating a certain angle B to 50Hz (Picture 2.b), after passing through the polarizer (1.8) the light beam is projected on a photomultiplier where it generates a alternating current signal.

Picture 2 - Light path

Note: 00 it is the polarizer's polarization axis 1.5 and PP is the polarizer's polarization axis 1.8.

Picture 2a: The polarized light generated by the polarizer 1.5 vibrates in the plane 00.

Picture 2b: After having passed through the rotation magnetic coil, the vibration plane of the polarizer's light, generates an opposite angle.

Picture 2c: After having passed through the sample, the polarized light's plane rotates by α_1 .

Picture 2d: Once the value indicated by the instrument has been balanced, the polarizer 1.5 rotates reversibly to position α_1 until the sample's optical rotation is compensated.

3. TECHNICAL SPECIFICATIONS

Scale's measurement range	$\pm 45^\circ$
Accuracy	$\pm(0.001^\circ + \text{measured value} \times 5/10000)$
Automatic digital screen for reading	
Minimum readable value	001 $^\circ$
Speed	1. $^\circ$ /Second
Monochromatic lamp's wavelength (sodium steam)	589,3nm
It includes tubes of 100 and 200 mm	
Electrical source	220v 50Hz
Weight	14 Kg
Outside dimensions	606x310x212 mm

4. SETTING UP

- ◆ Plug the equipment to the current intake 220V 50 Hz $\pm 10\%$, with the Schuko wire provided. The polarimeter 418 should be connected to a current wire provided with an earth wire. It should be handy to be disconnected in case of emergency.
- ◆ Once it is connected let the sodium lamp get warm for 5 minutes by pressing the switch on button.
- ◆ Then, press the *Light* button that feeds the lamp with direct current (D.C) and by pressing the *Measure* button the display turns on.
- ◆ When the indicating value switch is turned on, the tubes should indicate, the digits.
- ◆ Put the test tube containing distilled water or other bland solvent into the simple chamber, and close the cover. After the indicating value has got stabilized, press *Clear* button. If there are air bubbles in the test tube, the first thing should be done is to make the air bubbles float on the protruded tube-neck. The atomized waterdrops on the both ends of the light-through plane should be wiped dry. The screw nut of the test tube should not be screwed down too tightly, so as to avoid stress, otherwise, the readings may be affected. When placing the test tube, care should be taken for the marked position and direction.
- ◆ Press the knob *Repetition - Measurement* three times running. If the different is only 0.005, the parameter of polarimeter is correct.
- ◆ Empty the simple tube. Fill it with the simple to be measured. Put the tube into the sample chamber. The machine operates. When a red dot below the single digit lights up, take the reading.
- ◆ Press the *Repetition - Measurement* button for several times, and take the average value as the measuring result.
- ◆ If the polarized angle of the sample exceeds the measuring range, the instrument will be oscillating at $\pm 45^\circ$ back and forth: at this time, the test tube should be taken out, then the instrument will reset to zero automatically.