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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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### 3. Procedure

◆ In the case of young or sparkling wines, remove firstly the carbonic gas. Stir 250 mL of wine in a flask of 500 mL. The interior of the flask must be previously siliconated with 3 drops of a solution of water soluble silicone at 1%; when the flask is dry, add the wine.

#### 3. Distillation:

- ◆ Add 250 mL of wine into a volumetric flask and check the temperature.
- ◆ Pour the wine in a spherical distilling flask with some small pieces of porcelain and add 10 mL of calcium hydroxide to turn the wine alkaline. Add some drops of silicone to avoid froth formation.
- ◆ Collect the distilled product into the 250 mL volumetric flask; the flask must contain 10 mL of pure water in which the bent distillation adapter is submerged.
- ◆ Distil at least 200 mL of wine and add water up to 250 mL. Pour the volume into a 250 mL measuring cylinder. Wait 1 min and then check the temperature.
- ◆ Remove the thermometer, submerge the alcoholometer and leave to stand at least 1 min. Take note of at least 3 readings of the alcoholic grade. The scale of the alcoholometer must be read just under the meniscus.

**Note: Use a magnifying glass to read better the alcoholometer scale.**

### 4. Calculations and results

◆ The volumetric alcoholic grade is expressed in % of volume at 20 °C with 2 decimal numbers; the second decimal number must be approximated to 0 or 5. As the volumetric alcoholic grade depends on temperature, determinations must be always performed at 20 °C. If the temperature is not 20 °C, readings must be corrected by using the following chart:

- pH indicator paper

#### 4. Precautions

- ◆ Be extremely careful with the vapours; always wear safety glasses and avoid burns.
- ◆ Ethanol is very flammable. Keep always this product away from flames.
- ◆ The sulphuric acid is very corrosive.

#### 5. Procedure

- ◆ Put a round bottom flask in cold water and add 30 mL of ethanol, 50 g of ammonium acetate and 50 mL of sulphuric acid. Keep the flask in cold water during operation.
- ◆ Once reaction is finished, assemble the distillation equipment; connect the condenser with a tap to make the water flow through the jacket and warm gently the spherical flask, so as distillation is performed slowly.
- ◆ Collect the distillation produced until temperature reaches 80-85 °C and stop distillation.
- ◆ Look up in a chart of organic compounds boiling points and deduce which are the components of the initial mixture and which is the component that remains in the distilling flask.
- ◆ Prepare two saturated solutions of Na<sub>2</sub>CO<sub>3</sub> and CaCl<sub>2</sub>.
- ◆ Pour the distilled product into a separating funnel and wash it with the carbonate solution; mix gently and let the two phases (aqueous and organic) separate. Keep washing the solution until the pH of the upper phase is not acid (check the pH with indicator paper).
- ◆ To remove the rests of carbonate, wash the solution with calcium chloride and remove the aqueous phase.
- ◆ Pour the liquid that remains in the funnel into a porcelain dish with magnesium sulphate or any other dehydrating agent and let act for half an hour. Cover the dish with a watch glass.
- ◆ Wash, let dry and weight an Erlenmeyer flask with a cork stopper and filter the dish content.
- ◆ Weight the liquid obtained that will be mostly ethyl acetate.

#### 6. Calculations and results

- ◆ Calculate the performance of the process.

### Practice 4: Determination of volumetric alcoholic grade

#### 1. Principle

- ◆ The volumetric alcoholic grade is the number of litres of ethyl alcohol (ethanol) and equivalents (methanol, 2,3-butanediol) contained in 100 litres of wine; both volumes must be measured at 20 °C and expressed in volumetric alcoholic grades (%vol) with a precision of 0.1%.

#### 2. Material and reagents

- |                              |                              |
|------------------------------|------------------------------|
| - Distillation equipment     | - Hydrometer (alcoholometer) |
| - Heating mantle             | - Calcium hydroxide          |
| - Round bottom flask, 1 L    | - Phenolphthalein            |
| - Measuring flask, 250 mL    | - Silicone                   |
| - Measuring cylinder, 250 mL |                              |

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- ◆ Pour 35 mL of acetone and 35 mL of water into the distilling spherical flask and add some small pieces of porcelain.
  - ◆ Open the tap to make the water flow through the condenser jacket and warm uniformly the distilling flask.
  - ◆ The boiling point of the acetone is lower than the water's so the acetone will reach its boiling point first. The acetone vapour will move upward the fractionation column; a heat interchange will take place between the ascendant vapour and the descendant condensed liquid.
  - ◆ The acetone is collected in an Erlenmeyer flask.
  - ◆ When the boiling point of the water is reached, it will begin to distil and again, a heat interchange will take place.
  - ◆ It is very important to change the Erlenmeyer flask when the temperature is around 85 °C. Thus assuring that acetone distillation is finished and water distillation does not yet begin.
  - ◆ Measure the obtained acetone and water volume.
5. Calculations and results
- ◆ Make a chart and a graphic with the data of temperature and time.
  - ◆ Write the boiling points obtained for the acetone and water.


### Practice 3: Preparation of ethyl acetate

#### 1. Objectives

- ◆ To obtain an organic compound
- ◆ To combine different laboratory basic techniques

#### 2. Principle

◆ The action of a carboxylic acid over an alcohol produces an ester and water as shown by the following formula:

◆ This reaction is slow and the presence of an equilibrium makes impossible a good performance. The sulphuric acid is used as a catalyser and as a dehydrating agent, thus the equilibrium is displaced towards the right producing a large quantity of ester (ethyl acetate). Esters have a nice smell and are usually volatile.

◆ One of the techniques used in this practice is the distillation, which is used to separate mixtures of liquids. Separation is based on the different boiling points of the components of the mixture. When the mixture is warmed the liquids evaporate and condense later. The more volatile component is the first obtained.

#### 3. Material and reagents

- |  |  |
|--|--|
| - Distillation equipment                 | - Ammonium acetate   |
| - Beaker or Erlenmeyer flask             | - Ethanol 96°  |
| - Thermometer                            | - H <sub>2</sub> SO <sub>4</sub> - Concentrated sulphuric acid |
| - Separating funnel                      | - Na <sub>2</sub> CO <sub>3</sub> - Sodium carbonate           |
| - Watch glass                            | - MgSO <sub>4</sub> - Magnesium sulphate (anhydrous)           |
| - Two stands with retort ring and clamps | - CaCl <sub>2</sub> Calcium chloride                           |

**Fractional distillation****Practice 2: Separation of a mixture of acetone and water****1. Objectives**

- ◆ To understand the bases and applications of the fractional distillation
- ◆ To compare the efficiency of the fractional distillation against that of simple distillation

**2. Principle**

◆ It is a fractional distillation technique that consists of separating two liquids whose boiling points differ in less than 80 °C. This technique allows a better and more efficient separation of the liquids than successive distillations, thanks to a fractionation column.

◆ Warming the flask produces the evaporation of the solution. This vapour, during its travel along the column is condensed. Due to the ascendant vapours a heat interchange is produced between them and the descendant liquid and the more volatile component reaches the condenser and finally the collecting flask. The less volatile component, however, remains inside the distilling flask.


**3. Material and reagents**

- |  |                             |
|--|-----------------------------|
| - Round bottom flask   | - Thermometer               |
| - Vigreux column   | - Bent distillation adapter |
| - 75° bent connecting adapter                                      | - 2 Small Erlenmeyer flasks |
| - Condenser  | - Heating mantle            |
| - Measuring cylinder   | - Acetone                   |
| - Water  |                             |
| - Clamps, stand, rubber tubing, porcelain, grease, insulating tape |                             |

**4. Procedure**

Assembly the fractional distillation equipment as shown in **figure 4**:

figure 4

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#### 5. PRACTICES

##### Simple distillation

##### *Practice 1: Purification of methyl alcohol*

##### 1. Objectives

- ♦ To understand the bases and applications of the simple distillation
- ♦ To learn to assembly a simple distillation equipment

##### 2. Principle

♦ It is an easy technique that consists of separating the impurities of a liquid by heating it. When the liquid reaches its boiling point, distillation begins; the generated vapour is condensed and collected while impurities remain inside the distilling flask.

##### 3. Material and reagents

- Round bottom spherical flask
- 75° bent connecting adapter
- Condenser
- Thermometer
- Bent distillation adapter
- Clamps, stand, rubber tubing, porcelain, grease, insulating tape
- Erlenmeyer flask, 250 mL
- Heating mantle or Bunsen burner
- Measuring cylinder
- Methyl alcohol

##### 4. Procedure

♦ Assembly the simple distillation equipment according to figure 3. All the ground joints must be well greased.

♦ Pour 50 mL of alcohol inside the spherical or distilling flask.

♦ Warm the flask uniformly until the liquid reaches its boiling point.

♦ The generated vapour passes through a condenser where it cools and condenses; the liquid is collected into an Erlenmeyer flask.

♦ The non volatile particles or dissolved impurities do not boil and remain inside the spherical distilling flask; thus, the distilled product collected in the Erlenmeyer flask is pure alcohol.


♦ The thermometer situated at the bent connecting adapter indicates the temperature at which the boiling point of the alcohol is reached. This temperature must stay constant during all the distillation.

♦ At the beginning the temperature must be checked every 2 min and then every minute. When the liquid begins to boil, the temperature must be checked every 3 min.

##### 5. Calculation and results

- ♦ Write a chart with the data of time (min) and temperature (°C)
- ♦ Indicate the temperature at which the alcohol is distilled.

figure 3

nal region through which the vapour from the distilling flask passes and condenses, and the external region through which the cooling water flows.

◆ **(D) Bent distillation adapter:** It connects the condenser with the Erlenmeyer flask or the recipient to collect the distilled product.

◆ **(E) Vigreux column**

◆ **(F) Clips:** All the components are joined and fixed together by means of polyoxymethylene (POM) clips included with the equipment (*see fig. 2*).

#### 4. RECOMENDATIONS

Before starting with the distillation process make sure that:

◆ The capacity of the distilling flask is correct and enough to contain the volume to be distilled. The Nahita distillation equipments include a distilling flask of 250 mL, please consult our catalogue to find flasks of other capacities.

*Note: The flask must not be filled over 2/3 its capacity, neither must be excessively big.*

◆ All the ground joints are perfectly greased.

◆ Some small pieces of porcelain have been added to the distilling flask. Almost all liquids tend to overheat; this metastable state is interrupted periodically because vapours are emitted suddenly. To avoid this, 2 or 3 pieces of porcelain must be added to the distilling flask before starting heating. The pores of this substance provide places for groups of bubbles to form thus inducing a uniform boiling. If the temperature decreases under the boiling point, the liquid fills the pores of the substance losing all its effect and efficiency. Then, the liquid must be let to cool under a temperature below the boiling point before adding a new piece of porcelain; if it is added to an overheated liquid it can produce a violent boiling.

◆ All the components are joined together by means of the clips included with the equipment. Make sure that there are no tension forces that will cause breakage of the glass components.

◆ The condenser is connected to the water flow; the rubber tubing must be well connected to the water inlet/outlet.

◆ It is convenient to use a heating mantle for a uniform warming of the flask. This is especially necessary in case of flammable liquids.

◆ The heat source is adjusted so as the distillation speed is approximately 1 drop per second. If speed is higher there is no time to establish the liquid/vapour equilibrium and distillation is no good. If speed is too slow, the bulb of the thermometer will always be covered by drop of condensed liquid and temperature readings will not be correct.

◆ Do not heat the flask if it is empty; it could break.

◆ The water flows gently through the condenser. Using a rubber or latex tubing, connect the water inlet to a tap and the water outlet to a sink.

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

## 3. ASSEMBLY

**VERY IMPORTANT:** Before assembling the different components, all the ground joints must be greased.

Using the stands and clamps, assembly the different parts as indicated in the next figure:

- ◆ **(A) Distilling spherical flask:** With a ground joint at the neck; it contains the solution to be distilled. In case of the simple distillation equipment, the flask is fixed to the metallic stand (*G*) by means of a clamp (*H*). In the fractional distillation equipment, the clamp holds the Vigreux column (*E*). The distilling flask is put inside the heating mantle (*I*).
- ◆ **(B) Bent connecting adapter, 75°:** It is also known as distillation head and presents two ground joints, one of them is connected to the flask (*A*) and the other to the condenser (*C*). In the fractional distillation equipment, the bent connecting adapter is joined to the Vigreux column (*E*) and the condenser (*C*). A thermometer is introduced through the upper part of the bent connecting adapter; the bulb of the thermometer must stay under the lateral arm of the adapter so as it is in complete contact with the vapour flow and then, the temperature readings will be correct.
- ◆ **(C) Condenser:** It is joined to the bent connecting adapter (*B*) and fixed to a metallic stand (*G*) by means of a clamp (*H*). Along the condenser there are two regions: the inter-

The following accessories are not included but are necessary for a good assembly:

Letter	Description
G	Retort stand
H	Clamp for condensers
I	Heating mantle, 250 mL

Useful labware to perform a distillation:

- Erlenmeyer flask
- Bunsen burner
- Tripod stand
- Gauze iron wire
- Measuring cylinder
- Grease for ground joints

**Note:** All the accessories mentioned above are indicated for a standard configuration; to use a spherical flask with higher capacity or any other size for stands, heating mantles etc. please visit our web site ([www.auxilab.com](http://www.auxilab.com)).

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
the most volatile component.

Distillation techniques can be divided into different types:

#### Simple distillation

This is the simplest type. The assembly is composed of a spherical flask (distilling flask), a bent connecting adapter, a condenser, a bent distilling adapter and other material necessary for the assembly of the distillation equipment and later determinations.

The simple distillation technique is used to remove non volatile particles and separate a mixture of liquids with a high difference in their boiling points. As a general rule any mixture of liquids with a difference in their boiling points of at least 80 °C can be separated by simple distillation. The component with the lowest boiling point, it means the most volatile, boils faster and the first fraction of the distilled product is very rich in this component. The component with the highest boiling point, it means the less volatile, will be the main component of the last distilled fraction.

#### Fractional distillation

In this case, a fractionation column is connected between the distilling flask and the bent connecting adapter. In this type of distillation a series of minuscule separations are performed in a unique and continuous operation. The fractionation column provides a large surface for the ascendant vapours and the descendant condensed product to interchange heat, thus making possible that numerous partial evaporations and condensations take place all along the column. The component with the lowest boiling point keeps reaching its boiling temperature until the mixture is completely separated. Then, the temperature is increased up to the boiling point of the less volatile component so as it can reach the condenser. This change of temperature makes the distilled product richer in the most volatile component, having the same effect as multiple distillations.

## 2. DESCRIPTION

The Nahita distillation equipments are made of borosilicate glass with 29/32 ground joints:

- Simple distillation equipment >> Cod. 15000220.
- Fractional distillation equipment >> Cod. 15001220.

Letter	Code	Description
A	18015220	Spherical distilling flask, 250 mL
B	18117001	Bent connecting adapter, 75°
C	11215251	Liebig condenser, 250 mm
D	18127001	Bent distillation adapter, 100 mm
E	18230250	Vigreux column (only fractional distillation equipment)
F	63226229	POM clips for ground joints
J	72064100	Thermometer -10/110 °C

Thank you for choosing a Nahita distillation equipment made of borosilicate glass and with ground joints of 29/32. Its high quality will permit you to enjoy this equipment for years. We highly recommend to read thoroughly this manual to get the best results and a higher duration of the equipment.

## INDEX OF LANGUAGES

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## INDEX OF CONTENTS

1. THEORETICAL FOUNDATIONS .....	15
2. DESCRIPTION.....	16
3. ASSEMBLY.....	17
4. ADVICES .....	18
5. PRACTICES .....	19

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

## 1. THEORETICAL FOUNDATIONS

Distillation is a technique used to purify liquids by removing the non volatile impurities or by separating two liquids with different boiling points. In this technique the liquid is passed through different steps:

- Evaporation by heating
- Condensation of the vapour
- Collection of the condensed vapour in a separate recipient

The *boiling point of a liquid* is the temperature at which the vapour pressure of the liquid matches the environmental pressure. In the case of a mixture of liquids, the boiling point is reached when the sum of the partial vapour pressures of each liquid matches the environmental pressure.

In a process of distillation in the laboratory, the liquid is warmed up in a spherical flask (distilling flask). When the boiling point is reached, numerous bubbles of vapour are formed and move upward the flask. A part of the vapour is condensed while ascending and falls again into the distilling flask but the most of the vapour goes into the condenser and there condenses thanks to the cold water that flows through the external jacket of the condenser. The condensed liquid is transported to an Erlenmeyer flask or any other recipient through a bent distillation adapter. The non volatile particles remain as wastes in the distilling flask; if the initial solution is a mixture of liquids, then the vapour will be richer in

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