BALANZA ELECTRÓNICA ELECTRONICAL BALANCE BALANCE ÉLECTRONIQUE

Series/Series/Séries 5172, 5171

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita-Blue electronic balances, series 5171/5172. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- This manual is inseparable from the Nahita-Blue electronic balances 5171/5172, so it should be available for all the users of this equipment.
- You should carefully handle the balance avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the balance.
- To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler.
- This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- Overhaul is not covered by the equipment warranty.
- Operations made by non-qualified staff will automatically produce a loss of the warranty.
- Accessories (including their loss), are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- Manufacturer reserves the right to modify or improve the manual or equipment.

INDEX OF LANGUAGES

Spanish.....	2-11
English.....	12-22
French.....	23-33

INDEX OF CONTENTS

1. SPECIFICATIONS.....	14
2. DESCRIPTIO.N.....	15
3. INSTALLATION.....	16
4. USE OF THE BALANCE.....	19
5. MAINTENANCE AND CLEANING.....	20

SPECIFICATIONS

Precision balances series 5171 and 5172 work on high precision strain gauge load cell which implements high speed stabilization and high reliability.

Technical parameters

Code	55172230	55172250	55172310	55172330	55172350
Max. capacity. (g)	0-300	0-500	0-1000	0-3000	0-5000
Tare range (g)	300	500	1000	3000	5000
Readability	0.01				
Pan size (mm)	135 mm (round)			155x143 mm (square)	
Dimensions (LxWxH)	265x200x80mm				
Power supply	DC 9V				

Code	55171310	55171320	55171330	55171350
Max. capacity (g)	0-1000	0-2000	0-3000	0-5000
Tare range (g)	1000	2000	3000	5000
Readability	0.1			
Pan size (mm)	135 mm (round)		155*143mm (square)	
Dimensions (LxWxH)	265x200x80mm			
Power supply	DC 9V			

Working temperature: 5-35 °C

Humidity: 50-85%

Features


- High precision loadcell
- Full plastic housing
- Stainless steel pan
- Superbright easy to read LCD display with backlit
- Overload indicator on display
- Mains adapter supplied as standard
- Possibility of working with batteries or connected to the power supply
- Height adjustable feet
- Full capacity taring
- Simple 5 key menu for easy operation
- Stabilization time 1-1.5 seconds
- Automatic calibration
- Overload protection
- Below balance weighing facility
- Selectable measure units: g, oz, ct
- Piece counting function

2. DESCRIPTION

Structure

Control panel

1. **ON/OFF:** Turn the balance on/off
2. **UNIT:** Select the weighing unit
3. **PCS:** Piece counting function
4. **CAL:** Calibration function
5. **TARE:** Tare function

3. INSTALLATION

Choosing the installation place

The measuring performance of the balance is greatly influenced by the environment where it is installed.

Observe the following points to ensure a safe and accurate weighing.

Prohibitions

Do not use the balance anywhere exposed to explosive, combustible or corrosive gases.

This could cause fire or other kind of serious problems.

Instructions

Use the correct power supply and voltage with the balance.

Using an incorrect power supply or voltage with the balance will lead to fire or trouble with the balance.

Note also that if the power supply or voltage is unstable or if the power supply capacity is insufficient, it will not be possible to obtain satisfactory performance from the balance.

Prohibitions

Avoid locations where the balance will be exposed to any of the following

situations:

- ✘ Air flow from an air conditioner, ventilator, door or window.
- ✘ Extreme temperature changes
- ✘ Vibration from surroundings or nearby equipments
- ✘ Direct sunlight
- ✘ Dust, electromagnetic waves or a magnetic field

Instructions

Install the balance on a strong and stable flat table or floor.

Placing the balance in an unstable site could lead to troubles with the balance or injury to users.

When selecting the installation site, take into account the combined weight of the balance and the item to be weighed.

Unpacking and delivery inspection


Check that all of the items indicated below are included in the package, and that nothing has been damaged:

- Balance
- Pan
- Instruction manual
- Adapter

In case the balance or any component presents any damage due to transport tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Turning the power on

1. Connect the balance to the power supply or install 4 batteries (type AA 1.5 V) in its corresponding compartment. Press the power switch at the back to its On position (-).
2. Press On/Off key; display will automatically light on and go through the changes indicated below until finally displaying 0.0 g (series 5171) or 0.00 g (series 5172).

Note: The initializing sequence depends on the stability of the load cell, thus, the balance must not be located in an unstable surface or near wind blows.

Warming up

Before performing calibration on the balance or measuring its accuracy it is necessary to ensure that the balance is in a stable state.

For balance stabilization it is important that environmental temperature is stable.

Leave the balance on and in weighing mode for at least half an hour before performing calibration; this is called “warming up”.

Calibration

Perform balance calibration each time after moving it. Weights are required for external calibration balances.

Necessary calibration weight depending on maximum capacity of the balance

Max. capacity	Calibration weight
300 g	200 g
500 g	500 g
1000 g	500 g
2000 g	1000 g
3000 g	2000 g
5000 g	2000 g

Before performing calibration, warm up the balance in advance.

1. Remove any substance or object from the pan and press TARE key to change displayed value to 0.0g or 0.00g.
2. Keep pressed CAL key for 3 seconds until “---CAL---” is displayed.
3. Then the value of calibration weight will flash on display.
4. Put the calibration weight on the pan; after a few seconds “-----” will be displayed, wait until the value of calibration weight is shown on display.
5. Remove the weight; for a few seconds “-----” will be displayed until the balance enters stable weighing mode showing 0.0g or 0.00g.
6. It is recommended to make this operation twice to achieve the best calibration result.
7. **Lineal calibration:** keep pressed CAL key for a 3 seconds until “---CAL---” is displayed and then the value of calibration weight. Keep pressed PCS key to enter in Lineal calibration mode and follow the instructions on display.

Note: for lineal calibration is necessary the use of several calibration weights. It is recommended to perform lineal calibration only when big tolerance (error) is detected in repeatability or linearity of the balance. It is not recommended to make lineal calibration frequently.

4. USE OF THE BALANCE

Weighing function

Once the balance is warmed up and calibrated, it will stay in weighing mode showing on display 0.0g or 0.00g.

Put the sample to be weighed on the pan and when stable the weight of the sample will be displayed.

Tare function

1. Place a recipient on the pan; its weight will be displayed.
2. Press TARE key to change the value on display to 0.0g or 0.00g; this indicates that the weight of the recipient has been deducted.
3. Place the object or substance to be weighed in the recipient and when the balance stabilizes its weight will be displayed.

Piece counting function

1. Remove any substance or object from the pan and press TARE key to change displayed value to 0.0g or 0.00g.
2. Press PCS key to enter in counting mode; display will successively show “-COU-”, “- - -” and “10” (n^o reference n^o of pieces by default).
3. Press CAL key to select the desired reference number of pieces (10, 20, 50, 100, 150, 200, 250 or 500).
4. Place the corresponding number of pieces on the pan and press PCS key to save setting.
5. Remove the pieces from the pan and the balance is ready to be used in piece counting mode (weighing unit will change from g to pcs).

Note: In counting mode, the weight of the pieces must be even and the weight of an individual piece should not be less than the readability of the balance.

6. Press PCS key to return to weighing mode.

Weighing unit selection

Press UNIT key to select the desired weighing unit among g (grams), ct (carats), oz (ounces) y lb (pounds).

The weighing unit by default is g.

Overload

The weight of the sample cannot exceed the rated maximum capacity of the balance; when exceeded, the upper line “-----” will be displayed. Remove the sample immediately to avoid damages in the balance.

RECOMMENDATIONS OF USE

- Warming up is necessary before any operation is made to the balance.
- In TARE mode, the value of the tared object cannot exceed the rated maximum capacity of the balance.
- Calibration is necessary to ensure a reliable weighing.
- Remove batteries in case the balance will not be used for a long time.
- Do not mix new batteries to used batteries inside of the balance.
- Anode and cathode of the batteries need to be connected properly.
- Switch off the balance when it is not used.
- It is suggested to turn the round platter clockwise before take it off the balance.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the balance it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- Please follow the processes of use of this manual.
- This manual should be available for all users of this equipment.
- Prevent the balance from sudden movements and knocks, as well as from direct sunlight or air flows. The balance is a precision instrument, you must handle it carefully.
- The balance is supplied with an adapter. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- Never unplug the adapter by pulling the wire, do it from the base.
- Never use the balance in a wedged location as for example a shelf.
- Never use sharp objects as pens, etc... to press the buttons of the control panel; only use your fingers.
- Never place on the pan an object heavier than the maximum capacity of the balance, the sensor could be damaged.
- Neither submerge the balance nor spill liquids on it.
- When you are not using the balance for a long period of time, lock the rechargeable battery.
- If any liquid comes into contact with the electric parts of the balance please do immediately disconnect it from the net and send it to the technical service, as soon as possible, for its checking and adjustment.


- Always use original components and supplies. Other devices can be similar but they can damage the equipment.

Cleaning

- Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the balance and produce an early ageing of the equipment.
- Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

Troubleshooting

Problem	Cause	Solution
Upper line - - - -	Overload	Make sure the weight of the sample does not exceed the maximum capacity of the balance. Re-calibrate the balance
Under line _ _ _ _	Sample is not properly loaded on the balance	Re-calibrate the balance
Err-1	Too frequently turn on and off the balance	Turn off the balance, resume it after 3 seconds
Err-2	The balance is not stable	Wait for a few seconds for stabilization
	Low battery	Replace batteries or use the balance connected to the power supply

Note: According to the in force legislation regarding “Non-automatic weighting instruments” in which balances are included, by means of writ dating from 22nd October, 1994 (BOE 1/3rd/95), these balances must not be used for:

- *Commercial transactions*
- *Calculating taxes, tariffs, rates, indemnities and other similar canons*
- *Judicial surveys*
- *Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories*
- *Determining the price or total amount in sale price or in pre-packaged preparations*

Instructions on environmental protection

Nahita-Blue electronic balances series 5171/5172 can operate with rechargeable batteries. Batteries contain heavy metals (zinc, lithium, mercury, etc.) that could contaminate water, soils, etc.

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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