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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita-Blue analytical balance. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- This manual is inseparable from the Nahita-Blue analytical balance series 5134IN, 5134EX and 5133, so it should be available for all the users of this equipment.
- You should carefully handle the balance avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the balance.
- To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler.
- This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- Overhaul is not covered by the equipment warranty.
- Operations made by non-qualified staff will automatically produce a loss of the warranty.
- Accessories (including their loss), are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- Manufacturer reserves the right to modify or improve the manual or equipment.

INDEX OF LANGUAGES

Spanish.....	2-16
English.....	17-30
French.....	31-45

INDEX OF CONTENTS

1. SPECIFICATIONS.....	19
2. DESCRIPTION.....	21
3. INSTALLATION.....	22
4. USE OF THE BALANCE.....	23
5. MAINTENANCE AND CLEANING.....	28

1. SPECIFICATIONS

Analytical balances series 5134IN, 5134EX and 5133 are high performance equipments that offer fast and accurate measurements and use a high precision electromagnetic force balance sensor, improving the reliability of the balances still further.

Technical parameters


Model	External calibration	55134210	55134220
	Internal calibration	55134211	55134221
Accuracy class	I		
Max. capacity	120g	220g	
Min. capacity	10mg	10mg	
Division	0.0001g		
Verification scale value	0.001g		
Repeatability	±0.0002g		
±0.5e	0≤m≤50000		
±1.0e	50000<m≤200000		
±1.5e	200000<m		
Stabilization time	≤8 seconds		
Operating temperature	17.5°C~22.5°C, fluctuation range<1°C/h		
Relative humidity (RH)	50%~75%		
Pan size	Φ 80 mm		
Windshield	225x220x265 mm(L*W*H)		
Total dimensions	355x230x350 mm(L*W*H)		
Net weight	7.2kg		
Power	AC110-220V 50Hz / DC9V-2.2A		

Modelo	55133210	55133230	55133250
Accuracy class	Ⓜ	Ⓜ	Ⓛ
Max. capacity	100g	300g	500g
Min. capacity	20mg	20mg	20mg
Division	0.001g		
Verification scale value	0.01g		
Repeatability	±0.001g		
±0.5e	0≤m≤5000		
±1.0e	5000<m≤20000		
±1.5e	20000<m		
Stabilization time	≤5 seconds		
Operating temperature	15°C~30°C, fluctuation range<5°C/h		
Relative humidity (RH)	40%~80%		
Pan size	Φ90mm		
Windshield	225x220x265 mm (L*W*H)		
Total dimensions	355x230x350 mm(L*W*H)		
Net weight	7.2kg		
Power	AC110-220V 50Hz DC9V-2.2A		


- Electromagnetic force compensation technology
- Aluminum alloy die cast base and mains
- Stainless steel pan
- Easy to read LCD display with backlit
- Glass windshield with sliding doors
- Height adjustable feet
- RS232 and USB outputs
- Full capacity taring
- Below balance weighing facility
- Overload protection
- Selectable measure units: g, mg, ct, oz
- Piece counting and percentage functions

2. DESCRIPTION

Structure

Display and panel descriptions

-
1. **PRT:** Data print and menu.
 2. **OFF:** Switch off.
 3. **TAR:** Tare
 4. **ON:** Switch on.
 5. **CAL:** Calibration
 6. “±” Symbol of positive (over 0) and negative (below 0) weight.
 7. Stabilization symbol

3. INSTALLATION

Choosing the installation place

When analytical balances series 5134IN, 5134EX and 5133 are used under routine test lab or industrial measuring room, the weighing speed will be much quicker and the result more accurate. The conditions of installation place of the balance should be as below:

1. The working room should be clean and dry.
2. The balance should be placed on a solid, stable, horizontal and plain flat surface.
3. Avoid locations in which the balance can be exposed to any of the following situations:
 - b) Airflows coming from air conditioners, fans, doors or windows.
 - c) Vibration from surrounding or nearby equipments
 - d) Direct sunlight and radiation
 - e) Electromagnetic waves or fields
4. Do not use the balance anywhere exposed to explosive, combustible or corrosive gases.
5. Don't use the balance under the condition of high humidity or high dust.
6. When moved from a colder to a warmer place, the accuracy and reliability of the scale will be influenced by the moisture condensation inside the balance. In order to eliminate this influence, it is better to put the scale in the working place without power supply for 2 hours.
7. Avoid extreme temperatures or exposure to direct sunlight or air conditioners. The ideal working temperature is 20°C-25°C with a fluctuation temperature not higher than 1°C/h.
8. Keep the balance clean.
9. Do not stack material on the balance scale when not in use.
10. Use the correct power supply and voltage with the balance; the balance is supplied with a power adapter 110-220V/50-60Hz.

Unpacking and preliminary inspection

Check that all of the items indicated below are included in the package, and that

nothing has been damaged:

- Balance
- Pan
- Power cable (AC adapter)
- Instruction manual
- External calibration weight

In case the balance or any component presents any damage due to transport tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

4. USE OF THE BALANCE

Level adjustment

Once the balance is placed on its location, check the bubble level; if bubble is not well centered, turn the adjustable feet so as the bubble moves to the center and the balance is perfectly leveled.

Start

1. Connect the balance to the power supply and press ON key on the panel.
2. When the balance turns on, a countdown of 30 min will be displayed showing the necessary pre-warming time for a better performance and stable state of the balance.
3. Once pre-warming time has passed it is recommended to calibrate the balance. In internal calibration models once pre-warming time has passed, the balance will automatically make calibration.
4. After calibration, the balance will enter weighing mode.
5. To turn the balance off, press OFF key on the panel and display will light off. If the balance will not be used for a long time, disconnect it from the power supply.

Calibration

Internal calibration models (series 5134IN)

1. Press CAL key
2. Calibration using the internal weight will be made automatically.

If tolerance (error) after internal calibration is too high, then it is necessary to make external calibration. For this, proceed as follows:

1. Press simultaneously CAL and ON keys until F" appears at the upper side of

display; after releasing the keys CAL-200 (value of the necessary weight to make external calibration) will be displayed.

2. Put the calibration weight on the pan and wait until its weight is displayed (e.g. 200.0000g).
3. Remove the weight and the balance will enter weighing mode.

Note: Each time the external calibration is made it is necessary to make internal calibration first.

External calibration models (series 5134EX and 5133)

1. Press CAL key; the value of the necessary weight will flash on display.
2. Put the calibration weight on the pan; for this, open the glass door, put the weight on the pan and close the door again. Wait until the value of calibration weight stops flashing and is fixed on display (e.g. 200.0000g).
3. Remove the weight from the pan; for this, open the glass door, remove the weight from the pan and close the door again.
4. The balance will enter weighing mode and will display 0.000g or 0.0000g.

Note: Each balance includes the necessary weight to make external calibration.

Tare function

1. Put the recipient on the pan; the corresponding weight will be displayed.

+18.9001g

2. Press TAR key and Reading on display will automatically set to 0.000 g or 0.0000 g; tare is complete.

0.0000g

4. If the recipient is removed from the pan, the corresponding weight with negative value will be displayed.

-18.9001g

6. Press TAR key again and the value on display will be set to zero.

0.0000g

Weighing unit selection

Press PRT key to select the desired weighing unit among mg (UNT-0), g (UNT-1), ct (UNT-2) and oz (UNT-3).

The weighing unit by default is g.

Piece counting function

1. Keep pressed ON key until Set-0 is displayed.
2. Then, press PRT key to select the reference number of pieces (COU-10, COU-25 or COU-50).
3. Put the corresponding number of pieces on the pan, press CAL key and wait until the number of pieces is displayed.
4. Remove the pieces from the pan and the balance is ready to be used in piece counting mode (weighing unit will change from g to pcs).

Note: In counting mode, the weight of the pieces must be even and the weight of an individual piece should not be less than the readability of the balance.

To return to weighing mode:

1. Keep pressed ON key until Set-0 is displayed.
2. Press PRT key to select COU-00
3. Press TAR key and the balance will enter weighing mode.


Percentage function

1. Keep pressed ON key until Set-0 is displayed.
2. Then, press PRT key to select the parameter COU-100.0%.
3. Put on the pan the reference simple, press CAL and wait until 100.0% is displayed.
4. Remove the sample from the pan and the balance is ready to be used in percentage mode (weighing unit will change from g to %).

To return to weighing mode:


1. Keep pressed ON key until Set-0 is displayed.
2. Press PRT key to select COU-00
3. Press TAR key and the balance will enter weighing mode.

Ajuste del nivel de sensibilidad


To see the symbol displayed for each sensitive level please check *Display and panel description*

Weighing speed setting


To see the symbol displayed for each speed level please check *Display and panel description*

Data output setting

To see the symbol displayed for data output mode please check *Display and panel description*

Baudrate setting

RS232 interface

- Connection

	Balance (9 pins)	PC/Printer (9 pins)
RXD (Input)	2	3
TXD (Output)	3	2
GND (Ground)	5	5

- The baudrate by default is 1200 bps (see *Baudrate setting*)

- Data format: 10 bits, 0 as start bit, 1 as stop bit, 8 digits (ASCII code)

- No odd and even numbers adjusting

- Data output: by default is continuous mode. The data output mode can be changed into press output, timing output and continuous output (see *Data output setting*)

- Output data format

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
type	Space	Space or *	±	data	data	data	data or dot	data or dot	data	data	data	Unit1	Unit2	end	return

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the balance it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- Please follow the processes of use of this manual.
- This manual should be available for all users of this equipment.
- Prevent the balance from sudden movements and knocks, as well as from direct sunlight or air flows. The balance is a precision instrument, you must handle it carefully.
- The balance has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- Never unplug the balance by pulling the wire, do it from the base.
- Never use the balance in a wedged location as for example a shelf.
- Never use sharp objects as pens, etc... to press the buttons of the control panel; only use your fingers.
- Never place on the pan an object heavier than the maximum capacity of the balance, the sensor could be damaged.

- Neither submerge the balance nor spill liquids on it.
- When you are not using the balance for a long period of time, lock the rechargeable battery.
- If any liquid comes into contact with the electric parts of the balance please do immediately disconnect it from the net and send it to the technical service, as soon as possible, for its checking and adjustment.
- Always use original components and supplies. Other devices can be similar but they can damage the equipment.

Cleaning

- Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the balance and produce an early ageing of the equipment.
- Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

Troubleshooting

Problem	Cause	Solution
No display	No power supply; Fuse damaged; Power transformer is damaged;	Plug in adapter; Change the fuse; Change the power transformer; If problem persists, send the balance to the Technical Service for repair.
Unstable display	Bad working conditions; The windshield is open; Something between the scale pan and working table. The power exceeds its permissible value and is unstable; Unstable sample weight (moisture evaporation)	Improve the working condition, avoid vibration and airflow; Close the windshield; Remove the pan and clean well the balance surface; Connect the balance to power supply 110-220 V AC
Difference between displayed value and real value	The balance has not been calibration The weight of the recipient has not been tared The balance is not horizontal	Make calibration Make tare Adjust level of balance

Note: According to the in force legislation regarding “Non-automatic weighing instruments” in which balances are included, by means of writ dating from 22nd October, 1994 (BOE 1/3rd/95), these balances must not be used for:

- *Commercial transactions*
- *Calculating taxes, tariffs, rates, indemnities and other similar canons*
- *Judicial surveys*
- *Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories*
- *Determining the price or total amount in sale price or in pre-packaged preparations*

Instructions on environmental protection

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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