**IMANES –MAGNETS**

**REF:QLL1-QLL04**

**LIMADURA DE HIERRO**

**IRON FILLING**

**REF:QLL05**

**INDEX OF LANGUAGES- ÍNDICE DE IDIOMAS**

**ESPAÑOL .....2**

**ENGLISH .....4**

## IMANES Y LIMADURA DE HIERRO

### Descripción

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## MAGNETS AND IRON FILLING

Different models of magnets for performing experiments of magnetism and visualization of magnetic fields.

Code	Shape	Dimensions	Material	Quantity/pa
QLL001	Rectangular	75x15x10 mm	Alnico	2 uds/pack
QLL002	Cilyndrical	100x12 mm	Alnico	2 uds
QLL003	Ring	36x18x8 mm	Ceramics	1 ud
QLL004	Horeshoe	75x12x5 mm	Chromed steel	1ud

### Theory

Magnetism is a class of physical phenomena that are mediated by magnetic fields. Electric currents and magnetic moments of elementary particles give rise to a magnetic field which acts on other currents and magnetic moments. The most familiar effects occur in ferromagnetic materials, which are strongly attracted by magnetic fields. The most common ones are iron, nickel and cobalt and their alloys.

Each electron is a little magnet. Materials have electrons in all directions but a magnet has all electron in one direction, there is a magnetic effect depending on the number of electrons.

### Obsevatons

Each type of magnet has its code. In order to do experiments its normal to use iron filling. For the study of magnetic fields. Supplied in box of 500g.

REF:QLL005

### Other utilities of iron filling

Iron filling can be used to avoid experiments about homogeneous and heterogeneous mixtures. In case of heterogeneous mixture, there is a common example.

### Experiement

Mix sand, salt, water and iron filling to have a heterogeneous mixture. Then you have to separate all the components.

1. Mix all the components in a beaker.

Now separate

2. Filtering, you can separate water with salt and in the filter sand with iron filling.
3. Crystallization and magnetic separation, the sand and the iron filling can be separated.
4. Last step, evaporation of water in order to have the salt deposited.

Image: filtration, crystallization and evaporation of water.