INCLINATORIO
COMPASS INCLINOMETER
REF: QLL012

INDEX OF LANGUAGES-ÍNDICE DE IDIOMAS

ESPAÑOL.....	2
ENGLISH.....	3

Inclinatorio

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Inclinometer

The inclinometer is useful to measure the local inclination of the magnetic field and also for the representation of the magnetic field of a current.

To demonstrate the earth's magnetic field and to illustrate the thumb rule. It consists of an aluminum ring with a double pointed needle pivoting in the center, the ring is graduated from 0-360° in both clockwise and anticlockwise directions. The whole assembly is fitted on an aluminium pillar and mounted on a plastic base of approx. 142x80x40 mm

Use

Cautions

- Protect the device from humidity, dust and mechanical blows.
- Do not touch the magnetic needle.

The geometry of the lines of the Earth's magnetic field is strongly changed by static magnetic changes, steel frames on laboratory tables and installations, steel beams on the floor, ceilings and building walls. For this reason you can not avoid having deviations in the angles to expect.

Measure

The needle is oriented along the direction of the actual course of the Earth's magnetic field lines. With the scale plane horizontally, the blue side of the needle is oriented in a northerly direction, the needle is oriented so that it is at 0° (blue side of the needle shows north direction).

Then, the fork is turned 90° (vertical plane of the scale). The magnetic needle tilts with the blue side facing down.

The deviation of the magnetic needle from the horizontal is called a tilt. This is the different depending on place, in the northern parallel of Europe approx 50° is between 63° and 68°.

Magnetic effect of an electric current

For the realization of the experiment an adjustable DC source is additionally required.

Taking the scale plane horizontally, the blue side of the needle should be oriented so that is at 0° (blue side of the needle shows north).

The connecting bushes are connected to an adjustable DC source as the current increases the magnetic needle experiences an additional deviation. Changing the polarity of the source changes the direction of the deviation.