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VAN DE GRAAFF GENERATOR

Description

Very useful for electrostatic experiments where a continuous source of high voltage is needed. The aluminum base presents isolated sides made of plastic for a higher safety. An inner motor allows the continuous movement of a rubber belt between 2 pulleys, and this produces the charge transfer from the belt to the metal dome by means of charge collecting combs. The hollow dome is especially designed with a smooth, polished and free of any imperfection surface and without any sharp corner or edge to minimize charge leakage.

The motor operates on AC power supply 220-240 V, 50/60 Hz and is activated through a on/off switch.

The equipment presents a socket for 4 mm banana plug at the base for earth connection and another socket at the top for charge transfer or for attaching accessories. It is supplied complete with accessories and a discharge sphere with isolating handle.

When bringing earthed discharge sphere sufficiently close to the charge dome, transfer of the charge from dome to the ground takes place in the form of an electric spark jumping from dome to the discharge sphere; under favorable conditions, it can develop electric potential up to 200 kV with a spark of up to 70mm length.


Overview

1. Hollow metal globe over which the electric charge accumulates
2. Nylon belt which carriage the charge
3. Mains switch to power on /off.
4. Earthing socket of discharger
5. Discharger
6. Earthing socket of power supply/operation unit.

Handling instructions/ precautions

1. Equipment intended for use in educational establishments by student operators at least 11 years old under supervision of the responsible body.

2. Student operators in education with special needs must use equipment under strict supervision of responsible body with fully safety and precautions.
3. Installation cleaning and decontamination also must be done under the supervision of the responsible body.
4. Take care for the duty cycle of the apparatus, which is about 10 minutes aprox.

5. Read all stickers placed on the apparatus
6. Disconnect from mains before operating the apparatus

7. Ensure proper grounding of the apparatus

8. Equipment must be disconnect from the hazardous live voltage before Access.
9. Use discharge rod to discharge the charged sphere after every use of the equipment

10. Earthing symbol

11. Be cautious, static electricity exists

Features:

1. Hollow metal globe of 200mm diameter
2. Spark length of 40mm approx.
3. Duty cycle of the equipment is 10 minutes approx.
4. Easy replaceable nylon belt.

Technical specifications

Input Voltage	220VAC
Frequency	50-60Hz
Current Rating	0.3A
Power	66W(max.)
Fuse (Primary)	2A
Temperature Range	Room Temperature
Humidity	40% to 60%

Working procedure

1. Adjust the hollow metal globe on the cylindrical acrylic pipe, as


shown in the picture

2. Connect the earthing socket for discharger to the earthing socket of power supply unit using the connecting lead provided.
3. Plug-in the power cord in the mains in let of the apparatus and connect to the 220 V AC
4. Switch on the mains switch.
5. Motor inside the apparatus, will start rotating the nylon belt.
6. As the belt will start rotating it passes in front of the lower brush, it receives negative charge that escapes from its points due to the influence of the electric field around the lower roller which ionizes the air at the points.
7. As the belt touches the lower roller, it transfers some electrons, leaving the roller with a negative charge, which added to the negative charge in the belt generates enough electric field to ionize the air at the points of the upper brush.

8. Electrons then leak from the belt to the upper brush and to the terminal, leaving the belt positively charged as it returns down and the terminal negatively charged.
9. The hollow metal globe shields the upper roller and brush from the electric field generated by charges that accumulate at the outer surface of it, causing the discharge and change of polarity of the belt at the upper roller to our practically.
10. As the belt continues to move, a constant charging current travels via the belt, and the sphere continues to accumulate negative charge until the rate that charge is being los equals the charging current. The larger the sphere and the farther it is from ground, the higher will be its final potential.

Wiring diagram and PCB circuit connections

Maintenance

In case of fuse blown: disconnect it from mains then replace the fuse with spare fuse of 2A provided in the in built cartridge of mains inlet.

*Fuse not be charged by a student operator.

Verification after fuse replacement

Connect the mains to the apparatus


Switch on the rocker switch mounted on front panel.

In built indicator of rocker switch will glow on and apparatus will start working

In case of belt replacement:

Remove the hollow metal globe

Disassemble the upper roller tightened with butterfly nuts and pull out the belt.

Rotate the cylindrical acrylic pipe to the clockwise direction, as indicated on the upper panel of the power supply unit.

Pull out the belt from lower roller as well

Insert the new belt over the lower roller and assemble cylindrical acrylic pipe on the upper panel of power supply unit.

Now, insert another end of new belt over the upper roller and tightened that upper roller with the help of butterfly nuts as previous.

Verification after the belt replacement

Place the hollow metal globe

Connect the mains to the apparatus

Switch on the rocker switch mounted on front panel.

Belt will start moving on the lower and upper rollers.