

Rotary evaporator **RS 100 PRO**



model RS 100-PRO

We present the new rotary evaporator RSLAB model RS 100-PRO which expands our product line of evaporators. It includes some innovations: designed a detachable operating panel allows remote control with large LCD display for a digital control of lot of parameters as speed, temperature, timer.

Available in two different versions, with normal borosilicate glass set and with coated borosilicate glass set. This version includes a complete set of glassware with plastic coating to improve safety and to avoid contamination of samples with glass pieces in case of accidental breakage; this feature is necessary in laboratories of food and sanitary analysis.


Its range of applications is very large and provides reactions of distillation and separation liquid / liquid, liquid / solid, etc. This model can be used in many fields such as chemistry and biochemistry, food, industrial processes, research and education.

Features

- [01] Wide LCD display provides together informations of speed, temperature and operation time
- [02] Speed range from 20 up to 280 r.p.m.
- [03] Large 5L water-oil heating bath with heating temperature range of RT to 180°C
- [04] Digital temperature control and adjustable safety circuit
- [05] Designed a detachable operating panel allows remote control
- [06] Patented condenser (cooling surface 1,500 cm²) with excellent cooling effect
- [07] Motorized lift with quick-action, automatic release evaporating flask to top position in case of power failure
- [08] Adjustable final position recognition to protect operator and sample against breakage
- [09] Evaporating flask with an ejector, convenient to remove
- [10] Chemical-resistant double PTFE system and patented pressure spring provide excellent sealing
- [11] PC can control and document all parameters via USB



rotary evaporator RS100 PRO



1

Patented condenser(cooling surface 1,500 cm²) with excellent cooling effect

3

Chemical-resistant double PTFE system and patented pressure spring provide excellent sealing

5

Wide LCD display provides together informations of speed, temperature and operation time

4

Large 5L water-oil heating bath with heating temperature range of RT to 180°C

technical specifications

Code	PGB006	PGB007
Type of glass	Borosilicate	Borosilicate coated
Power [W]	1400	
Motor	DC brushless motor	
Speed range [r.p.m.]	20 - 280	
Reversible direction of rotation	Yes	
Temperature [°C]	RT - 180	
Heat control accuracy [°C]	+/- 1	
Temperature display	LCD	
Lift	Motorized	
Stroke [mm]	150	
Timer	Yes	
Time setting range [min]	1 a 999	
Dimensions [D x P x H mm]	465x457x583	
Weight [Kg]	15	
USB	Yes	
Protection Class acc. DIN EN60529	IP20	

Display



rotary evaporator RS100 PRO

modelo RS 100-PRO accessories / spare parts

[01] Complete glass set



Including condenser, evaporating flask 1 L (29/32)
and receiving flask 1 L (35/20)


Code	Type
PGG014	Borosilicate glass
PGG015	Borosilicate glass, plastic-coated

[02] Evaporating flasks 29/32



Code	Volume	
Borosilicate glass	Borosilicate glass, coated	
PBJ001	PBJ007	50 mL
PBJ002	PBJ008	100 mL
PBJ003	PBJ009	250 mL
PBJ004	PBJ010	500 mL
PBJ005	PBJ011	1000 mL
PBJ006	PBJ012	2000 mL

[04] Vapor tube 29/32



Code	NS
PGG020	NS 29/32

[03] Receiving flasks 35/20



Code	Volume	
Borosilicate glass	Borosilicate glass, coated	
PGG024	PGG029	100 mL
PGG025	PGG030	250 mL
PGG026	PGG031	500 mL
PGG027	PGG032	1000 mL
PGG028	PGG033	2000 mL

model RS 100-PRO accessories / spare parts

rotary evaporator RS100 PRO

[05] Connectors



Code	NS
PGG016	NS 29/32, NS24/29
PGG017	NS 29/32, NS19/26
PGG018	NS 29/32, NS14/23

[06] Spider with 5 flasks



Code	Volume*
PGG023	5x50 mL (29/32)*

* not included

[07] Foam brake 29/32



Code	Volume
PGG019	250 mL

[08] Seal



Code
PGG021

vacuum pump



Referencia	PXB005
Vacío (máx.)	13 mbar
Velocidad	0 mmHg 25L/min 100 mmHg 22L/min 200 mmHg 20L/min 300 mmHg 12L/min 400 mmHg 7.5L/min 550 mmHg 5L/min
Velocidad motor	1450 rpm
Peso [Kg]	7

rotary evaporator RS100 PRO

modelo RS 100-PRO accessories

digital vacuum regulator



Code	PXB006
Model	DV380
Display	The default display is the pressure in the attached piece of equipment in units of mm of Hg (torr).

Modes of regulation

Meter: Evacuates the attached equipment to the set point pressure entered in the controller.

Off: Isolates the vacuum pump from the attached equipment. In this position the equipment can be opened to the atmosphere while the vacuum pump is running.

Full vacuum (vacío total): Opens the regulator continuously to the vacuum pump (independent of set point pressure). The pressure in the attached equipment is lowered to the limit achievable by the vacuum pump..

recirculating chiller



Code	PXB007
Model	F250
Temperature range	-10°C à +40°C
Temperature stability	+/- 0.5°C
Temperature display	LED
Pump capacity flow rate	15 L / min.
Pump capacity flow pressure	0.35 bar
Cooling capacity	20°C 0.25 KW
	10°C 0.22 KW
	0°C 0.21 KW
	-10°C 0.18 KW
	-20°C 0.09 KW
Filling volume	1.7 - 2.6 Litres
Dimensions [H x L x P cm]	52 x 42 x 24 cm
Weight [Kg]	27 Kg