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- Yellow (red + green) + blue = white
- Cyan (green + blue) + red = white
- Magenta (red + blue) + green = white

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NEWTON CONLOR DISK MOUNTED ON MOTOR (Ref. QLG002)

USER INSTRUCTIONS

Purpose

The Newton colour disk is used to demonstrate that white light is a resultant of additive mixing of different spectral colours.

Description

The instrument is composed of a multi-coloured disk having sectors of various spectral colours in proper proportions. The disc is connected to a motor and mounted on a stand. When the disc turns by the action of the motor, all colours will be observed to disappear by merging into a single white colour.

Necessary additional elements

- Power supply, ref. QLN001
- Red/Black connection cables, ref. QLN005 (black) y QLN006 (red)

Procedure

- Place the Newton colour disk on a stable surface.
- Connect the instrument to the power supply and select a voltage of 6 V; the disk will start rotating.
- At this rotating speed, the eye is not able to follow the movement of the colour sectors and produces their

additive mixing. Due to technical causes, a grey colour, not white, is observed.

- It is advisable to apply a frontal white light to the colour disk.

Some theory about primary and secondary colours

When white light is dispersed by a prism, seven colours (red, orange, yellow, green, blue, indigo and violet) can be observed on a white screen. However, the combination of only three colours (red, Green and blue in their proper intensities) is necessary to give white light.

By using three light boxes with filters red, green and blue respectively and making the three beams overlap, white light is observed where all the three colours overlap. When just red and green lights overlap, yellow light is seen. On mixing green and blue light cyan colour is seen. Similarly, on mixing blue and red light magenta is seen.

Red, green and blue colours are called primary colours because they cannot be made by the combination of the other colours, they can be combined to form all the other colours and all three combined together give white light. Yellow, cyan and magenta, the three colours formed by the combination of various pairs of primary colours are known as secondary colours. When a secondary colour is combined with the remaining primary colour white light is seen.

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Colours which add together to give white light are called complementary colours.