**BLOQUES  
CALORIMETRICOS**

**CALORIMETER  
BLOCKS**

**ALUMINIO-QLJ004**

**ACERO- QLLJ005**


**COBRE- QLJ006**

**LATÓN- QLJ007**

**RESISTENCIA**

**HEATING RESISTANCE**

**QLJ008**

**INDEX OF LANGUAGE-ÍNDICE DE IDIOMAS**

ESPAÑOL..... 2

ALUMINIO ..... 2

ACERO ..... 3

COBRE..... 4

LATÓN ..... 5

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
## CALORIMETER TYPES

### CALORIMETER BLOCK, ALUMINIUM

It is a simple calorimeter facilitating quick experimental determination of the specific heat capacity of Aluminium. It comprises cylindrical metal block of mass adjusted to  $1\text{kg} \pm 2\%$ . The metal block is drilled with two holes, a large central hole, about 12.5 mm diameter to accept special immersion heater and a smaller hole about 7.5 mm diameter, to accept thermometer or temperature sensor. The diameter of Aluminium block is 76 mm and height is 84 mm.


#### **Additional required components (not included):**

- Power supply 1-12 V AC/DC (ref. QLN001)
- Immersion heater (ref. QLJ008)
- Thermometer
- Stopwatch

#### **How to use**

1. Weigh the calorimeter block and record its mass.
2. Place the calorimeter block on a heat proof mat surrounded by insulating material, so that the heat losses are kept to the minimum.
3. Insert the immersion heater and the thermometer into the appropriate hole. Drop before some oil or water into the thermometer hole to ensure good thermal contact between the thermometer and the block.
4. Set up the circuitry according Fig.1.
5. Switch on the power supply and adjust it to give a current of about 4 A. Switch again the power supply off.
6. Before starting the experiment, wait for a few minutes. Then, take note of the initial temperature of the block.
7. Switch the power supply on and start measuring time with stopwatch.
8. Wait until the temperature has risen about  $20^{\circ}\text{C}$ .
9. Record the time and final temperature

The specific heat capacity can then be calculated using the relation:

$$I \cdot U \cdot t = m \cdot c \cdot (\theta_2 - \theta_1)$$

With  $I$ : Current,  $U$ : Voltage,  $t$ : Time,  $m$ : Mass of the calorimeter block,  $c$ : Specific heat capacity,  $\theta_1$ : Initial temperature,  $\theta_2$ : Final temperature

**Precautions:** There is a risk of burns from heater or calorimeter, therefore, allow the apparatus to cool before moving it.

## STEEL

It is a simple calorimeter facilitating quick experimental determination of the specific heat capacity of Steel. It comprises cylindrical metal block of mass adjusted to  $1\text{kg} \pm 2\%$ . The metal block is drilled with two holes, a large central hole, about 12.5 mm diameter to accept special immersion heater and a smaller hole about 7.5 mm diameter, to accept thermometer or

temperature sensor. The diameter of Steel block is 44 mm and height is 89 mm.


### Additional required components (not included):

- Power supply 1-12 V AC/DC (ref. QLN001)
- Immersion heater (ref. QLJ008)
- Thermometer
- Stopwatch

### How to use

1. Weigh the calorimeter block and record its mass.
2. Place the calorimeter block on a heat proof mat surrounded by insulating material, so that the heat losses are kept to the minimum.
3. Insert the immersion heater and the thermometer into the appropriate hole. Drop before some oil or water into the thermometer hole to ensure good thermal contact between the thermometer and the block.
4. Set up the circuitry according Fig.1.
5. Switch on the power supply and adjust it to give a current of about 4 A. Switch again the power supply off.

6. Before starting the experiment, wait for a few minutes. Then, take note of the initial temperature of the block.
7. Switch the power supply on and start measuring time with stopwatch.
8. Wait until the temperature has risen about 20°C.
9. Record the time and final temperature

The specific heat capacity can then be calculated using the relation:

$$I \cdot U \cdot t = m \cdot c \cdot (\theta_2 - \theta_1)$$

With  $I$ : Current,  $U$ : Voltage,  $t$ : Time,  $m$ : Mass of the calorimeter block,  $c$ : Specific heat capacity,  $\theta_1$ : Initial temperature,  $\theta_2$ : Final temperature

**Precautions:** There is a risk of burns from heater or calorimeter, therefore, allow the apparatus to cool before moving it.

## COPPER

It is a simple calorimeter facilitating quick experimental determination of the specific heat capacity of Copper. It comprises cylindrical metal block of mass adjusted to  $1\text{kg} \pm 2\%$ . The metal block is drilled with two holes, a large central hole, about 12.5 mm diameter to accept special immersion heater and a smaller hole about 7.5 mm diameter, to accept thermometer or temperature sensor. The diameter of Copper block is 44 mm and height is 79 mm.


### Additional required components (not included):

- Power supply 1-12 V AC/DC (ref. QLN001)
- Immersion heater (ref. QLJ008)
- Thermometer
- Stopwatch

### How to use

1. Weigh the calorimeter block and record its mass.
2. Place the calorimeter block on a heat proof mat surrounded by insulating material, so that the heat losses are kept to the minimum.

3. Insert the immersion heater and the thermometer into the appropriate hole. Drop before some oil or water into the thermometer hole to ensure good thermal contact between the thermometer and the block.
4. Set up the circuitry according Fig.1.
5. Switch on the power supply and adjust it to give a current of about 4 A. Switch again the power supply off.
6. Before starting the experiment, wait for a few minutes. Then, take note of the initial temperature of the block.
7. Switch the power supply on and start measuring time with stopwatch.
8. Wait until the temperature has risen about 20°C.
9. Record the time and final temperature

The specific heat capacity can then be calculated using the relation:

$$I \cdot U \cdot t = m \cdot c \cdot (\theta_2 - \theta_1)$$

With  $I$ : Current,  $U$ : Voltage,  $t$ : Time,  $m$ : Mass of the calorimeter block,  $c$ : Specific heat capacity,  $\theta_1$ : Initial temperature,  $\theta_2$ : Final temperature

**Precautions:** There is a risk of burns from heater or calorimeter, therefore, allow the apparatus to cool before moving it.

### BRASS

It is a simple calorimeter facilitating quick experimental determination of the specific heat capacity of Brass. It comprises cylindrical metal block of mass adjusted to  $1\text{kg} \pm 2\%$ . The metal block is drilled with two holes, a large central hole, about 12.5 mm diameter to accept special immersion heater and a smaller hole about 7.5 mm diameter, to accept thermometer or temperature sensor. The diameter of Brass block is 44 mm and height is 85 mm.


### Additional required components (not included):

- Power supply 1-12 V AC/DC (ref. QLN001)
- Immersion heater (ref. QLJ008)
- Thermometer

- Stopwatch

### How to use

1. Weigh the calorimeter block and record its mass.
2. Place the calorimeter block on a heat proof mat surrounded by insulating material, so that the heat losses are kept to the minimum.
3. Insert the immersion heater and the thermometer into the appropriate hole. Drop before some oil or water into the thermometer hole to ensure good thermal contact between the thermometer and the block.
4. Set up the circuitry according Fig.1.
5. Switch on the power supply and adjust it to give a current of about 4 A. Switch again the power supply off.
6. Before starting the experiment, wait for a few minutes. Then, take note of the initial temperature of the block.
7. Switch the power supply on and start measuring time with stopwatch.
8. Wait until the temperature has risen about 20°C.
9. Record the time and final temperature

The specific heat capacity can then be calculated using the relation:

$$I \cdot U \cdot t = m \cdot c \cdot (\theta_2 - \theta_1)$$

With  $I$ : Current,  $U$ : Voltage,  $t$ : Time,  $m$ : Mass of the calorimeter block,  $c$ : Specific heat capacity,  $\theta_1$ : Initial temperature,  $\theta_2$ : Final temperature

**Precautions:** There is a risk of burns from heater or calorimeter, therefore, allow the apparatus to cool before moving it.

## HEATING RESISTANCE

It is a especial heater designer for use with metal block calorimeters. The heating element is enclosed in a stainless steel tube with two flying leads coming out of the sealed tube body for electrical connections. It operates on 12V, 50W.