**BOMBA DE VACÍO POR RECIRCULACIÓN DE AGUA  
CIRCULATING WATER VACUUM PUMP  
POMPE À VIDE Á RECIRCULATION D'EAU**

**Nahita**

Ref. PXC001

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**Vacuum gauge needle does not move or jumps**

- If the needle does not move, this may be due to blockage of the vent hole at the rear of the vacuum gauge. Remove the vacuum gauge for cleaning.
- The needle jumps because the check valve or the vacuum tube is leaking. The check valve must be cleaned and the two sections of the vacuum tube must be tightened.

**Motor start-up anomaly**

- If the motor does not start, it may be due to a power failure or a blown fuse.
- After starting the motor, it beeps without starting. The power must be cut off immediately, then turn the ventilation wheel inside the engine cover a few times by hand.

**Attention**

- Do not allow the pump head to idle away from the water for long periods of time to prevent wear of the mechanical seal. The motor must be well grounded and the circulating water must be kept clean to prevent the vacuum pump from blocking.

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Spanish ..... 2-4  
 English ..... 5-7  
 French ..... 8-10

**INDEX OF CONTENTS**

1 Introduction.....5    4 Instructions of use .....6  
 2 Features .....5    5 Troubleshooting .....6  
 3 Specifications .....6

**1. INTRODUCTION**

This vacuum pump uses the water jet method to generate vacuum; a centrifugal pump propels the water contained in the tank at high speed through two vacuum tubes. It is a multi-purpose device which, in addition to providing vacuum, can be used to supply cooling water for applications such as evaporation, distillation, crystallisation, filtration, sublimation, etc. It is a very useful apparatus in laboratories of universities, pharmaceutical, food and other sectors.

**2. FEATURES**

- With this equipment it is possible to save a lot of water. The water contained in the tank is recycled due to water recirculation, avoiding the wasteful generation of vacuum using tap water and a vacuum tube. In addition, the same water can be used for a long time, which is especially important in laboratories where availability is low or pressure is insufficient. In addition to providing vacuum conditions, this device also allows the establishment of an external cooling water circulation circuit for various applications where cooling water is required. In these cases, to prevent water temperature rise, a continuous supply of tap water entering the tank may be necessary.
- Ease of use. The suction of the centrifugal pump is directly submerged in the water. The equipment is small in size and light in weight, making it easy to move from place to place. The tank lid can be opened to observe the water level and refill if necessary.
- Equipped with two vacuum ports, which can work separately or simultaneously. Each port has its own vacuum gauge to observe the vacuum pressure. A check valve is installed in the suction port to prevent circulating water from flowing into the vacuumed equipment when pump operation is stopped.

- Corrosion resistant and non-polluting. The internal parts of the centrifugal pump are made of stainless steel, the water tank of ABS plastic, materials that are not affected by corrosive substances such as acids and alkalis, also there is no contamination in the laboratory, such as that produced by oils.
- Stable and reliable operation, with low noise.

**3. SPECIFICATIONS**

<b>Vacuum level (gauge pressure)</b>	-0,09 MPa
<b>Vacuum flow rate (at each port)</b>	10 L/min
<b>Tank capacity</b>	15 L aprox.
<b>Centrifugal pump flow rate</b>	60 L/min
<b>Centrifugal pump head</b>	8 metros
<b>Alimentación eléctrica</b>	220V, 50 Hz
<b>Consumo</b>	280 W
<b>Nivel de ruido</b>	< 45 dB
<b>Dimensiones (largo x ancho x altura)</b>	395 x 290 x 450 mm
<b>Peso</b>	10 kg

**4. INSTRUCTIONS FOR USE**

- Preparation for work. Place the equipment in a horizontal position on the workbench. When using it for the first time, open the tank lid and add clean cold water (you can also add water through the hose connected to INLET). When the water level is about to reach the height of the overflow nozzle behind the water tank, stop adding water. Switch on the equipment; add more water if necessary.
- Connect the suction nozzle of the equipment requiring vacuum to one of the two suction nozzles of the vacuum pump, and check that the circulating water valve, located at the rear, is closed. Turn on the main switch to start suction. The vacuum pressure can be observed on the vacuum gauge corresponding to the suction nozzle used.
- When the pump works continuously for a long time, the temperature of the water in the tank increases, affecting the vacuum level. To avoid this effect, the hose from the water source (tap water) can be attached to the INLET connector, the overflow nozzle serves as an outlet to the drain, and the tap water flow rate can be properly controlled, so that the water temperature in the tank does not rise and the vacuum level remains stable.
- When cooling water supply is required, the water inlet and outlet hoses of the equipment to be cooled are connected to the OUTLET nozzle and the INLET nozzle respectively. If the water contained in the tank increases in temperature to such an extent that it affects the cooling, then tap water must be continuously added to the tank. Excess water will flow out of the overflow hole to the drain. For the cooling water supply the valve at the rear must be open.

**5. TROUBLESHOOTING**

- Low vacuum**
- There are foreign objects inside the vacuum tubes. Remove foreign objects.
  - There are foreign objects in the check valve. Remove foreign objects.
  - The water inlet of the centrifugal pump is blocked. Clean the water inlet.