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
Index

1. General Information	12
2. Packing list and accessories	12
3. Technical specification	12
4. Information about safety in the use of the equipment	13
4.1 Risks the operator is subjected to	13
4.2 Staff qualification	13
4.3 Use with hazardous substances	13
4.4 Maintenance	13
4.5 Warranty	13
5. Equipment Description	14
6. Installation	14
6.1 Bucket/Tank	14
6.2 Connection to the mains	14
7. Operation	14
7.1 Quick operation	16
7.2 Setting screen	17
7.3 Programming screen	17
7.4 Selecting a program for operation	18
7.5 Operation screen	18
7.6 Circulation pump	18
7.7 External probe	19
7.8 Communications	19
7.9 Alarms	19

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ENGLISH VERSION

1. General Information

- Handle the parcel with care.
- Check that contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list and accesories

Units are supplied with its owns manuals.

The immersion thermostat 3000544 needs some accesories to operate.

Stainless steel buckets:

Code	Liter	High (cm)	Wide (cm)	Length (cm)
6000390	9	20	34	28
6000391	12	20	37	35
6000392	20	20	55	35
6000393	27	25	55	35
6001093	45	23	70	56

Metacrilate buckets: (operate up to 60°C)

Code	Liter	High (cm)	Wide (cm)	Length (cm)
1000394	8	15	18	31
1000544	9.5	14.5	16.5	39
1000395	13	15	23	40
1000396	30	20	30	50

Extending support:

Code: 6001094

3. Technical specification

	3000544	3000550	3000551
Dimensions: High/wide/depth (mm)	280/180/190	360/550/350	410/550/350
Power supply voltage:	230 VAC 50Hz		
Electrical power (W):	2060		
Weigth: (Kg)	4.5	10	12
Fuse:	Replaceable (Type: 10x5mm)		
Temperature range:	From... ambient T ^a + 5°C ... to ... 200°C		
Temperature resolution:	0.1°C		
Timed operation:	Hours, minutes		
Pump flow	10L/min		

4. Information about safety in the use of the equipment

**Hazard
risk**

**Electric
risk**

**Risk of high
temperature contact**

4.1 Risks the operator is subjected to

- Risk of burns
- Electrical risk

4.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

4.3 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

ATTENTION!

Burning risk

Take proper precautions
(use of thermal gloves, clamps,...)
to extract samples from inside the equipment.

It can be at a high temperature.

4.4 Maintenance

Unplug the equipment from the mains before you handle inside it.

**The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage.
Take it to one of J.P.SELECTA's authorized technical services.**




It is recommended to clean the accessible parts of the equipment with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

4.5 Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

5. Equipment Description

Immersion thermostats are widely used on laboratories for thermal treatment of samples, or heating/cooling other equipments.

The equipment includes a circulation pump that can be used to improve the homogeneity within a bath or to perform a closed liquid circulation circuit for cooling or heating.

The unit could work in a continuous or timed operation.

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

6.1 Bucket/Tank

A bucket or tank will be needed to operate the thermostat. See J.P. SELECTA, s.a.u. product range or paragraph 2 of this manual.

Place the equipment on a flat level stable surface, suitable for the equipment weight.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

7. Operation

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.






When starting the computer, the starting screen will be shown (Fig 7.1).

- **SETTINGS:**
 - Time and date: menu to set the date and time.
 - Temperature correction.
 - Temperature alarm limits.
- **PROGRAMMING**
 - Allows you to program the working temperature & time.
- **TASK**
 - Allows you to select and launch the operating temperature, the working time and the pump flow.

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
Fig 7.1 Initial screen


SETTINGS SCREEN**INITIAL SCREEN****PROGRAMMING SCREEN****OPERATION SCREEN**


7.1 Quick operation


1 Press over «TASK».

2 Press over temperature display
(Frame changes to yellow color)

3 Set the temperature with  and 
(SET 40°C on the picture)

4 Press over .

Notes:

In a similar way, hours and minutes of working time could be programmed and also the pump flow rate.

The temperature set from this screen does not remain recorded once equipments are switched off.
(See paragraph 7.2)


Fig 7.2 Setting screen

7.2 Setting screen

7.2.1 Set time and date

Press on Date/time button. Press successively on the box to select the day, month, year, hour, minute and second settings.

Use Up and Down icons to entry a value.

7.2.2 Set temperature, language and screen contrast

In a similar way, set the temperature, language and screen contrast by pressing the icons.

7.2.3 Working parameters protected under password

There are several parameters under password "0000" which affect the equipment operation.

. Press "Insert password" and follow the instructions:

password: 0000

. Press CAL to access these parameters:

log interval:	1 to 60min
external sensor:	yes/no
calibration:	-100 to +100
lower temp. limit:	-10
upper temp. limit:	-10

a) Log interval

It is the frequency of data collection for the logging.

b) External sensor

Select the value **YES** if you want to work with an external probe connected to the connector on the back side of the equipment

c) Calibration (OFFSET)

Use only this setting if you notice a deviation between the temperature indicated by the equipment and the one indicated by a standard thermometer.

For example: the equipment indicates 40.0°C, the pattern indicates 40.5°C. Enter a value of -0.5°C to correct the deviation.

d) The lower and the upper limit of temperature

These are the maximum and the minimum limits of temperature adjustments. This allows you to adjust the working temperature range of the equipment.

e) To update the firmware of the device, click on "Info" and the UPGRD button will appear at the bottom of the screen.

Contact the technical service to obtain the updated file.

7.3 Programming screen

Up to 10 programmes of temperature and time can be stored.

Each program has the following parameters:

Temperature:	from ambient +5°C ... 200°C
Time:	HH:MM
Undertemp:	Alarm value for temperature drop.
Overtemp:	Alarm value for overtemperature.

From main menu, press over «PROGRAMMING».


Fig 7.3. Programming screen


Fig 7.4. Selecting a program


Fig 7.5. Operation screen


Fig 7.6.1 INTERNAL circulation.


Fig 7.6.2 EXTERNAL circulation.

Fig 7.3 screen is displayed.

Select a program (from 00 to 09) with keys:

Set working temperature and the other parameters:

Press: TEMP: 16.0°C

Keys will be shown to select temperature.

In a similar way press:

TIME: 00:00

Sets de working time.
A time setting value of 00.00 means «continuous operation».

UNDERTMP: 0°C

Sets the undertemperature alarm value.

OVERTMP: 0°C

Sets the overtemperature alarm value.

Press **SAVE** to save .

7.4 Selecting a program for operation

To select a stored program press over SEL (Fig 7.4).
(SEL indication changes its color to yellow).

Once selected, this program in placed on the TASK screen.

7.5 Operation screen

On the operation screen «TASK», the program selected (PROG 04) on the programming screen is shown (Fig 7.5) and it starts when pressing over .

This screen shows:

Actual temperature
Temperature set
Working time set (SET)
Elapsed time

7.6 Circulation pump

The liquid circulation pump could operate in two configuration modes:

- INTERNAL circulation:

Water (or other) will be pumped from horizontal pump output (See fig 7.6.1).

Use this configuration to made water circulation into the bath vessel in order to homogenize the temperature from sample load.

- EXTERNAL circulation:

Water (or other) will be pumped from vertical pump output (See fig 7.6.2).

Use this configuration to made water circulation out of the bath vessel in order to temperate an external equipment. This is the factory setting.

7.7 External probe

The process can be controlled by means of an external probe.

Connect the probe to the rear connector and place it in the media you want to control the temperature. To adjust it, follow the steps in section 7.2.3.

7.8 Communications

7.8.1 USB

Access the SETTINGS > LOGS DOWNLOAD screen to download the logs of the processes carried out. Then follow the instructions that appear on the screen.

7.8.2 RS-232

RS232 port configuration. RS232 communication protocol based on technical specification SPT-0004-B.

Transmission speed	56.000bps
Data	8 bits
Stop	1 bit
Parity	No parity

7.9 Alarms

Error	Screen	Description	Type
01	Level	Safety or level thermostat	A
02	Overtemp	High temperature level	A
03	Undertemp	Low temperature level	A
04	External RTD opened	Open external probe	A
05	External RTD shortcircuited	Shortcircuited external probe	A
06	Internal RTD opened	Open internal probe	A
07	Internal RTD shortcircuited	Shortcircuited internal probe	A
08	Power Fail	Without net	B

J.P. SELECTA, s.a.u. 2017