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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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8- CLEANING

- ◆ The equipment has no lateral ventilation grilles and front panel is plane so, all surfaces are easy to clean. All surfaces can be cleaned with soft and standard detergents used in the laboratory.
- ◆ Before cleaning, always disconnect the equipment from the net and let it cool.

Attention!: Never use organic solves as xylol for cleaning the equipment.

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At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

- ◆ Put the slices obtained with the microtome on the surface of the water of the bath. The slices, at the beginning joined by the ends, will be extended softly thus removing the wrinkles and folds produced during slicing; slices will finally stay floating on the water.
- ◆ Then, by using a microscope slide, take the slices. This way it is possible to have one or several slices, obtained from the same paraffin block, in one microscope slide.
- ◆ Remove the excess of water, let dry and if necessary use an oven at a temperature similar to that of the bath to evaporate completely the water and add the sample to the slide, at least for 48 hours. Finally stain the sample for a later observation with the microscope.
- ◆ Once operation with the bath has finished press ON/OFF (stop). The equipment will stop heating and by indicating “_XX”, shows the cooling temperature at any moment.
- ◆ Turn the equipment off by pressing general switch.

Note: Never turn the equipment on with the tank empty.

7- SAFETY RECOMMENDATIONS

- ◆ Read carefully the instruction manual before using the equipment.
- ◆ During operation, the equipment should be placed on a horizontal, plane and stable surface and correctly ventilated. Keep it away from heat sources.
- ◆ Before turning the equipment on make sure it is enough water in the tank.
- ◆ While the instrument is working, the surfaces can be very hot. Danger of burns!
- ◆ Use only non-flammable liquids, preferably distilled water.
- ◆ Never leave or store flammable substances near the equipment.
- ◆ Though it is properly isolated, avoid water get inside the equipment.
- ◆ Make sure the equipment is well grounded and power supplies meets with parameter specified.
- ◆ Press gently the keys of the front panel; never use sharp objects to press them since front panel could be damaged.

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


5- CONTROL PANEL**6- OPERATION MODE**

Fill the tank of the bath up to 2/3 its maximum capacity and proceed as follows:

- ◆ Turn the equipment on by pressing the general On/Off switch. Display will show the room temperature by indicating “_XX” (e.g.: “_22”)
- ◆ Set the working temperature:
 - Press Set.; display will show “PXX” (e.g.: “P22”)
 - By pressing ▲ and ▼ select the desire working temperature (the temperature of water must be always lower than melting temperature of paraffin).

Note: The selected nominal value is kept in the memory of the instrument after turning it off.

- Press ON/OFF. The system starts working and display, by indicating “CXX”, shows the temperature in the tank at any moment (e.g.: “C35”).

Note: For visualizing the set temperature value at any moment, as well as for setting a new value it is necessary to press “ON/OFF” (stop) previously and proceed according to the above steps.

- ◆ Wait until the system reaches the selected working temperature and stabilizes completely.

TAKE ALL NECESSARY PRECAUTIONS TO AVOID ANY ACCIDENT DUE TO A CONTACT WITH HOT WATER OR SURFACES

2- FUNCTIONS AND FEATURES


- ◆ The PTC heating material used allows a fast and uniform heating with low inertia for energy saving.
- ◆ The high quality temperature survey integration block has been used to make it precise in temperature survey and credible in performance.
- ◆ The microprocessor control and LED digital display allows an easy and intuitive working temperature adjustment.
- ◆ With function of memorization to keep the set temperature automatically after startup.

3- MAIN TECHNICAL SPECIFICATIONS

Capacity	~2L
Tank volumen	240x180x54 mm
Temperature range	RT up to 90°C
Range for constant temp. control.	-3°C a +5°C
Power supply	AC 220±10% (50Hz)
Consumption	550 w
Fuse	6A
Dimensions	340x350x105 mm

4- PARTS OF THE EQUIPMENT

1. Tank
2. Power On/Off switch
3. Control panel

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Spanish..... 2-6
 English..... 7-11

INDEX OF CONTENTS

1. Introduction 7
 2. Functions and features 8
 3. Main technical specifications 8
 4. Parts of the equipment 8
 5. Control panel 9
 6. Operation mode 9
 7. Safety recommendations 10
 8. Cleaning 11

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

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1- INTRODUCTION

The bath for paraffin sections is one of the necessary and main used equipment's in Pathological Anatomy laboratories and for tissue processing in hospitals, universities and research laboratories. The advanced microprocessor control and heating element used, together with the compact, light and stable design make of this equipment an advance, reliable and easy-to-use instrument.