

UVStudio

For UV/Visible Spectrophotometers

Instruction Manual

Contents

Preface	1
Introduction	1
Installation	1
PC requirements	1
Installation of UVStudio software	1
Installation of USB driver	4
Connecting the spectrophotometer to the computer	5
Run UVStudio software	6
Start running	6
Basic Operation	7
Enter an application module	7
Back to main interface	9
Return to previous interface	9
Do blank/Scan baseline	9
Measuring sample	9
Enter data /curve list	9
Measurement Results Operation	10
Calibration and System Settings	10
Calibration	11
Bandwidth selection	12
Display resolution setting	12
Light source management	12
Language selection	13
Printer setting	14
General settings	14
Measurement	15
Multi-wavelength	15
Kinetics	17
Time scanning	19
Quantitation	21
Nucleic Acid/Protein	28
Spectrum	31
File management	33
Formula editor	34
Performance verification	35
Establishing performance validation methods	36
Verify performance	42
Copyrights and Warranties	43

Preface

Introduction

Controlling the spectrophotometers from a computer running UVStudio software offers additional data display, manipulation and storage capabilities beyond those available in the firmware and build-in software.

UVStudio software includes fixed wavelength (single or multiple wavelengths), quantitative analysis, spectral scanning, time scanning, and biological methods (DNA/RNA/protein). The software also features a comprehensive validation application that enables you to verify the continued performance of the instrument.

Installation

PC requirements

Hardware

- | | |
|------------------------|--|
| - Processor | Intel or AMD at 1.0GHz or above |
| - RAM | 4GB or above |
| - Hard disk | 100GB free space |
| - Display | 1280x800 or higher, color |
| - USB ports | USB 2.0 or USB 3.0x2 |
| - Input/Output Devices | Keyboard/Mouse or Touch Screen |
| - Operating System | Microsoft Windows 8/10/11 Pro (32-bit or 64-Bit) |

Operating System

- | | |
|--------------------------------|-------------------------------|
| - Operating System | Microsoft Windows 8/10/11 Pro |
| - Automatic standby/sleep | Disable |
| - Hard disk automatic shutdown | Disable |
| - Local computer time change | Disable |

PC security


- | | |
|----------------------|-----------------------------|
| - Automatic Updates | Disable |
| - Firewall | On |
| - Antivirus software | Installation is not allowed |

Installation of UVStudio software


1. Insert the UVStudio software installation USB disk into the USB port of your computer. Open the USB disk and double-click "Setup.exe" to start the installation. Click "Next";




2. Enter the "User Name" and "Company Information" and click "Next";




3. Click "Change..." Select the installation directory and click "Next";




4. Select the shortcut folder and the range of users to which the shortcut applies, and click "Next";




5. Display the installation information that has been set or selected and click "Next";



6. The installation program copies the required files and database to the installation directory, and makes the corresponding settings, and displays the "Installation Success" interface after completion. Click "Finish" to end the installation of UVStudio software.






Installation of USB driver



1. Open the "CP210x_VCP_Win7_8_10" in the USB flash drive, select the appropriate installation file according to the different operating systems, and double-click to start the installation. Click "Next";
 - CP210xVCPInstaller_x64.exe : Windows 8/10/11 64-bit
 - CP210xVCPInstaller_x86.exe : Windows 8/10/11 32-bit



2. The software license agreement is displayed, select "I accept this agreement" and click "Next";




3. The installation program copies the required files to the installation directory, and displays the "Installation Success" interface after completion. Click "Finish" to end the installation of USB Driver.




Connecting the spectrophotometer to the computer

1. Connect the USB cable to the USB-A port on the computer and USB-B port on the back of the instrument. Connect the

USB key to the USB-A port on the computer;



2. Press the start button of the instrument, wait for the instrument startup to complete, and enter the "Device Manager" → "Ports (COM and LPT)" to check the port number on the PC;



Run UVStudio software

Start running

1. Power on the spectrophotometer.
2. Open the sample room lid, check and remove anything in the light path that might block the light, and close the sample room lid.
3. Double-click the "UVStudio" icon on the desktop to run the software, initialize and system calibration in turn, and then enter the main interface after completion.

Initialization	
System	Passed
Light Source Initialization	Passed
Grating Initialization	Passed
Adjustable Slit Initialization	Passed
Automatic Sample Holder Initialization	Not Configured
Filter Wheel Initialization	Passed
Dark Current Correction	Passed
Wavelength Calibration	Passed
Energy Checking	Passed
System Baseline Checking	Passed
Warm-up	Remaining 14 Minute

Retry Skip

Note

If the pop-up window "Do you want run offline versions?" , it may be


- 1. The instrument has not been started or the USB communication cable is not connected.*
- 2. The default port number is not the currently connected device.*

After checking and confirming, click the "x" button to select the port number of the connected device and click "OK" to connect.





Basic Operation

Enter an application module

Main interface, press an application icon to enter the module interface.




	<p>Photometry</p> <p>Measure the photometric value of a sample at a single wavelength.</p>
	<p>Multi-wavelength</p> <p>Measure the photometric value of a sample at multiple wavelengths.</p>
	<p>Kinetics</p> <p>Measure the change in absorbance or absorbance change rate over time at a specified wavelength.</p>
	<p>Time Scanning</p> <p>Measure the change of photometric value with time at a single wavelength.</p>
	<p>Quantitation</p> <p>Establish a standard curve and measure the concentration of the sample using a standard curve.</p>
	<p>Nucleic Acid/Protein</p> <p>Measure DNA, RNA and protein concentrations using built-in methods or new methods.</p>
	<p>Spectrum</p> <p>Measure the photometric curve of a sample over a range of spectra.</p>
	<p>Account Management (only for UVStudio Audit Trail edition)</p> <p>Assign user rights, operate log management.</p>

	<p>Performance Verification</p> <p>Verify the technical performance of the instrument.</p>
	<p>File Management</p> <p>Efficient management of user files, browsing, copying, renaming, and deleting operations.</p>
	<p>Formula Editor</p> <p>Users can add calculation formula for special applications or exploratory studies according to their needs.</p>
	<p>System</p> <p>Set system parameters and system calibration.</p>

Back to main interface

Measurement interface, press  back to main interface.


Return to previous interface

Method/Settings/Data list/Curve list interface, press  back to previous interface.

Do blank/Scan baseline

Measurement interface,


- Split beam models: Put the "Reference" in the measurement channel.
- Double beam models: Put the "Reference" in the reference channel and in the measurement channel.

Press  to do blank/scan baseline.

Measuring the sample

Measurement interface,








- Split beam models: Put the "Sample" in the measurement channel.
- Double beam models: Put the "Reference" in the reference channel and "Sample" in the measurement channel.

Press  to measure the sample.

Enter data/curve list


Measurement interface, press  to enter data/curve list.


Measurement Results Operation

	Open Open a stored file, load data or parameters.
	Save Save data, parameters to storage.
	Print Print test report.
	Export to MS Word Export file to MS Word format.
	Export to MS Excel file format Export file to MS Excel format.
	Export to PDF file format Export file to PDF format.
	Delete Delete the selected results.

Calibration and System Settings




Select the icon  in the main interface. Display options to calibrate the system and configure the basic instrument settings.



Calibration

Select the tab **Calibration** in the System interface. Remove something from the measurement channel, close the sample chamber cover, and select the item **Wavelength calibration**, **System baseline calibration** or **Dark current calibration** to do calibration.

Important Before performing the calibration, you must remove something from the measurement channel, close the sample chamber cover, and maintain this state throughout the calibration process.



Bandwidth selection

Select the tab **Bandwidth** in the System interface. Select the bandwidth which you want.

Important After changed the bandwidth, you must do a calibration (wavelength calibration, system baseline calibration and dark current calibration).

The screenshot shows a dialog box titled "Spectral Bandwidth" with a teal header. It contains a list of numerical options: 0.5, 1.0, 2.0, 4.0, and 5.0, each followed by a horizontal line for selection. At the bottom, there is a teal bar with a white "X" icon in the center.

Display resolution setting


Press the tab **Display resolution**. Select the resolution of the required display digits according to different measurement modes.

The screenshot shows a dialog box titled "Display Resolution" with a teal header. It contains a table with measurement modes and their corresponding resolutions. At the bottom, there is a teal bar with a white "X" icon on the left and a white checkmark icon on the right.

Measurement Mode	Resolution
Abs	0.000
%T	0.00
Conc.	0.0
r ²	0.0000

Light source management

Press the tab **Light source**, which displays the usage timing of the light source and the switching wavelength of the light source.




- **Turn on/off lamp**

Click  to turn on/off Tungsten lamp/Deuterium lamp.

Note


When using a single light source for a long time, you can turn off another unused light source to save light source life and energy.

- **Set the light source switching wavelength**

Press the current switching wavelength, pop up the numeric keyboard, input the switching wavelength, and press the button  to complete the setting.

Note

After modifying the light source switching wavelength, you need to calibrate system baseline.


Click  to reset the Tungsten lamp/Deuterium lamp timing.

Note

To ensure the accuracy of the light source timing, perform the corresponding light source clear after replacing the light source. Do not operate at other times.


Language selection

Press the tab **Language** and select the desired language.




Printer setting

Press the tab **Printer** and select the installed printer.



General settings


Click the tab **General** to enter.



·**Warm up**

Click "Warm up After Start-up (min)", pop up numeric keyboard, input value 15~120 minutes.

·**Calibrate system baseline when startup**

Click the icon  to turn on/off the "Calibrate System Baseline After Start-up" option.

·**Set sample holder type**


If you change to an auto sample holder or micro volume test base, you will need to set the sample holder to the matching type.



Measurement


Multi-wavelength

Multi-wavelength mode is used to measure the absorbance or transmissivity of the sample at multiple wavelengths.




Main interface, click the icon  to start a **Multi-wavelength** application.



- 1 **Multi-wavelength measurement** interface, click the button  to set the measurement parameters. The method can be saved or directly called from the memory. Click  to accept the new parameters and return to the measurement interface.



Photometric mode	2 photometric modes: absorbance and transmittance.
Measurement Wavelength Number	1~20 wavelengths are available, wavelength range: 190~1100nm.
Cycles	1, 2, 3, 5, 10, 20, 30, 50 times can be selected, the instrument will calculate the average as the final output.
Calculation	The instrument can choose to embed the calculation formula to directly calculate the result. The formula can support user customization (in the custom method module).
Path Length	The width of the cuvette used for the measurement.
Equivalence Transformation	Activation automatically converts measurements from different path cuvettes to values in the 10mm light path length.

- 2 Put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to do blank.
- 3 Put the sample in the measurement channel, close the sample chamber cover, and click the button  to measure and calculation results.
- 4 Repeat step 3 to measure more samples.
- 5 Click  to switch to list mode to browse the list of measurement results.


No.	Calculated	Measured	Date & Time	Memo	Select
0001	0.159	A(635.0)=1.118 A(465.0)=...	2024/7/23 14:48:34		<input checked="" type="checkbox"/>
0002	0.107	A(635.0)=0.763 A(465.0)=...	2024/7/23 14:48:53		<input checked="" type="checkbox"/>
0003	0.119	A(635.0)=0.609 A(465.0)=...	2024/7/23 14:49:10		<input checked="" type="checkbox"/>

Kinetics

Kinetics mode is used to measure the rate of change of the sample.





Main interface, click the icon to start a **Kinetics** application.



- 1 **Kinetics measurement** interface, click the button to set the measurement parameters. The method can be saved or directly called from the memory. Click the button to accept the new parameters and return to the measurement interface.

Measurement Mode	Kinetic rate >	Path Length (mm)	10 >
Wavelength (nm)	546.0 >	Equivalence Transformation	<input type="checkbox"/> OFF
Background Correction	<input type="checkbox"/> OFF	Coordinate Auto Scaling	<input checked="" type="checkbox"/> ON
Correction Wavelength (nm)	>	Y Maximum	>
Total Time (s)	120 >	Y Minimum	>
Sampling Interval (s)	1 >	Auto-printing	<input type="checkbox"/> OFF
Differential Interval(s)	1 >	Autosave	<input type="checkbox"/> OFF
Criteria for Judging Linearity (%)	10 >		

Measurement mode	2 measurement modes: kinetics and kinetic rate.
Wavelength	Measured wavelength, range: 190~1100nm.
Background correction	Background correction switch, can be set according to actual needs.
Correction wavelength	Background corrected wavelength, range: 190~1100nm.
Total Time	Measure the total sampling time required.
Sampling Interval	Sampling interval.
Delay/Derivative	Waiting time before starting sampling/time of participating activity calculation.
Coefficient	Coefficient of activity calculation equation.
Path Length	The width of the cuvette used for the measurement.
Equivalence Transformation	Activation automatically converts measurements from different path cuvettes to values in the 10mm light path length.
Coordinate Auto Scaling	Whether to automatically adjust the coordinates based on the data.
Y Maximum	The maximum value of the ordinate (valid only when the coordinates are fixed).
Y minimum	Minimum value of the ordinate (valid only when the coordinates are fixed).
Auto-printing	Automatically print curves and results after measurement is complete.
Autosave	Automatically save curves and results after measurement is complete.



- 2 Put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to do blank.
- 3 Put the sample in the measurement channel, close the sample chamber cover, and click the button  to measure and get sampled data and draw the curve.
- 4 Repeat step 3 to measure more samples.

Time scanning



Time scanning mode is used to measure photometric value changes of the sample.

Main interface, click the icon  to start a **Time scanning** application.



- 1 **Time scanning measurement** interface, click the button  to set the measurement parameters. The method can be saved or directly called from the memory. Click  to accept the new parameters and return to the measurement interface.

Photometric mode	2 photometric modes: absorbance and transmittance.
Wavelength	Measurement wavelength, range: 190~1100nm.
Total	Measure the total sampling time required.
Interval	Sampling interval.
Path Length	The width of the cuvette used for the measurement.
Equivalence Transformation	Activation automatically converts measurements from different path cuvettes to values in the 10mm light path length.
Coordinate Auto Scaling	Whether to automatically adjust the coordinates based on the data.
Y Maximum	The maximum value of the ordinate (valid only when the coordinates are fixed).
Y minimum	Minimum value of the ordinate (valid only when the coordinates are fixed).
Auto-printing	Automatically print curves and results after measurement is complete.
Autosave	Automatically save curves and results after measurement is complete.

- 2 Put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to do blank.
- 3 Put the sample in the measurement channel, close the sample chamber cover, and click the button  to measure and get sampled data and draw the curve.
- 4 Repeat step 3 to measure more samples.

Quantitation

Quantitation mode is used to measure sample concentration by establishing and using a standard curve.

Main interface, click the icon  to start a **Quantitation** application.

ppm

5.20

A1	546.0								
	0.5229								

Calculation	0.5229	Date & Time	2024/7/23 15:27:01	Memo	
-------------	--------	-------------	--------------------	------	--

Operations Method Blank Calibration slot difference Measure

Establish quantitative methods

- 1 **Quantitation measurement** interface, click the button to enter into **Method** interface.

Quantitation - Method

12.02

ppm

-0.04

0.0000

Abs

1.2041

Standard Curve Equation: $C = 10.0203 \cdot A - 0.0412$ $R^2 = 0.998661$

Standard samples

Conc.	Abs	Conc.	Abs	Conc.	Abs
0.00	0.0000				
5.00	0.5230				
7.00	0.6856				
10.00	1.0034				

Measure method: Single Wavelength

Wavelength (nm): 546.0

Cycles: 1

Equivalence Transformator: Closed

Path Length (mm): 10

Equation Form: $C = f(\text{Abs})$

Fitting: Linear

Zero Intercept: Closed

Calibration method: Measuring Standard Samples

Number of Standard Samples: 4

Unit: ppm

Threshold:

New Method Open Method Save Method Export Method Print Method OK

- 2 **Method** interface, click the button to start a new measurement method and enter the parameter setting interface.

Quantitation - New Method ↩


Measurement Method					Single Wavelength A=A1 >	Equation Form	Abs=f(C) >
λ1	λ2	λ3	λ4	λ5			
546.0	650.0	-----	-----	-----		Fitting	Linear >
λ6	λ7	λ8	λ9	λ10		Zero Intercept	<input type="checkbox"/> OFF
-----	-----	-----	-----	-----		Calibration method	Measuring Standard Samples >
λ11	λ12	λ13	λ14	λ15		Number of Standard Samples	4 >
-----	-----	-----	-----	-----		Unit	ppm >
λ16	λ17	λ18	λ19	λ20		Threshold	>
-----	-----	-----	-----	-----		Cell Deviation Correction	<input type="checkbox"/> OFF
Cycles					1 >		
Path Length (mm)					10 >		
Equivalence Transformation					<input type="checkbox"/> OFF		


Open Parameters
Save Parameters
Extract Parameters
Next →

Measurement Method	Built-in single wavelength, dual wavelength difference, dual wavelength ratio, three wavelengths, area 5 ways, and support custom formula.
Wavelength	Measuring wavelength, range: 190~1100nm.
Cycles	1, 2, 3, 5, 10, 20, 30, 50 times can be selected, the instrument will calculate the average as the final output.
Path Length	The width of the cuvette used for the measurement.
Equivalence Transformation	Activation automatically converts measurements from different path cuvettes to values in the 10mm light path length.
Equation Form	Equation form: $C=F(Abs)$ and $Abs=F(C)$.
Fitting	Three ways of fitting are provided: Linear, Quadratic, Cubic.
Zero Intercept	Turning on the representative fitting curve will directly cross the zero point, and closing will represent the fitting curve but the zero point.
Calibration Method	Three ways to generate a standard curve: inputting equation coefficient, measuring standard samples and inputting standard sample values.
Unit	Built-in 19 commonly used concentration units: -, %, ppm, ppb, g/l, mg/l, µg/l, ng/l, g/dl, mg/dl, µg/dl, mg/ml, µg/ml, ng/ml, µg/µl, ng/µl, mol/l, mmol/l, IU, and support input custom units.
Number of standard samples	The number of standard samples can be selected (only valid for standard sample calibration and standard sample input), quantity: 2~20.
Threshold	Upper and lower limits of measurement results.

3 Establish standard curve

3.1 Establish standard curve by inputting the equation coefficient

- 1) Set **Calibration** to **Input Equation Coefficient**, set other measurement parameters according to measurement requirements, and click the button  to start.
- 2) **Input Equation Coefficient** interface, click the coefficient K0–Kn to pop up the keyboard and enter the coefficient.



Quantitation - New Method - Inputting Equation Coefficients

Standard Curve Equation

$$C = K1 \cdot A + K0$$

Coefficient K0 0.012 >

Coefficient K1 3.44 >

Coefficient K2 >

Coefficient K3 >

Standard Curve

340.57


ppm

0.00

0.0000 Abs 100.0000


$$C = 3.4400 \cdot A + 0.0120$$




Previous Blank Correct Cell Deviation Measure Finish


- 3) After the method is completed, the standard curve and related information are displayed. Click the button  to complete and return to **Method** interface.

3.2 Establish standard curve by measuring standard samples

- 1) Set **Calibration** to **Measure standard samples**, set other measurement parameters according to measurement requirements, and click the button  to start.





- 2) **Measure standard sample** interface, put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to do blank.
- 3) Put the standard sample in the measurement channel, close the sample chamber cover, and click the button  to measure.
- 4) Repeat step 3 to measure all standard samples.
- 5) Click the concentration cell to enter the corresponding concentration value.
- 6) After the method is completed, the standard curve and related information are displayed. Click the button  to complete and return to **Method** interface.




3.3 Establish standard curve by inputting standard samples






- 1) Set **Calibration** to **Input standard samples**, set other measurement parameters according to measurement requirements, and click the button  to start.



- 2) Click the **Abs** cell to input Abs value and click **Concentration** cell to input the corresponding concentration value.
- 3) After the method is completed, the standard curve and related information are displayed. Click the button  to complete and return to **Method** interface.



Measure sample


- 1 **Quantitation measurement** interface, click the button  to enter into **Method** interface.
- 2 **Method** interface, load or new a method, click the button  to accept and return to measurement interface.
- 3 Put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to do blank.
- 4 Put the sample in the measurement channel, close the sample chamber cover, and click the button  to measure and calculation results.
- 5 Repeat step 4 to measure more samples.
- 6 Click  to switch to list mode to browse the list of measurement results.

Quantitation - Results

No.	Memo	Abs	ppm	Date & Time	Select
1		0.5229	5.20	2024/7/23 15:27:01	✓
2		0.5229	5.20	2024/7/23 15:27:26	✓
3		0.6856	6.83	2024/7/23 15:27:28	✓
4		0.6856	6.83	2024/7/23 15:27:29	✓
5		0.6856	6.83	2024/7/23 15:27:31	✓
6		1.0038	10.02	2024/7/23 15:27:33	✓
7		1.0038	10.02	2024/7/23 15:27:34	✓
8		1.0039	10.02	2024/7/23 15:27:35	✓
9		1.0038	10.02	2024/7/23 15:27:36	✓
10		1.0038	10.02	2024/7/23 15:27:37	✓

Nucleic Acid/Protein

Nucleic Acid/Protein measurement mode is used to measure DNA, RNA and protein concentrations using built-in methods or new methods.

Main interface, click the icon  to start a **Nucleic Acid/Protein** measurement application.


Nucleic Acid/Protein


No.	Calculation
0001/0001	C(DNA)=A1*50


197.35


A1 260.0	A2 280.0	A3 230.0	A1/A2	A1/A3					
3.9469	3.8297	4.0000	1.03	0.99					
Date & Time	2024/7/23 15:32:13		Memo						

Operations
Method
Blank
Measure





1 **Nucleic Acid/Protein** measurement interface, click the button  to enter into **Method** interface.



2 **Method** interface, click the button  to start a new measurement method and enter the parameter setting interface.



Measurement Method	Built-in Nucleic Acid (dsDNA, ssDNA, RNA), Protein A-280 (1Abs=1mg/mL, BSA, IgG, Lysozyme), DNA Quantitation 1 (260/280), DNA Quantitation 2 (260/230), Protein-Lowry, Protein-BCA, Protein-CBB, Protein-Biuret measurement method.
Wavelength	Measure wavelength, range: 190~1100nm.
Samples	dsDNA, ssDNA, RNA, Protein, etc.
Cycles	1, 2, 3, 5, 10, 20, 30, 50 times can be selected, the instrument will calculate the average as the final output.
Path Length	The width of the cuvette used for the measurement.
Equivalence Transformation	Activation automatically converts measurements from different path cuvettes to values in the 10mm light path length.
Equation Form	Equation form: $C=F(\text{Abs})$ and $\text{Abs}=F(C)$.
Fitting	Three ways of fitting are provided: first order, second order, third order.
Zero Intercept	Turning on the representative fitting curve will directly cross the zero point, and closing will represent the fitting curve but the zero point.
Calibration Method	Three ways to generate a standard curve: inputting equation coefficient, measuring standard samples and inputting standard sample values.
Unit	Built-in 19 commonly used concentration units: -, %, ppm, ppb, g/l, mg/l, $\mu\text{g/l}$, ng/l, g/dl, mg/dl, $\mu\text{g/dl}$, mg/ml, $\mu\text{g/ml}$, ng/ml, $\mu\text{g}/\mu\text{l}$, ng/ μl , mol/l, mmol/l, IU, and support input custom units.
Threshold	Upper and lower limits of measurement results.


- 3 **Method** interface, load or new a method, click the button  to accept and return to measurement interface.
- 4 Put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to do blank.
- 5 Put the sample in the measurement channel, close the sample chamber cover, and click the button  to measure and calculation results.
- 6 Repeat step 5 to measure more samples.
- 7 Click  to switch to list mode to browse the list of measurement results.


Nucleic Acid/Protein - Results							
No.	µg/ml	A1/A3	A1/A2	Detail data	Date & Time	Memo	Select
1	25.52	0.69	1.65	A1(260.0nm)=0.5104 A2(280.0nm)=0.3099 ...	2024/7/23 15:33:...		✓
2	31.67	0.74	1.47	A1(260.0nm)=0.6334 A2(280.0nm)=0.4322 ...	2024/7/23 15:33:...		✓
3	70.00	0.86	1.17	A1(260.0nm)=1.4000 A2(280.0nm)=1.2000 ...	2024/7/23 15:33:...		✓
4	78.49	0.87	1.15	A1(260.0nm)=1.5698 A2(280.0nm)=1.3696 ...	2024/7/23 15:33:...		✓

Spectrum

Spectrum mode is used to obtain the photometric curve of a sample over a range of wavelengths.

Main interface, click the icon  to start a **Spectrum** application.





1 **Spectrum scanning** interface, click the button  to enter into **Method** interface.

Spectrum - Method ↩



Photometric Mode	Abs >	Equivalence Transformation	<input checked="" type="checkbox"/>
Start Wavelength (nm) (≤ 1100.0)	600.0 >	Coordinate Auto Scaling	<input checked="" type="checkbox"/>
Stop Wavelength (nm) (≥ 190.0)	500.0 >	Y Maximum	1. >
Scan Step (nm)	1.0 >	Y Minimum	0. >
Scan Speed	Medium >	Auto-printing	<input checked="" type="checkbox"/>
Cycles (1~99)	1 >	Autosave	<input checked="" type="checkbox"/>
Cycle Interval	10 >		
Path Length (mm)	10 >		

New Method
Open Method
Save Method
Finish

Photometric mode	2 photometric mode: Abs and %T.
Start wavelength	Scan start wavelength, range: 190~1100nm.
Stop wavelength	Scan stop wavelength, range: 190~1100nm.
Scan Step	7 wavelength intervals selectable 0.1, 0.2, 0.5, 1.0, 2.0, 5.0, 10.0nm.
Scan Speed	3 scan speed selectable: Quick, Medium, Slow.
Cycles	Scan times.
Cycle Interval	Interval between 2 scans.
Path Length	The width of the cuvette used for the measurement.
Equivalence Transformation	Activation automatically converts measurements from different path cuvettes to values in the 10mm light path length.
Coordinate Auto Scaling	Whether to automatically adjust the coordinates based on the data.
Y Maximum	The maximum value of the ordinate (valid only when the coordinates are fixed).
Y minimum	Minimum value of the ordinate (valid only when the coordinates are fixed).
Auto-printing	Automatically print curves and results after measurement is complete.
Autosave	Automatically save curves and results after measurement is complete.

- Put the reference in the measurement channel and also in the reference channel in double beam models, close the sample chamber cover, and click the button  to scan baseline.
- Put the sample in the measurement channel, close the sample chamber cover, and click the button  to scan


sample and draw the curve.

- 4 After the scan is complete, click the button  to zoom the graph as needed.
- 5 Click the button  to retrieve the value of each point (peak) of the curve and mark it for a specific point.

File management

File management is used to manage saved parameters, methods, and measurement files.

Main interface, click the icon  to start File management.



The screenshot shows the 'File Management' window with a table of files. The table has columns for Selected, File Name, Directory, Measurement Module, FileType, Operator, and Date&Time. Below the table are buttons for Copy, Rename, Select All/Unselect, and Delete.

Selected	File Name	Directory	Measurement Module	FileType	Operator	Date&Time
<input checked="" type="checkbox"/>	SpectrumScan-1	C:\TestFile\	Spectrum	Measurement Results	Guest	2024/7/23 15:39:26
<input checked="" type="checkbox"/>	DNA-1	C:\TestFile\	Nucleic Acid/Protein	Measurement Results	Guest	2024/7/23 15:34:21
<input checked="" type="checkbox"/>	QuaMethod-1	C:\TestFile\	Quantitation	Measurement Results	Guest	2024/7/23 15:27:59
<input checked="" type="checkbox"/>	QuaMethod-1	C:\TestFile\	Quantitation	Method	Guest	2024/7/23 15:26:48
<input checked="" type="checkbox"/>	TimeScan-1	C:\TestFile\	Time Scanning	Measurement Results	Guest	2024/7/23 15:17:53
<input checked="" type="checkbox"/>	Kinetics-1	C:\TestFile\	Kinetics	Measurement Results	Guest	2024/7/23 14:55:25
<input checked="" type="checkbox"/>	MultiWave-1	C:\TestFile\	Multi-wavelength	Measurement Results	Guest	2024/7/23 14:49:57

Browse file list

File management interface displays the data files stored in the storage, including data, methods and parameters. Click **Measurement Module** and **File Type** to set the filter conditions and display the corresponding file types.

Delete files


Select the file you want to delete to make it selected, then click the button **Delete** to delete it.

Copy files

Select the file you want to copy to make it selected, then click the button **Copy**, pop up the **Path selection** form, select target path and click the button



Rename a file


Click the file name you want to rename, pop up the keyboard, enter a new file name, and click the button  to complete the modification.

Formula editor


Formula Editor is used to write formulas for photometric calculations.




Main interface, click the icon to start **Formula editor**.



Add calculation formula

Formula list interface, click the button **Add formula** to enter to **Formula editor** interface. Enter the formula name as required by clicking the label **Formula Name**, and select the applied module, and enter the calculation formula and click the button  to complete.



Delete calculation formula


Formula list interface, select the method you want to delete and click the button **Delete formula** to delete the method.

Performance verification

Performance Verification is used to verify that the instrument's performance indicators are good.





Main interface, click the icon to start **Performance Verification**.




The performance validation module has two built-in validation templates to meet the United States Pharmacopoeia (USP) and European Pharmacopoeia (EP) for direct use by the user. Users can modify based on the templates or create new test methods and save them for use during testing.


Establishing performance validation methods

Performance Verification interface, click the button  to enter the **Method** interface.




1 Selecting Validation Items

Click **New Method** to create an empty method template. In the **Validation Items** pop-up window, select the items that need to be validated and click .




2 Edit wavelength accuracy verification

Click the **Wavelength accuracy** tab in the verification items column on the left. Enter the parameter editing interface for wavelength accuracy.



- 1) Click the button **Add** to add a test point;
- 2) Click the cell in the **Standard** column to select a standard;




Information If you are using a Standard that is not in the list, you can click the Add button to add it to the list, or you can remove the selected Standard from the list.

- 3) Click the cell in the **Certified** column to enter a certified wavelength test point;


- 4) Click the cell in the **Accuracy** column to enter the allowable deviation for wavelength accuracy;
- 5) Click the cell in the **Repeatability** column to enter the allowable deviation for wavelength repeatability;
- 6) Click the cell in the **Cycles** column to enter the number of tests.

3 Edit photometric accuracy verification

Click the **Photometric accuracy** tab in the Verification Items column on the left. Enter the parameter editing interface for photometric accuracy.




- 1) Click the button **Add** to add a test point;
- 2) Click the cell in the **Standard** column to select a standard;



- 3) Click the cell in the **Wavelength** column to enter a wavelength test point;

- 4) Click the cell in the **Certified** column to enter a certified photometric test point;
- 5) Click the cell in the **Accuracy** column to enter the allowable deviation for photometric accuracy;
- 6) Click the cell in the **Repeatability** column to select the repeatability calculation method and enter the allowable deviation for photometric repeatability;



Repeatability Calculation

SD (Standard Deviation)

RSD (Relative Standard Deviation)


NA (Not Calculated)

✕

- 7) Click the cell in the **Cycles** column to enter the number of tests.

4 Edit photometric linearity verification


Click the **Photometric linearity** tab in the Verification Items column on the left. Enter the parameter editing interface for photometric linearity.



Performance Verification - Method

Wavelength accuracy	Standard	Linearity (%)	Wavelength(nm)														
	PDC	0.999	235.0														
Photometric accuracy	Number of Standard Samples	Cycles	Unit														
	6	6	ppm														
Photometric linearity	<table border="1"> <thead> <tr> <th>Conc-1</th> <th>Conc-2</th> <th>Conc-3</th> <th>Conc-4</th> <th>Conc-5</th> <th>Conc-6</th> </tr> </thead> <tbody> <tr> <td>20</td> <td>40</td> <td>60</td> <td>80</td> <td>100</td> <td>120</td> </tr> </tbody> </table>					Conc-1	Conc-2	Conc-3	Conc-4	Conc-5	Conc-6	20	40	60	80	100	120
Conc-1	Conc-2	Conc-3	Conc-4	Conc-5	Conc-6												
20	40	60	80	100	120												
Stray light																	
Resolution																	


- 1) Click the cell in the **Standard** to select a standard;




- 2) Click the parameter item to set the parameter;
- 3) Enter the concentration of the standard sample in the table.

5 Edit stray light verification

Click the **Stray light** tab in the Verification Items column on the left. Enter the parameter editing interface for stray light.




- 1) Click the button **Add** to add a test point;
- 2) Click the cell in the **Standard** column to select a standard;



- 3) Click the cell in the **Photometric Mode** column to select Abs or %T;
- 4) Click the cell in the **Wavelength** column to enter the test wavelength;
- 5) Click the cell in the **Stray Light** column to enter the allowable deviation for stray light.

6 Edit resolution verification

Click the **Resolution** tab in the Verification Items column on the left. Enter the parameter editing interface for resolution.



- 1) Click the cell in the **Standard** to enter standard name;
- 2) Click **Scanning Range** to set the scanning wavelength range;
- 3) Click **Resolution** to enter the allowable deviation for resolution.

7 Finished

Click the button **Save method** to save the verification method. Establishment of the validation method is complete.


Verify performance

1 Preparing

- The instrument must be warmed up for more than 30 minutes before performing performance verification.
- The instrument must first calibrate the dark current and the system baseline before verification.
- The standard material used for verification must be within the validity period of the test.

Note *Standard filters used to verify the performance of the instrument are not supplied with the device.*

2 Open the performance verification method

Performance Verification interface, Click the button **Method, Open method** to open a performance verification method, click the button .

Performance Verification 🏠

Wavelength Accuracy and Repeatability


Measure

Standard	Certified	Accuracy Allowed	Repeatability Allowed	Measured-1	Measured-2	Measured-3	Measured-4	Measured-5	Measured-6	Accuracy	Repeatability	Result
Holmium	241.2	0.3	0.1									
Holmium	278.1	0.3	0.1									
Holmium	287.6	0.3	0.1									
Holmium	333.5	0.3	0.1									
Holmium	361.1	0.3	0.1									
Holmium	416.6	0.3	0.1									
Holmium	451.3	0.3	0.1									
Holmium	485.3	0.3	0.1									
Holmium	536.9	0.3	0.1									


Photometric Accuracy and Repeatability

Measure


Standard	Wavelength (nm)	Certified(Abs)	Accuracy Allowed	Repeatability Allowed	measured-1(Ab)	measured-2(Ab)	measured-3(Ab)	measured-4(Ab)	measured-5(Ab)	measured-6(Ab)	Accuracy	Repeatability	Result
Neutral fil...	465.0	0.9547	0.008	SD/0.005									
Neutral fil...	546.0	1.0044	0.0081	RSD/0.5%									
Neutral fil...	590.0	1.1079	0.0081	RSD/0.5%									
Neutral fil...	635.0	1.1192	0.0081	RSD/0.5%									



Operations




Method



System Calibration

3 System calibration

Remove anything from the measurement and reference channels and close the sample chamber lid. **Performance Verification** interface, click the button  to do system calibration.

4 Verify

Click the button **Measure** and follow the prompts of the software to put the **reference** or **standard** to complete the test.

Copyrights and Warranties

We warrant the UVStudio software storage disk, USB key and cabling against defects in materials and workmanship for a period of one year from the date of original retail purchase. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, has been modified without the written permission, or if the serial number has been removed or defaced.

In no event will We are liable for direct, indirect, special, incidental, or consequential damages resulting from any breach of warranty, or under any legal theory, including lost profits, downtime, goodwill, damage to or replacement of equipment and property, and any costs of recovering, reprogramming or reproducing any program or data used with UVStudio.

This manual and the software described in it are copyrighted with all rights reserved. Under the copyright laws, this manual or the software may not be copied, in whole or part without written consent of us, except in the normal use of the software or to make a backup copy. The same proprietary and copyright notices must be affixed to any permitted copies as were affixed to the original. This exception does not allow copies to be made for others, whether or not sold, but all material purchased with all (with all backup copies) may be sold, given or loaned to another person.

