

**ESPECTROFOTÓMETRO 4211/50**  
**4211/50 SPECTROPHOTOMETER**  
**SPECTROPHOTOMÈTRE 4211/50**

REF. / CODE / RÉF. HJD002



Este manual es parte inseparable del aparato por lo que debe estar disponible a todos los usuarios del equipo. Le recomendamos leer atentamente el presente manual y seguir rigurosamente los procedimientos de uso para obtener las máximas prestaciones y una mayor duración del mismo.

*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

*Ce manuel est une partie indissociable de l'appareil et doit être mis à la disposition de tous les utilisateurs de l'équipement. Nous vous recommandons de lire attentivement ce manuel et de suivre scrupuleusement les procédures d'utilisation afin d'obtenir des performances maximales et une plus longue durée de vie de l'appareil.*

**LANGUAGE INDEX**

Spanish .....	1-18
English .....	19-35
French .....	36-52

**TABLE OF CONTENTS**

Safety .....	20
Package Contents .....	20
Unpacking .....	20
Installation .....	20
1. Environment Required .....	20
2. Install Spectrophotometer .....	21
Overview .....	21
Symbols .....	21
Main Specifications .....	22
Description of Appearance and Keys .....	22
1. Appearance .....	22
2. Keypad .....	23
3. Description of Keys .....	23
Functions .....	23
Getting Started .....	24
Important Guidelines .....	24
General Operating .....	25
Measuring .....	26
1. Photometry .....	26
2. Quantitation .....	26
3. Utility .....	29
Troubleshooting .....	31
Repair and Maintenance .....	32
1. Daily Maintain .....	32
2. Spare Parts Replacement .....	32
Warranty .....	35
Equipment Disposal .....	35

## SAFETY

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.



- Do not open the device.
- Disconnect the device from the mains supply before carrying out maintenance work or changing the fuses.
- The inside of the device is a high-voltage area Danger!
- Do not use the device if it is damaged, especially if the main power cable is in any way damaged or defective.
- Repairs may only be carried out by the authorized service technicians.
- The device must be connected to a power outlet that has a protective ground connection.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- Do not allow any liquid to enter into the device.
- Do not operate the device in a hazardous location or potentially explosive environment.

## PACKAGE CONTENTS

Description	Quantity
■ Spectrophotometer	1PC
■ 10mm Glass Cuvette	4PCS
■ 10mm Quartz Cuvette	2PCS
■ Power Cord	1PC
■ User Manual	1PC
■ Dust Cover	1PC

## UNPACKING

Open the package, according to carefully check the packaging packing list items, if found inside the packaging are missing or damaged items please contact us and authorized contractual partners.

## INSTALLATION

### 1. Environment Required

To ensure the best performance, the following conditions are required:

- The best working temperature range is 16—35°C and the humidity is 45—80%.
- Keep it as far as possible away from the strong magnetic or electrical fields or any electrical device that may generate high-frequency fields.
- Set the unit up in an area that is free of dust, corrosive gases and strong vibrations.

- Remove any obstructions or materials that could hinder the flow of air under and around the instrument.
- The power requirement is  $110\pm 11V/60\pm 1Hz$  or  $220\pm 22V/50\pm 1Hz$ .
- Use the appropriate power cord and plug into a grounded outlet.
- If the local voltage is not stable, a voltage regulator is required.
- Be away from direct sunlight.

## 2. Install Spectrophotometer

### Placement

Place the instrument on the stable table carefully.

### Install Printer (Printer is Optional Accessories)

Check to confirm instrument power switch is turned off. Connect the printer's data cable to the instrument's parallel port.

### Connect the Power Cords

Check to confirm instrument power switch is turned off. Connect the power cables of the instrument and printer to two separate power outlets.





## OVERVIEW

4211/50 spectrophotometer is an electrical measure instrument which is widely used in the laboratories.

- |                               |               |
|-------------------------------|---------------|
| ■ Use Frequency:              | Intermittence |
| ■ Excessive Voltage(Current): | No            |
| ■ Pollution Class:            | Class 1       |

## SYMBOLS

The following chart is an illustrated glossary of the symbols that are used in this manual.

	Caution, Danger!
	Caution, High Voltage!
	Caution, Hot!
	This equipment will be recycled at the end of its useful life.

## MAIN SPECIFICATIONS

Optical System	Single beam
Wavelength Range	200-1000 nm
Wavelength Accuracy	$\pm 2$ nm
Wavelength Repeatability	0,8 nm
Photometric Range	-0,3-3A, 0-200%T
Photometric Accuracy	$\pm 0.5\%T$
Photometric Repeatability	0.3%T
Spectral Bandwidth	4 nm
Stray Light	0,3%T@220nm&360nm
Stability	$\pm 0,002A/h@500nm$
Work Mode	Photometry, Quantitation
Interface	USB, Parallel(printer)
Power Requirement	110 V CA/60 Hz o 220 V CA/50 Hz
Dimensions	490 x 360 x 210
Weight	14 kg

## DESCRIPTION OF APPEARANCE AND KEYS

### 1. Appearance

#### Front View



## Back View









- 1 LCD Display
- 2 Keypad
- 3 Lid of Sample Room
- 4 Rod
- 5 LCD Contrast Adjust
- 6 Printer port
- 7 USB port
- 8 Cover of Fan
- 9 Power Socket
- 10 Power Switch
- 11 Cover of Cooling Vents

## 2. Keypad



## 3. Description of Keys

	<b>SET Key:</b> Set Parameters
	<b>GOTO Key:</b> Set Wavelength
	<b>ZERO Key:</b> Blank
	<b>PRINT Key:</b> Print measuring result
	<b>Function Key:</b> Functions according to the screen
	<b>UP, DOWN Keys:</b> Scroll menu/data and set Y scale

## FUNCTIONS

**Photometry**

Display results as Abs, %T or Energy.

**Quantitation**

By using a Standard Sample to establish the Standard Curve.

## GETTING STARTED

The following chart describes the basic operation of the instrument.

### Turn On and Self-check

Switch on the power. Then the instrument begins to self-check and 20 minutes' warm up. Self-check includes the following steps: Turn on lamps - Check Sensor - Initialize AD - System position - Get Dark Current - Warm up.

Self-test . . .

Warm up 20 minutes,  
Any key to skip

After warm up, instrument displays Main Interface.

WL: 500.0nm		0.000A	
100.0%T			
Basic		Quantitative	

## IMPORTANT GUIDELINES

- Reagents and dilution buffers can cause cauterization and other damage to health.
- Samples (nucleic acids, proteins, bacteria cultures) can be infectious and cause serious damage to health.
- During sample preparation, measuring procedures and maintenance and cleaning work, observe all local laboratory safety precautions (e.g. wear protective clothing and gloves, use of disinfectant) regarding the handling of sample material.
- Dispose of measuring solutions and cleaning and disinfectant materials in accordance with the relevant local laboratory regulations.


## GENERAL OPERATING

### Select Application

Main interface, press the key  (left) to enter into.

### Set Wavelength

Test interface, press key  to set wavelength,   to modify wavelength value, then press key  (left) to go to wavelength and blank.



### Set Parameters

Press  enter into setup interface,   to select items or input parameters,  (left) to confirm.

### Delete the test result and stored data

Test Interface, press the key , then press key   to select "Clear Data, not Print",  (left) to delete.

### Blank

Put the Reference in the light path, press  to do blank.

### Measure Samples

Put the samples in the light path, press  (left) to measure.

### Print the test results

Test Interface, press the key    select "Print, clear data", press the key  (left) to print.

### Store the Standard Curve

After got the Standard Curve, press   input the file name and press  (left) to save.

### Load the Standard Curve


"Quantitative" interface, press   to select "Load Curve", press   to choose the curve you want, press  (left) to open.

## MEASURING

### 1. Photometry

#### Step 1. Start Photometry

Main Interface, press key  (left) to choose “Basic”.



#### Step 2. Set Wavelength

Press  to set wavelength, press   to input wavelength value, press  (left) to go to wavelength.

#### Step 3. Blank

Put the Reference in the light path and press  to do blank.

#### Step 4. Measure samples


Put the sample in the light path, and then the result displays on the screen automatically, press  (left) to record.

WL: 500.0nm 0.000A			
Nb.	WL	%T	Abs.
1	500.0	100.0	0.000
2	500.0	100.0	0.000
Test			Cancel

### 2. Quantitation


#### Step 1. Start Quantitation

Main Interface, press key  (right) to choose “Quantitative”.



## Step 2. Establish or call Standard Curve

<input checked="" type="radio"/> Create Curve	WL: 500.0nm
<input type="radio"/> Load Curve	<input checked="" type="radio"/> Coefficient
<input type="radio"/> Clear Curve	<input type="radio"/> Standard Curve
OK	Cancel




WL: 500.0nm	<input checked="" type="radio"/> Coefficient
	<input type="radio"/> Standard Curve
OK	Cancel

2 methods to establish Standard Curve:

### Method 1: Establish Standard Curve by inputting coefficients

1) **Establish.** Press  , to select "Coefficient", then press  (left) to confirm.

WL: 500.0nm	<input checked="" type="radio"/> Coefficient
	<input type="radio"/> Standard Curve
OK	Cancel




WL: 500.0nm	Work WL.:
	500.0 nm
OK	Cancel

2) **Set wavelength.** Press   to input wavelength value, press  (left) to confirm.

3) **Set coefficients K and B.** Press   to input coefficient K, press  (left), same way set B.

WL: 500.0nm	Work WL.:
	500.0 nm
OK	Cancel




WL: 500.0nm	Coefficient K:
	00000
OK	Cancel

### Method 2: Establish Standard Curve by using Standard Samples

1) **Establish.** Pulse   para seleccionar "Standard Curve" y, a continuación,  (izqu.).

WL: 500.0nm	<input type="radio"/> Coefficient
	<input checked="" type="radio"/> Standard Curve
OK	Cancel




WL: 500.0nm	0.000A
	Please insert Blank:
OK	Cancel

2) **Set Wavelength.** Press  to set wavelength, press   to input wavelength value, press  (left) to go the setting value.

3) **Blank.** Put the Reference in the light path, press  (left) to do blank.

4) **Setup number of Standard Samples.** Press   , to input the quantity of standard samples(no more than 9), press  (left) to confirm.


WL: 500.0nm	
Coefficient K: 00000	
OK	Cancel



WL: 500.0nm	
Number: 3	
OK	Cancel

5) **Calibrate Standard Samples.** Put the corresponding standard samples in the light path as the screen indicates, press   , to input the concentration, press  (left) to confirm, complete finish all the standard samples.

WL: 500.0nm	
Number: 3	
OK	Cancel




WL: 500.0nm	
Insert 1# Standard.	
Input 1# Conc: 00000	
OK	Cancel

#### Load the Stored Curve

Press   to choose "Load Curve", press   to choose the curve, press  (left) to confirm.


<input type="radio"/> Create Curve <input checked="" type="radio"/> Load Curve <input type="radio"/> Clear Curve	
OK	Cancel



$C=K \times A+B$	WL.
$C=1.000 \times A+0.2$	500.0 <
$C=0.02 \times A+0.32$	470.0
OK	Cancel

Press  (left) to enter the test mode after building or loading standard curve.

**Step 3. Enter into Measuring Interface.** Press  (left) to enter into the Quantitation Measuring Interface.




**Step 4. Blank.** Put the Reference in the light path, press  to do blank.

**Step 5. Measure Samples.** Put the sample to be tested in the light path, press  (left) to measure. Then the test result will display in the data sheet. Repeat this step to finish measuring all the samples.


500.0nm		0.000A	
Nb.	WL	Abs	Conc.
1	500.0	0.039	0.078
2	500.0	0.042	0.084
3	500.0	0.041	0.082

### 3. Utility

Main Interface, press  to go into utility setting.




**Test Mode** Press   to choose "Test Mode", press  (left) to enter, press   to choose "Abs", "%T", "Energy", press  (left) to confirm.



**Turn On/Off D2 Lamp** Press   to choose "D2 Lamp", then press  (left) to enter.  
Press   to choose "On" or "Off", press  (left) to turn on/off.


<input type="radio"/> Test Mde		
<input checked="" type="radio"/> D2 Lamp	On	
<input type="radio"/> WLamp	On	
OK		Cancel



√ D2 Lamp	On	
<input checked="" type="radio"/> On		√
<input type="radio"/> Off		
OK		Cancel

**Turn On/Off W Lamp** Press   to choose "W Lamp", then press  (left) to enter.  
Press   to choose "On" or "Off", press  (left) to turn on/off.

<input type="radio"/> Test Mde		
<input type="radio"/> D2 Lamp	On	
<input checked="" type="radio"/> WLamp	On	
OK		Cancel




√ WLamp	On	
<input checked="" type="radio"/> On		√
<input type="radio"/> Off		
OK		Cancel

**Get Dark Current** Keep the light path without anything blocking, press   to choose "Dark Current", then press  (left) to resample Dark Current.

**Note:** During the course, open the lid of the compartment is prohibited.


<input checked="" type="radio"/> Dark Current		
<input type="radio"/> Reset WL		
<input type="radio"/> Version		
OK		Cancel



Get Dark Current!		
-------------------	--	--

**Reset Wavelength** Keep the light path without anything blocking, press   to choose "Reset WL", then press  (left) to reset wavelength.

<input type="radio"/> Dark Current		
<input checked="" type="radio"/> Reset WL		
<input type="radio"/> Version		
OK		Cancel



Calibrating WL!		
-----------------	--	--

**About Version** Press   to choose “Version”, press  (left) to view version information, press any key to return.



## TROUBLESHOOTING

Review the information in the table below to troubleshoot operating problems.

Problem	Cause	Solution
Power on, no response	<ol style="list-style-type: none"> <li>1. Power cord connection is not reliable</li> <li>2. Fuse burning</li> </ol>	<ol style="list-style-type: none"> <li>1. Improve connection</li> <li>2. Replace fuse</li> </ol>
Measurement uncertainty	<ol style="list-style-type: none"> <li>1. Warm up is not enough</li> <li>2. Sample is not Stable</li> <li>3. The concentration of sample is too high</li> <li>4. Power Supply Voltage Low or not Stable</li> <li>5. Lamp damage or lamp life maturity</li> </ol>	<ol style="list-style-type: none"> <li>1. Warm up more time</li> <li>2. Improve the sample</li> <li>3. Dilute the sample</li> <li>4. Improve the Power Supply</li> <li>5. Replace lamp</li> </ol>
Dark Current Error when self-check	<ol style="list-style-type: none"> <li>1. The lid of the compartment is open during self-check</li> </ol>	<ol style="list-style-type: none"> <li>1. Close the lid, restart</li> </ol>
System Calibrate Failed	<ol style="list-style-type: none"> <li>1. Something block the Light path</li> </ol>	<ol style="list-style-type: none"> <li>1. Remove it, calibrate again</li> </ol>
Power on, back light is OK, but nothing display on the screen or display is not clear	<ol style="list-style-type: none"> <li>1. Display Contrast problem</li> </ol>	<ol style="list-style-type: none"> <li>1. Adjust the contrast potentiometer</li> </ol>
Measurements inaccurate	<ol style="list-style-type: none"> <li>1. Cuvettes were contaminated</li> <li>2. Samples were contaminated</li> <li>3. Bad matching of the cuvettes</li> <li>4. Dark current error</li> </ol>	<ol style="list-style-type: none"> <li>1. Clean the cuvettes</li> <li>2. Improve the samples</li> <li>3. Improve the matching of the cuvettes</li> <li>4. Resample dark current</li> </ol>

## REPAIR AND MAINTENANCE

### 1. Daily Maintain

---

#### Check the Compartment

After measurement, the cuvettes with sample solutions should be taken out of the compartment in time. Or the volatilization of the solution would make the mirror go moldy. Users must pay more attention to the corrosive sample and liquid easy to volatilize. Any solution remains in the compartment should be wiped off immediately.

#### Surface Cleaning

If paint drops fall on the instrument cover, wipe them off immediately with a damp towel. Organic solution is forbidden to be used to clean the cover. Please wipe off the dirt on the cover timely.

#### Cleaning the Cuvettes

After every test or after a solution change, the cuvettes should be cleaned carefully, or the remains on the surface would cause measuring error.

### 2. Spare Parts Replacement

---

#### Fuse replacement



**Danger! Be sure to switch off the power and unplug the socket before replacement!**

#### Step 1. Tools preparation

Prepare a 3×75 Flat Blade screwdriver.

#### Step 2. Switch Off the power supply

Switch off the power supply, and unplug the socket.


#### Step 3. Take out the Fuse Seat

Take out the fuse seat by the screwdriver.



#### Step 4. Put the new fuse

Pick out the spare fuse (3.15A/250V) and replace it to the working position.



### Step 5. Reset the fuse seat

Replace the fuse seat in the power socket.

### Step 6. Switch on the power

Plug the socket and switch on the power.

### Lamps replacement



**Hot! Wait 20 minutes before open the lamp chamber after power off to avoid scald!**

### Step 1. Tools preparation

Prepare a 6×150mm Cross Blade screwdriver and a pair of glove.

### Step 2. Power Off

Switch off the power supply and unplug the socket.


### Step 3. Open the cover

Unscrew the 4 screws indicated (each side with 2 screws) and remove the cover.




**Step 4. Open the cover of the light chamber**



Unscrew the 2 screws on the light chamber cover and remove it.

**Step 5. D2 lamp replacement**

Unscrew the 2 screws on the D2 Flange (No.1), unplug the connector in the Power Board (No. 2) and remove the D2 lamp. Draw on the cotton gloves and replace a new lamp. Fix the 2 screws and plug the connector again.


**Step 6. W lamp replacement**

Pull out the defected W lamp and draw on the cotton gloves. Insert the new W lamp as deep as possible on the lamp seat. Be sure to keep the filament in the same direction as the old one face.



### Adjust the position of the W lamp

Switch on the power(the Switch Mirror should be placed to the position as indicates). Observe the light spot, it should in the center of the entrance hole. If the light spot deviates to Left or Right, then loosen the No.1 screws and move the lamp seat to Left or Right until it focus on the center of the slot. Then fix the screws. If the light spot deviates to Up and Down, then loosen the No.2 screws and move the lamp seat Up and Down until the light spot focus on the center of the slot. Then fix the No. 2 screws again.



### Step 7. Finish

Reset the cover of the light chamber and fix the screws. Reset the cover of the instrument and fix the screws. Install the rod of the sample compartment, then the work is finished.

## WARRANTY AND EQUIPMENT DISPOSAL

AUXILAB S.L. warrant that this product will be free from defects in material and workmanship for a period of 2 years from date of purchase except the lamps. Lamps have a warranty of 1000 hours usage time or 6 months max. If a defect is present, we will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.



It's your responsibility to correctly dispose of your equipment at lifecycle -end by handing it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you!