PÉNDULO DE HILO

PENDULUM

REF. QLD005

INDEX OF LANGUAGES-ÍNDICE DE IDIOMAS

ESPAÑOL2

ENGLISH4

PÉNDULO DE HILO

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PENDULUM

Designed for performing a number of pendulum related experiments quickly and accurately. It comprises a meter rule with zero at the bottom, graduated in millimeters and reading every centimeter, attached vertically on a support rod. Towards the upper end of the rod a hanger is fixed, from which 3 pendulum bobs, each of 25mm diameter and of different material (wood, steel and aluminium) are suspended. A sliding index mounted on the meter rule can be adjusted for vertical position and has a projection through which first pendulum cord passes thus controlling the effective length of suspension of first pendulum. The complete assembly is mounted on a stable metal stand. The hanger has provision for controlling the length of suspension of each pendulum.

Aim

To study the variation of the time period for simple pendulum with different lengths and bobs of different masses.

Theory

The time period of a simple pendulum is given by:

$$T = 2\pi \sqrt{\frac{l}{g}} \quad (1)$$

Where l is the length of simple pendulum and g is acceleration due to gravity. The length l of the pendulum is the distance between the point of suspension and the center of the bob.

From equation (1) the acceleration due to gravity is:

$$g = 4\pi^2 \frac{l}{T^2} \quad (2)$$

From equation (1) it can also be obtained that period T is proportional to the square root of pendulum length l and thus the graph between l and T will be an exponential curve. Also, the square of period (T^2) is proportional to pendulum length l and in this case the graph will be a straight line. By determining the mean value of l/T^2 from the slope of the graph, the acceleration due to gravity g can be calculated using equation (2).

Procedure

1. Place the apparatus on a stable table
2. Determine the diameter of each bob using vernier calipers and then radius of each bob.
3. Displace the bobs from mean positions and release. Ensure that the bobs do not vibrate.
4. Note the time for 20 oscillations for each bob.
5. Repeat steps 3 and 4 at least for three times and take the mean.
6. Calculate the time period of oscillation of each bob.
7. Plot a graph between l taking it along x-axis and T along y-axis

8. Find the value of T^2 for various observations
9. Plot the graph between l taking it along x-axis and T^2 along y-axis.

Observation table

Sr. nº	Mass of the bob (g)	Time for 20 oscillations (t) (sec)				Period $T=t/20$
		(i)	(ii)	(iii)	Mean	
1.	Wooden bob mass, $m_1 =$					
2.	Steel bob mass, $m_2 =$					
3.	Aluminum bob mass, $m_3 =$					

Conclusions

1. Time period T of a simple pendulum increases with the increase in the length l of the bob pendulum.
2. Time period T of the pendulum does not depend upon the mass and type of the bobs.
3. The graph between l and T is an exponential curve and the graph between l and T^2 is a straight line.

Precautions

1. The bobs should oscillate without spinning.
2. The experiment must be done in a room free of air currents that could disturb the oscillation of the pendulum.