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SET OF 8 TUNING FORKS, STEEL

REF :QLD007

INDEX OF LANGUAGES-ÍNDICE DE IDIOMAS

ESPAÑOL.....	2
ENGLISH.....	4

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DESCRIPCIÓN

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SET OF 8 TUNING FORKS, STEEL

Description

With the best quality, plain shanks and with frequency clearly marked. Each set is composed of 8 forks, one each of diatonic scale frequencies: C1-(256), D-(288), e- (320), f (341.3), G-(384), A-(426.6), B-(480) Y C2-(512).

The complete set is supplied in a box.

The sound

The speed of sound in the atmosphere is 343,2m/s (a 20°C temperature, 50% humidity and at the level of the sea. The speed depends on the medium. The tuning forks is a metal rod U-shaped. The sound is emitted by the tuning forks with a frequency. We can use the tuning forks to calculate the speed of the sound.

If we know the frequency of the tuning fork we can calculate the speed of the sound in the air, as we can see in the figure.

We have a graduated cylinder with water. Put the tuning fork on top very closed and make it vibrate. We go down the water until start resonance, a maximum intensity.

The distance of the cylinder L is that we need to calculate.

$$f = \frac{2n + 1}{4} \frac{v_s}{L} \quad n = 0, 1, 2, 3, \dots$$


Where:

F = frequency

N= vibration mode

Vs= sound speed m/s

L= distance in meters

Materials

- Graduated cylinder
- Tuning forks
- Water
- Meter if it is necessary

Practice

1. Select a tuning forks with a specific frequency in this case f= 426.6 Hz.
2. Go down the water until resonance
3. Note the distance L in meters
4. Calculate
5. Vs is the speed

Example

We have a tuning forks with a frequency f=440Hz. Then remove the water until l= 58 when vibrates. Calculate.

$$440 = \frac{2 \cdot 1 + 1}{4} \frac{v_s}{0.58} \quad v_s = 340.3 \text{ m/s}$$

Storm distance

Since we see a ray we can count the time until we see the thunder and calculate the distance with the formula. The light speed is 300000 km/s and the sound 0.340km/s.

X=v.t

X is the distance we want to know

V speed of sound

T time

Example

We see a ray and in 3 seconds we hear the thunder.

V=340 m/s and t=3.

V= 340 m/s x 3s = 1020m