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## ***USER'S MANUAL***

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*This manual should be available for all users of this equipment. It is advisable to read this manual and follow the use instructions carefully in order to get the best performance and extend the working life of the equipment.*

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## 1. Configuration and features

- The single annular phase contrast kit is composed of:

- 3 achromatic phase contrast objectives (10x, 40x y 100x)
- 3 ring diaphragms
- 1 condenser, with centring screws for the ring diaphragms
- 1 centring eyepiece

- In each ring diaphragm is indicated the magnification of the contrast phase objective to be used with:

- **10**: observations with the 10x phase contrast objective
- **40**: observations with the 40x phase contrast objective
- **100**: observations with the 100x phase contrast objective

## 2. Assembly and use

- Remove the objectives from the nosepiece of the microscope and screw the phase contrast objectives in their place, following an increasing magnification order clockwise.

- Put the ring diaphragm that matches the magnification of the contrast phase objective you are going to use (it is recommended to start the observations with the lowest magnification objective, 10x) in the phase contrast condenser. To do so, introduce the diaphragm by the bottom of the condenser and hold it with the smaller screw. Remove the condenser of the microscope and place the phase contrast condenser with the ring diaphragm in its place. Hold the condenser tightly with the microscope holding screw.

- Switch on the microscope and turn the nosepiece until the 10x phase contrast objective is placed in the path of the light beam.

- Remove one of the eyepieces of the microscope and place the centring eyepiece in its place

- Looking through the centring eyepiece, focus sharply the light ring (clear ring) and the phase plate (dark ring). To do so, lighten the screw of the centring eyepiece and adjust the height of the upper part until the image is correctly focused. If you don't see the illumination ring properly, adjust the height of the condenser

- Using the centring screws of the condenser to move the ring diaphragm, match the illumination ring and the phase plate (Fig. 1).


Figure 1

- Remove the centring eyepiece and place the eyepiece of the microscope; the equipment is ready to make phase contrast observations.

- When you change to an objective with higher magnification, you must remove the phase contrast condenser, place the corresponding ring diaphragm and repeat the procedure of adjustment and centring previously described.