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POWER SUPPLY

Manual de uso / User's manual

Referencia/Code

ZFD014

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1. Safety

Electrophoresis power supplies use high output voltages that are electrically isolated from earth ground through a protective impedance to minimize the risk of electrical shock to the user. The following guidelines should be observed and followed when using a power supply.

Electrophoresis power supplies have passed test for operation at temperatures between 0° and 40°C, with relative humidity between 0 and 95% non-condensing.

Operating the power supply outside these conditions is not recommended by our company and will void the warranty.

1. To ensure adequate cooling of the power supply, be sure that there is at least 6 cm clearance around the power supply. Do not block the fan vents at the rear of the unit.
2. Always connect the power supply to a 3-prong, grounded AC outlet, using the 3-prong AC power cord provided with the power supply
3. Do not operate the power supply in extreme humidity (less than 95%) or where condensation can short the internal electrical circuits of the power supply.
4. When taking the power supply into a cold room, the unit can be operated immediately. However, when removing the power supply from the cold room, let the unit equilibrate to room temperature for a minimum of 2 hours before using it.
5. Never connect a high voltage output lead to earth ground. This defeats the floating electrical isolation of the power supply and exposes the user to potentially lethal high voltages.

This instrument is intended for laboratory use only.

This product conforms to the class A standards for Electromagnetic Emissions, intended for laboratory equipment applications. It is possible that emissions from this product may interfere with some sensitive appliances when placed nearby or on the same circuit as those appliances. The user should be aware of this potential and take appropriate measures to avoid interference.

2. Introduction

The C series Power supply has the following Technical specifications and Characteristic features:

You can switch on the equipment with the bottom located in the right face of the equipment. There are five bottoms on the frontal face: run/stop, edit, enter and two arrows to increase or decrease the numbers.

Technical specifications

- Type of Output: Constant-Voltage, Current or Power
- Increment: 1V, 1mA, 1W
- Time Range: 1min~99hr, 59min
- Display: LCD display with backlight
- Output jacks: Four sets of output jacks
- Output voltage: 10-600V
- Output current: 1-500mA
- Energy: 1-300W
- Supply: 220+-10V, 50Hz

Characteristic features

- Molding shell, touch keys, dual core microprocessor
- intelligent control;
- Indicate the preset value and the actual output value at the same time;
- It can store 10 electrophoresis methods;
- With automatic memory function;
- With standard, timing operation function;
- With constant voltage, constant current, constant power, misoperation, fault intelligent prompt functions;
- Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating.

3. Setup and Operation

- a. Turn on the power switch, enter into Setting State.

EDIT[0]	
U= 100V	0v
I= 500mA	0mA
P= 300W	0w

Possible values:

U: 10-600V

I:1-500mA

P:1-300W

Time: 1seg- 99hr 59min

Increment 1 unit.

- b. The cursor will flash at "U=□□□V". Set the value of constant voltage by pressing "▲" or "▼" button.
- c. The cursor will flash at "I = □□□mA" when press "ENT" button. Set the value of constant current by pressing "▲" or "▼" button.
- d. The cursor will flash at "P = □□□W" when press "ENT" button. Set the value of constant power by pressing "▲" or "▼" button. The range of power can be set from 1 to 300W.

The cursor will return and flash at "U = □□□V" when you press "ENT" button. Set the value of constant voltage again by pressing "▲" or "▼" button.

- e. Power supply in Setting State now:
Press "RUN/STOP" button, it will enter into Operating State.

EDIT[0]	RUN! 00:01
U= 100V	▶ 100v
I= 500mA	350mA
P= 300W	35w

- f. Power supply in Operating State now:
Press "RUN/STOP" button, it will stop output and return to Setting State.

- g. In Setting State, continue to press "EDIT" button, it will enter into Editing State and display:
"SAVE [□]" → "LOAD [□]" → "T= □□:□□" → "QUIT" → "SAVE [□]"...

EDIT[0]	
U= 100V	0v
I= 500mA	0mA
P= 300W	0w

- h. At "SAVE □", press "▲" or "▼" button to select number of memory program. The range of number is from 1 to 9. Then press "ENT" button to save current method into this program.
- i. At "LOAD [□]", press "▲" or "▼" button to select number that needed. The range of number is from 1 to 9. Then press "ENT" button to load the current method from program.
- j. At "T= □□:□□", press "▲" or "▼" button to set operating time.
- k. The range of timer is from 1min to 99hs and 59mins. Then press "ENT" button to save it.

EDIT[0]	
U= 100V	0v
I= 500mA	0mA
P= 300W	0w

- l. At "QUIT", press "ENT" button to quit Editing state, and return to Setting State.

4. Error Alarm

Type 1: When the power supply run with Over-load, it will stop automatically and alarm with buzzer, then display "ERROR1!".Turn off the switch immediately.

Type 2: When the power supply run with No-load, It will stop automatically and alarm with buzzer, then display "ERROR 2!".Turn off the switch immediately.