ANALIZADOR DE HUMEDAD HALÓGENO
HALOGEN MOISTURE ANALYZER
ANALYSEUR D'HUMIDITÉ-THERMOBALANCE

Ref. KCD001

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Sample weight is overweight	1. Exceed the maximum load.	1. Reduce the load immediately.
Overweight or ultralight during calibration	1. The calibration weight is incorrect 2. Before calibrating the moisture meter for weighing, there is an object on the weighing pan.	1. Take the object to zero and calibrate 2. Calibrate the moisture meter with a calibration weight
The weighing result is unstable (the data fluctuates)	1. Too much airflow in the working environment 2. The workbench where the moisture meter is located is unstable.	1. Airflow should be avoided as much as possible 2. The moisture meter should be placed on a stable workbench.
Incorrect weighing result	1. It is not cleared before the scale 2. The moisture meter is used without calibration, or the calibration code is inaccurate.	1. Touch the tare button 2. Recalibrate the moisture meter.
There is no response when touching a function button	1. Possible instantaneous interference 2. The power supply voltage is incorrect, or the power cable is not connected properly.	1. It can be turned off/on again or plugged in again 2. Switch to a normal power supply or connect well the power cable.

INDEX OF LANGUAGES

Spanish 2-9

English 10-17

French 18-25

INDEX OF CONTENT

Preface 10

1. Overview 11

2. Structure 12

3. Features 12

4. Technical parameters 12

5. Installation 13

6. Operation 14

7. Maintenance 16

8. Troubleshooting 16

PREFACE

In order to get a better experience, please read this manual carefully and follow the safety operation regulations!

Please keep this instruction manual for reference when needed!

NOTES

- Please ensure that only personnel with relevant training can operate and use this instrument.
- Please comply with safety regulations, personal safety, accident prevention and other related regulations.
- Please fully insert the power plug of the moisture analyzer into the power socket, please do not use a power source other than the specified one.
- Do not plug or unplug the power plug with wet hands.
- Do not damage, modify, pull, excessively bend, or twist the power cord, and do not place heavy objects on the power cord.
- Please place the instrument on a stable, clean, non-slip, dry and fireproof surface.
- Ensure that the moisture meter should be placed around enough space to prevent heat accumulation and overheating (at least 1 meter above the moisture meter).

- Please pay attention to check the instrument and accessories before each use to ensure no damage.
- Note! Before operating the instrument, please warm up for at least 30 minutes to ensure the accuracy of the moisture analyzer□
- Note! During the heating process of the moisture meter, the upper cover part will heat up. Please do not place anything, let alone touch it directly with your hands. Switch the upper cover with the handle of the front cover.
- Note! Do not cover the black round holes of the upper cover during the heating process to maintain good heat dissipation. Because the temperature of the drying unit will rise, do not place flammable materials on or around the instrument. Please be careful when removing the sample. Or other sample containers make it hotter and easy to cause burns.
- Please select the appropriate protective device according to the type of measuring medium.
- When handling toxic and volatile media, please pay attention to it in a suitable fume hood.
- When handling flammable and explosive media, please keep it cool and dry.
- When handling corrosive media, please take a small amount of material for testing.
- Please turn off the power after the operation is completed and unplug the power plug.
- Do not disassemble and adjust the parts of the instrument at will. When the spare parts are damaged, please only use the original spare parts for replacement.
- When the equipment fails, the power supply should be cut off in time, and contact our company for maintenance guidance or return to the factory for maintenance as soon as possible.
- The weighing pan and shell should be gently scrubbed with a soft cloth and toothpaste frequently and must not be scrubbed with a strong dissolving agent.
- When cleaning the instrument, please turn off the instrument and cut off the power before proceeding!

1. OVERVIEW

The electronic halogen moisture tester adopts the principle of halogen lamp pyrolysis weight and is a fast moisture tester. While the moisture analyzer measures the weight of the sample, the halogen lamp heating unit and moisture evaporation channel quickly dry the sample. During the drying process, the moisture meter continuously measures and instantly displays the missing moisture content “%” of the sample. After the drying process is completed, the final determination is made; the moisture content value of is locked and displayed.

Compared with the international oven heating method, the halogen heating can reach the maximum heating power in the shortest time, and the sample is quickly dried at high temperature. The test results are in good agreement with standard oven method, and the test efficiency is far higher. Generally, it only takes a few minutes to complete the measurement. The operation is simple, the test is accurate, and the moisture value, dry value, measurement time, temperature value, and other data can be displayed separately. Besides it has the possibility of connecting with computers and printers.

1. Nine-pin connection mode: 5 pins grounded, 3 pin connected to send.
2. There are 9 kinds of printing parameters.

■ Look at the test data

The user can recall and view the stored test data. Proceed as follows:

1. Select the set parameters or the samples that have been tested and saved (with changes)
2. Touch to return to the homepage
3. Touch the saved record to view different test moisture parameters.

■ Data interface

Serial interface parameter setting:

1. Baud rate 9600
2. No surprises

7. MAINTENANCE

1. When the moisture meter is used frequently, the moisture meter should be continuously energized to reduce the warm-up time and keep the moisture meter in a relatively stable state. If the moisture meter is not used for a long time, the power should be turned off.
2. The moisture meter should be kept clean. Beware of dust and other things entering the moisture meter. The weighing pan and housing should always be wiped clean with a soft cleaning agent and soft cloth. Do not scrub with strong solvents such as alcohol or gasoline. The moisture meter should not be placed in an environment with corrosive gas.
3. According to the degree of use of the moisture meter, periodic inspections should be made.

8. TROUBLESHOOTING

Failure	Reason	Solution
The display does not light up	<ol style="list-style-type: none"> 1. The moisture meter is not connected to the power supply normally 2. The moisture meter shows that the switch is not turned on 3. The fuse wire is damaged 	<ol style="list-style-type: none"> 1. Try to connect to the power supply 2. Adjust the fuse wire, if it burns out again, it needs to be repaired.
Ultra-light or zero ultra-light	<ol style="list-style-type: none"> 1. The load is under light load without the upper tray 2. The pan support is not installed properly 3. The internal calibration number may be damaged. 	<ol style="list-style-type: none"> 1. Reinstall the pan support 2. Touch any key to recalibrate the moisture meter to weigh.
Overweight or zero overweight	<ol style="list-style-type: none"> 1. There are heavy objects on the weighing pan when starting up 2. Exceed the maximum load 3. The internal memory calibration number may be destroyed. 	<ol style="list-style-type: none"> 1. Remove the heavy objects on the weighing pan 2. Reduce the load immediately 3. Touch any key to recalibrate the moisture meter to weigh.

1. When you need to modify the drying temperature, touch the display box of the drying temperature, and set the drying temperature in the pop-up interface for setting the temperature, and touch the confirm key to return to the timing drying interface.
2. When you need to modify the drying time, touch the display box of the drying time, and set the drying time in the pop-up interface for setting the drying time, and touch the OK key to return to the timing drying interface.
3. After placing the empty tray, if the weight is not zero, touch the tare button to reset the weight to zero. After the weight is zero, place 3-5 grams of the sample. Close the cover. After the weight is stable, touch the dry start button and the lamp will light up. Start to dry and enter the drying interface.
4. In the drying interface, you can see the changes in the temperature in the cavity of the moisture meter, the weight of the sample, and the heating time. Touch the conversion key to cycle through the water loss rate, drying rate, moisture regain rate, wet weight rate, and touch the emergency stop button to immediately stop drying.

When the drying time reaches the set drying time, the moisture meter will stop drying. Enter the drying time result view interface, which can display 9 water measurement parameters:

- Sample weight (G)
- Dry weight (g)
- Heating time (min; sec)
- Drying temperature
- End way
- Water loss rate (i.e., water content) = $(G-g)/G$
- Drying rate = g/G
- Moisture regain = $(G-g)/G$
- Wet weight rate = G/g

Automatic drying

Automatic drying means that when the weight of the sample is no longer decreasing, the drying ends. Except that the automatic drying does not require drying time, other operations are the same as timing drying.

Current record

The current record is to view the last drying parameters of the sample of the moisture analyzer.

1. On the main interface, touch the current record key to enter the current record interface.
2. In the current record interface, you can view the drying parameters. Print, save or delete. If the instrument has not been dried after startup, the touch screen will display no current record and return to the main interface.

Keep records

Save records means that after the moisture meter completes a drying test, press the Save button to save the test results in the instrument (up to six groups can be stored).


1. On the main interface, touch the save record key to enter the save record interface.
2. In the save record interface, touch the sample number button to view the drying parameters of the sample. Print or delete.

6.4. Print settings

This moisture meter can print test data or the parameters of the sample to be tested through an external micro printer.

1. Select the parameter that has been set or the sample that has been tested and saved and touch the print button to print the parameters and test data of the test sample.
2. Print via RS232 interface.

2. STRUCTURE

3. FEATURES

1. Fully automatic heating cabin door double-drive ultra-quiet opening to prevent the influence of artificial lifting of the lid and the weight, and at the same time meet the requirements for automatic detection.
2. 5-inch high-definition smart touch screen, real-time dynamic analysis of the curve, direct display of moisture percentage/dry weight percentage/dry weight.
3. The extreme heating temperature of 180°C meets the test of most samples.
4. 9 groups of recent test records can be retrieved at any time and the test results can be printed.
5. The stainless-steel heating inner chamber improves the efficiency of sample testing.

4. TECHNICAL PARAMETERS

1. Surrounding environmental conditions

The technical parameters are valid in the following environments:

- Recommended ambient temperature: 15°C-30°C; the analyzer can be operated in an environment between 5°C and 40°C, and the measurement accuracy of the moisture rapid tester can be guaranteed under extreme temperatures.
- Relative humidity: 50%-75%, no condensation at 30°C.
- Warm-up time: wait at least 30 minutes after connecting the power supply to the moisture rapid tester. When it comes back from the standby state, the rapid moisture tester can be used immediately.

2. Power supply

- Input voltage: 220VAC, 3A, 50Hz±1Hz
- Voltage fluctuation: 220VAC±10%
- Power load: the maximum power is 350W during heating.
- Protection: protect from dust and moisture
- Pollution degree: Level 2

Code	KCD001
Readability	0.001g
Maximum Weighing Value	100g
Readability Of Moisture Content Determination	0.01%
Accuracy Of Moisture Measurement	±0.5% (sample≥3g)
Temperature Control Adjustment Tolerance	±1°C
Heating Temperature Range Setting	50°C- 180°C
Heating Time Range Setting	0-99min
End Point Control	Timed automatic
Sample Drying Parameters	9 types (see later print for details)
Moisture Content Measurement Range	0-100%
Communication	Serial port RS232
Weighing Pan Size	Φ100mm
Dimensions (L×W×H)	350mm×185mm×140mm
Net Weight	3.0kg
Gross Weight	4.5kg
Heating Source	High-efficiency halogen lamp

5. INSTALLATION

1. The choice of placement

■ Keep the surface of the console stable and level when operating the rapid moisture analyzer. Choose a safe and adequately ventilated location. Samples with corrosive or toxic fumes and other hazardous materials need to be specially prepared for placement. Please try to avoid the rapid temperature fluctuations, excessive humidity, air flow, vibration, electromagnetic fields, heat, or direct sunlight.

2. Installation of components

- Place a windshield
- Place the sample tray holder and rotate until it is positioned
- Place the stainless-steel sample pan flat on the weighing pan support

3. Connect power

■ Plug the correct end of the provided power cord into the power input slot on the back of the moisture rapid tester, and then connect the other end to the power output socket.

Note: In order to obtain the best test results, please use it after power-on for at least 30 minutes.

1. Adjustment level

■ The device is equipped with a level meter and two leveling feet to compensate for the slight unevenness on the weighing operation table. When the blister is placed in the center of the level meter, the instrument is completely level.

Note: Every time the position changes, the level must be readjusted

6. OPERATION

6.1. Matters needing attention (related to measurement accuracy)

- As a precision analysis instrument, the moisture meter should be operated carefully by the user, and the sample to be tested should be spread evenly on the weighing pan.
- Typical samples should be 3g-5g, and the tare weight should be removed before measurement.
- Please set up the heating parameters (heating temperature, end method) before weighing the sample.
- In order to ensure the accuracy of the test, the sample should be in powder state, and the large particle sample should be crushed with a grinder first.
- In order to measure correctly, wait for the moisture meter to cool down to 50□ before using it for the second time.

6.2. Switch of moisture meter

In the shutdown state, touch the screen to start the moisture meter, and perform the self-test inside the instrument, and enter the main interface after the end.

6.3. Operation of the moisture meter

Weight calibration

■ The moisture meter can be calibrated with external weights. In order to improve and ensure the accuracy of the weighing data, the weight calibration should be carried out before the first weighing.

The weight calibration steps are as follows:

1. Please open the top cover when calibrating.
2. Clean the sample pan, the sample pan should be placed in place, touch the weight calibration key on the main interface to enter the calibration interface.
3. Touch the calibration key to start weight calibration.
4. Put the touch screen prompt to load the weight on the sample pan and close the cover.
5. It shows that the calibration is finished, please remove the weight, and take off the weight. (If it does not return to zero, then repeat the above operation and weight calibration).

Clock calibration

■ When the time and date are deviated and need to be reset, the steps for setting the time and date are as follows:

1. Touch the clock calibration key on the main interface to enter the clock interface.
2. After setting the current date and time, touch the confirm button to return to the clock interface. The setting is successful

Timing drying

Timed drying means that the drying ends when the preset drying time is reached.

1. On the main interface, touch the timing drying key to enter the timing drying interface.