# DISCO CROMÁTICO DE NEWTON

## NEWTON COLOUR DISC

REF. QLG001

### INDEX OF LANGUAGES-ÍNDICE DE IDIOMAS

ESPAÑOL.....2

ENGLISH.....3

## Objetivo

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## Purpose

It is used to demonstrate that white light is a resultant of additive mixing of seven different colours

## Description

It comprises of a multi-coloured circular disc having sectors of various spectral colours in proper proportions for demonstrating that white light is composed of all the spectral colours through the additive mixing of the colours. The disc has a toothed wheel at the back and gear with the driving toothed wheel with a handle. On rotating the disc at sufficient speed, all the colours will be observed to disappear by merging into a single white colour. Complete assembly mounted on a stable base

## Primary and secondary colours

When white light is dispersed by a prism seven colours (red, orange, yellow, green, blue, indigo and violet) are seen on white screen. However, only three coloured lights need to be combined to give white light. These are red, green and blue mixed in the correct intensities.

On using three ray boxes fitted with red, green and blue filter, the coloured rays can be made to overlap. Where all the three colours overlap, white light is seen. When just red and green lights overlap, yellow light is seen. On mixing green and blue light cyan colour is seen. Similarly, on mixing blue and red light magenta colour is seen.

Red, green and blue colours are called **Primary colours** because they cannot be made by the combination of the other colours, they can be combined to form all the other colours and all three combined together give white light.

Yellow, cyan and magenta, the three colours formed by the combination of various pairs of primary colours are known as **Secondary colours**. When a secondary colour is combined with the remaining primary colour white light is seen.

Yellow (red + green) + blue = white

Cyan (green + blue) + red = white

Magenta (red + blue) + green = white

Colours which add together to give white light are called **Complementary colours**.

## Procedure

1. Place the Newton colour disc on a balanced table.
2. Hold the apparatus by just pressing the base of the apparatus with one hand.
3. Rotate the driving wheel in clockwise direction with the help of the handle fitted with the wheel.
4. On rotating the disc at a sufficient speed, all colours will be observed to disappear by merging into de single white colour.