**MICROPIPETAS AUTOCLAVABLES
AUTOCLAVABLE MICROPIPETTES
MICROPIPETTES AUTOCLAVABLES**

DIGIPETTE

ES	MANUAL DE INSTRUCCIONES	02
EN	INSTRUCTION MANUAL	08
FR	MANUEL D'UTILISATION	14

DESCRIPCIÓN

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DESCRIPTION

Digipette pipettes are available as fixed or variable, general purpose micropipettes for sampling and dispensing accurate amounts of liquid in micro volumes.

These pipettes operate on the air displacement principle and use disposable tips at any standard make.

All pipettes are equipped with a built-in tip ejector.*Except 10 mL pipette MGB009

They cover the volume range at 0.1ul-10000ul.

variable volume pipettes

<i>Reference</i>	<i>Volume</i>
MGB001	0,1-2.5 µL
MGB002	0.5-10 µL
MGB003	2-20 µL
MGB004	5-50 µL
MGB005	10-100 µL
MGB006	20-200 µL
MGB007	100-1000 µL
MGB008	500-5000 µL
MGB009	500-10000 µL*

fix volume pipettes

<i>Reference</i>	<i>Volume</i>
MGD001	10 µL
MGD002	20 µL
MGD003	25 µL
MGD004	50 µL
MGD005	100 µL
MGD006	200 µL
MGD007	500 µL
MGD008	1000 µL

multi (8) channel pipettes

<i>Reference</i>	<i>Volume</i>
MGG001	0.5-10 µL
MGG002	5-50 LµL
MGG003	10-100 µL
MGG004	20-200 µL
MGG005	30-300 µL

multi (12) channel pipettes


<i>Reference</i>	<i>Volume</i>
MGG006	0.5-10 µL
MGG007	5-50 µL
MGG008	10-100 µL
MGG009	20-200 µL
MGG010	30-300 µL

Note: special volume pipettes other than specified above may be supplied on demand.

** Without tip ejector; manual tip ejection*


PIPETTING TECHNIQUES

A. Forward Pipetting

1. Press the operating button to the first stop.
2. Dip the tip attached to the pipette into the solution to a depth of about 2-3 mm, and slowly release the operating button. Wait for a while, and then withdraw it from the liquid touching it against the edge at the reservoir to remove excess liquid adhering to the outer surface of the tip.
3. Dispense the liquid into the receiving vessel by gently pressing the operating button to the first stop. After a second, press the operating button to the second stop. This will empty the tip completely. Remove the pipette from the vessel sliding it up the wall of other vessel.
4. Release the operating button to the ready position.

B. Reverse Pipetting

1. Press the operating button to the second stop.
2. Dip the tip attached to the pipette into the solution to a depth of about 2-3 mm, and slowly release the operating button. This action will fill the tip with a volume that is larger than the set volume. Wait 1-2 seconds and withdraw the tip from the liquid, touching it against the edge of the reservoir to remove excess liquid.
3. Dispense the liquid into the receiving vessel by pressing the operating button gently and steadily to the first stop. This volume is equal to the set volume. Hold button in this position. Some liquid will remain in the tip, which should not be dispensed.
4. The liquid remaining in the tip can be dispensed back into the original solution by pressing the button to the second stop or disposed together with the tip.
5. Release the operating button to the ready position.

Note : *Reverse pipetting technique is recommended for viscous solutions, solutions having tendency to foam or for dispensing very small volumes.*

PIPETTING RECOMMENDATIONS

- ◆ Aspirate liquid into the pipette only when a tip is attached to its tip cone.
- ◆ Allow liquids, tips, and pipettes to equilibrate to the ambient temperature.
- ◆ Pre-rinsing of tip 5 times with the liquid to be dispensed is recommended. This is important especially when dispensing liquids which have a viscosity and density different from water.
- ◆ Pre-rinse the tip several times before use when pipetting liquids at temperature different from ambient.
- ◆ While pipetting, the pipette should be vertically straight and tip should be dipped only a few millimeters into the liquid.
- ◆ Always control the push button movements with the thumb for consistency.
- ◆ Wipe the tip only if there is liquid on the outside of the tip, being careful to avoid touching the tip's orifice.
- ◆ Don't keep pipette in your hand while not working, to avoid transferring of body heat
- ◆ Use the correct pipette tip designed for use with the particular pipette.
- ◆ Select the correct pipetting technique (e.g. Reverse, forward etc.) depending on the nature of the liquid.
- ◆ Using excessive force to turn the push button outside the range specified for it may jam the mechanism and damage the pipette.
- ◆ Use always filter tips when handling aggressive reagents

MAINTENANCE

- ◆ To maintain the best results from your pipette, each unit should be checked every day for cleanliness. Particular attention should be paid to the tip cone(s).
- ◆ This pipette has been designed for easy in-house service. However, we also provide complete repair service. Please return your pipette to your local distributor for repair by Digipette Technical Service. Before returning, please make sure that it is free from any contamination.
- ◆ Check the performance of your pipette regularly e.g every 3 months and after every in-house service or maintenance.

Cleaning Your Pipette

- ◆ To clean your pipette, use ethanol and a soft cloth or lint-free tissue. It is recommended to clean the tip cone regularly.

Perform cleaning and lubrication of internal parts at regular intervals to keep smooth pipetting operation.

In-House maintenance

1. Hold down the tip ejector.
2. Place the tooth of the opening tool between the tip ejector and the tip ejector collar to release the locking mechanism.

3. Carefully release the tip ejector and remove the ejector collar.

4. Place the wrench end of the opening tool over the tip cone and turn it anti-clockwise

5. After removing the tip cone, wipe the piston, the O-ring and the tip cone with ethanol and a lint-free cloth.

6. Before replacing tip cone, it is recommended to grease the piston slightly using the silicone grease provided.

7. After reassembling, use the pipette (without liquid) several times to make sure that the grease is spread evenly.

AUTOCLAVING

This entire pipette can be autoclaved (sterilized) by steam at 121 °C for 20 minutes at 1 atm pressure. After autoclaving, let the pipette cool down and dry for 6 hours before use.

Verify pipette calibration after each autoclaving cycle.

CALIBRATION

Each pipette has been checked & calibrated at factory with procedure conforming to EN-ISO 8656 standards. It is recommended to check the calibration at least once a year, for regularly used pipette.

Can consult the verification procedure in our technical support <https://www.auxilab.es/en/technical-information/>

Checking calibration

1. Fit new tip onto the tip cone.

2. Pipette distilled water into pre-weighed beaker at least five times & record the weight each time.

Check the pipette at maximum volume, at 50% of maximum volume and at minimum value or 10% of maximum volume, whichever is bigger.

3. Compare the results with the permitted limits given in the calibration report or with the variation limits set by your own laboratory.

Important Notes

1. Procedure should take place at 20° (± 0.5 C) constant temperature.
2. The weighing beaker, distilled water, pipette & tips must be at the same temperature.
3. Use an analytical balance with 0.01 mgs readability.
4. Pre-rinse the tip 3 to 5 times before pipetting.
5. Divide the weight of the water by its density (at 20°C , 0.9982) to get the volume.

Recalibration

1. Place the service tool into grooves at the base of the push button as shown below. Turn it clockwise to increase & anticlockwise to decrease the volume.
2. Repeat the 'Checking calibration' procedure.

STORAGE

- ◆ When not in use, it is recommended that your pipette is stored in a vertical position.
- ◆ Leaving the pipette on its side can cause liquids to leak into the body of the pipette and cause corrosion.

TROUBLE SHOOTING

1. Droplets left inside the tip:
 - Unsuitable tip, non-uniform wetting of the plastic: **use new tip.**
2. Leakage or Pipetted volume too small:
 - Tip holder (cone) scratched or damaged: **change the tip cone.**
 - Organic solvent as liquid: **aspirate and discard the organic solvent several times before actual pipetting by the same tip.**
 - Tip incorrectly attached: **attach firmly.**
 - Unsuitable tip: **use new tip.**
 - Foreign particles between tip and tip cone: **clean the tip cone.**

—Insufficient amount of grease on piston and O-ring: **clean and grease O-ring.**

—O-Ring not correctly positioned or damaged: **change the O-ring.**

3. Inaccuracies:

—Incorrect operation: **follow instructions carefully.**

—Calibration altered: **recalibrate according to instructions.**

—Unsuitable for the particular liquid pipetting technique: **use correct pipetting technique.**

—Instrument damaged: **send for repair.**

4. Push button jammed or moves erratically:

—Piston contaminated or penetration of solvent vapors: **clean and grease O-ring and piston.**

5. Tip ejector jammed or moves erratically:

—Tip cone contaminated from outside: **Remove ejector collar and clean tip cone's outer surface with ethanol.**

6. Volume setting is not properly click stopped:

—Click stop mechanism damaged: **Send for repair.**

7. Push button does not turn for volume setting:

—Use of excessive force beyond the range of pipette : **Send for repair.**

WARRANTY 3 YEARS

These pipettes are warranted to be free from defects in material and workmanship for a period of three years from the date of purchase. Your product will be duly repaired upon prompt notification in compliance with the following conditions:

◆ Defects or damages caused by physical & chemical abuse or normal wear and tear or resulting from the pipette being used in any manner other than as instructed in the manual are not covered by the warranty.

◆ The warranty is invalidated by non-factory modification, which will immediately terminate all liabilities on us for the product or damages caused by its use.

◆ The buyer & its customer shall be responsible for the product or use of products as well as any supervision required for safety.

◆ If requested, the product must be returned to us or the distributor in well packed & insured manner and all shipping charges must be paid.

◆ Routine cleaning and recalibration are not covered under the terms of warranty.

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