PLATINA TERMOSTATIZABLE HEATING STAGE

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment. Please inquire at the community administration for the authorized disposal location.

- ◆ When you are not using the equipment for a long period of time please make sure it is unplugged from the net in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the equipment by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the equipment (see Maintenance).

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The equipment is supplied with a standard Schuko wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the equipment by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the stage and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

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
Figura 2

- ◆ Connect the heating stage to the temperature controller by inserting the stage wire (1.2) in the corresponding plug (2.4) of the controller.
- ◆ The equipment is supplied with a standard Schuko wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz \pm 10% provided with earth wire and the other end to the temperature controller connector (2.8).

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz \pm 10%.

Note: The tension can be changed to 110 V, 60 Hz \pm 10%. For this purpose change the position of the tension selector placed at the back of the temperature controller (2.6).

Turn the temperature controller on by pressing the On/Off switch (2.5) at the back. When turning it on, the upper display (2.1) will show the temperature of stage at that moment, while lower display (2.2) will show the last programmed working temperature.

- ◆ Select the working temperature by pressing buttons for temperature selection (2.3). Please, bear in mind that the working temperature range is from RT up to 50 °C. The selected working temperature will be shown on the corresponding display (2.2).
- ◆ Once selected the working temperature, the stage will start heating. Real temperature of the stage will be shown on the upper display (2.1).
- ◆ When real temperature of the stage reaches programmed temperature, then the observation of samples can start.
- ◆ Once observation is finished, turn the controller off and disconnect it from the power supply.
- ◆ Before removing the heating stage wait until it cools.
- ◆ When you are not using the heating stage for a long period of time, please make sure it is disconnect it from the net and protected from dust (avoiding possible accidents and prolonging its life time).

Security

- ◆ The equipment must be used by previously qualified staff that knows how to operate it thanks to the user manual.
- ◆ Do not put the equipment near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid the presence of flammable or toxic products in the working area.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the working area.

3. TECHNICAL SPECIFICATIONS

| | |
|-------------------------------|------------------------------------|
| <i>Temperature range</i> | RT up to 50 °C |
| <i>Consumption</i> | 50 W |
| <i>Power</i> | AC, 220 V / 50 Hz or 110 V / 60 Hz |
| <i>Stage dimensions (LxW)</i> | 130x150 mm |

4. INSTALLATION / SETTING UP

Preliminary inspection

◆ Unwrap the equipment, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all the protective items and, without connecting the equipment to the net, make sure that it does not present any damage because of the shipment. In case the equipment presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Heating stage with connection wire to temperature controller
- Temperature controller
- Schuko wire
- Fuse 3 A
- Roll of adhesive tape
- User manual

We will only accept the equipment return within 15 days after delivery provided it comes in its original wrapping.

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ◆ Previously clear the microscope stage by removing all accessories as clamps, etc.
- ◆ Put the heating stage on the microscope stage, so as central holes of both stages match (the thicker side of heating stage must be placed at the front side). At the bottom, the heating stage has two strips of adhesive tape that will help to fix it on the microscope stage.

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2. DESCRIPTION

- 1.1 Heating stage
- 1.2 Connection wire of stage to temperature controller.
- 1.3 Temperature controller


Figure 1

- 2.1 Display of stage real temperature
- 2.2 Display of programmed working temperature
- 2.3 Buttons to select working temperature
- 2.4 Connection plug
- 2.5 On/Off switch
- 2.6 Tension selector
- 2.7 Fuse holder
- 2.8 Socket


Figure 2

INDEX OF LANGUAGES

| | |
|------------------|--------|
| 1. Spanish | 2 - 8 |
| 2. English | 9 - 15 |

INDEX OF CONTENTS

| | |
|------------------------------------|----|
| 1. USES OF THE INSTRUMENT | 10 |
| 2. DESCRIPTION | 11 |
| 3. TECHNICAL SPECIFICATIONS | 12 |
| 4. INSTALLATION / SETTING UP | 12 |
| 5. MAINTENANCE AND CLEANING | 14 |

1. USES OF THE INSTRUMENT

Zuzi heating stage allows keeping the samples to be observed at a constant and specific temperature. Thus, it is very useful for observation of biological samples (usually at 37°C) as for example semen, blood, etc. in both medical and veterinary laboratories.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi Heating stage. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Heating stage, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the stage avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the heating stage.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler.

This equipment is protected under the Warranties and consumer goods regulation (10/2003).

- ◆ Overhaul is not covered by the product's warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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