BALANZAS DE BOLSILLO POCKET SCALE BALANCE DE POCHE

ES

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EN

Thanks for choosing our Nahita-Blue pocket scale which is easy to operate, please read the user manual carefully before you use it.

1. Specification parameter

Power control: Auto power off after 3 minutes without operation (Adjustable)

Working temperature: 0-30 °C

LCD Display: 5 digital LCD display

Backlight: Blue

Power supply: 2xAAA batteries

Scale Size: 127X77.5X18.5mm

Net weight: 118g

Salver Size: 69X64mm

Tare

Counting

Lower power alert

Unit switch: g/oz/ozt/dwt/ct/l/gn

2.Operation

Put 2xAAA battery into the battery room of the scale, place the scale on a solid and flat board without wind or vibrancy.

1. Press the button "⏻", the screen displays "0.0" or "0.00" in 3 seconds

2. Put your good on the salver. The unit code is twinkling when it has not been steady; You can read the weight when the unit code stop blinking.

3. Press the button "M" to switch unit: g/ oz/ ozt/ dwt/ ct/l/gn;

4. Press the button "T" to tare the scale.

5. Press the button "T" for 3 seconds to operate counting (See the detail under 8.How to operate Counting Function)

6. Press the button "⏻" to turn off the scale at working status.

7. Calibrating Operation

Press the button "⏻" for 4 seconds till the screen displays "0", then release till screen displays a weight (for example:500.0; if you do not have right calibration weights, you can press the button "M" to select a weight). Please put right calibration weights on the salver, you can see the data flashing on the screen and then shows the value of the present calibration weight. The calibration is completed. (You should calibrate the scale when the temperature changes or other reasons).

8. How to operate Counting Function

8.1 Please prepare 25, 50, 75 or 100 pieces of samples. When the scale is ready to weight goods

8.2 Place the samples on the salver, then you can see the weight. Please press the button "T" for 3 seconds, the screen shows the quantity P25PCS (P50, P75, P100), you can press the button "M" to select a right quantity-25, 50, 75 or 100 pieces,

8.3 Finally, press the button "T" for confirmation. At this time, the screen displays the quantity. You can press the button "T" to come to weighing status, then you can weight more goods. Press the button "T" to read quantity when you need. Go to 8.1 if need to do re-counting.

9. "⏻" is a Multi-function button
- Press "⏻" to switch on the scale.
- Press "⏻" again to power off the scale.
- Press "⏻" 4 seconds and display 0 to do the calibration.

3. Notice

1. Place the scale on a solid and flat board without wind and vibrancy.
2. Do not over load, otherwise the sensor will be damaged; Do not bump the salver when you put the goods on the salver, put it carefully.
3. Keep away from mobile phones, $\geq 50\text{Cm}$, The mobiles will affect the accuracy of the scale.
4. Read the weight after the unit code stop blinking.
5. Replace batteries when it shows low power.
6. Working temperature: 10-30 °C. Please calibrate the scale if the temperature is too low/high. Please place the scale for 30 minutes first, let the temperate is even, then make the calibration.
7. The screen display "FULL" when it is over weight.
8. Please do not disassemble! Careful storage!

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