**MANÓMETRO EN U**  
**MANOMETER U**

**REF : QLB007**

**INDEX OF LANGUAGES-ÍNDICE DE IDIOMAS**

ESPAÑOL ..... 2  
ENGLISH..... 4

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## Descripción

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## MANOMETER U

### Description

Crystal manometer with U form and a key of step on the top of one way. It has a support. The scale is 50-0-50 with 1mm divisions.

Detail of the key of step

### Theory

Manometer is an instrument useful to measure the pressure in fluids, it is often used to measure the difference between two pressures taking as reference the local pressure.

Pressure is defined as the press per unit of area that makes a liquid or gas on a surface. Pressure can be measured as atmosphere (atm), newton per metre, pascal (Pa).

$$1 \text{ atm} = 101.325 \text{ Pa} = 760 \text{ mm Hg.}$$

We have to consider that the majority of manometer measure the difference between a liquid and the atmospheric pressure, thus we have to add the result to the atmospheric pressure to achieve our result.

This is because the manometer has two columns, one is connected to our gas that we want to measure the pressure and the other is opened. In both cases there is a calibration to help us in the process.


In the figure, P1 is the pressure that a gas makes over the water because this one is connected, the other column is P1 this is the atmospheric pressure, the difference Ah between the two columns is that we have to add to the atmospheric pressure to measure the pressure of the left column.

The unit of this manometer is the meter of water column. Meter of water column is the unit of pressure equal to the pressure that is exerted by a column of water on earth gravity. It symbol is mca and is a multiple of mmca.

$$1 \text{ mca} = 100 \text{ cmca} = 1000 \text{ mmca}$$

$$1 \text{ mca} = 9896,65 \text{ Pa}$$

$$1 \text{ atm} = 10,33 \text{ mca}$$

### Process

1. Put water into the manometer to measure the pressure
2. Connect the left column to a gas or vacuum pump with a specific pressure to measure
3. Add the difference  $\Delta h$  to the atmospheric pressure
4. Change and study some pressure units