Refractive index correction chart with distilled water temperature

Temperature °C	Refraction Index	n _F -n _C Mean dispersion
10	1.33369	0.00600
11	1.33364	0.00600
12	1.33358	0.00599
13	1.33352	0.00599
14	1.33346	0.00599
15	1.33339	0.00599
16	1.33331	0.00598
17	1.33324	0.00598
18	1.33316	0.00598
19	1.33307	0.00597
20	1.33299	0.00597
21	1.33290	0.00597
22	1.33280	0.00597
23	1.33271	0.00596
24	1.33261	0.00596
25	1.33250	0.00596
26	1.33240	0.00596
27	1.33229	0.00595
28	1.33217	0.00595
29	1.33206	0.00594
30	1.33194	0.00594
31	1.33182	0.00594
32	1.33170	0.00593
33	1.33157	0.00593
34	1.33144	0.00593
35	1.33131	0.00592
36	1.33117	0.00592
37	1.33104	0.00591
38	1.33090	0.00591
39	1.33075	0.00591
40	1.33061	0.00590

**INSTRUCTIONS ON ENVIRONMENT PROTECTION**

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non-adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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
- ◆ When observing through the eyepiece, search for the divisory line between light and dark colored fields and do as the instructions in the previous case.

c) To test translucent solids

- ◆ The object to be measured must be smooth and polished faced.
- ◆ Put one or two naphthalene solid and clean solid surface to place it on the refracting prism. Open the reflecting mirror (1.6) and adjust to most appropriate angle, using it as incidence surface for measuring.
- ◆ Reading meted is the same as previous cases.

Picture 4

d) To measure sugar concentration in sugary liquids

- ◆ The procedure is similar to liquid samples, but the reading can be done on a Brix scale (larger scale 0- 95%) reading directly the sugar concentration in the liquid.

e) Refractive index estimation at different temperatures

- ◆ Insert the thermometer into its socket (1.7), connect the thermostat to the prepared connectors (1.8) and (1.13) for bottom prism and (1.12) for top prism, thanks to the thermostat temperature can be read in the ranges required. Before each reading, wait a couple of minutes to stable the temperature.

6. MAINTENANCE

- ◆ To ensure accuracy and to avoid instrument damage, please, pay attention to your refractometer maintenance.
- ◆ Follow the instructions and warnings relating to this manual.
- ◆ You should always have this manual by hand, so as any user can read it.
- ◆ Preserve the refractometer from heavy shocks. Look for a dry place away from dust, direct sunlight, or air flows.
- ◆ Since it is a precision instrument it must be handled carefully.
- ◆ It must be kept in a dry and well-ventilated room to preserve optical parts from humidity and mold.
- ◆ If you deal with corrosive liquids, you must clean all the instruments immediately to prevent them from rust damages you should include in this list optical and mechanical parties and also painted surfaces.
- ◆ After each use, you must clean the refractive prisms with a non-fuzz gauze, use alcohol or water.
- ◆ Liquid samples must not contain solid impurities. Solid samples must have polished and flat faces to avoid scratching the prism.
- ◆ When not using the refractometer for long periods, keep it in a box with a draining substance (Silica Gel).
- ◆ Never remove any of its mechanical components. In addition to losing the warranty, it could also lead to a faulty performance of the equipment.
- ◆ For any inquiry, please contact Zuzi official service.

Installation

- ◆ Please read carefully this manual before operating with this equipment.

BEFORE CONNECTING OR OPERATING WITH THIS INSTRUMENT, READ THESE INSTRUCTIONS THOROUGHLY.

- ◆ Place the refractometer in a flat horizontal table, avoiding humidity, heat, dust, and vibrations.

Calibration

- ◆ Before using the first time, you must check if the refractometer is well calibrated. For doing this, a known refraction index ($n_D = 1.5168$) test piece is included.
- ◆ If you put two naphthalene bromide drops on refraction prism and you place then the smooth test piece item over them, you can see through the eyefield, the line that separates light and dark colored parts. Once the dispersion is corrected, you agree inferior side scale with the amount in the test piece. If the line that divides light and dark parties do not meet in the center of the crossed lines seen in the top, you need to turn the calibration screw until you have them met.
- ◆ Once the refractometer is calibrated, the screw does not be manipulated while measuring. If you doubt about the refraction index estimation when testing, please, calibrate it again from the beginning.
- ◆ The ridged incident prism face, the smooth refraction prism face and test piece must be cleaned with the aid of an absorbent cotton gauze wet with a 1:1 absolute alcohol and ethermixture in order to avoid from the dirtiness that have an effect on image brightness and measurement accuracy.

Testing work***a) To test transparent and translucent liquid***

- ◆ Drop the tested liquid onto prism surface using a dropping tube. Cover the incident prism and lock with the knob (1.3). The sample must be homogeneous and spread prism surface wide without bubbles.
- ◆ Open the shield (1.4) and check the reflecting mirror (1.6) is closed, then adjust the ocular (1.1) based on the dioptric correction required until you see the crossed lines sharply.
- ◆ Turning the refraction index adjustment knob (1.11) look for light and dark field divisory line, then turn knob (1.10) to adjust the dispersion in a way this line is seen without any color and as clear as possible. Once you have done this, turn the knob (1.11) again to move divisory line up to the middle of both crossed lines.
- ◆ The value in the bottom scale tells us the accurate liquid refraction index. If the scale is not lighted enough, turn the condenser (1.9) towards light focus.

b) To test transparent solids

- ◆ Solid to be measured must be smooth and polished faced.
- ◆ Put one or two naphthalene bromide drops on the refraction prism surface. Clean all the tested object surface properly to get the best contact.

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

Prism refractive angle ϕ and prism refractive index n_2 , are calculated when knowing the critical angle i . Sample refraction index can be known.

4. TECHNICAL SPECIFICATIONS

Measurement range	
Refractive index	1.3000 to 1.7000 nD
^a Brix	0 to 95%
Measurement accuracy	
Refractive index	± 0.0005 nD
^a Brix	± 0.25%
Illumination	By natural light
Thermometer range	-50/70°C
Working temperature	0 – 80°C
Outside dimensions	270x170x395 mm
Weight	4.3 Kg

5. OPERATING INSTRUCTIONS

Preliminary examination

- ◆ Remove refractometer from the packing material and examine it carefully to make sure no damage has occurred during shipping. If there is any damage, notify it to your dealer/forwarder agency.
- ◆ Check all the accessories that must be included with the equipment:
 - Digital thermometer up to 70°C
 - Screwdriver
 - Test piece with known refractive index
 - Instruction manual
 - Solution of monobromide naphthalene
- ◆ Keep all packing material; any return of the equipment must be done in its original packing and with all the accessories supplied.

3. WORKING PRINCIPLES

Refractometry

◆ When a light beam arrives to the divisory surface between two different mediums, one part is reflected and the other is bended (it comes into the second medium), being spread with an angle and speed different to first medium. Both values are measured by Refractive index. Each substance has its own refractive index value.

◆ Refractive index is related to mass, charge and the number of particles in the substance through which the light beam is transmitted. There are different kinds of refractometers, but they all are based on the same principle.

◆ The basic refractometer principle is based on the following refraction law:

$$n_1 \times \text{sen} \sigma' = n_2 \times \text{sen} \sigma''$$

Being n_1 and n_2 medium refractive index at both limit sides between different sides medium, where σ' is the incident angle and σ'' refractive angle.

◆ When a light beam enters an optically lighter medium from a thicker one, the incident angle is smaller than refractive angle and if we change the first one until the second one reaches 90° , the incident angle is named critical angle.



◆ As shown in picture 3, when a light beam falls on surface AB from any angle, its refractive angle is always bigger than i , if we observe the emergent light beam through a microscope, it appears the view field divided in two halves, one dark and the other one light colored. If we place the two halves separation line just in the middle, this is the critical angle position.

◆ $ABDC$ is the refracting prism whose refractive index is n_2 , and the sample to be measured (liquid or transparent solid) is placed on AB surface. Its refractive index is n_1 .

◆ From refraction law, we obtain the following equation:

$$\begin{aligned} n_1 \times \text{sen} 90^\circ &= n_2 \times \text{sen} \sigma \\ n_2 \times \text{sen} \beta &= \text{sen} i \\ \phi &= \sigma + \beta \\ \sigma &= \phi - \beta \end{aligned}$$

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1. APPLICATIONS OF THE INSTRUMENT


The refractometer is used to measure the refractive index n_D and mean dispersions $n_F - n_C$ of transparent and translucent liquid or solid, mainly the transparent liquid. Zuzi 325 Refractometer includes a thermostat to measure refractive indexes n_D within the range of 0°C and 70°C. See joint chart. Refractive indexes are important optical constants in substance and can be used to determine optical performance, purity and concentration, etc.

It is possible to get sugar percentage in saccharose solutions within a range of 0 - 95% Brix, equal to refractive index 1.333 - 1.531.

This instrument can be widely used as indispensable tool in petrological, oil and fat, pharmaceutical, food, cosmetic or chemical industry, as well as universities and research institutions.

2. DESCRIPTION

Picture 1

- | | |
|---|--|
| 1.1. Eyepiece | 1.9. Scale illuminator condenser |
| 1.2. Calibration screw | 1.10. Scale knob for dispersion graduation |
| 1.3. Close prism control | 1.11. Refraction index adjustment knob |
| 1.4. Incident illumination window | 1.12. Top prism circuitation thermostat connections |
| 1.5. Prisms | 1.13. Bottom prism circuitation thermostat connections |
| 1.6. Refraction prism mirror | |
| 1.7. Thermometer socket | |
| 1.8. Bottom prism circuitation thermostat | |

Thanks for choosing a Abbe Refractometer ZUZI. It is advisable to read carefully this manual and follow the processes of use.

ZUZI develops its products putting more emphasis on Ergonomics and Security for users. Its high mechanical and optical quality will allow you to have this equipment available for many years.

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In order to use this equipment, bear in mind the following points:

- ◆ This manual should be available for all the users of this equipment.
- ◆ As an accuracy instrument, it should be used carefully avoiding the sudden movements and knocks.
- ◆ Never dismantle the different pieces of the instrument to repair it by yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler.

ATTENTION!!

IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

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