Table 2

S. No.	Magnitude $\vec{A}$ (N)	Magnitude $\vec{B}$ (N)	Angle $\vec{B}$ ( $^{\circ}$ )	Magnitude $\vec{C}$ (N)	Angle $\vec{C}$ ( $^{\circ}$ )

Table 3

S. No.	Experimental magnitude $\vec{C}$ (N)	Analytical magnitude $\vec{C}$ (N)	Difference (%)	Experimental angle $\vec{C}$ ( $^{\circ}$ )	Analytical angle $\vec{C}$ ( $^{\circ}$ )	Difference (%)

Precaution: The pulleys do not have to be attached extremely tight, only tight enough to keep them from falling off. If they are too tight, they can damage the force table. Do not put the ring on the post or the string on the pulleys at this time!

**Note:** Similarly we can work with a system of four forces in equilibrium.

**Instructions on environment protection**

- At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances.
- It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.
- The materials are recyclable as mentioned in its marking.
- By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.
- Please inquire at the community administration for the authorized disposal location.

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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**Procedure**

**Balance Forces**

1. On the force table, set up two pulleys and hanger with same weigh, facing opposite directions (180° difference in angle).
2. The force of each will be equal to  $F_{grav} = m \times g$ . Check whether the two forces are equal and opposite by examining the ring at the center of the force table, which should not move.
3. Notice that if the components of the vectors associated with these forces are added, the resultant vector will have zero magnitude. This is how to determine that all forces are in equilibrium.
4. Now add little mass on one hanger and check equilibrium.
5. Find the maximum mass one can place in hanger 2 and still maintain experimental equilibrium.

**Analytical Calculation**

1. This experiment will consist of three forces in equilibrium. Two forces will be known, while the third will be found first analytically, using the theory of vectors, and then experimentally. For this, keep  $\vec{A}$  at 0° for the duration.
2. Note that if  $\vec{A}$  and  $\vec{B}$  are known and  $\vec{C}$  when added to the system, causes the two forces to be in equilibrium, then  $\vec{C}$  is of equal magnitude but in the opposite direction to the sum ( $\vec{A} + \vec{B}$ ).
3. Calculate the magnitude of  $\vec{A}$  and  $\vec{B}$ . Use the fact that  $F_{grav}=m \times g$  and that 1 Newton is a unit of force equal to  $kg \times m \times s^{-2}$
4. Using the theory of vectors, calculate what magnitude  $\vec{C}$  would be if it was the sum ( $\vec{A} + \vec{B}$ ).
5. Using the theory of vectors, calculate what angle  $\vec{C}$  would be if it was the sum ( $\vec{A} + \vec{B}$ ).

**Experiment**


1. Set up the two forces on the force table keeping  $\vec{A}$  at 0° and  $\vec{B}$  at 120° and place equal mass on both hangers.
2. Place a third pulley at some angle at which you expect the system to balance
3. Set up the third force  $\vec{C}$  by adding weights and changing the angle until equilibrium is reached.
4. Repeat the step for different masses.
5. Determine the percent difference from the analytical result by calculating

$$\frac{\text{experimental result}}{\text{theoretical result}} \times 100$$

**5. Calculation**

**Table 1**

S. No.	A		B	
	Mass	Angle	Mass	Angle

## 4. Experiment to study equilibrium of forces

### Components required

- Force Table
- Pulley
- Weight hanger
- Weight set
- String
- Ring

### Other components required


- Scissors
- Level
- Graph paper
- Pencil
- Ruler
- Protractor
- Calculator

### Set up

1. Attach the three legs into the brackets on the underside of the tabletop by screwing.
2. Push the pin into the hole in the center of the table until it is flush with the bottom of the table.
3. Check the center ring to make sure it is smooth and free of sharp edge (the string should be able to slide smoothly over the ring without catching).
4. Tied the string with loops over the ring to allow the string to slide freely around the ring.
5. Attach the pulleys to the disk at the desired locations (may need to move them during the experiment).
6. Make sure the string is short enough so that when the ring is placed on the center post, the hangers do not touch the floor.
7. Straighten the string out so they are not tangled or crossed over each other. Lay the string over the pulleys, as shown in the diagram.
8. Place the level on the tabletop and adjust the leveling feet on the legs until the table is level.

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




The components  $F_x$  and  $F_y$  in the x and y directions of the vector  $F$  are related to the magnitude  $F$  and angle  $\theta$  by:

$$F_x = F \cos \theta \text{ and } F_y = F \sin \theta$$

Conversely

$$F = \sqrt{F_x^2 + F_y^2} \text{ and } \theta = \arctan \left( \frac{F_y}{F_x} \right)$$

When several forces act on a point, their sum can be obtained according to the rules of vector algebra. Graphically, the sum of two forces  $F = F_1 + F_2$  can be found by using the parallelogram rule illustrated in Fig. 2 or, equivalently, by the head-to-tail method illustrated in Fig. 3.

The sum of the vectors can also be derived analytically by adding their components:

### 3. Condition of equilibrium

An object is in translational equilibrium when the vector sum of all the forces acting on it is zero. In this experiment we shall study the translational equilibrium of a small ring acted on by several forces on an apparatus known as a force table. This apparatus enables one to cause the forces of gravity acting on several masses ( $F = mg$ ) to be brought to bear on the small ring. These forces are adjusted until equilibrium of the ring is achieved. You will then add the forces analytically by adding their components and graphically by drawing the vectors and determining if they add to zero using the rules for the addition of force vectors listed above.

**INDEX OF LANGUAGES**

Spanish ..... 2-7  
 English ..... 8-12

**INDEX OF CONTENTS**

1. Description .....8  
 2. Addition of forces .....8  
 3. Condition of equilibrium .....9  
 4. Experiment to study equilibrium of forces 10  
 5. Calculation .....11

**1. Description**

Used for verifying the laws of composition and resolution of forces utilizing vector diagrams of several concurrent forces. Comprising a thick aluminium circular table, about 40cm in diameter and finished with scratch resistant black epoxy coating. Raised outer rim of the table has circular scale graduated 0 to 360°×1° and marked every 10°. Table mounted on a heavy, stable cast metal tripod base through a metal pipe, with base having three leveling screws. Also included are four pulleys in clamp that can be attached anywhere along the rim. A detachable pin, positioned at the center of the table indicates the balancing for forces through the centering of ring, tied to the strings carrying masses and suspended through the pulleys, around it. Supplied complete with set of strings tied to the rings, center pin and 4 sets of slotted masses with hanger, each set having one hanger 100g and masses 1 each of 10, 20, 50 and 100g.

**2. Addition of forces**

Forces are one of a group of quantities known as vectors, which are distinguished from regular numbers (known as scalars) by the fact that a vector has two quantities associated with it, a magnitude and a direction (related to a coordinate axis of the system you are dealing with). These properties completely characterize a vector. A vector may alternatively be described by specifying its vector components. In the case of the Cartesian coordinate system (the system we will be primarily dealing with) there are two components, the x-component and y-component. These two properties also completely characterize a vector. Vectors, and in the case of this lab, force vectors, can be represented pictorially by an arrow pointing in the direction of action of the force, with a length proportional to the strength (magnitude) of the force.

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