# ESPECTROFOTÓMETROS SPECTROPHOTOMETERS SPECTROPHOTOMÈTRE

## Modelos / Models/ Modèle 4251/50

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi spectrophotometer 4251/50. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual is inseparable from the Zuzi spectrophotometer model 4251/50, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the spectrophotometer avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the spectrophotometer to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the spectrophotometer.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to [asistencia@auxilab.es](mailto:asistencia@auxilab.es).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

**ENGLISH**

 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

## INDEX OF LANGUAGES

Spanish .....	2-18
English .....	19-34
Français .....	35-52

## INDEX OF CONTAINS

1. USES OF THE INSTRUMENT .....	20
2. DESCRIPTION .....	20
3. TECHNICAL SPECIFICATIONS .....	21
4. INSTALLATION / SETTING UP .....	21
5. MAINTENANCE AND CLEANING .....	33
6. TROUBLESHOOTING .....	34
ANNEX I: CE CERTIFICATE .....	52




## 1. USES OF THE INSTRUMENT

Zuzi spectrophotometer model 4251/50 is a versatile instrument that performs precise quantitative and qualitative analysis in the visible and ultraviolet spectral range. It is very useful in conventional laboratories and in a wide range of fields such as petrochemical and mining industry, agricultural engineering, environmental preservation, biochemistry, hospitals and research laboratories.


## 2. DESCRIPTION

- 1.1 LCD display and control panel
- 1.2 Sample compartment
- 1.3 Cell changer manual knob
- 1.4 Power On/Off

- 1.5 Fuse holder
- 1.6 Socket
- 1.7 USB and parallel ports
- 1.8 Manual power switch 110/220V

- 2.1 SET: selection of the different functions
- 2.2 GOTO  $\lambda$ : selection of wavelength
- 2.3 ZERO: blank
- 2.4 PRINT: printing of results
- 2.5 START/STOP: start/stop testing
- 2.6 RETURN: exit to previous screen or cancel operation
- 2.7 ENTER: confirm the input data or selected function
- 2.8 CLEAR: clear all characters when inputting data or clear test results on display
- 2.9  $\blacktriangledown/\blacktriangle$ : browse functions for selection
- 2.10 [0-9]: numeric keypad for data input or menu selection
- 2.11 [-]: input "minus" symbol
- 2.12 [.] : input dot

### 3. TECHNICAL SPECIFICATIONS

<i>Code</i>	54251050
<i>Wavelength range</i>	190-1100 nm
<i>Band width</i>	2 nm
<i>Optical system</i>	Single beam, grating 1200 lines/mm
<i>Wavelength accuracy</i>	± 0.5 nm
<i>Wavelength repeatability</i>	0.3 nm
<i>Wavelength resolution</i>	± 0.1 nm
<i>Photometric range</i>	-0.3/3 A, 0/200% T
<i>Photometric accuracy</i>	± 0.3 %T
<i>Photometric repeatability</i>	± 0.2 %T
<i>Spray light</i>	0.05% T A 220 nm, 340 nm
<i>Stability</i>	± 0.002 A/h to 500 nm
<i>Detector</i>	Silicon photodiode
<i>Sample compartment</i>	4 cells of 10 mm
<i>Lamps</i>	Tungsten and deuterium
<i>Outputs</i>	USB and parallel (printer) ports
<i>Power</i>	AC 220V/50Hz or AC110V/50HZ
<i>Dimensions</i>	480x360x160 mm
<i>Weight</i>	14 Kg

### 4. INSTALLATION / SETTING UP

#### Preliminary inspection

Unwrap the spectrophotometer, take off the involving plastic and take off the protection in which it comes fitted. Take off all the protective items and, without connecting the spectrophotometer to the net, make sure that it does not present any damage because of the shipment. In case the spectrophotometer presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- 4 glass cells
- 2 quartz cells
- Schuko wire
- User's manual
- Warranty

Compatible with MWave Professional software (Ref. 94251000), not included

**We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.**

#### Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.**

Please put the spectrophotometer on top of a horizontal, plane and stable table making a free space of at least 30 cm per side. Do not put the spectrophotometer near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid inflammable or toxic substances in the working area.

The spectrophotometer is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz  $\pm$  10% provided with earth wire and to the other end to the spectrophotometer connector.

**Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz  $\pm$  10%.**

### Setting up

- ◆ Press On/Off (1.4).
- ◆ The equipment initializes and begins an auto-test with the following steps:
  - First of all the equipment checks the position of filters and then the position of the lamps and AD converter; then, the equipment initializes the W/D2 lamp warm up, wavelength reset and Dark current; finally, after checking printer, battery and system the screen displays the main menu.
  - If data in memory have been lost, the instrument will directly calibrate the system without any user's choice.

**Note: Do not open the sample compartment during auto-test.**

### BASIC OPERATION

#### 1. Blank


- ◆ Put the cell with the blank sample into the light path.
- ◆ Press ZERO (2.3); the equipment will make the blank automatically.

**Note: 1. Blanking is automatic after every wavelength change.**

**2. Do not open the sample compartment during blanking.**

#### 2. Set wavelength (Fig. 3)

Figure 3

To set the wavelength value at any measurement mode, proceed as follows:

- ◆ Press GOTO  $\lambda$  (2.2).
- ◆ Use the numeric keypad (2.10) to input the desired wavelength value.
- ◆ Press ENTER (2.7) to save the input value and make the blank; after blanking, the instrument will return to the measurement mode main menu.

**Note: 1. If changing the wavelength is not necessary, press RETURN (2.6).**

**2. Wavelength range is 190-1100, any value out of this range will be invalid.**

### SAMPLE ANALYSIS

Zuzi spectrophotometer 4251/50 presents 3 different measurement modes: basic, quantitative and kinetics.

#### 1. Basic mode


- ◆ Put the cell with the blank sample in the light path.
- ◆ Using  $\nabla/\blacktriangle$  (2.9) select Basic mode on the main menu displayed after finishing auto-test.
- ◆ Press ENTER (2.7) to enter the basic mode display. The instrument will make the blank automatically.
- ◆ Press RETURN (2.6) to return to the main menu.

#### Select test mode (Fig. 4)

The equipment presents three measurement modes (Transmittance, Absorbance and Energy). To select the desired measurement mode proceed as follows:

- ◆ Press SET (2.1) and then select the desired measurement mode by using  $\nabla/\blacktriangle$  (2.9).
- ◆ Press ENTER (2.7) and the symbol  $\surd$  will be displayed beside the selected mode.
- ◆ Press RETURN (2.6) to return to the basic mode display and make the blank at the same time.


Figure 4

Measuring (Fig. 5)

- ◆ Press START/STOP (2.5) to enter the basic mode results display.
- ◆ Put the sample into the light path and press START/STOP (2.5) again. The test result will be displayed.
- ◆ Whenever changing the sample press START/STOP (2.5) to perform measurement and visualize the result on display.

Figure 5

**Note: 1. Up to 5 results can be displayed at the same time; a total of 200 results can be saved.**

**2. On the Basic mode results display it is possible to change wavelength [GOTOΛ(2.2)], select test mode [SET (2.1)], print results [PRINT (2.4)] and make the blank [ZERO (2.3)]. To clear all displayed data press CLEAR (2.8).**

2. Quantitative mode

To determine the concentration of different samples from a curve equation in which the absorbance and concentration are related:

$$C=K \cdot A+B$$

C: concentration

A: absorbance

K and B: curve coefficients

To create a curve, two methods can be followed:


- Coefficient: this method is used when the equation of the curve is known. First of all, the equation coefficients must be introduced and then, proceed to the analysis of the samples of unknown concentration.

- Standard curve: this method is used when the equation of the curve is unknown. First of all, the curve must be created from standard samples of known concentration and then, proceed to the analysis of the samples of unknown concentration.

Proceed as follows:

On the main menu display, select Quantitative with ▼/▲ (2.9) and press ENTER (2.7); the instrument will show the method choosing display (Fig. 6).


Figure 6

2.1 Standard curve method

- ◆ Select Standard curve and press ENTER (2.7); the pre-test interface will be displayed (Fig. 7)


Figure 7

Set wavelength

- ◆ On the pre-test interface, press GOTO λ (2.2) and input the desired wavelength value by using the numeric keypad (2.10) (Fig. 8).
- ◆ Press ENTER (2.7) to save the input value.

Figure 8

Blank


- ◆ Put the blank sample into the light path and press ZERO (2.3)

Set parameters

Before performing the quantitative analysis several parameters must be set: concentration unit, number of standard samples to obtain the standard curve and concentration of those standard samples.


- ◆ On the pre-test interface, press SET (2.1); the instrument will enter the parameter setting display (Fig. 9).

Figure 9


- ◆ Select the concentration unit:
  - On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Unit and press ENTER (2.7); 8 different units will be displayed (Fig. 10).
  - Select the appropriate measurement unit by using the numeric keypad (2.10) or ▼/▲ (2.9) and press ENTER (2.7).
  - If it is not necessary to change the unit press RETURN (2.6).

Figure 10

- ◆ Input the number of standard samples:
  - On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Number and press ENTER (2.7) (Fig. 11).
  - By using the numeric keypad (2.10) input the number of standard samples of known concentration that will be used and press ENTER (2.7). Note: the equipment allows a minimum of 1 standard samples and a maximum of 9.

Figure 11


- ◆ Input the concentration of standard samples:
  - First of all, blank must be performed. Put the blank sample into the light path and, on the pre-test interface, press ZERO (2.3).
  - On the parameter setting interface, move the cursor with ▼/▲ (2.9) to select Concentration and press

ENTER (2.7) (Fig. 12).

- Put the sample number 1 into the light path, input the concentration of that sample by using the keypad (2.10) and press ENTER (2.7).
- Repeat this procedure with each Standard sample. Once the last concentration has been input, the equipment will return to the parameter setting display.

Figure 12


Note: 1. The absorbance or transmittance value of Standard samples will not be displayed; they will be directly saved in the RAM.

- 2. Concentration range is 0-9999; any value out of this range will be invalid.
- 3. Standard samples must be input from lower to higher concentration.

◆ Visualization of standard curve:


- On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Display curve and press ENTER (2.7); the just established standard curve will be displayed (Fig. 13). The equation of the curve will also be displayed at the lower side of the screen and will be saved on the equipment memory.

Figure 13

- If an error message appears (Fig. 14) check samples and input data; press any key and repeat the procedure.


Figure 14

Selection of a curve:

- On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Select curve and press ENTER (2.7) (Fig. 15).

Figure 15

- The equations of the previously established curves (the newly established curve will be displayed at the end of the group). A total of 200 curve equations can be saved.
- Move the cursor with ▼/▲ (2.9) to select the desired equation and press ENTER (2.7); the equipment will return to the parameter setting display.
- To delete a equation, select it and then press CLEAR (2.8).

Test and print

Once the concentration curve has been selected or created proceed as follows:


- ◆ On the pre-test interface, press START/STOP (2.5) to enter the results display.
- ◆ After blanking, put the sample of unknown concentration into the light path and press START/STOP (2.5). The

result will be displayed (Fig. 16). The procedure is the same for all samples.

- ◆ To print results, press PRINT (2.4).


Figure 16

## 2.2 Coefficient method


- ◆ Select Coefficient and press ENTER (2.7); the equipment will display the formula selecting interface (Fig. 17).

Figure 17

- ◆ Select the desired type of formula by using the keypad (2.10) or with ▼/▲ (2.9) and press ENTER (2.7); the Coefficient pre-test interface will be displayed (Fig. 18).

Figure 18

Note: Operation is the same for both formulas since they can be converted one into each other; select one or other only depends on user's preferences.


### Set wavelength

- ◆ On the pre-test interface, press GOTO λ (2.2) and input the desired wavelength value by using the keypad (2.10).
- ◆ Press ENTER (2.7) to set the input value.

### Set parameters

- ◆ Before performing the quantitative analysis it is necessary to set the concentration unit and the value of coefficients K and B of the equation.
- ◆ On the pre-test interface, press SET (2.1); the equipment will enter the parameter setting display (Fig. 19).


Figure 19

- ◆ Set coefficient K:

- On parameter setting display, move the cursor with ▼/▲ (2.9) to select Coefficient K and press ENTER (2.7) (Fig. 20).

Figure 20


- Input the value of coefficient K by using the numeric keypad (2.10) and press ENTER (2.7). Note: the value of coefficient K must be in the range 0-9999, any value out of this range will be invalid.

◆ Set coefficient B:

- On parameter setting display, move the cursor with ▼/▲ (2.9) to select Coefficient B and press ENTER (2.7) (Fig. 21).

Figure 21


- Input the value of coefficient B by using the numeric keypad (2.10) and press ENTER (2.7). Note: the value of coefficient B must be in the range -9999/9999, any value out of this range will be invalid.

◆ Select the concentration unit:

- On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Unit and press ENTER (2.7); 8 different units will be displayed (Fig. 22).

Figure 22

- Select the appropriate measurement unit by using the numeric keypad (2.10) or ▼/▲ (2.9) and press ENTER (2.7).

- If it is not necessary to change the unit press RETURN (2.6).


Blank

◆ Put the blank sample into the light and press ZERO (2.3)

Test and print

- ◆ Once the curve coefficients have been set, proceed as follows:
- ◆ On the pre-test interface, press START/STOP (2.5) to enter the results display.
- ◆ Put the sample into the light path and press START/STOP (2.5). The result will be displayed (Fig. 23). The procedure is the same for all samples.
- ◆ To print results, press PRINT (2.4).
- ◆ To delete the displayed results press CLEAR (2.8).


Figure 23

**3. Kinetics mode**

◆ On the main menu display, select Quantitative with ▼/▲ (2.9) and press ENTER (2.7); the equipment will show the pre-test interface of the kinetics mode (Fig. 24).

Figura 24


### Set wavelength

- ◆ On the pre-test interface, press GOTO  $\lambda$  (2.2) and input the desired wavelength value by using the numeric keypad (2.10).
- ◆ Press ENTER (2.7) to save the input value.

### Set parameters

- ◆ Before performing the kinetics analysis several parameters must be set:
- ◆ On the pre-test interface, press SET (2.1); the equipment will enter the parameter setting display (Fig. 25).


Figure 25

#### ◆ Set time interval:

- On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Time interval and press ENTER (2.7) (Fig. 26).
- Input the time interval (in seconds) by using the keypad (2.10) and press ENTER (2.7). Note: the time interval range is 1-200; any value out of this range will be invalid.
- If no changes are necessary press RETURN (2.6).


Figure 26

#### ◆ Set the total time:

- On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Total time and press ENTER (2.7) (Fig. 27).

Figure 27


**Note: 1. The total time of the analysis must be in the range 1-120000; any value out of this range will be invalid.**  
**2. The maximum value of total time is correlative to the time interval since the system only allows a maximum of 1000 measurements per analysis. Thus, if the time interval is 1 s, the maximum total time will be 1000 s, however, if the time interval is 200 s, the maximum total time can not be 200000 s but 120000 s (600 measurements altogether).**


#### ◆ Set test mode:

- On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Test mode and press ENTER (2.7) (Fig. 28).


Figure 28

- Select the desired test mode ▼/▲ (2.9) and press ENTER (2.7); the symbol √ will be displayed beside the selected mode.
- Press RETURN (2.6) to return to the parameter setting display.
- ◆ Set the upper limit:
  - On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Upper limit and press ENTER (2.7) (Fig. 29).
  - Input the value of the upper limit by using the numeric keypad (2.10) and press ENTER (2.7). Note: the Absorbance measurement range is -0.3/3 and the Transmittance's 0-200.
  - If no changes are necessary press RETURN (2.6).


- ◆ Set the lower limit:
  - On the parameter setting display, move the cursor with ▼/▲ (2.9) to select Lower limit and press ENTER (2.7) (Fig. 30).
  - Input the value of the lower limit by using the numeric keypad (2.10) and press ENTER (2.7). *Note: the Absorbance measurement range is -0.3/3 and the Transmittance's 0-200.*
  - If no changes are necessary press RETURN (2.6).

### Sample analysis

- ◆ After setting the time interval, the total test time and the upper and lower limits, proceed as follows:
- ◆ Press RETURN (2.6) to return to the pre-test interface.
- ◆ Put the blank sample into the light path and press ZERO (2.3)
- ◆ Put the sample into the light path and press START/STOP (2.5) to enter the results display.
- ◆ Press START/STOP (2.5) again to start testing (Fig. 31). As test goes on, the kinetics curve and the time of analysis will be displayed.
- ◆ To interrupt the test press START/STOP (2.5); by pressing again this button, the test will start again from the beginning.


**Note: the curve can not be printed; to check data of each measurement, select Curve data on the parameter setting display.**

To display data of each measurement:

- ◆ Press SET (2.1) to go into the parameter setting interface.
- ◆ Move the cursor with ▼/▲ (2.9) to select Curve data and press ENTER (2.7).
- ◆ Data of the different measurements will be displayed in groups of 6 (Fig. 32).
- ◆ Press PRINT (2.4) to print data.


Figure 32

UTILITIES

- ◆ On the main menu display select Utility with ▼/▲ (2.9) and press ENTER (2.7). A list with all the utilities of the equipment will be displayed (Fig. 33).
- ◆ Use ▼/▲ (2.9) to browse all utilities and press ENTER (2.7) to go into the selected function.
- ◆ Press RETURN (2.6) to go back to the main menu.

Figure 33


1. Deuterium lamp on/off (D2 Lamp On/Off) (Fig. 34)

Note: Zuzi recommends to have the deuterium lamp switched off when working out of the ultraviolet range.

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select D2 Lamp and press ENTER (2.7) to enter the deuterium lamp state selection display.
- ◆ By using ▲/▼ (2.9) select On or Off and press ENTER (2.7). Beside the selected state the symbol √ will be displayed.
- ◆ Press RETURN (2.6) to go back to the function selection display.

Note: When the deuterium lamp is switched on, it will take a few seconds to warm up.

Figure 34


2. Tungsten lamp on/off (W Lamp On/Off) (Fig. 35)

Note: Zuzi recommends to have the tungsten lamp switched off when working out of the visible range.

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select W Lamp and press ENTER (2.7) to enter the tungsten lamp state selection display.
- ◆ By using ▲/▼ (2.9) select On or Off and press ENTER (2.7). Beside the selected state the symbol √ will be displayed.
- ◆ Press RETURN (2.6) to go back to the function selection display.


Figure 35

### 3. Set printer

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Printer and press ENTER (2.7) to enter the printer setting display (Fig. 36)

Figure 36


- ◆ Move the cursor with ▼/▲ (2.9) to select Printer style and press ENTER (2.7):
  - If printer is outside the equipment select Table style.
  - If printer is inside the equipment select Panel style.
- ◆ Whenever changing the roll of paper, move the cursor with ▼/▲ (2.9) to select Scroll paper and press ENTER (2.7). Repeat the operation until the new paper spreads out.
- ◆ Move the cursor with ▼/▲ (2.9) to select Reset and press ENTER (2.7); the printer will be reset.

### 4. Set time and datec

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Time and date and press ENTER (2.7) to enter the time and date setting display (Fig. 37).
- ◆ Input the data of hour and date by using the keypad (2.10). Use ▼/▲ (2.9) to move the cursor.

Figure 37


### 5. Dark Current (0%T) (Fig. 38)

Zuzi recommends to carry out this function whenever changing the environmental conditions such as temperature, voltage or environmental light. For this purpose:

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Dark current and press ENTER (2.7).
- ◆ After a few seconds, the equipment will automatically return to the function selection display and will be ready to operate.

Figure 38


### 6. Wavelength calibration (Reset WL) (Fig. 39)

Zuzi recommends calibrating the wavelength whenever suspecting of it. For this purpose:

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Reset WL and press ENTER (2.7).
- ◆ After a few minutes, the equipment will automatically return to the function selection display and will be ready to operate.


Figure 39

### 7. Lamp life (Fig. 40)

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Lamp life and press ENTER (2.7).
  - Move the cursor with ▼/▲ (2.9) to select Display Lamp life and press ENTER (2.7). The time of use of both lamps will be shown. Note: It is recommended to replace the lamp when the time of use is higher than 1000 hours.
  - Whenever replacing the lamp the time of use must reset to 0. For this, move the cursor with ▼/▲ (2.9) to select Reset D2 life (deuterium lamp) or Reset W life (tungsten lamp) and press ENTER (2.7)


Figure 40

### 8. Default parameters (Fig. 41)

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Load default and press ENTER (2.7). To establish the parameters by default select Yes and press ENTER (2.7).


Figure 41

### 9. Lamp change point (Fig. 42)

- ◆ This function allows setting the wavelength at which tungsten and deuterium lamp exchange. The wavelength of change should be between 300 and 400 nm, any value out of this range will be invalid.
- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Lamp Change and press ENTER (2.7).
- ◆ Input the wavelength value by using the keypad (2.10) and press ENTER (2.7).

Figure 42

### 10. Version

To display the hardware and software version of the equipment:

- ◆ On the function selection display, move the cursor with ▼/▲ (2.9) to select Version and press ENTER (2.7); the corresponding version of hardware and software will be displayed.

### Security

- ◆ The spectrophotometer must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the spectrophotometer on top of a horizontal, plane, stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the spectrophotometer near any warm supply (burners, blowlamps, etc), nor expose it directly to the

sun. Avoid vibrations, dust and dry environments.

- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the spectrophotometer for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the spectrophotometer by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the spectrophotometer (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

## 5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The spectrophotometer is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the spectrophotometer by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the spectrophotometer (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.

### Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the spectrophotometer and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.**

**The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.**

**Please inquire at the community administration for the authorized disposal location.**

## 6. TROUBLESHOOTING

<i>TROUBLE</i>	<i>CAUSE</i>	<i>SOLUTION</i>
The equipment does not turn on	The wire is not connected	Connect the wire
	The power source is broken	Change the power source
	Voltage error	Check the voltage of the electric installation
	Fuse blown or defect in an electric component	Contact the Zuzi Technical Assistance Department
100% T (0.000A) can not be selected	Light beam blocked: holder misaligned	Check the sample holder
	The lamp is worn out or defective	Contact the Zuzi Technical Assistance Department for lamp replacement
	Lamp is not aligned	Contact the Zuzi Technical Assistance Department
	Defect in an electric component	Contact the Zuzi Technical Assistance Department
The equipment strays and there is noise	Not enough period of warm up	Leave 20 minutes to warm up
	Use of glass cell in the UV range	Use a quartz cell
	The lamp is worn out or defective	Contact the Zuzi Technical Assistance Department for lamp replacement.
	Unstable power source	Stabilize the voltage
	Significant change of temperature The lamp is not adjusted	Adjust correctly the lamp Contact the Zuzi Technical Assistance Department
	Defect of detector or defect in an electric component	Contact the Zuzi Technical Assistance Department
Incorrect readings	Not enough volume of sample	Fill the cell with sample
	Equipment out of calibration	Contact the Zuzi Technical Assistance Department
	Bubbles or particles in sample	Check sample preparation, analytical procedure
	Lost of vapours of sample	Adapt the sample to the equipment Use the adequate ventilation
	Wrong blanking	Check the wavelength and the analytical procedure Make Dark current and calibrate the equipment

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