AGITADOR MAGNÉTICO CON CALEFACCIÓN MULTIPLACA
MAGNETIC STIRRER WITH HEATING MULTIPLAQUES
AGITATEUR MAGNÉTIQUE AVEC CHAUFFAGE MULTIPLAQUES

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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TABLE OF CONTENTS

1. INTRODUCTION	4
2. TECHNICAL SPECIFICATIONS	4
3. OPERATION	5

1. INTRODUCTION

The Nahita-Blue multiposition magnetic stirrer with heating is designed to stir and heat low-viscosity liquids or solid-liquid mixtures, with synchronized or independent control of stirring, temperature, and time. It offers easy operation, stable performance, and continuous speed regulation.

2. TECHNICAL SPECIFICATIONS

Code	LBD019
Model	670/4
Stations	4
Max stirring volume	4x2 L (H□O)
Motor	Brushless motor
Control	Display and knob
Station material	Stainless Steel ceramic coated
Station dimension	Ø135 mm
Speed	50-1500 rpm (increments 10 rpm)
Plate temperature	RT-310°C
Stirring bar size	30 mm
Stirring rotation	Clockwise
Timing	1min-24h
Power	1200 W
Voltage	110-220V 50/60 Hz
Permissible temperature and humidity	5-40°C / 80% HR
Dimensions	410x320x92 mm

3. OPERATION

■ Using steps:

- Connect the device to the power supply and switch it on.
- Press the right knob to toggle between speed/temperature/time settings and turn the right knob to adjust the desired parameters.
- Turn the left knob to switch between the different stations; press and hold it for 2 seconds to synchronize the current station's parameters with the others (when synchronizing stations, if the number of probes in the four positions is inconsistent, synchronization will not be possible).
- Press and hold the right knob for 2 seconds to start/stop.

■ Instrument use environment:

- Indoor use only
- Altitude: ≤ 2000 m
- Working temperature range: $+5$ °C a $+40$ °C
- Humidity range: $\leq 80\%$
- No vibrations or airflow that could affect performance.
- No conductive dust, explosive gases, or corrosive gases in the surrounding air.
- Maintain a safety distance of 10 cm around the instrument while it is in operation.

■ Safety precautions:

- Do not unplug or plug in the power connector or toggle the power button with wet hands.
- Do not unplug the power plug while the instrument is powered on.
- Do not perform maintenance or cleaning while the instrument is powered on.
- Do not install the instrument on an uneven, swaying, or vibrating work surface.