MINICENTRÍFUGA DIGITAL MINIFUGE-GD06 DIGITAL MINICENTRIFUGE MINIFUGE-GD06

Ref. GDC006

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6. Troubleshooting

The following table lists the faults of the instrument and the causes of the related faults and the treatment methods. You can follow the prompts to eliminate them. If the user still cannot solve the faults or the faults are not included in the list below, the user should contact the production immediately. Unit maintenance staff.

After an abnormality occurs, turn off the power, unplug the plug, and then start the fault after troubleshooting.

Table 1: Fault phenomenon, cause and troubleshooting

Fault phenomenon	Causes and troubleshooting
The instrument does not rotate	1. Check if the power socket and wiring are good, and whether the power socket has power. 2. The power switch is not in good contact. 3. If it cannot be excluded, please contact the manufacturer's maintenance personnel.
Suddenly stopped while the instrument is running	1. Is the power connection good? 2. Determine the time limit for use, whether the carbon brush is within the service life (carbon brush life is about 1600 hours). 3. The motor is overheated or overloaded. Turn off the power and wait for 10-20 minutes before restarting. 4. Does the instrument have a burnt smell? If it is, please unplug the power immediately and contact the manufacturer's maintenance personnel. 5. The voltage range may be too large. Check if the power supply voltage meets the specified AC100 250V.
Transparent cover can't be opened or closed.	1. The transparent cover cannot be opened, and there is any foreign matter under the cover plate. 2. The transparent cover cannot be closed, and the triangular shaped button above the cover plate is worn out.
Instrument noise is relatively large	1. The rotor has a critical speed and some vibrations of the machine are normal. 2. Check if the rotor is locked. 3. Check the symmetry balance of the rotor load and check the level of the instrument. 4. Check if the rotor is installed. 5. Check the drive shaft and rotate it by hand. If it cannot rotate smoothly, there may be a problem with the drive shaft or the motor. If it cannot be excluded, please contact the manufacturer's service personnel.
Panel LCD does not light after turning on the power switch	1. The power supply is not connected, check the power distribution. 2. Switch the remaining indicators to check if they are normally lit. If not, please contact the manufacturer. 3. If it cannot be excluded, please contact the manufacturer's maintenance personnel.
E3 code	1. Compruebe si la tapa transparente está bien cerrada o si el cierre de la tapa no se aprieta automáticamente. 2. Si no se soluciona póngase en contacto con el fabricante.
Beyond the above fault content	1. Whether the transparent cover is tightly closed or the buckle wear cover does not tighten automatically. 2. Exclude the first reason, please contact the manufacturer's maintenance personnel! Solved by the manufacturer's professionals!

Company website:
www.auxilab.es
www.rogo-sampaic.com

Stop instrument

Centrifuge at the end of the process can be directly pressed START/STOP, the key stops, the standby device stops running, the transparent upper cover is opened → the sample is taken out → the power is turned off (the transparent upper transparent upper cover is closed without the need to lock when the instrument is no longer used).

4.6.2 [About calculation of Centrifugal Force]

The relative centrifugal force is typically thousands of times the Earth's gravity (g), which is a measure of the efficiency of various instruments to separate or precipitate objects. The calculation of the centrifugal force is related to the centrifugal speed and the centrifugal radius, and is calculated according to the following equation:

$$RCF = 11.18 \times n / 1000^2 \times r$$

r: is the radius of the centrifuge, in cm;

n: centrifugal speed, unit rpm (revolutions per minute)

Note: "The maximum centrifugal force value is related to the maximum centrifugal radius";

5. Maintenance matters

Do not attempt to clean the centrifuge when the power cord is plugged in or the power switch is turned on. If the equipment or accessories are contaminated by pathogenicity, toxicity or radioactive materials, it is the responsibility of the user to perform proper purification. This should be explained when developing the "Laboratory Biosafety Manual". Cleaning and Decontamination In addition to the other requirements mentioned in this chapter, please contact our company to check if they are likely to damage the equipment. If you plan to return equipment or accessories to the company or request repair, you must ensure that the equipment or accessories are clean and harmless to the human body! Do not use organic solvents because it can decompose the grease in the motor bearings; during the cleaning process, liquids, especially organic solvents, cannot be contacted to the motor shaft and the bearings.

5.1 Do not use sharp objects to collide with the rotor. Prevent bumps during handling and disassembly. Prevent cracks in the rotor during use due to scratches or trauma.

5.2 Regularly check the rotor assembly (especially the bottom of the test tube hole) for corrosion spots, grooves, and small cracks. If any of the above conditions are found, stop using the rotor and contact the production unit.

5.3 Normally, the rotor is cleaned once a week. If it is separated from salt solution or other corrosive samples, please wash it immediately after use. If the sample is found to be spilled during use, dipped or dripped on the rotor, it should be immediately drained and partially cleaned.

5.4 When cleaning the rotor, clean it with a mild detergent dampened with a sponge or cotton cloth, then wash off the detergent with distilled water. Do not sprinkle or spray the rotor with water as the liquid may be left somewhere and cause corrosion. Allow to invert and dry after washing.

5.5 Check the locking of the transparent cover and the main body casing, and whether the cover can be opened and closed smoothly.

5.6 Carefully wipe the surface of the instrument with a soft cloth dampened with a mild detergent, then remove the detergent with a clean, damp cloth.

5.7 If there is a major problem with the instrument, if the power supply is not working properly, there is a burnt smell inside the machine, the rotor is broken, and the instrument is shaking sharply during operation. Please contact the manufacturer or the authorized dealer of the manufacturer in time. Repair by yourself under the guidance of maintenance personnel!

INDEX OF LANGUAGES

Spanish	2-11
English	12-20
French	21-30

INDEX OF CONTENTS

Instrument use environment	12	4.2 Turn On/Off Power	16
Safety warning	13	4.3 [Installation / Disassembly]	16
Safety warning sign description	13	4.4 Test tube placement method	17
Meaning of the safety instructions	13	4.5 Operation Panel Description	18
Safety precautions	14	4.6 Start/Stop Instrument	18
1. Product description	14	4.6.1 Example of use	18
2. Product main technical parameters	14	4.6.2 Centrifuge force	19
2.1 Rotor description	15	5. Maintenance matters	19
3. Appearance function introduction	15	6. Troubleshooting	20
4. Operation instructions	16	Table fault phenomenon	20
4.1 Open/Close	16		

This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT

Instrument use environment

To ensure the continued safety of the machine, the following factors may damage the centrifuge:

- Chemical effects;
- Environmental impacts, including natural UV radiation;
- Corrosion and wear of hood parts and other safety parts.

- Indoor use
- Altitude: ≤ 2000 m
- The applicable ambient temperature range of the instrument is +5 °C ~ +40 °C
- The working temperature range of the instrument is ≤ 80%
- The working model of the instrument working power supply and its main technical parameters
- There must be places in well ventilation room
- There are no vibrations and airflows around the performance
- There is no conductive dust in the surrounding air, explosive gas and corrosive gas

Safety warning

- Use this machine for the first time, please read this manual carefully!
- MiniFuge GD06 mini centrifuges can only be operated by trained and authorized personnel
- Equipment maintenance can only be done by Auxilab Group or Auxilab Group authorized agents
- Do not use the following materials in the centrifuge:
 - Flammable and explosive materials.
 - Strong chemical action materials.
 - Toxic or radioactive substances, or pathogenic microorganisms.
- Only qualified maintenance personnel can use the appropriate tools to perform maintenance on the Mini Fuge GD06 mini centrifuge system.
- If the operator encounters a situation not mentioned in this manual, please contact the authorized agent of Auxilab Group or Auxilab Group to inquire about the correct treatment.
- Try to use the accessories provided by Auxilab Group. If the user wants to use other accessories, Auxilab Group will not be responsible for the adverse consequences. However, the user can apply to Auxilab Group to verify that the accessory meets Auxilab Group's
- The Mini Fuge GD06 mini centrifuge must be inspected and maintained at regular intervals.

Safety warning sign description

Note: Please read the instruction manual carefully before using the instrument!

Attention: Danger of high voltage! Risk of electric shock!

Meaning of the safety instructions


In order to avoid damage to surrounding people or damage to surrounding objects and the environment, please be sure to follow all safety instructions in this user manual.

In addition to accident prevention, environmental protection, and recognized professional rules in safety and occupation, including the country where the user is located and the local laws and regulations for the installation of the centrifuge, they must be carefully observed.

Neglecting the consequences of safe operating procedures

Any violation of safe operating procedures, laws, regulations and rules will lead to harm to people, objects and the environment.

4.5 [Operation Panel Description]

1. MENU button: switch speed / time (min. / sec.);
2. speed display area;
3. parameter increment key;
4. door cover open indicator;
5. time display area;
6. parameter decrement key;
7. Start/stop button;
8. Operation indicator;
9. Door cover closing indicator;
10. Time unit indicator character;

4.6 [Start/Stop Instrument]


4.6.1 Example of use

Start the instrument

Open the transparent cover → install the required rotor → symmetrically place the sample → turn on the power → close the transparent cover;

For example: If the speed required is 1200 rpm and the time is 10min. Press MENU, the speed indicator on the display panel will flash, press key 3 to increase the speed to 1200. Next press MENU again, the time indicator will flash on the display, press Key 3 or Key 6 to adjust time up or down to 10min. The display will flash twice when set.

The settings are as follows:


Press STAR/STOP when the key is activated, the machine will stop the operation according to the set time countdown, and the buzzer alarm will prompt the completion of the centrifugation; The first time the instrument is turned on, due to the completion of the centrifugation;

The first time the instrument is turned on, due to the vibration generated during transportation, the static brushing of the carbon brush for a long time will cause a slight carbon brush friction noise at the start up. If the instrument is idling for 2-3 minutes, the motor will operate normally without worrying about the instrument quality!

[Installation]

A" and "B" directions in the "mounting rotor" (installation diagram 3), that is, the position indicated by the arrow "A" of the rotor and the position of the rotor elastic buckle "B" B", to ensure that the positioning holes of the two positions are the same. Hold the rotor symmetrically with both hands and press it vertically. When you hear the sound of "Carbazole", it indicates that the rotor has been locked with the elastic buckle of the rotor. It should be further verified whether the installation is in place. At this time, grasp the rotor with both hands and slowly lift it up. The head is firmly fixed on the elastic buckle of the rotor, indicating that it has been correctly installed in place; (Note: The configuration of the adapter is the same!)

[Desinstalación]

The picture above shows the universal rotor of the MiniL mini centrifuge (disassembly diagram 1)

According to the direction of the arrow in (disassembly diagram 1), hold the rotor with both hands and squeeze the two thumbs inward, slowly lift the rotor up and take it out; (Note: The disassembly method is the same!)


Never use a strong pull out of the rotor! Will cause the rotor to deform or damage the motor! The face with the Arabic numerals printed on the head is front. Should be facing up, the reverse may cause the machine to be damaged and the personnel to be hurt!

- **DO NOT USE FOR MORE THAN 4-5 MINUTES AT A TIME.**
- **WAIT 4-5 MINUTES BETWEEN CENTRIFUGATION CYCLES.**

RISK OF HEATING OF PLASTIC COMPONENTS!

4.4 [Test tube placement method]

Symmetrically balanced placement of 12 --, 8 --, 6 --, 4 and 2 tube angles (Figure 1)

Safety Precautions

- Do not plug or unplug the power plug and toggle the power button when the hands are in possession of liquid!
- Do not plug or unplug the power plug when the instrument is powered!
- It is strictly forbidden to maintain and scrub the instrument under power on state!
- Do not allow the sample volume to exceed 70% of the tube capacity for centrifugation!
- It is strictly forbidden to install the instrument on the work surface with unevenness and shaking vibration!
- When the rotor is not completely stopped, it is strictly forbidden to put the hand or other objects on the rotor, in an attempt to replace the instrument brake!

1. Product description

MiniFuge GD06 adopts international advanced design concept and manufacturing technology, the shape is smooth and beautiful, the structure is compact and firm, the high versatility, the simple operation and so on, it is a rare instrument for scientific research experiments. This handheld centrifuge is ideal for microtube filtration and rapid centrifugation, micro blood cell separation, microbial sample processing, PCR experiment zone centrifugation, and prevention of liquid wall hanging in the centrifuge tube.

2. Product main technical parameters and rotor parameter description


Product number	MiniFuge-GD06
Input power	AC100-250V/50/60Hz
Power fuse	PPTC /self recovery fuse / no need to replace
Input power [W]	45W
Control method	ST16 bit MCU/PWM/speed regulation
Display method	Black and white LCD
Motor	DC24V/DC permanent magnet motor
Effective centrifuge time	1-99 min. /1-59 seg.
Speed [r/min.]	500~12000 rpm±5%
Speed increase	500 rpm
Maximum relative centrifugal force [×g]	9660×g
Misoperation/fault alarm	Voice prompt + display code
Size: L×W×H [mm]	194×229×120
Weight (no rotor)	2Kg
Noise at maximum speed	≤47dB(A)

Fastest acceleration time	≤12 s
Fastest deceleration time	≤16 s
Protection level	IP21
Allowable ambient T ^a /HR	5-40 °C / 80 %
Testing standard	EN 61010-1 EN 61010-2-020 EN 61326-1 EN 61010-3-2/A2
Note: Ambient temperature 23 °C / at 23 °C ambient temperature	

2.1 MiniL rotor parameter description MiniFuge-GD06

Rotary number	Capacity (ml × n° tubes)	Max speed. (r/min)	Max Rcf. (g)	Test type tube
#1 GDF008	1.5-2ml ×12 tubes + PCR 0.2ml×8 tubes×4 row	12000	9500	Conical bottom microtube
#2 GDF012	PCR 0.2ml×12strips×4 row	12000	5960	PCR strips
#3 GDF011	5ml × 4tubes	12000	9660	Conical bottom
#4 GDF005	Blood tubes capillaires × 12	12000	8371	Tubes ø1.5×40 mm
#5 GDF006	5ml 4 tubes / 1.8ml×4 tubes (using adapter GDF009)	7000	3180	Cryotubes with stopper

3. MiniFuge-GD06 centrifuge components

- | | | |
|----------------------------|------------------------------|-----------------------------|
| 1. Transparent cover | 5. Lower casing | 9. Power switch |
| 2. Top casing | 6. Non-slip feet | 10. AC100 250V power socket |
| 3. Rotor | 7. Closure | |
| 4. Operation display panel | 8. Cover release lever (tab) | |

4. Operation instructions

Note: Before using the instrument, please ensure that the instrument is transported without damage!

- Appearance of the instrument is brand new!
- The release tab opens the cover easily and without unusual noises.
- The rotor elastic buckle is tightly fitted on the motor spindle and well bonded to the casing.
- The power switch and power cord are securely fastened to the casing.

The equipment can be used safely if all the above conditions are met

If there is any of the above, please contact the manufacturer or the authorized dealer of the manufacturer in time, and do not repair it yourself without the guidance of professional maintenance personnel!

4.1 [open/close]


[Open cover] The cover lever (8) in diagram 1 is pressed in the direction of the operation panel and lifts upward to open the cover.

[Close] Slowly lower the transparent top cover until you hear a closing sound (click) from the lever (tab).

4.2 [Turn On/Off Power]

The power switch button in the appearance diagram 1 is set to the “I” position to turn the power on, and the “O” position to turn the power off;

4.3 [Installation / Disassembly of the rotor]

The picture above shows the universal rotor of MiniL mini centrifuge (installation diagram 1)