ANILLO Y ESFERA DE GRAVENSANDE

GRAVESANDE'S RING AND BALL

REF: QLJ001

INDEX OF LANGUAGE-ÍNDICE DE IDIOMAS

SPANISH.....2

ENGLISH4

ANILLO Y ESFERA DE GRAVESANDE

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GRAVESANDE'S RING AND BALL

A simple and economical device for demonstrating the thermal expansion of solids. It is composed of a brass ring and ball each one of them mounted on separate insulated handles. When cold, the ball of 19 mm diameter easily passes through the ring but does not pass when heated. If the ring is heated on the same extent than the ball, then it passes through the ring.

Theory

Dilatation

Metals are materials that have high dilatation, in part because of their conductivity. Sometimes the change of a material can be perceived to undergo the changes of ambient temperature. Thermal expansion is the tendency of matter to change in shape, area, and volume in response to a change in temperature. The degree of expansion divided by the change in temperature is called the material's coefficient of thermal expansion and generally varies with temperature.

Theoretical basic

The volume V_2 of a body with a coefficient of thermal expansion δ , that at the beginning has V_1 volume at a T_1 temperature, through by a second temperature T_2 $V_2 = V_1(1 + \delta\Delta t)$. If we talk about density is the same $D_2 = D_1(1 + \delta\Delta t)$.

The knowledge of the coefficient of thermal expansion acquires a high importance in industry and construction. For example the railroad tracks are changing their temperature constantly due to the atmosphere. If the temperature rises the railroad tracks can increase in volume and thus move from its original position being very dangerous. In order to minimize this

problem tracks are heated and cut before use it to avoid natural dilatation. This process is called neutralization of stress.

In comparison with anomalous dilatation of water. In this case, water is dilatated if the temperature rises and is contracted when decreases, but in a room temperature. In the freezing point, at 0° the behavior is the opposite, it is very important to save the life.

Testing


Materials

- Gravesande's ring and ball
- Bunsen lighter
- Beaker
- Beaker holder
- Termometer

Steps

1. Check that both the ring and the ball have the same temperature.

2. Check that the ball can pass through the ring easily.
3. Put water in the beaker.
4. Switch on the bunsen lighter to heat the water.
5. Introduce the Gravesande ball into the water to increase its temperature.
6. Try to introduce the ball into the ring once heated
7. Introduce the Gravesande ring into the water to heat the ring at the same temperature as the ball.
8. Pass the ball through the ring when both have the same temperature.

What happened when is heated only the ball ?

The ball can't pass through the ring.

Explanation

The ball has increased in volume.

Cubic dilatation is the rise in volume that a body experiments when is heated.

Thus the buildings, bridges... are designed carefully taking into account the dilatation and contraction that they can experiment with the variation of temperature.

Other example:

Some bridges are 1m more long aprox. in summer than in winter. This fact is important to avoid the bridges bend.