Módulo de transferencia / Transfer module
Referencia / Code: ZFD006
Manual de Instrucciones / Instruction Manual

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SAFETY PRECAUTION

WHEN USED CORRECTLY, THESE UNITS POSE NO HEALTH RISK. HOWEVER, THESE UNITS CAN DELIVER DANGEROUS LEVELS OF ELECTRICITY AND ARE TO BE OPERATED ONLY BY QUALIFIED PERSONNEL FOLLOWING THE GUIDELINES LAID OUT IN THIS INSTRUCTION MANUAL.

ANYONE INTENDING TO USE THIS EQUIPMENT SHOULD READ THE COMPLETE MANUAL THOROUGHLY.

THE UNIT MUST NEVER BE USED WITHOUT THE SAFETY LID CORRECTLY IN POSITION.

THE UNIT SHOULD NOT BE USED IF THERE IS ANY SIGN OF DAMAGE TO THE EXTERNAL TANK OR LID.

INTRODUCTION


Electro-blotting is a very useful technique in research laboratories for transferring macromolecules (proteins or nucleic acids) from the gel in which they are embedded to a membrane (normally nitrocellulose membrane), followed by detection methods for specific molecules.

Transfer module is mainly used for transferring proteins from a polyacrylamide gel when carrying out the Western Blot (WB) technique.

PACKING LIST

Part Name	Quantity
Trans core	1
Transfer splint	2
Spongy cushion	4

COMPONENTS

CARE AND MAINTENANCE

Units are best cleaned using warm water and a mild detergent. **Water at temperatures above 60° C can cause damage to the unit and components.**

The tank should be thoroughly rinsed with warm water or distilled water to prevent buildup of salts, but care should be taken not to damage the enclosed electrode and vigorous cleaning is not necessary or advised.

Air drying is preferably before use.

The units should only be cleaned with the following: RNase Decontamination

This can be performed using the following protocol:

Clean the units with a mild detergent as described above.

Wash with 3% hydrogen peroxide (H₂O₂) for 10 minutes.

Rinse with 0.1% DEPC-(diethyl pyrocarbonate) treated distilled water,

Caution: DEPC is a suspected carcinogen. Always take the necessary precautions when using.

RNaseZAP™ (Ambion) can also be used. Please consult the instructions for use with acrylic gel tanks.

OPERATION METHOD

1. Add a small amount of transfer buffer (39 mmol/L glycine, 48 mmol/L Tris, 0.037% SDS, 20% methanol) in a shallow tray, place the spongy cushion and filter paper soak in the buffer.
2. Cut the filter paper to the same size with the gel, and put the filter paper on the spongy cushion, put the gel on the filter paper, and then put the cellulose nitrate membrane and filter paper on the gel, finally cover with a spongy.
3. Carefully move the filter paper, spongy cushion, membrane and the gel to the transfer splint and fasten the bayonet, insert the transfer splint to the core of the tank, add the buffer until the transfer splint immerse in the buffer.

Note: 1. Please pay attention of the direction of positive and negative.

2. Before the experiment, please cool the electrophoresis buffer in the refrigerator.

4. Carefully place the lid on the tank and connect to a power supply, start the transfer electrophoretic.

Note: In order to ensure a good result, please keep the temperature at 4° C.

5. After the experiment, take out the transfer splint, remove the cellulose nitrate membrane, according to the experimental purpose to proceed the next step detection.