CENTRÍFUGA ANGULAR DIGITAL RS-0408

RS-0408 DIGITAL ANGLE CENTRIFUGE

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Spanish 2-15
 English 16-29
 French 30-43

INDEX OF CONTENTS

Safety Instructions..... 16
 1. Intended use 18
 2. Specifications 18
 3. Required Operational Conditions 19
 3.1 Basic operational conditions 19
 3.2 Transport and storage conditions 19
 4. Installation 19
 4.1 Location 19
 4.2 Connection of the power cord and grounding 19
 5. Structure 20
 6. Control panel 20
 7. Rotor Preparation 21
 7.1 Prepare the samples 21
 7.2 Inject the samples into tubes 21
 7.3 Keep the tubes balanced 21
 7.4 Inspect the rotor 21
 7.5 Symmetrically load centrifuge tubes into rotor 21
 8. Operation 22
 8.1 Normal operation 22
 8.2 RCF operation 24
 8.3 Pulse operation 24
 9. Maintenance 24
 9.1 Cleaning 24
 9.2 Rotor Installation 25
 10. Troubleshooting 26
 10.1 Possible problems and solutions 26
 10.2 How to open the door 26
 11. Instructions for the rotor and tubes 27
 11.1 Rotor 27
 11.2 Tubes 28
 12. Calculate RCF 29

SAFETY REMINDER

Common safety precautions

Carefully read the following safety precautions for a thorough understanding.

- Follow the instructions and procedures described in this manual to operate this centrifuge safely.
- Carefully read all safety messages in this manual and the safety instructions on the instrument.
- Safety messages are labeled as indicated below. They are in combination with signal words of “WARNING” and “CAUTION” with the safety alert symbol to call your attention to items or operations that could be dangerous to you or other persons using this instrument. The definitions of signal words are as follows:

3. Autoclaving of PA, PC and PP centrifuge tube

PA begins softening at the temperature of 120°C, while PC and PP begin softening at 130°C. Generally, PA may be sterilized for 30 min at 115°C (0.7 kg/cm²), while PC and PP may be sterilized for 20 min at 121°C (1.0 kg/cm²). Too high temperature would result in deformation of centrifuge tube.

When autoclave is used, take the following steps:

- Place the centrifuge tube upright with opening facing upward. If the centrifuge tube is placed in an inclined or horizontal manner, it will deform due to gravity.
- Remove the threaded cover and inner cover to prevent deformation or crack of the centrifuge tube.
- Take the centrifuge tube only when the autoclave cools down to the room temperature.

4. Service life of centrifuge tube

The service life of plastic centrifuge tube depends upon the nature of sample, rotor speed and centrifugation temperature. When the plastic centrifuge tube is used for centrifugation of conventional neutral samples (pH5-pH9), its estimated service life at the max. rotation speed is as follows:

- When used at the max. rotation speed:
- High-quality centrifuge tube (PA, PC, PP): 30-50 times
- Conventional centrifuge tube: about 10 times.
- The service life of centrifuge tube is also related to the cleaning and sterilization conditions and might decrease somewhat.
- Note: Never use any centrifuge tube with cracks on it.

12. RCF CALCULATION

Relative centrifugal force (RCF) can be calculated using the following formula:

$$RCF = 1,118 \times r \times n^2 \times 10^{-6}$$

r-rotation radius, unit: cm

n-rotation speed, unit: rpm

3. Precautions

The density of sample that centrifuge rotor can separate is less than 2.0g/ml. If the density of sample to be separated exceeds 2.0g/ml, please calculate the permissible rotation speed using the following formula:

$$\text{Allow speed (rpm)} = \text{Max. Speed} \times (2,0 \text{ (g/ml)} / \text{sample density (g/ml)})^{1/2}$$

4. Sterilization

Rotor body is made of high-strength plastics and cannot be autoclaved or radiated by ultraviolet beam.

11.2 Centrifuge tubes

1. Cleaning and sterilization

Conditions			Materials		
			PA	PC	PP
Cleaning	Fluid cleaning	Acid cleaner (pH5 or lower)	X	X	X
		Acid cleaner (higher than pH5)	O	O	O
		Alkaline cleaner (higher than pH9)	O	X	O
		Alkaline cleaner (pH9 or lower)	O	O	O
		Neutral cleaner (pH7)	O	O	O
		70°C hot water	O	O	O
	Ultrasonic	Neutral cleaner (pH7)	O	O	O
Sterilization	Autoclaving	115°C (0.7kg/cm2) 30min	O	O	O
		121°C (1.0kg/cm2) 20min	X	O	O
		126°C (1.4kg/cm2) 15min	X	X	X
	Boiling	15min to 30min	O	O	O
	Ultraviolet	200-300nm	X	X	X
	Gas	Ethylene oxide	O	X	O
Formaldehyde		O	O	O	

PA: Polialómero PC: Policarbonato PP: Polipropileno

2. PC centrifuge tube cleaning

PC material has relatively low chemical stability to alkaline solvent, therefore any use of cleaning agent with pH value of over 9 should be avoided. Some neutral cleaning agents still have pH value of over 9 after being diluted as recommended by the vendor, therefore use of cleaning agent with pH value of 7-9 is recommended.

WARNING: Personal Danger

Warning notes indicate any condition or practice, which if not strictly observed, could result in personal injury or possible death.

CAUTION: Possible damage to instrument

Caution notes indicate any condition or practice, which if not strictly observed or remedied, could result in damage or destruction of the instrument

NOTE: Notes indicate an area or subject of special merit, emphasizing either the product's capability or common errors in operation or maintenance.

- Do not operate this centrifuge in any manner not described in this User manual. When in doubt or have any troubles with this centrifuge, ASK FOR HELP.
- The precautions described in this User manual are carefully developed in an attempt to cover all the possible risks. However, it is also important that you are alert for unexpected incidents. Be carefully operating this centrifuge.

WARNING

- This centrifuge is not explosion-proof. Never use explosive or flammable samples.
- Do not install the centrifuge in or near places where inflammable gases are generated or chemicals are stored.
- Do not place dangerous material within 30cm around the centrifuge.
- Make sure to prepare necessary safety measures before using samples that are toxic, radioactive or contaminated with pathogenic micro-organisms at your own responsibility.
- If the instrument, rotor and/or accessories that has been contaminated by solutions with toxic, radioactive or pathogenic materials, clean it according to the decontamination procedure that you are specified.
- If you require services at site, please sterilize and decontaminate it in advance, and then notice the service center involved in the details of the particular materials.
- Do not handle the power cord or turn on or off the POWER switch with wet hands to void electrical shocks.
- For safety purposes, do not enter within 30cm around this centrifuge while it is in operation.
- While the rotor is rotating, never forcedly release the door lock.
- Unauthorized repairs, disassembly, and other services to the centrifuge except by our service center are strictly prohibited.

- This centrifuge must be located on one firm and level table.
- Make sure the centrifuge is horizontal before running.
- Make sure the angle between the door and cover is greater than 70 degrees when open the door.
- Be careful not put your fingers or hands between the door and cover when the door off.
- Do not move or relocate this centrifuge while it is running.
- If fluid spills in the rotor chamber, please promptly clean and dry with a dry cloth to avoid sample contamination.
- Ensure to remove any objects and fragments of the tubes dropped inside the rotor chamber before running this centrifuge.

■ Cautions on rotors:

- 1) Always check for damages on the rotor surface before using it. Do not use the rotor in question.
- 2) Do not set the centrifuge speed beyond the speed limitation of the rotor kits (the combination of rotor and adapters). Make sure to run it below the lowest permissible speed.
- 3) Do not exceed the imbalance tolerance.
- 4) Do not over fill centrifuge tubes.
- 5) Ensure the lid is fastened onto the rotor if there is one before running.

■ In case of anomalies such as strange noises during operation, please stop it immediately and contact our service center. Keep a record of the error code if displayed and inform the service personnel of the code.

■ Intense vibrations could possibly damage the centrifuge. Contact our service center in case of such suspicion.

1. INTENDED USE

This device is a medical device (laboratory centrifuge) and this centrifuge is mainly used to separate blood, urine or feces samples. Before using this instrument, please undergo specialized training and operate it strictly according to the instructions.

2. SPECIFICATIONS

Maximum speed	4500rpm (Range: 300-4000rpm, increment: 100rpm)
Maximum RCF	1900×g
Maximum capacity	15ml×8 / 10ml×12
Timer	1 min-99 min-HOLD (continuous operation)
Noise	≤60dB(A)
Driving Motor	Brushless DC motor
Safety devices	Single door locks, over speed and status diagnosis system
Power requirements	Single phase 100V-240V, 50/60Hz
Ambient conditions: -Set-up site -Altitude -Ambient temperature -Max. humidity -Excess-voltage category -Pollution degree	Indoor No more than 2000m 2°C ~ 40°C 80% II 2
Device protection class	I
EMC	EN/IEC 61326-1 / Class A - FCC / Class A
Dimensions	286 (L) × 367 (D) × 227 (H)
Weight	9,2 kg
Additional features	LED display, sound signal at the end of operation

11. ROTOR AND CENTRIFUGE TUBES


PRECAUCIÓN

Carefully read the user manual and correctly install and use the rotor.

Don't exceed the max. speed allowed by test tubes and adaptors. The max. speed allowed by certain adaptors is lower than the max. speed of the rotor.

11.1 Rotor

1. Rotor structure

2. Rotor, tubes and adapters

Rotor	Tube	n° of Tubes	Dimensions (Ø×L) mm	Use of adapter
Angular	Falcon type 15ml	8	17×120	No
	Vacutainer type	12	13×75	No
	Vacutainer type	12	13×100	Yes (with pad)
	Vacutainer type	12	16×107	Yes (without pad)
	Collection tube (EU)	12	15×92	Yes (without pad)

10. COMMON FAILURES AND SOLUTIONS

10.1 List of common failures

This centrifuge is capable of self-diagnosis. When the centrifuge fails, the time display area will show the error code to identify the cause of failure

Failure	Possible cause	Solution
No display after powering on	No power supply to power socket	Check the power and reconnect the power supply
Abnormal vibration of centrifuge	The rotor does not match with the drive shaft The sample installed asymmetrically	Reinstall the rotor Weigh using a balance and install it symmetrically
Error code	E-02 Upper cover door failed	Door opened during operation Press the key START/STOP when the door is open
	E-10~86	See the service manual

10.2 How to open the upper cover

1. When the centrifuge is powered on

Open the upper cover only when the rotor is not running.

- When the rotor fully stops, press the uncovering button to open the upper cover

2. When the centrifuge is powered off

If a sudden power outage occurs, the upper cover may be opened manually as well.

- Check whether the rotor is running.
Listen carefully to confirm that the rotor fully stops.

- Upper cover opening

Confirm the rotor has stopped and press the opening button to open the upper cover

3. REQUIRED OPERATIONAL CONDITION

3.1 Basic operational conditions

- Power supply: 100V-240V, 50/60Hz
- Ambient temperature: 2°C~40°C
- Relative humidity: ≤80%
- No vibration or airflow around
- No airborne dust with charge, explosive and corrosive gases around

3.2 Storage temperature

- Temperatura de almacenamiento: -40°C~55°C
- Relative humidity: ≤93%

4. INSTALLATION

Users must strictly comply with the installation instructions contained in this chapter.

WARNING

Improper power supply may damage centrifuge.
Make sure the power source conforms to the requirement before connecting

4.1 Location

- Place the centrifuge on a firm, flat and level bench top, ensure the four feet of this centrifuge stand on the table firmly. Avoid slippery surface or surface prone to vibration.
- The recommended ambient temperature is 20°C±5°C. Avoid temperature over 30°C or direct exposure to sunlight.
- Keep clear of the centrifuge at least 10cm on both sides and at least 30cm behind it to guarantee the cooling efficiency.
- Keep away from heat source or liquid leakage to avoid elevated sample temperature or centrifuge failures.

4.2 Connection of power cord and grounding


WARNING

To avoid electrical shocks, ensure your hands are dry when touching the power cord.
This centrifuge must be properly grounded.

Proper earth protection in compliance with local safety requirements is compulsory.

5. STRUCTURE

Drive shaft

We recommend to clean the drive shaft with a piece of soft cloth and apply a thin layer of silicone grease on it periodically

Upper cover

Clean or sterilize the upper cover in the same manner as (1) centrifuge.


Rotor

If the rotor is put out of service for a prolonged period, please remove the rotor from the centrifugal chamber and place the rotor upside down to empty the rotor hole in order to prevent corrosion. If any sample leaks into the rotor hole, flush the rotor hole with clean water and apply a thin layer of silicone grease on the rotor surface after it dries up.


The rotor needs periodical maintenance. We recommend to clean the rotor every three months by cleaning the tube hole and spindle hole and applying a thin layer of silicone grease thereafter.

9.2 Rotor installation

Installation

Adjustment

Before tightening the rotor, rotate the rotor with hand while carefully observing for any visible vibration. If yes, please remove the rotor and rotate it to a certain angle before remounting it until the rotor has no visible vibration. Then tighten it.

Error display

■ It will automatically stop if any failure occurs in the running process, with the error code indicated on the time display window. By looking up it from table 11-1, the cause of error can be found, and appropriate action may be taken

3. Stop running

When the running time is out or by pressing the key , the centrifuge begins deceleration to stop

- When the rotor fully stops, the buzzer will beep to remind the user that the operation is over.
- After the rotor fully stops, the red LED and green LED indicators will blink alternately, and you will see a blinking text "OPEN" in the time display area to remind you that the operation is over. If there is no operation within 1 minute, the red LED indicator goes out and the green LED indicator stays on. If there is any operation within 1 minute, both the red and the green LED indicators will go out.

The upper cover lock is opened

- After the end of operation, the upper cover may be opened by pressing the door-opening button.
- After the end of operation, the program will automatically save the set parameters of this operation.
- When the device is started up again, the program will automatically retrieve the set parameters of the last operation.

Open the upper cover and take out the samples

9. MAINTENANCE

9.1 Cleaning

CAUTION

Cleaning or sterilizing the centrifuge without complying with the instructions contained herein might cause damage to the centrifuge.

Centrifuge


The color of housing might change and the label thereon might fall off if the centrifuge is exposed to ultraviolet for a prolonged period of time. After use, please cover the centrifuge with cloth to avoid exposure to light.

Please clean the centrifuge using cloth or sponge soaked with neutral cleaning agent. The centrifuge may be sterilized using cloth soaked with 70% alcohol.

Centrifugal chamber

- Never pour water or other solvents directly into the centrifugal chamber, otherwise these solvents might enter into the drive unit and cause corrosion or damage to the bearing.
- Clean the chamber using cloth or sponge soaked with neutral cleaning agent and sterilize the centrifuge using cloth soaked with 70% alcohol

6. OPERATION PANEL

Markings	Name	Function
	Blood program key	Press this key to load the parameter settings for blood separation.
	Urine program key	Press this key to load the parameter settings for urine separation.
	Fecal program key	Press this key to load the parameter settings for fecal separation.
	SET key	Press this key to choose the parameters to be input.
	Parameter input key	Parameter increase
	Parameter input key	Parameter decrease
	ENTER key	Parameter decrease
	Start/Stop key	When the speed is zero, it will start running. While the centrifuge is running, press this key to stop.

Speed display area

Time display area

The screen display is shown in Figure 6.2, in which the speed is set at 3000rpm and the set running time is 88 minutes.

7. ROTOR PREPARATION

7.1 Prepare the samples

7.2 Transfer the samples into tubes

CAUTION

Do not over fill the tube with sample and avoid liquid leakage.
Do not over load the rotor.

7.3 Keep the tubes in balance

■ Although this centrifuge may be used with visually confirmed balance, we recommend that samples be weighed using a balance to ensure balanced centrifuge tube in order to prolong the service life of centrifuge, with the deviation of no more than $\pm 1.5g$.

■ Never intentionally run the centrifuge with rotor imbalance, no matter within the tolerance or not.

7.4 Inspect the rotor

■ Check the rotor for damages before using

■ If any damage is found on the rotor, please stop using it.

■ Only use rotor from the manufacturer

7.5 Symmetrically load centrifuge tubes onto the rotor

■ Make sure the rotor and motor shaft are tightened. Otherwise, the rotor may be moved off while rotating and cause damage to the centrifuge and rotor.

■ Never press the unlocking button while the rotor is in operation

8. OPERATION

PRECAUCIÓN

While the centrifuge is in operation, don't move the centrifuge or lean against it.

Before running the centrifuge, keep the centrifugal chamber clean and be sure to remove any foreign object from the centrifugal chamber, such as fragments of centrifuge tube.

If any abnormal phenomenon such as strange noise occurs during the operation, please immediately stop the centrifuge, contact the service center and inform them of the error code on the LCD display

8.1 Normal operation

Turn on the power switch and the self-test will launch automatically. Then the most recent program setting is loaded and displayed as shown in Figure 8.1

Speed display area

Time display area

■ Speed settings: 3000 rpm

■ Time settings: 88 min.

1. Set the operating parameters

Press the SET key gently to switch between different parameters. The parameter will blink on the LCD display when it is chosen. Press the parameter input keys to increase or decrease the selected parameter. Press and hold the parameter input keys, the value of the selected parameter will change faster. Press the parameter input keys briefly to change the parameter value slowly. The increment is 100rpm for the rotation speed, and 1 minute for the timer. Press the ENTER key to confirm the parameters.

Set the rotation speed

■ Press the SET key to choose the rotation speed in rpm.

■ When the speed parameter is selected, it will blink on the LCD screen.

■ Increase the speed from 300rpm by 100rpm each time.

■ Press the parameter input key to increase the parameter; press the parameter input key to decrease the parameter.

■ Press and hold the parameter input key or to accelerate the parameter change.

Set the running time

■ Press the SET key to make the time value blink and enter into the time setting mode.

■ Press the parameter input key or input the set time within a range of 1-99 min.

■ HD refers to unlimited time for continuous running.

Preset program setting

■ Press the blood program BLOOD key to load parameter settings for the blood sample.

■ Press the urine program URINE key to load the parameter settings for the urine sample.

■ Press the fecal program FECAL key to load the parameter settings of the fecal sample.

Modify the preset mode parameters

■ After setting the speed and running time, press and hold the SET key + blood program BLOOD key for over 5s to save the current parameter settings as the blood presetting program.

■ After setting the speed and running time, press and hold the SET key + urine program URINE key for over 5s to save the current parameter settings as the urine presetting program.

■ After setting the speed and running time, press the SET key + fecal program FECAL key for over 5s to save the current parameter settings as the fecal presetting program.

2. Start running

Press the running key START/STOP to start operation

■ The upper cover is locked, the rotor starts rotating and the red LED indicator is lit.

■ Begin timing after startup, with the time display indicating the remaining running time.