CENTRÍFUGA ANGULAR DIGITAL RS-0412
RS-0412 DIGITAL ANGLE CENTRIFUGE
CENTRIFUGEUSE ANGULAIRE NUMÉRIQUE
RS-0412

RS Lab

Ref. GBF013

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Spanish 2-15
 English 16-29
 French 30-43

INDEX OF CONTENTS

Safety Instructions..... 16
 1. Intended use 18
 2. Specifications 18
 3. Required Operational Conditions 19
 3.1 Basic operational conditions 19
 3.2 Transport and storage conditions 19
 4. Installation 19
 4.1 Location 19
 4.2 Connection of the power cord and grounding 19
 5. Structure 20
 6. Control panel 20
 7. Rotor Preparation 21
 7.1 Prepare the samples 21
 7.2 Inject the samples into tubes 21
 7.3 Keep the tubes balanced 21
 7.4 Inspect the rotor 21
 7.5 Symmetrically load centrifuge tubes into rotor 21
 8. Operation 22
 8.1 Normal operation 22
 8.2 RCF operation 24
 8.3 Pulse operation 24
 9. Maintenance 24
 9.1 Cleaning 24
 9.2 Rotor Installation 25
 10. Troubleshooting 26
 10.1 Possible problems and solutions 26
 10.2 How to open the door 26
 11. Instructions for the rotor and tubes 27
 11.1 Rotor 27
 11.2 Tubes 28
 12. Calculate RCF 29

SAFETY REMINDER

Common safety precautions
 Carefully read the following safety precautions for a thorough understanding.
 ■ Follow the instructions and procedures described in this manual to operate this centrifuge safely.
 ■ Carefully read all safety messages in this manual and the safety instructions on the instrument.
 ■ Safety messages are labeled as indicated below. They are in combination with signal words of “WARNING” and “CAUTION” with the safety alert symbol to call your attention to items or operations that could be dangerous to you or other persons using this instrument. The definitions of signal words are as follows:

3. Autoclaving PA, PC and PP tubes

PA begins softening at about 120°C, PC and PP at about 130°C. Autoclave PA tubes at 115°C (0.7kg/cm²) for 30 minutes, PC and PP tubes at 121°C (1.0kg/cm²) for 20 minutes. If a certain temperature is exceeded, the tubes may be deformed.

When use a sterilizing chamber, please operate as follows:

- Place tubes in vertical position, mouths upward. If tubes are placed sideways, they may deform into an oval shape due to gravity.
- Remove locking nut and lid to prevent from deformation or rupture.
- Wait until the sterilizing chamber cools down to the room temperature before removing tubes.

4. Conditions and life expectancy of tubes

The life expectancy of plastic tubes depends on the characteristics of samples, speed of the rotor used, temperature applied and so on. When the plastic tubes are used for ordinary aqueous samples (pH between

5.0 and 9.0), their life expectancies are defined as follows. Be operated at the maximum speed:

- High quality tubes (PA, PC, PP): 30-50 operations
- Ordinary tubes (PA, PC, PP): around 10 operations
- Life expectancy of tubes also depends on the pretreatment conditions such as cleaning and sterilization, lifetime can be cut down.

Notice: Do not use damaged or cracked tubes.

12. CALCULATE RCF

The Relative Centrifugal Force (RCF) can be determined with the following calculation formula:

$$RCF = 1,118 \times r \times n^2 \times 10^{-5}$$

r- rotating radius, unit: cm
 n- rotating speed, unit: rpm

3. Notice

The centrifuge can separate samples with a density lower than 2.0g/ml. If the sample density is over 2.0g/ml, please calculate allowable speed depending on the following formula.

Allow speed (rpm) = Max. Speed × (2,0 (g/ml) /densidad muestra (g/ml)) ¹/₂

4. Autoclaving

The rotor is made of plastic, cannot be high-pressure sterilization and UV irradiation, only ordinary sterilization can be used.

11.2 Tubes

1. Cleaning and sterilizing tubes

Condition			Material		
			PA	PC	PP
Cleaning	Cleaning fluids	Acidic (pH5 or lower)	X	X	X
		Acidic (higher than pH5)	O	O	O
		Alkaline (higher than pH9)	O	X	O
		Alkaline (pH9 or lower)	O	O	O
		Neutral (pH7)	O	O	O
		Warm water(up to 70°C)	O	O	O
	Ultrasonic	Neutral detergent (pH7)	O	O	O
Sterilization	Autoclaving	115°C (0,7kg/cm²) 30min	O	O	O
		121°C (1,0kg/cm²) 20min	X	O	O
		126°C (1,4kg/cm²) 15min	X	X	X
	Boiling	15 to 30 minutes	O	O	O
	Ultraviolet sterilization	200-300nm	X	X	X
	Gas sterilization		Ethylene oxide	O	X
Formaldehyde			O	O	O

PA: Polyallomer PC: Polycarbonate PP: Polypropylene

2. Cleaning PC tubes

PC material is low in chemical resistance against alkaline solutions. Avoid using neutral detergents with pH higher than 9. Note that pH of some neutral detergents is still higher than 9 even if diluted according to the manufacturer’s instructions. Use detergent with its pH between 7 and 9.

WARNING: Personal Danger

Warning notes indicate any condition or practice, which if not strictly observed, could result in personal injury or possible death.

CAUTION: Possible damage to instrument

Caution notes indicate any condition or practice, which if not strictly observed or remedied, could result in damage or destruction of the instrument

NOTE: Notes indicate an area or subject of special merit, emphasizing either the product’s capability or common errors in operation or maintenance.

- Do not operate this centrifuge in any manner not described in this User manual. When in doubt or have any troubles with this centrifuge, ASK FOR HELP.
- The precautions described in this User manual are carefully developed in an attempt to cover all the possible risks. However, it is also important that you are alert for unexpected incidents. Be carefully operating this centrifuge.

WARNING

- This centrifuge is not explosion-proof. Never use explosive or flammable samples.
- Do not install the centrifuge in or near places where inflammable gases are generated or chemicals are stored.
- Do not place dangerous material within 30cm around the centrifuge.
- Make sure to prepare necessary safety measures before using samples that are toxic, radioactive or contaminated with pathogenic micro-organisms at your own responsibility.
- If the instrument, rotor and/or accessories that has been contaminated by solutions with toxic, radioactive or pathogenic materials, clean it according to the decontamination procedure that you are specified.
- If you require services at site, please sterilize and decontaminate it in advance, and then notice the service center involved in the details of the particular materials.
- Do not handle the power cord or turn on or off the POWER switch with wet hands to void electrical shocks.
- For safety purposes, do not enter within 30cm around this centrifuge while it is in operation.
- While the rotor is rotating, never forcedly release the door lock.
- Unauthorized repairs, disassembly, and other services to the centrifuge except by our service center are strictly prohibited.

- This centrifuge must be located on one firm and level table.
- Make sure the centrifuge is horizontal before running.
- Make sure the angle between the door and cover is greater than 70 degrees when open the door.
- Be careful not put your fingers or hands between the door and cover when the door off.
- Do not move or relocate this centrifuge while it is running.
- If fluid spills in the rotor chamber, please promptly clean and dry with a dry cloth to avoid sample contamination.
- Ensure to remove any objects and fragments of the tubes dropped inside the rotor chamber before running this centrifuge.

■ Cautions on rotor

- (1) Always check for damages on the rotor surface before using it. Do not use the rotor if an abnormality is found.
- (2) Do not set the centrifuge speed beyond the allowable speed of the rotor kits (rotor or adapters).
- (3) Do not exceed the allowable imbalance.
- (4) Use the rotor and tubes within their actual capacities.
- If any abnormal condition occurs during operation, please stop it immediately and contact your distributor to request technical service.
- Vibrations are likely to damage the centrifuge, contact the service center if abnormality observed.

1. INTENDED USE

The centrifuge is used for separation of human blood and urine samples. Operator should be trained before using the centrifuge by reading this manual.

2. SPECIFICATIONS

Maximum speed	4500rpm (Range: 300-4500rpm, increase: 100rpm)
Maximum RCF	2400×g, increase: 100×g
Maximum capacity	12×10ml, 8×15ml
Timer	30-59 seconds / 1-99 minutes / Continuous operation (HOLD)
Noise	56dB
Driving Motor	Brushless DC motor
Safety devices	Lid locking, overspeed detector, automatic internal diagnostics
Power requirements	Single-phase current 110V-240V, 50/60Hz
Ambient conditions: -Set-up site -Altitude -Ambient temperature -Max. humidity -Excess-voltage category -Pollution degree	Indoor only Up to 2000 m above sea level 2°C ~ 40°C 80% II 2
Device protection class	I
EMC	EN/IEC 61326-1 / Class A - FCC / Class A
Dimensions	(D)301 ×(W)354 ×(H)217
Weight	6 kg
Additional features	Speed/RCF switch, short-time run function (pulse operation), sound-alert function

11. INSTRUCTIONS FOR THE ROTOR AND TUBES


CAUTION

Read the instructions thoroughly, to properly load and use rotor.
Do not exceed the allowable maximum speed of tubes and adapters.

11.1 Rotor

1. Rotor structure

2. Rotor, tubes and adapters

Rotor	Rotor No.	Tubes	Dimension (φ×L mm)	Use of adapter	Maximum speed (rpm)	Maximum RCF (×g)
Angular	Falcon type 15ml	8	17×120	No	4500	2400
	Vacutainer type 5ml	12	13×75	Yes, with mat	4500	2218
	Vacutainer type 7ml	12	13×100	Yes, without mat	4500	2400
	Vacutainer type 10ml	12	16×107	No	4500	2400
	Collection tube (EU) 7,5-8,2ml	12	15×92	Yes, without mat	4500	2400

10. TROUBLESHOOTING

10.1 Possible problems and solutions

This centrifuge has a self-diagnostic function. If a problem occurs, an error/warning code will be displayed on the time display screen, and the operator can determine the malfunction with the alarm code below.

Symptom	Causes	Solutions
Nothing appears on the screen when the POWER is turned on.	Building power circuit breaker trips.	Remove the trouble and turn on the POWER.
Abnormal vibration	Rotor do not match with spindle ·Samples are imbalance	Install again the rotor · Weighting scales, install symmetrically
Alarm code appeared on the time display screen	E-02 Door fault	The door opened in running. · The button is pressed while the door opening.
	E-06 Set wrong speed	The setting speed exceed the allowable range.
	E-10~86	·Read the service manual

Alarm codes E-1~E-9 are related to incorrect operation/programming. You can continue running the centrifuge after implementing corrective procedures.

10.2 How to open the door

1. In the case of power on

The door just can be opened while the power is on and rotor stops rotating

- Turn on the power switch, release the door automatically.
- The door will be released automatically once the operation is finished.
- It is available to release the door by press DOOR button once the rotor stops.

2. In the case of power outage

The door cannot be opened automatically if there is a power outage. It is available to be opened manually as follows.

- Ensure if the rotor has stopped rotating.
 - Listen carefully to ensure no rotating sound can be heard.
- Insert a screwdriver into the hole to open the door.
 - Hole is located on the right side of the unit.
 - Insert a screwdriver into the hole and push forward to release the door

3. REQUIRED OPERATIONAL CONDITIONS

3.1 Basic operational conditions

- Power: 110V-240V, 50/60Hz
- Ambient temperature: 2°C~40°C
- Relative humidity: ≤80%
- No vibration and airflow around
- No electric dust, explosive and corrosive gases around

3.2 Transport and storage conditions

- Storage temperature: -40°C~55°C
- Relative humidity: ≤93%

4. INSTALLATION

This section describes the instructions that you should abide when install the centrifuge to ensure your safety and the optimum performance. Before moving the centrifuge, the rotor must be removed.

WARNING

Improper power supply may damage centrifuge.
Make sure the power source conforms to the required power supply before connecting

4.1 Location

- Place this centrifuge on a firm, flat and level surface, ensure the four feet of this centrifuge stand on the counter firmly. Avoid installing on a slippery surface or surface prone to vibration.
- Ideal ambient temperature is 20°C±5°C, avoid placing the centrifuge in direct sunlight if temperature exceeds 30°C.
- Keep clear of the centrifuge at least 10cm on both sides and at least 30cm behind it to guarantee the cooling efficiency.
- Keep away from heat or water to avoid sample temperature issues or centrifuge failures.

4.2 Connection of the power cord and grounding


WARNING


To avoid electrical shocks, ensure your hands are dry when touching the power cord.
This centrifuge must be grounded properly.

The outlet must be grounded according to local safety requirements.

5. STRUCTURE

6. CONTROL PANEL

Button	Name	Function
	Select button	Press the button to choose the program which you want to modify.
	Pulse button	The speed can be accelerated and held at the speed when pressing Pulse on.
	Door opening button	Press the button to open the door The button is not available when the centrifuge is running.
	Start/Stop button	Press the button to start running. The centrifuge will brake to stop running if pressed during centrifugation.
	Parameter button	Clockwise rotate to increase program values. Rotate anti-clockwise to decrease parameter values. Press the button, shift between speed and RCF display

Drive shaft

We recommend regular maintenance for drive shaft. You can wipe the drive shaft with soft cloth, and then apply a thin coat of silicon grease.

Door


Clean and sterilize the door using the same method as the section (1) above

Rotor


If not in use for a long term, then detach the rotor, clean it if necessary and turn upside down to dry the tube holes and keep clean.
 For sample leaks in the rotor, rinse the rotor with water. Apply a thin coat of silicon grease to the rotor when it is completely dry.
 The rotor should be checked every 3 months to ensure the tube and rotor holes keep are clean and apply a thin coat of silicon grease.

9.2 Rotor installation

Installation

Adjustment

Observe here

Before lock the rotor, rotate it, and observe carefully if there is obvious vibration, if so, please take off the rotor, turn some angle and install it again, until the rotor rotates smoothly, then, lock it firmly.

8.2 RCF operation

Turn on the power switch.

Set a RCF (Relative Centrifugal Force) value

- Press the select button  and choose speed unit $\times g$, the speed symbol will flash into RCF value input status.
- If no button is pressed after the speed value has flashed after 5 seconds, the input mode will be shut down.
- Rotate program button  to input a RCF value, RCF increment is $100 \times g$.

Set operation time

Please refer to the section 8.1.

8.3 Pulse operation

This function is used to remove the residual samples adhered to the interior of the tubes or for quick spins.

Note: The button works only while the rotor stopped and the door is locked.

- Turn on the power switch and load the rotor to the shaft, tighten the rotor lid and make sure it is in secured position, and then close the door.
- The centrifuge goes into preparation mode and displays last running values.
- Press  knob and hold, the centrifuge will speed up to the setting speed. While releasing the  knob during acceleration, the centrifuge will start to decelerate and stop

9. MAINTENANCE

9.1 Cleaning

CAUTION

If do not follow the recommended instructions for cleaning or disinfecting this may damage the centrifuge.

Centrifuge

If the centrifuge is exposed to ultraviolet rays for a long time, the color of the door may be changed. After using, cover the centrifuge with a piece of cloth to protect it from direct exposure.


If the centrifuge needs cleaning, clean it with a cloth or sponge moistened with a neutral detergent solution.

Sterilize the centrifuge by wiping with a cloth moistened with 70% ethanol solution.

Rotor chamber

Do not directly pour water, neutral detergent or disinfectant solution into the rotor chamber, otherwise fluids may leak into the drive units and cause corrosion or deterioration to the bearings

- If the rotor needs cleaning, clean with cloth or sponge moistened with a neutral detergent solution
- Sterilize the centrifuge by wiping with a cloth moistened with 70% ethanol solution.

Main interface is as figure 6-2. The speed is set to be 4500rpm, the door lock is released, and the running time is 12 minutes. When speed symbol  is rotating, this indicates the centrifuge is running. Time symbol  displays the ratio of working to time setting.

7. ROTOR PREPARATION

7.1 Prepare the samples

7.2 Inject the samples into tubes

CAUTION

Do not overload samples into the centrifuge which will cause leaking.
Do not exceed the actual capacity allowed in the User's Manual.

7.3 Keep the tubes balanced

- Although the centrifuge can accept sample balancing by eye, we recommend that you keep this centrifuge in a well-balanced condition to extend its life expectancy.
- Never intentionally run the centrifuge under an unbalanced condition even though the allowable imbalance is not exceeded

7.4 Inspect the rotor

- Check the rotor for corrosion or scratches before using
- If any abnormality such as corrosion or scratches are found, stop using the rotor and contact our service center.
- Only manufacturer's rotors must be used with the unit.

7.5 Symmetrically load centrifuge tubes into rotor

- Make sure the rotor and motor shaft are tightened. Otherwise, the rotor may be moved off while rotating and cause damage to the centrifuge and rotor.

8. OPERATION


CAUTION

Do not push or lean against the centrifuge while it is running.
Do not run the centrifuge when fragments or sample solutions are left in the centrifuge chamber. Always keep the centrifugal chamber clean.
If the centrifuge makes strange noise during operation, stop it immediately and contact our service center. Notify them of the warning code if displayed

8.1 Normal operation

Turn on the power switch, centrifuge will display the running interface last time after passing the self-diagnostic checks, see figure 9-1 below:

- Speed: 4500rpm.
- Running time: 12 minutes.
- The door lock is released.

1. Set the operation programs

Press the **SELECT** button to select required program. The parameter can be modified when the program is flashing. Rotate the program button clockwise to increase parameter value. Rotate the program button clockwise to decrease parameter value. Rotate the program button faster, and the parameter value will increase faster. The minimum speed increment is 100 rpm, the minimum time increment is 1 second.

Set the speed

- Press the select button **SELECT** until the speed rpm is displayed.
- When the speed button is selected, the speed symbol will flash the speed value.
- The minimum speed value you can set 500rpm, the minimum increment is 100 rpm.
- Rotate program button clockwise to increase speed value. Rotate the program button anti-clockwise to decrease speed value.
- You can speed-up set the speed value by rotating program button quickly.
- There is a circulating function to increase/decrease the speed values. Rotate the program button clockwise to change settings from small → large → maximum → minimum. Rotate the program button anti-clockwise to change settings from large → small → minimum → maximum.

Set the time

- Press select button **SELECT**, time value flashes in the time setting mode.
- Rotate the program button to set running time from 30 seconds to 99 minutes.
- When the time displays HD, this is a continuous running mode.

2. Start the operation

Press button **START STOP** to start running

- Timer will operate once the speed setting value is reached; the screen displays the remaining run time.

View and modify the operation programs

- Operation programs can be modified after the centrifuge reaches the set speed.
- Pressing the select button **SELECT**, returns the display to the program interface and displays setting programs. Press the select button **SELECT** to the desired program. When flashing, rotate parameter button to modify values. Release the button after 5 seconds, and the centrifuge will return to normal operation mode and run according to the new value.

Warning display

- If an error occurs during the operation, the centrifuge will brake to stop automatically, and display the error code on the time/display area. The error code can be checked in the table 10-1, and corrective actions can be applied accordingly.

3. End the operation

The centrifuge will brake when it reaches the set time or **START STOP** button is pressed.

- When the rotor stops rotating, the centrifuge will start beeping to alert the operation has finished

Open the door

- The door can be released automatically when the operation has stopped.
- With the door closed, you are able to press the **DOOR** button to open it.
- After ending the operation, the program will store the setting parameters of this operation and will recall these parameters when restarting the program.

Open the door and take out the rotor and samples.