MANUAL CUBETA DE ELECTROFORESIS VERTICAL
VERTICAL ELECTROPHORESIS CELL
CUVE ÉLECTROPHORÈSE VERTICALE

Nahita

Ref. ZFD005

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
- Put the core with installed plates in the gel casting frame, hang over the hooks on both sides of the core, press down the wrenches to make the bottom of the glass plates be sealed by the rubber. See image below:

6. GEL POURING


- Pour the gel into the gel blocks carefully so as not to generate bubbles.
- Insert the comb carefully and make sure there is no air bubbles under the ends of the comb teeth.
- Leave the gel by itself and wait it to concrete.
- Pull the comb(s) out carefully and move the core with the gel to the main tank.

7. RUNNING THE GEL

- Put the core into the bottom tank.
- Fill 1x buffer into the upper tank, as well as the bottom tank.
- Load the samples into the wells with pipettes and take care not to damage the wells or induce bubbles.
- Carefully cover the tank with lid and connect it with a power supply.
- Typically gels run under 150V- 200V. Be noted that, generally higher voltage enables faster electrophoresis but poorer quality of sample resolution.
- Run electrophoresis.

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INDEX OF CONTENTS

1. SPECIFICATIONS	5	5. SETTING UP	6
2. ACCESSORIES INCLUDED	5	6. GEL POURING	7
3. SAFETY INSTRUCTIONS.....	6	7. RUNNING THE GEL.....	7
4. CARE AND MAINTENANCE	6		

1. SPECIFICATIONS

Code	ZFD005
Number of gels	1 ó 2
Gel size (WxL)	82x88 mm
Glass plate size (WxL)	100x100 mm
Gel thickness	0,75 mm / 1,0 mm / 1,5 mm
Number of samples	11 ó 15 per gel
Buffer volume	750 mL
Dimensions (LxWxH)	150x120x115 mm
Net weight	1,5 kg

2. ACCESSORIES INCLUDED

- Combs: 12
- Wedged-shaped plates: 2
- Gel shovel: 1
- Dummy plate: 1
- Spacer plates with 0.75 mm spacer: 2
- Spacer plates with 1.0 mm spacer: 2
- Spacer plates with 1.5 mm spacer: 2
- Concave plates: 4

3. SAFETY INSTRUCTIONS

Safety precaution when used correctly, these units pose not health risk. However, these units can deliver dangerous levels of electricity and are to be operates only by qualified personnel following the guidelines laid out in this instruction manual. Anyone intending to use this equipment should read the complete manual thoroughly. The units must never be used if there is any sign of damage to external tank or lid.

4. CARE AND MAINTENANCE

Cleaning vertical units are best cleaned using warm water and a mild detergent. Water at temperatures above 60°C can cause damage to the unit and components. The tank should be thoroughly rinsed with warm water or distilled water to prevent buildup of salts but care should be taken not to damage the enclosed electrode and vigorous cleaning is not necessary or advised. Air drying is preferably before use. The units should only be cleaned with warm water with a mild concentration of soap or other mild detergent. Compatible detergents include dishwashing liquid, Hexane and Aliphatic hydrocarbons. The units should not be left in detergents for more than 30 minutes. The unit should never come into contact with the following cleaning agents, these will cause irreversible and accumulative damage: Acetone, Phenol, Chloroform, Carbon tetrachloride, Methanol, Ethanol, Isopropyl alcohol. RNase decontamination can be performed using the following protocol: Clean the units with a mild detergent as described above. Wash with 3% hydrogen peroxide (H₂O₂) for 10 minutes. Rinsed with 0.1% diethyl pyrocarbonate (DEPC) treated distilled water.

Caution: DEPC is a suspected carcinogen. Always take the necessary precautions when using.

5. SETTING UP

5.1. Instructions for fitting electrode cables

- Note the position of the lid on the unit to keep the correct polarity and the correct orientation of the cables; black is negative and red positive.
- Remove the lid from the unit. Note if the lid is not removed, fitting the cables may result in un-tightening of the gold plug and damage to the electrode.
- Screw the cables into the tapped holes as fully as possible to that there is no gap between the lid and the leading edge of the cable fitting.
- Refit the lid.

5.2. Vertical Gel Casting

- Clean glass plates, flush with distilled water, let them dry in the air.
- Overlap a spacer plate with a notched one, put the two glasses into the core, keep their bottoms contacting the table and insert a cuneal plate closely outside the two plates.
- Use the same method to complete the other side.

Notice:

- **Spacer plate should be put outside and the notched one inside, otherwise the experiment can't be carried on.**
- **The operation must be done on a level surface.**
- **The bottom ends of glass plates must be perfectly aligned.**