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FREE FALL APPARATUS WITH DIGITAL COUNTER

Introduction

Free fall apparatus with digital counter has been designed to calculate the value of acceleration due to gravity. The digital counter device can be also used as a stopwatch.

The free fall apparatus consists of an electromagnet mounted in a plastic case to hold/release a metallic ball. As soon as the ball gets released from the electromagnet, the timer starts counting; when the ball hits the platform with an incorporated impact sensor, the timer stops. The time taken by the sphere from top to bottom will be displayed on the LED display of digital counter to calculate the acceleration due to gravity.

Product overview

11. 5 digit LED display
12. Indicators for timer in μ second, millisecond and second mode
13. Output sockets for the free fall base
14. Push button to disconnect the electromagnet
15. Output sockets for connecting the digital counter with the electromagnet
16. Input sockets for the free fall pulse from impact sensor

17. Selector to select the mode edge triggered or level triggered when digital counter is used as a stopwatch
18. Push button to reset the timer
19. Start/Stop switch
20. Platform with impact sensor

Features

- Indicators for the timer result in μ seconds, millisecond and second mode.
- 5 digit LED display
- SMPS based, universal input voltage
- High accuracy impact sensor

Technical specifications

Input voltage	90-240 V AC
Frequency	50-60 Hz
Current rating	20 mA
Power	5 W (max.)
Temperature range	Room temperature
Humidity	40-70%
Protection grade	IP20

Working procedure

Calculation of the acceleration due to gravity

- Connect the “9 V DC” output socket of digital counter to “9 V DC” input sockets of free fall stand using patch cords.
- Connect the “Electromagnet” sockets of the digital counter to the electromagnet sockets of the free fall stand, using patch cords.
- Connect the “Sensor plate” sockets of the digital counter to the “To Timer” sockets of the free fall stand, using patch cords.
- Plug-in the power cord with an input voltage of 90-240 V AC.
- Put into contact the metallic ball with the electromagnet (solenoid) mounted on the free fall stand. Due to magnetization, the electromagnet will hold the metallic ball.
- Press the corresponding button in the digital counter to release the ball from the electromagnet.
- Once the ball is released, the timer will start and the ball will fall on the sensor plate of free fall stand.
- Impact sensor placed in the platform will send a pulse from free fall stand to the digital counter to stop the timer.
- Read out the timer taken by the ball between releasing from electromagnet and fall on sensor plate.
- Fill in below table with obtained data

Seri es nº	Inicial distance between ball and sensor plate (S)	Time taken by the ball to reach at the sensor plate (t)	$g = 2S/t^2$
1.			
2.			
3.			

Formula used

$$S = ut + \frac{1}{2} at^2$$

Where S is the initial distance between the metallic ball and the sensor plate

t , time taken by the ball to reach the sensor plate once it is released from the electromagnet

u , initial speed of the ball

a , acceleration of the metallic ball; in this case $a = g$ since it is a free fall vertical movement

As the initial speed of metallic ball is zero, therefore

$$S = \frac{1}{2} gt^2 \rightarrow g = 2S/t^2$$

Press reset switch after every observation to obtain accurate results.

Use digital timer as a stopwatch

- Do not connect the free fall stand with the digital counter
- Plug in the power cord with an input voltage of 90-240 V AC

- Press the Start/Stop button on the front panel to start and stop the timer.
- Below Start/Stop button there is a slider switch to select the timer mode. According to the symbols there are two interruption modes:
 - o Positive edge triggered (slider on the upper side): in this mode, pressing once Start/Stop button, the timer will start and continues itself until pressing the button again to stop the timer.
 - o Level triggered (slider switch on the lower side): in this mode, timer will continue itself while keeping pressed Start/Stop button.

Precaution and maintenance

- In case of fuse blown, replace the 0.5 A fuse in the built-in cartridge of inlet socket.
-  This symbol is mentioned for AC MAINS voltage. Do not touch.
-  This symbol is mentioned for the earthing points in the apparatus.
- In case of not working properly do not open the apparatus to repair it by yourself; contact with technical service department.