**APARATO DE CAÍDA LIBRE CON CONTADOR
DIGITAL**
FREE FALL APPARATUS WITH DIGITAL COUNTER

Ref. QLB003

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
This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

Equipment

- Digital timer for g by free fall
- Secondary device (any controlled device) (not included)
- Connecting leads

Procedure

1. Connect the connection as shown in the figure.
2. Slide the slider switch to the level triggering mode .
3. The timer will take +5V to +9V as a high (1).
4. When the pulse will goes to high then the timer will starts to count the time until the pulse/signal will goes to low or zero level position.
5. The reading on the display shows the time taken by the pulse at high position.

Precaution and Maintenance:

- In case of fuse blown, replace the 0.5A fuse, in the inbuilt cartridge of inlet.
- Symbol is mentioned for AC MAINS voltage. Do not touch.
- Symbol is mentioned for the earthing points in the apparatus.
- Do not open the apparatus, in case of not working properly, as this apparatus is non repairable at user end.

Instructions on environment protection

- At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances.
- It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.
- The materials are recyclable as mentioned in its marking.
- By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.
- Please inquire at the community administration for the authorized disposal location.

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1. Press the push button mounted on the device (9) to release the ball from electromagnet.
2. Once the ball will release, the timer will start and it stops when the ball hit to the sensor plate, the mercury tilt switch mounted under the sensor plate will release a pulse from free fall stand to the device to stop the timer.
3. Read out the timer taken by the ball between releasing from solenoid and fall on sensor plate.
4. Press Reset Switch after every observation to obtain the accurate results.
5. Repeat the experiment 2-3 time.

Put the data in the below table:

S. No.	Distance in between ball & sensor plate (S)	Time taken by ball to reach at the sensor plate (t)	$g = \frac{2h}{t^2}$
1			
2			
3			

Formula Used

$$s = ut + \frac{1}{2}at^2$$

s = Distance between metallic sphere and sensor plate

t = time taken by the ball to reach at the sensor plate

u = initial speed of the ball

a = acceleration of the metallic sphere, Now a = g in vertical direction

$$s = \frac{1}{2}gt^2$$

As the initial speed of the metallic sphere is zero

Put all the value in the formula and calculate the acceleration due to gravity.

Working in the level triggering mode (for other device)

1. Indicators for timer in μ second, millisecond, second mode.
2. 5 digit LED display.
3. Push button to reset the timer.
4. Start-Stop switch.
5. Selector to select the mode edge triggered or level triggered.
 - Edge trigger for free fall apparatus
 - Level trigger (high) for other compatible devices.
6. Input from the sensor plate
7. DC output supply for electromagnet.
8. 9V DC output supply for biasing the free fall apparatus.
9. Push button to release the ball from electromagnet.
10. Electromagnet input terminals.
11. 9V DC input.
12. Steel ball
13. Connecting leads
14. Output from sensor plate.
15. LED indicator of output pulse.

6. Specification

Input Voltage	90-240 AC
Frequency	50-60Hz
Current Rating	20mA
Power	5W (max)
Temperature	Room temperature
Humidity	40% to 70%
Ingress protection code	IP 20

7. Working procedure (for edge mode)

To count the value of "g".

1. Conecte las tomas "9V DC" del aparato a las tomas "9V DC input" del soporte de caída libre 1- Connect the "9V DC" sockets of the device to the "9V DC input" sockets of free fall stand using connecting leads.
2. Connect the "Electromagnet" sockets of the device to the electromagnet terminals of the free fall stand.
3. Connect the "Sensor Plate" sockets of the device to the "Sensor Output" sockets of the free fall stand.
4. Plug-in the power cord with an input voltage of 90-240V AC mains.
5. Switch On the device with the help of rocker switch provided on the back panel.
6. Slide the slider switch on upper side i.e edge triggering mode.
7. Touch the MS ball to the solenoid mounted on the free fall stand. Due to magnetization, that solenoid will hold the ball with itself.
8. Press the reset button to initialize the operations.

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
Level Triggering: In level triggering the circuit will become active when the gating or clock pulse is on a particular level. This level is decided by the designer. We can have a negative level triggering in which the circuit is active when the clock signal is low or a positive level triggering in which the circuit is active when the clock signal is high.


Edge Triggering: In edge triggering the circuit becomes active at negative or positive edge of the clock signal. For example if the circuit is positive edge triggered, it will take input at exactly the time in which the clock signal goes from low to high. Similarly input is taken at exactly the time in which the clock signal goes from high to low in negative edge triggering. But keep in mind after the input, it can be processed in all the time till the next input is taken.

4. Experiment

In this experiment we will determine the acceleration of gravity 'g' by measuring the time of flight for balls dropped from a known height. The acceleration due to gravity does not depend on the mass of the ball.

The falling body in this experiment will be a metal ball, which falls freely from the catch at the top of the apparatus to the pad below it. The apparatus is designed to record the time of flight for the ball. When the ball is in the catch, it forms part of an electrical circuit. When the ball is released, the circuit breaks and the timer start. The pad acts as a switch, such that when the ball hits it the timer is stopped.

5. Product overview

3. Theory

A free falling object is an object that is falling under the sole influence of gravity. Any object that is being acted upon only by the force of gravity is said to be in a state of free fall. There are two important motion characteristics that are true of free-falling objects:

- Free-falling objects do not encounter air resistance.
- All free-falling objects (on Earth) accelerate downwards at a rate of 9.8 m/s/s (often approximated as 10 m/s/s).

This object has a downward acceleration, denoted by (g), toward the center of earth. In order to calculate the value of the gravitational acceleration, we will use the free fall experimental setup illustrated in Figure 1. The steel ball is fixed to the releasing mechanism of the free fall apparatus. Allowing the ball to fall a fixed distance h toward the Receptor Plate, the free fall timer will compute the time elapsed for the ball to fall that distance.

The position of the steel ball, starting from rest at time $t = 0$ and undergoing constant downward acceleration along the y -direction, can be understood using the following equation

$$s = ut + \frac{1}{2}at^2$$

Where y_0 is the initial position and g is the acceleration due to gravity. If h is the falling distance the body has traveled from its starting point during time t , then above equation can be written

$$h = \frac{1}{2}gt^2$$

Where

s = height / displacement m (meters)

u = initial velocity ms^{-1}

v = final velocity ms^{-1}

g = acceleration due to gravity ms^{-2}

t = time s (seconds)

h = distance the ball falls through (m)

This equation shows that g can be easily determined by measuring free fall time's t as a function of height h . Therefore, we can determine the acceleration of gravity.

Triggering: This means making a circuit active. Making a circuit active means allowing the circuit to take input and give output. For example, supposed we have a flip-flop. When the circuit is not triggered, even if you give some input data, it will not change the data stored inside the flip-flop nor will it change the output Q or Q' . Now there are basically two types of triggering. The triggering is given in form of a clock pulse or gating signal. Depending upon the type of triggering mechanism used, the circuit will become active at specific states of the clock pulse.

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1. Description

Digital Counter “g” By Free Fall can be used to calculate the value of gravity. This apparatus can also be used as a timer and stopwatch. “g” By Free Fall apparatus is consisting an electromagnet (9 V DC) on the upper plate of the free fall stand. 4 mm color coded sockets are also mounted to supply the input voltage. Sensor plate (base plate) is consisting the mercury tilt switch for high accuracy.

2. Features

- SMPS based, universal input voltage.
- Timing in μ -second, millisecond, second with indicators (Auto-Change).
- 5-digit seven segment display.
- Positive edge trigger & level trigger mode.
- Sturdy push switches for stopwatch & reset.
- 4mm color coded safety sockets.