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

Index

1. General information	8
2. Packing list and accesories	8
3. Technical specification	8
4. Information about safety in the use of the equipment.....	8
5. Equipment description	9
6. Installation	9
7. Operation.....	9
8. Error messanges.....	12
9. Maintenance	12

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Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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1. General Information

- Handle the parcel with care.
- Check that contents coincide with the packing list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list and accessories

Units are supplied with their own manuals.

Immersion thermostats 3000540 & 3000543 need some accessories to operate.

Stainless steel buckets:

Code	Liters	Height(cm)	Width(cm)	Length(cm)
6000390	9	20	34	28
6000391	12	20	37	35
6000392	20	20	55	35
6000393	27	25	55	35
6001093	45	23	70	56

3. Technical specification

Part Number	3000543	3000548	3000549	3000540
Dimensions: High/wide/depth (mm)	280/180/190	360/550/350	410/550/350	280/180/190
Power supply voltage:	230 VAC 50Hz / 60Hz			
Electrical power (W):	2000W			
Weight: (Kg)	4	9,2	10,8	4
Fuse:	Replaceable (10x5mm)			
Temperature range*:	+ 5°C ... to 150°C			+100°C
Display resolution:	0.1°C			
Stability:	± 0.1°C			
H ₂ O pump flow:	12L/min			

* Always 5°C above room temperature

ATTENTION! Risk of burns

Take precautions
(use of thermal gloves, clamps, etc...)
to take samples from the bath.

They could be at a high temperature.

Electric Risk

**Risk of high
temperature contact**

4. Information about safety in the use of the equipment

4.1 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

4.2 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

4.3. Maintenance

Unplug the equipment from the mains before removing the oven's cover to handle inside.

The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.

4.4. Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.


Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

5. Equipment Description

Immersion thermostats are widely used on laboratories for thermal treatment of samples and bath's heating.

The equipment includes a circulation pump that can be used to improve the homogeneity within a bath or to perform a closed liquid circulation circuit for heating.

The unit could work in a continuous or timed operation. Once the operating time has elapsed, the unit stops (only **DIGITERM** models)

The level switch inhibits the thermostat operation if no liquid level is detected (only **DIGITERM** models)


Fig 7.1. Fixing clamp

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

6.1 Liquid tank

A bucket or tank will be needed to operate the thermostat. See J.P. SELECTA, s.a.u. products range or paragraph 2 of this manual.

Fix the thermostat in the tank with the fixing clamp.

Place the equipment on a flat level stable surface, suitable for the equipment weight.

Select a suitable fluid for working temperatures, water, silicone oil, antifreeze, etc.

To work in a comfortable way, please leave a free space of 20cm around the equipment.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied.


For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency match the one indicated in the equipment plate.

7. Operation


All the equipment functions are set with the touch screen keys. Gently press with your finger in the center of each key.

Use the fixing clamp (see Fig 7.1) in the rear part of the equipment.

7.1 Control panel

1. Stop / Start cycle Key
2. Reduce value key
3. Increase value key
4. Adjustment (change between temperature / time) (only **DIGITERM** models).

7.2 Starting up

1. Press the start switch and at the same time the screen will light up.
 2. Press the Up/Down keys (3/2) until the desired temperature value is selected in the green display (the above one).
 3. Press key (4) if you want to select an operating time (only **DIGITERM** models).
 4. Press key (1) to start the cycle.
 5. The heating symbol will light up and flash.
 6. Press key (1) again to stop the incubation process.
- Note: during operation, temperature or time adjustment value can be changed.

REMARKS:

- * When using the time function, the green display (the above one) will indicate the remaining time. If you want to see the target temperature value, press key (4) (only **DIGITERM** models).

7.3 Safety thermostat

The equipment has a safety element to prevent overheating.

In case the thermostat triggers, press the rear button to reset the device.


Fig 7.2 Rear part view


Fig 7.3 Cooler coil

7.6a Thermostat bulb for working at temperatures below 120°C

7.6b Thermostat bulb for working at temperatures above 120°C

7.4 Operation under ambient temperature (only DIGITERM models)

The unit allows to operate at temperatures under ambient (to a minimum of +25°C)

Coolant liquid circulation is needed, through the coolant coil (See figure 7.3).

For instance, to get the bath at 25°C, made circulate tap water through the coolant coil, so you can get an operation temperature on 25°C.

If tap water is under 18°C the unit could be set the water tank to 25°C.

7.5 External probe (only DIGITERM models)

Connect a Pt100 type probe to control the process from an external probe. The equipment automatically detects the probe change.


NOTE: the external probe must be located at a point where the heating element takes effect, otherwise the equipment will heat without control.

7.6 Bulb clamp change

If you are going to work at high temperatures (between 120°C and 150°C) with oils or other liquids with a higher density than water, it is recommended to change the position of the safety thermostat bulb to prevent the thermostat from jumping unnecessarily.

Use the attached clamping flange to modify the position of the bulb, as indicated in figure 7.6b.

- Remove the bulb from the retaining clips.
- Screw it to the resistance using the attached flange.


Fig 7.4 **INTERNAL** circulation

Circulation pump

The water circulation pump (or other liquid) could operate in two configuration modes:

- **INTERNAL circulation:** Water will be pumped from horizontal pump output (See fig 7.4).

Use this configuration to achieve a good temperature homogenization in the liquid in the vessel.

Fig 7.5 **EXTERNAL** circulation

- **EXTERNAL circulation:** Water will be pumped from vertical pump output (See fig 7.5).

Use this configuration to send the liquid to an external circulation circuit for tempering, cooling, etc ... an external equipment.

8. Error messages

P.FAIL Power failure. The oven continues with the process.

Err will appear in green on the top line of the display.

TSEG/SAFE

The safety thermostat is open.
(in the DIGTERM models, the buoy detects lack of water)

SENS

Open temperature probe.
The control does not detect the temperature probe.

HEAT

The resistance does not work.
The equipment detects that the temperature does not increase.

9. Maintenance

Before removing the oven lid to handle inside, disconnect it from the mains socket.

The oven internal electronic circuits manipulation by unauthorized personnel can cause damages that are difficult to repair. Be sure to take the equipment to one of the technical services authorized by J.P. SELECTA, s.a.u.

CLEANING:

We recommend the following products for cleaning the different pieces of the equipment:

Stainless steel cleaning: Alcohol.
Covers and plastics cleaning: Alcohol with a cotton cloth.