



GW2145



GW1160, GW4060



GW4190



GW2145 - GW1160
GW4060 - GW4190

GLASSWARE WASHER

USER MANUAL

READ THE USER MANUAL

This manual is an integral part of the appliance.

Take good care of it and keep it to hand throughout the appliance's life cycle.

This manual and all the information it contains must be read carefully before using the appliance.

Failure to read, misunderstanding or incorrect interpretation of the instructions provided in this manual may lead to misuse of the appliance, put the operator at risk and considerably reduce the appliance's performance.



Installation, maintenance and any repairs must be carried out by authorised technical staff.

Repairs of the appliance performed by unauthorised staff may put the user's safety at risk, and void the warranty cover.

Any components must always be replaced with genuine Smeg spare parts.

The use of the device in breach of the instructions provided by the manufacturer may jeopardise the specified protection level (safety of the appliance) and the warranty cover for it (see point 5.4.4 of IEC 61010-1:2001).

The manufacturer declines any liability for uses other than those described in this manual.



Consumables (detergents, air filters, thermal printer paper, etc.) are not covered by the warranty, except for any manufacturing defects.

The Warranty does not cover any parts found to be faulty due to negligent or careless use, improper use, failure to comply with the appliance's operating instructions, incorrect installation or maintenance, maintenance or repairs performed by unauthorised staff or made with non-genuine parts, damage during transport, or any circumstances in which the appliance's defects cannot be traced back to manufacturing faults. Any work relating to installation and connection to the intake and drain systems, and the maintenance work described in the operator's manuals, are also excluded from the Warranty.



Any accessories may not be installed on the appliance by the user; this must be done by Smeg authorised technical staff.

To request technical documentation relating to accessories from Smeg: instruments@smeg.it (ref. 5.4.4.c IEC61010-2-040:2005)

The information in this manual is provided for guidance only. The contents and the equipment described may be subject to change without notice.

19 390 4630 04	08/02/2016
Manual Ed.	Date

INFORMATION AND AFTER-SALES SERVICE FOR SMEG INSTRUMENTS DIVISION PRODUCTS

Our Sales Department staff will be able to provide you with information about prices and special offers. Our After-Sales Department will be able to provide you with guidance about keeping your appliance functioning correctly and put you in touch with your nearest authorised Service Centre.

**INTERNATIONAL CUSTOMERS
PLEASE CONTACT YOUR LOCAL DISTRIBUTOR**

instruments@smeg.it

service.instruments@smeg.it

Our whole product offering can be viewed at:

www.smeg-instruments.com

CONTENTS

1	INTRODUCTION.....	6
1.1	GW FAMILY - INTENDED USE	7
1.2	DEFINITION: “RESPONSIBLE AUTHORITY” IN RELATION TO THE DEVICE.....	7
1.3	WASHING TROLLEYS AND ACCESSORIES.....	8
1.3.1	Smeg WD-PRINTE printer – GW2145, GW1160, GW4060 series	8
1.3.2	Smeg WD-PRINT9 Printer – GW4190 Series	8
1.4	PERISTALTIC PUMP P3 (NOT INSTALLABLE ON GW2145).....	9
1.5	LIST OF MODELS.....	10
2	KEY TO THE SYMBOLS USED IN THE MANUAL AND ON THE APPLIANCE	11
3	GENERAL WARNINGS	12
3.1	ACCESSING AND REUSING THE DEVICE AFTER AN INCOMPLETE CYCLE	13
3.2	DOOR OPENING	14
3.3	MANUAL DOOR RELEASE PROCEDURE	15
3.4	GW4190 OPENING THE DOOR OF THE SIDE CABINET - SERIES GW4190 ONLY	15
4	DESCRIPTION OF CONTROLS AND PROGRAMS	17
4.1	CONTROLS.....	17
4.2	CONTROLS - DETAILED DESCRIPTION	18
4.2.1	A0 PARAMETER DISPLAY	19
4.3	INDICATOR LEDS - INTRODUCTION	20
4.4	INDICATOR LEDS - DETAILED EXPLANATION	21
4.5	SELECTING THE WASHING PROGRAMME	23
4.6	BRIEF TABLE OF THE PROGRAMS INSTALLED.....	25
4.6.1	RESINS - SOFTENER, REGENERATION AND WASHING PHASE	26
4.6.2	DRYING PHASE – Only for GW4060, GW4190 series	26
4.6.3	CUSTOM PROGRAMS	28
4.7	A0 THERMAL DISINFECTION PARAMETER	29
4.8	PREPARING THE LOAD FOR WASHING AND DISINFECTION	30
4.9	END OF CYCLE	31
4.10	INTERRUPTING A PROGRAM IN EXECUTION	32
4.11	PROGRAM SPECIAL OPTIONS.....	33
4.11.1	DELAYED START - Setting the hours.....	33
4.11.2	DELAYED START - Starting the program and display.....	33
4.11.3	TEMPORARY DISABLING OF DEMI WATER.....	34
4.11.4	Deleting “Delayed start” - “No Demi” settings	35
4.12	RESET PROCEDURE.....	36
4.13	PRINTING THE CYCLE AND APPLIANCE PARAMETERS.....	37
5	SETUP MODE - “FUNCTIONS”	37
5.1	ACCESSING THE SETUP MODE / ENTERING THE PASSWORD	38
5.2	ACCESS AND PARAMETER MODIFICATION	39
5.3	“SET” PARAMETERS	41
5.3.1	APPLIANCE DATA STORAGE, OVERWRITING, “Set -b”	44
5.4	“CLOC” – DATE AND TIME PARAMETERS.....	46
5.5	SETTING THE “PrEn” “PROGRAM ENABLING” FUNCTION.....	46
6	INSTALLATION INSTRUCTIONS	49
6.1	USE OF THE WATER SOFTENER	49
6.2	USING DETERGENT AND NEUTRALISING AGENT	50
6.2.1	LIQUID DETERGENT INTAKE SYSTEM	50
6.2.2	DETERGENT LEVEL SENSORS - OPTIONAL	52

6.3	USING THE DETERGENT DISPENSER – GW SERIES ONLY	53
6.4	RECOMMENDED DETERGENTS	54
7	ALARMS AND WARNING MESSAGES	56
8	CLEANING AND MAINTENANCE	64
8.1	CLEANING THE DEVICE AND ITS PARTS.....	64
8.1.1	GW2145 - filter elements.....	66
8.1.2	GW1160 - GW4060 – GW4190 - filter elements.....	66
8.2	IF THE DEVICE IS TO BE OUT OF USE.....	67
8.3	REUSING THE DEVICE AFTER A PERIOD OUT OF USE	67
8.4	TROUBLESHOOTING.....	67
9	ROUTINE INSPECTIONS AND MAINTENANCE	69
9.1	DAILY	69
9.2	WEEKLY	69
9.3	EVERY SIX MONTHS	69
9.4	EVERY YEAR.....	69
10	TECHNICAL DATA	70
11	PRODUCT DIMENSIONS - MEASUREMENTS IN MM	72
11.1	GW2145	72
11.2	GW1160, GW4060 SERIES.....	72
11.3	GW4190 SERIES.....	73
12	POSITIONING THE APPLIANCE	73
13	ELECTRICAL SYSTEM REQUIREMENTS.....	77
13.1	GW2145 - POWER SUPPLY CABLE AND DISCONNECTION DEVICE	77
13.2	GW1160, GW4060, GW4190 - POWER SUPPLY CABLE AND DISCONNECTION DEVICE.....	77
14	WATER CONNECTION FITTINGS.....	79
14.1	WATER INTAKE.....	79
14.1.1	WATER INTAKE CONNECTION	79
14.1.2	WATER SUPPLY TAP POSITION.....	79
14.1.3	WATER SUPPLY SYSTEM REQUIREMENTS.....	80
14.1.4	PAD ACCESSORY FOR UNPRESSURISED DEMI WATER	81
14.2	WATER DRAIN	82
15	DATA CONNECTION.....	84

1 INTRODUCTION

GW2145, GW1160, GW4060, GW4190 – professional glassware washing machines, laboratory glassware washers (GW – glassware washer).

Smeg GW1160, GW4060, GW4190 family products comply with all the prerequisites of the relevant safety regulations and the standards in force for electrical equipment.



This manual describes products series:

- **GW2145 Series**, laboratory glassware washers, 45cm, without drying system.
- **GW1160 Series**, laboratory glassware washers, 60cm, without drying system.
- **GW4060 Series**: laboratory glassware washers, 60cm, equipped with forced air drying system.
- **GW4190 series**: laboratory glassware washers, 90cm, with forced-air drying system and side cabinet to take detergent cans.

The diagrams and images contained in this document refer to **GW4060** series. For the **GW2145, GW1160, GW4190** series the same considerations and procedures apply, unless otherwise specified.

The table below summarises the product's main characteristics. (For a full list of models and their configurations refer to the "LIST OF MODELS" point below).

	GW – basic model
CONTROL	ELECTRONIC with MICROPROCESSOR (with Time and Date functions)
DOOR OPENING	Automatic
DETERGENT DISPENSING	Powder detergent dispenser DD (integral in inner door) and peristaltic pump P2
OPTIONAL DETERGENT DISPENSING	peristaltic pump P1 peristaltic pump P3 (not installable on GW2145)
CONTROL OF DETERGENT DISPENSING	Optional
DETERGENT CAN LEVEL CONTROL	Optional
TEMPERATURE CONTROL	Single Platinum PT probe (TL1)
STEAM CONDENSER	On special models with suffix "C" (e.g. GW1160 Cxx, GW4060 Cxx, GW4190 Cxx)
DRYING SYSTEM	Yes, forced air (only on GW4060, GW4190 series)
ABSOLUTE DRYING FILTER	Optional (and only on GW4060, GW4190 series)
RAILS IN CHAMBER FOR TOP BASKET	Yes
WATER INTAKE CONNECTIONS	Mains and Demineralised
PRESET PROGRAMS	Specific LABORATORY programs, TOTAL OF 16 PROGRAMS of which 6 are CUSTOM (customisable using SMEG TRACELOG software ¹ . The software is not included in the product's basic outfit)
COMMUNICATION	RS232 serial (For PC or Printer)

¹The TRACELOG software enables the user to communicate with SMEG "WD" and "GW" instrument and glassware washers.

Once the connection has been made, the user can:

- **Keep a record of the disinfection cycles performed.**
- Monitor the appliance's status.
- Download new versions of the appliance's software.
- Modify programs.
- Save the appliance parameter setup.

1.1 GW FAMILY - INTENDED USE

Laboratory glassware washer.

The appliance is built to provide the following function:

- **Washing of laboratory glassware with chemical or thermal disinfection;**
- **The appliance cannot be used to sterilize instruments or any other device.**



Any use other than that described in this manual constitutes misuse.

SMEG declines all liability for uses other than those stated above.

SMEG declines all liability for any damage caused by the use of the appliance to wash glassware or instruments not approved by their manufacturers for automatic decontamination by means of thermal disinfection.

The appliance may only be used by specifically trained staff. Smeg will train the staff assigned to use the appliance at the time of installation.

Smeg declines all liability in the event of malfunctions or accidents caused by use of the appliance by untrained staff.

The training of the staff responsible includes specific information on the possible risks involved in the use of the appliance, and training in the safest possible way of carrying out the operating procedures.

The Smeg installation engineer is also responsible for notifying the responsible authority of the USER and SUPERUSER passwords to access the setup parameters. The responsible authority must keep these passwords in a safe place.

*It is the duty of the appliance's **RESPONSIBLE AUTHORITY** to ensure that those using the equipment have been suitably trained in its operation, its safe use and routine checks, and that this training is suitably maintained.*

Staff training should be checked regularly.

The installation engineer is responsible for ensuring that it operates correctly after commissioning.
Safety information supplied in compliance with 5.4.101.1 IEC61010-2-040:2005



Instrument manufacturers' instructions should always be followed when choosing the most appropriate disinfection treatment.

In particular, it is important to check the compatibility of the load for treatment with the specific washing cycle chosen, in terms of the maximum temperatures reached and the chemicals used. Information supplied in compliance with point 5.4.4.r IEC61010-2-040:2005.

1.2 DEFINITION: "RESPONSIBLE AUTHORITY" IN RELATION TO THE DEVICE

Responsible authority: "Person or group responsible for the use and ordinary maintenance of the unit and for operator training.

Definition taken from point 3.5.13 of the IEC61010-1:2001 standard.

The responsible authority should be clearly identified within the facility where the appliance is used, (for example by recording the relative names and responsibilities on corporate forms).

1.3 WASHING TROLLEYS AND ACCESSORIES



N.B: *the device is supplied without the washing trolleys and instrument holder baskets. Contact Smeg for details of the range of racks and baskets best suited to your washing requirements.*

Visit www.smeg-instruments.com for information about the standard range of products available.



*Any accessories may not be installed on the appliance by the user; this must be done by **Smeg authorised technical staff**.*

Any components must always be replaced with genuine Smeg spare parts.

1.3.1 Smeg WD-PRINTE printer – GW2145, GW1160, GW4060 series

The appliance's RS232 port can be used to connect the external Smeg WD-PRINTE printer, supplied with its own power transformer.

This accessory should be used **only with thermal paper, 57.5mm ± 1mm**, recommended weight 55 g/m².

To operate this accessory, simply connect the transformer to the power socket and connect the RS232 port to the RS232 socket on the rear of the glassware washer.

Refer to the accessory's manual for detailed instructions.



fig. 1 – Smeg desktop printer.




fig. 2 – RS232 cable supplied with the printer and position of the RS232 socket on the device.

1.3.2 Smeg WD-PRINT9 Printer – GW4190 Series

The Smeg "WD-PRINT9" is a panel printer that can be installed inside the side cabinet on GW4190 series devices.



fig. 3 – WD-PRINT9 accessory.

1.4 PERISTALTIC PUMP P3 (NOT INSTALLABLE ON GW2145)

Peristaltic pump P3 is always optional: contact Smeg for advice about the model best suited to your needs.



Note for GW4190 series: the device's cabinet is designed to take two 5 litre cans of Smeg detergent. To take a third detergent can, for use in combination with P3, the cabinet must be fitted with the Smeg "TANKBSK" (Tank Basket) accessory.

1.5 LIST OF MODELS

GW1160, GW4060 model configuration chart.

Symbol used:

SYMBOL	MEANING
●	feature present
○	optional accessory, installable on the model.
-	feature not present and not installable on the model.

GW MODELS	GW2145	GW1160	GW4060	GW4190	GW1160S	GW4060S	GW4190S	GW1160SC	GW4060SC	GW4190SC	GW1160C	GW4060C	GW4190C
Features													
DD - Powder detergent dispenser in inner door	●	●	●	●	-	-	-	-	-	-	●	●	●
P1 - Peristaltic detergent pump	○	○	○	○	●	●	●	●	●	●	○	○	○
P2 - Neutralising agent peristaltic pump	●	●	●	●	●	●	●	●	●	●	●	●	●
P3 - Optional peristaltic pump	-	○	○	○	○	○	○	○	○	○	○	○	○
Drying System	-	-	●	●	-	●	●	-	●	●	-	●	●
Absolute filter (dryer air)	-	-	○	○	-	○	○	-	○	○	-	○	○
Steam condenser	-	-	-	-	-	-	-	●	●	●	●	●	●
Control of detergent dispensing	○	○	○	○	○	○	○	○	○	○	○	○	○
Detergent can level control	○	○	○	○	○	○	○	○	○	○	○	○	○
Acquastop	-	○	○	○	○	○	○	○	○	○	○	○	○
Side cabinet to accommodate 3 detergent cans	-	-	-	●	-	-	●	-	-	●	-	-	●
Electrical main switch	-	-	-	●	-	-	●	-	-	●	-	-	●

ELECTRICAL CONFIGURATION

A number of different electrical connection options are available for each 60cm and 90cm model referred to in the table.

The names of the models have different suffixes to identify the power supply rating. Example:

- "GW4060-1", "GW1160-1" and "GW4190-1" refer to 230V single-phase models: "1N~"

- "GW4060-3", "GW1160-3" and "GW4190-3" refer to three-phase models, 230V between phases, without neutral: "3~"

- "GW4060", "GW1160" and "GW4190". If the model does not have a suffix it will be three-phase, with neutral, 400V between phases: "3N~".

- "GW4060U", "GW1160U" and "GW4190U". The "U" suffix refers to "North American models", convertible 208V models: 208V 3~ 60Hz * / 208 V ~ 60Hz (* default connection).



For the 45 cm model: GW2145, only the single-phase model is available for 230V ~ 50Hz electrical connection with Schuko plug

Contact Smeg for advice about the model best suited to your needs instruments@smeg.it

2 KEY TO THE SYMBOLS USED IN THE MANUAL AND ON THE APPLIANCE

The following is the key to the symbols used on the appliance and in this manual, as required by point 5.4.4.e of IEC61010-2-040:2005.



Read with particular attention.

(Symbol which appears beside particularly important instructions or warnings).



Warning, danger: refer to the manual.

(The symbol appears on the appliance's technical dataplate to emphasise that staff must read the manual before using the device. The symbol appears in the manual next to safety instructions)




Warning, hot surface.

(the symbol is affixed to the appliance close to parts which may heat to high temperatures and constitute a burn hazard - do not touch parts near this symbol. The symbol appears in the manual to emphasise safety instructions relating to the burn hazard).



Warning, electric shock hazard.

(The symbol is affixed close to live parts - the device's power supply must be disconnected before doing any work on these parts. Never touch live parts unless the electricity supply is disconnected).



Transport and unpacking precautions.



Flammability hazard.



Biohazard.



At the end of its working life, the product must be consigned to a disposal plant for recovery and recycling in accordance with the relevant legislation in the country of installation. Contact the specialist disposal consortia. This appliance is marked according to the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.



CE marking; the manufacturer guarantees compliance with the EU directives applicable to the product. This symbol appears on the machine's technical data label and in this manual.

3 GENERAL WARNINGS

	<p>Leaning or sitting on the open door of the device might cause it to tip over, putting people in danger. The door is not designed primarily to support loads.</p> <p><i>The maximum weight which can be loaded on the door, including the weight of the instrument trolley, must never exceed</i> <i>GW1160, GW4060, GW4190: 37kg</i> <i>GW2145: 18kg</i></p> <p><i>- if the weight of the instrument trolley and baskets is excluded, the maximum load must never exceed:</i> <i>GW1160, GW4060, GW4190: 23kg.</i> <i>GW2145: 15Kg</i></p> <p>For optimal DRYING, the load for processing must never exceed 15 kg. (only for GW4060, GW4190 Series).</p> <p><i>With use of the glassware washer, localised or general discolouring of the heating element may occur. This is normal since it is due to the operating mode and does not reduce the appliance's effectiveness.</i></p> <p><i>In the event of malfunction, disconnect the appliance from the electricity supply and turn off the water tap. Then contact your nearest authorised Service Centre.</i></p> <p><i>Open the door carefully, first waiting for the wash cycle to end.</i></p> <p><i>The appliance has an automatic door opening system; do not force the door open while a program is in progress.</i></p> <p><i>The appliance must only be used by staff suitably trained in its operation.</i></p> <p><i>The chamber of the appliance is not designed for users to climb into it. The user must never climb into the chamber; this might put his safety at risk (ref. 7.102 IEC61010-2-040:2005).</i></p>
	<p>Do not place flammable substances inside the device. Do not use flammable detergents.</p> <p><i>Never place alcohol or solvents such as turpentine, which might cause an explosion, inside the appliance. Do not place materials dirty with ash, wax or paint inside the appliance.</i></p>
	<p>DEMI WATER ABSENCE</p> <p><i>If demineralised water is not available, the user is responsible for ensuring that the quality of the water supplied to the medical device does not cause mineral salts or other substances to be deposited on the treated instruments, rendering their subsequent use unsafe.</i></p>
	<p>Never touch the heating element immediately after the end of a wash programme.</p> <p><i>Do not touch any residual liquids left inside the washing chamber; scalding hazard (ref. 7.102.c IEC61010-2-040:2005).</i></p>
	<p>When moving the appliance around, a forklift truck or pallet truck must be used.</p> <p><i>Before leaving the factory, the base of the appliance is secured to a pallet, which is used for lifting and transporting it.</i></p> <p>Do not use appliances which have been damaged in transit!</p> <p><i>If in doubt, contact your dealer.</i></p>
	<p>Once decommissioned, the appliance must be rendered unusable. Cut the power supply cable after removing the plug / disconnecting the cable from the power socket.</p>

3.1 ACCESSING AND REUSING THE DEVICE AFTER AN INCOMPLETE CYCLE

The instructions relating to the device's safety in the event of an incomplete operating cycle are provided in compliance with points 5.4.4.g and 13.1.102 of IEC61010-2-040:2005.

WARNING

Warning applicable if the device is used for processing of contaminated materials.

If a disinfection cycle is interrupted (by the user or due to an alarm generated by the appliance itself): take care when handling the elements and the load in general in the washing chamber.

The load and the internal parts of the appliance might be biologically contaminated/infected.

Before handling tools or performing any maintenance: Perform a complete disinfection cycle or, if you cannot perform a complete Thermo-disinfection cycle, handle tools with caution (using protection devices required for the handling of infected instruments, eg. gloves, lab coat).




DANGER, HOT SURFACES

The appliance performs a thermal disinfection cycle using water at high temperatures (up to 93°C) and detergents. If, in the event of a failure, there is water in the chamber when the door is open: avoid contact with the skin, burning hazard and risk of irritation due to the toxicity of the chemicals used.

Never touch the heating elements inside the chamber.

Contact Smeg authorised technical staff in the event of a failure.



Do not touch
the heating
elements

fig. 4 – View of inside of chamber, detail of heating elements.



3.2 DOOR OPENING

The following is a brief preview of the procedure for opening the door of the device for easier access to the chamber. The points which follow provide a complete description of all the controls.

The appliance is equipped with an **automatic door locking system**.

Do not force it open; proceed as follows:

- Connect the appliance to the electricity supply.
- For GW4190 series only: open the hatch of the side cabinet and turn the main switch on.
- Hold down the button with the **On/Off** symbol for 2 seconds to switch the appliance on.
- Press the **Door opening** button: the appliance gives a **beep** when it receives the signal. Release the button and the door will open automatically after a few moments.

BUTTONS TO BE PRESSED TO OPEN THE DOOR			
	On/Off		Door opening




fig. 5 – Selecting automatic door opening

To open the door without connecting the appliance to the electricity supply (useful for emergencies or in case of power blackouts), the manual release system has to be used (see next point).

3.3 MANUAL DOOR RELEASE PROCEDURE

If necessary, due to a malfunction or power blackout, the appliance can be opened manually by releasing the lock using a screwdriver with stem $\varnothing 4\text{mm}$; be gentle, taking care not to damage the device.

1. This procedure is only possible with the appliance in standby status, with no cycle running.
2. Warning: before opening the appliance manually, **disconnect it from the electricity supply.**
3. Insert a screwdriver or a rod 4 mm in diameter into the hole in the bottom of the centre of the front panel.
4. **Push the pin upward until the lock is released.**




fig. 6 – Diagram showing how to release the lock by hand in any emergency.

3.4 GW4190 OPENING THE DOOR OF THE SIDE CABINET - SERIES GW4190 ONLY

GW4190 series models have a side cabinet with a hatch which opens like a drawer. Once the hatch is open, users have full access to the components inside:

- **Main switch** (circuit-breaker),
- **Panel printer** (if installed),
- **Detergent cans.**

To open the cabinet hatch, simply pull it towards you using the stainless steel handle provided.

The following are located above the panel that carries the master switch and printer:

- **RS232 port:** for data transmission – connection by means of Smeg TRACELOG software.
- **“Data output” switch:** used to activate the panel printer or data transmission on the RS232 port. These two devices can never operate simultaneously. The data output switch activates the RS-232 port when turned to “I” or the printer when turned to “II”.




fig. 6 – Opening the GW4190 side cabinet hatch.



fig. 6 – GW4190: Once the side cabinet hatch is open, the user has access to the device's master switch (and the printer, if installed).




fig. 6 – GW4190: the RS232 data port and the data output switch (which directs data either to the transmission port or to the printer) are located above the master switch panel.

4 DESCRIPTION OF CONTROLS AND PROGRAMS

4.1 CONTROLS

All the device’s control and monitoring devices are placed together on the front panel.










fig. 7 – Control panel. The display and indicator LEDs are in the middle and the control buttons at the sides and bottom.



CONTROL PANEL SYMBOLS			
SYMBOL	MEANING	SYMBOL	MEANING
	On/Off		Functions
	Drying Decrease (only on GW4060, GW4190 series)		Drying Increase (only on GW4060, GW4190 series)
	Start / Stop		Door opening
	Selection		Selection



In the figures mentioned in this document, we use only the GW4060 front panel layout. The same considerations and procedures apply to the GW1160, GW4190 models, unless otherwise specified.

4.2 CONTROLS - DETAILED DESCRIPTION

BUTTONS	DESCRIPTION		
	<p>On/Off The button is operational with the door opened or closed. The button turns the device’s interface on and off, allowing its use. Hold it down for 2 seconds for both switch-on and switch-off. Once the appliance has been switched on, the code of the selected program appears on the display (e.g. “Pr 03”). When the interface is switched off, the display shows a flashing “OFF” message. When the door is open, the display shows 4 horizontal dashes and none of the other interface buttons are active. The device is still powered up even when the interface is off.</p>		
	<p>Drying +/- (only on GW4060, GW4190 series) The buttons increase/decrease the intensity (and duration) of the drying phase. Holding down the “Drying–” button selects a separate drying program; the display shows the message “Pdry”. If the cycle set has a drying phase, these buttons are enabled and the drying LED illuminates. If the cycle set does not include a drying phase, the buttons are not enabled and the drying LED remains off. These are “first touch” buttons; there is no need to hold them down. The drying duration set is displayed as follows: - 3 horizontal bars on the display; maximum time setting, default time - 2 horizontal bars, medium duration - 1 horizontal bar, short duration</p>		
	<p>Start/Stop This button starts the selected cycle or interrupts a cycle in progress. It must be held down for 2 seconds. Two beeps sound when a cycle is started or stopped. The button is also used (again holding it down) to confirm the value of an appliance parameter during modification.</p>		
	<p>Functions Hold down: access to the appliance setup parameters, password required. Refer to the SETUP MODE - “FUNCTIONS” section for a detailed explanation. Short press: whenever the Functions button is pressed, the display screen changes; the various screens show the parameters of the cycle selected in combination with the LED which explains their meaning. For a full understanding of the values shown on the display, note that they depend on the appliance’s status: 1. Appliance in standby 2. Washing cycle running</p>		
	LED	ON DISPLAY	Meaning of LED on
		1 : 16	<p>With appliance in standby The value is the nominal duration of the cycle. The time LED also illuminates. The duration is guideline; it applies in test conditions and with correct electrical connection. With cycle running Time left to end of current cycle</p>
		-- : 90	<p>With appliance in standby Maximum temperature reached in the cycle set With cycle running - During washing: Tank temperature - During drying: Inlet air temperature (function introduced from firmware version 4.6.xx.yy).</p>

BUTTONS	DESCRIPTION		
	A ₀	600	<p>With appliance in standby A0 value associated to the selected cycle, if the cycle includes a thermal disinfection phase.</p> <p>With cycle running Display function which starts from the thermal disinfection phase: displays the A0 value reached. The value reached during the thermal disinfection phase is displayed during the phases which follow. [The parameter is displayed on two consecutive display screens from firmware version 1.11.0.33. The value in thousands is indicated by the suffix "t", in units with the suffix "u"]</p>
	-	F5	<p>With appliance in standby Number of phases in the program</p> <p>With cycle running Current phase N.B.: Washing cycles always start with a short draining phase to bring the device to the standard starting condition: this phase is called "phase 0".</p>
	-	PR 03	<p>If the user does not give any commands for 5 seconds, the display shows the program selected.</p>
	<p>Door opening The door is opened automatically by means of an electric lock; opening is only permitted in conditions of safety for the user. The user simply presses the button; there is no need to hold it down. The button also provides an "Escape" function; it is pressed to quit the procedure for displaying or modifying an appliance parameter. In this case, it has to be held down (2 seconds).</p>		
	<p>Selection Program selection buttons. The letter "P" and the name of the selected program appear on the display. These buttons can be used to increase or decrease the current value during modification of an appliance parameter.</p>		

4.2.1 A0 PARAMETER DISPLAY

The parameter is displayed on two consecutive display screens from firmware version 1.11.0.33. The A0 parameter can reach high values (eg typical A0 = 12000) and will appear on two consecutive screens; thousands display screen (t-Thousands) and units screen (u-units), the value in thousands is indicated by the suffix "t" units with the suffix " u ". The transition from one screen to the next occurs automatically after 3 seconds or by short press of the "Functions" key.




fig. 8 – The A0 value is displayed split into thousands (t) and units (u)

4.3 INDICATOR LEDS - INTRODUCTION









fig. 9 – In the front panel display area: digit area in the centre with indicator LEDs around the perimeter

LED SYMBOLS			
SYMBOL	MEANING	SYMBOL	MEANING
	No detergent	A ₀	A0, thermal disinfection
	Maintenance: filter replacement (Only for GW4060, GW4190)		Temperature
	Appliance maintenance		Time: time left on display
	Lack of salt	T _d	Thermal disinfection, identifies cycles with thermal disinfection
	Alarm		Drying (Only for GW4060, GW4190)

4.4 INDICATOR LEDS - DETAILED EXPLANATION

1 - Symbols associated to a warning	
SYMBOL	MEANING
	<p>No detergent</p> <p>Lights up if the detergent level sensors are installed and one of them is at "minimum" level. The warning is displayed at the end of the cycle and, in combination with the LED, the segment display shows the code of the detergent to which the warning refers:</p> <p>P1: code "A-:68" P2: code "A-:69" P3: code "A-:70"</p> <p>Activated at the end of the cycle and when the user attempts to start a new program. The washing cycle can be started with this light on by pressing the Start/Stop button again. If both alarms are present (cans P1 and P2 both empty), the Start/Stop button has to be pressed twice to start the cycle. (The device's internal memory records the event).</p> <p>Replace the can of detergent which is running out with a new can as soon as possible.</p>
	<p>Maintenance: filter replacement</p> <p>Only for GW4060, GW4190 series. Active if the dryer is present and the absolute filter is installed. (The absolute filter is optional on the GW series).</p> <p>The LED indicates that the absolute filter needs to be replaced.</p> <p>Operation of the light is triggered by the number of drying operating hours, set at 500h when the filter is new. The filter must be replaced by Smeg authorised staff.</p>
	<p>Appliance maintenance</p> <p>The appliance counts the number of cycles performed and the LED comes on to alert the user that maintenance is required.</p> <p>These are inspection and maintenance operations, scheduled every 1000 cycles, which must be performed by Smeg authorised staff to keep the appliance safe and in good working order.</p>
	<p>Lack of Salt</p> <p>Salt must be added to the water softener reservoir in the chamber.</p> <p>Salt is necessary to enable the softener to reduce the hardness of the intake water.</p>
	<p>Alarm</p> <p>The machine displays an anomaly, which can be generated by a fault or by the detection of abnormal conditions. A numerical code on the segment display indicates the alarm which has been triggered. The alarms table in this manual details the components which may have caused the event for each alarm code, to simplify diagnostics and troubleshooting.</p> <p>A record should be kept of any alarms to allow the After-Sales Service to be provided with a detailed description of the malfunction.</p>

2 - Symbols associated to a washing cycle parameter

When the LEDs listed below light up, the value shown on the display is a parameter of the washing cycle in progress or being selected.

When the cycle has not yet started: pressing the **functions** button displays the cycle parameters.

When the operating cycle is in progress: pressing the **functions** button displays the value reached by the parameter associated to the LED.

SYMBOL	MEANING
	<p>A0, thermal disinfection</p> <p>With cycle running: The LED is on when the digits in the middle of the screen show the A0 value reached.</p> <p>With appliance in standby: when the LED is on the digits in the middle of the screen show the A0 value associated to the cycle being selected.</p> <p>Parameter calculation and display are only enabled if the cycle has a thermal disinfection phase with temperature of at least 80°C.</p>
	<p>Temperature</p> <p>With cycle running: The LED lights up when the display shows the temperature in the chamber.</p> <p>With drying in progress: The LED lights up when the display shows the temperature in the dryer pipe (function introduced from firmware version 4.6.xx.yy).</p> <p>With appliance in standby: the value on the display shows the max. temperature associated to the cycle being selected.</p>
	<p>Thermal disinfection</p> <p>Identifies a cycle with thermal disinfection; the LED flashes when a thermal disinfection phase is in progress.</p>
	<p>Drying</p> <p>Only for GW4060, GW4190 series. The LED comes on to identify a cycle with drying phase, even during selection. The LED flashes when the drying phase is in progress.</p>
	<p>Remaining time</p> <p>With cycle running: The LED lights up when the display shows the time remaining.</p> <p>With appliance in standby: the value on the display shows the estimated time for performance of the program.</p>




fig. 10 – **Functions** button. With the cycle in progress or cycle selected, pressing the **Functions** button the screen displays the parameters of the washing cycle: A0 value, Temperature and Time Remaining.

4.5 SELECTING THE WASHING PROGRAMME

The device has an interface, with segment display and 10 LED-lit symbols, for user-device communications.

To run a program:

1. Close the door by shutting it against the appliance until the lock clicks shut. **The buttons are only enabled with the door closed.**
2. If the appliance is off, switch it on by pressing the **On/Off** key (hold down for 2 seconds)
3. Use the **Selection** buttons to select the program.
4. Once the selection has been made, the **Functions** button can be used to scroll through the program parameters.
5. Only for GW4060, GW4190 series: If the program includes a drying phase, its duration can be adjusted using the **Drying Decrease** and **Drying Increase** buttons.
6. Press the **Start/Stop** button to start the cycle (hold it down for 2 seconds).
7. Two beeps sound to indicate the start of the cycle.
8. Once the selected program has been successfully completed, the display shows **End**".
9. **Interrupting the cycle:** A cycle in progress can be interrupted by pressing the **Start/Stop** button. **A** After the cycle has been interrupted, the options available vary depending on conditions in the chamber:
 - i. If the temperature in the chamber is below **40°C**: the door can be opened. If the door is closed again within 1 minute, the cycle can be restarted from the phase where it was interrupted by holding down the **Start/Stop** button.
 - ii. If the temperature in the chamber is above 40°C: a reset cycle is required since the door cannot be opened.




fig. 11 – Standard procedure for starting a cycle. 1 Select the program using the Selection buttons. 2 Check the program parameters: Functions button. 3 Start the program: hold down the Start/Stop button.

CONTROL PANEL SYMBOLS			
SYMBOL	MEANING	SYMBOL	MEANING
	On/Off		Functions
	Drying Decrease (only for GW4060, GW4190)		Drying Increase (only for GW4060, GW4190)
	Start / Stop		Door opening
	Selection		Selection



fig. 12 – Example of a **Display screen**. Program 5 selected
The program includes a **drying phase** – **Drying LED on**.
The program includes a **thermal disinfection phase**– **"Td" LED on**.

NOTE - CYCLE COMPLETION TIMES

The completion times shown on the display are guideline: they may vary due to factors such as the intake water temperature and pressure.

At the end of each program, a **resin regeneration** phase may be automatically tripped, followed by **resin washing**, to ensure **that the integral softener operates correctly**. These procedures do not take place at every cycle but are triggered as needed, depending on the intake water hardness and the amount of water used.

The procedures are not included in the theoretical cycle time shown initially on the display.

In order for the real times to be close to the theoretical stated times, the appliance parameters must be set appropriately, so that the time passed is calculated on the basis of the electrical connection as made.



4.6 BRIEF TABLE OF THE PROGRAMS INSTALLED

For details of the various washing program phases, always refer to the enclosed “[GW2145-GW1160-GW4060-GW4190 Program Table](#)”; the brief table below is purely guideline. The programs listed below may be subject to modification for the purposes of improvement.



PROGRAMS with THERMAL DISINFECTION

In programs normally used in **GW** Glassware Washers, **thermal disinfection** takes place at the **beginning of the washing cycle**, for the most effective removal of contamination and grease from the material for processing.

This process is different from the one used and recommended in Hospital Instrument Washers, where this phase concludes the cycle so that the instruments are maintained in a disinfected state.



CUSTOM PROGRAMS

The **Custom programs** are programs which can be modified to meet specific requirements.

Programs can only be modified using the Smeg TRACELOG software; contact Smeg for information concerning the use and installation of the program.

ID	GW - GLASSWARE WASHER
PR 1	Prewash
PR 2	Plastics 75°C [+ Drying]
PR 3	Standard glassware 80°C [+ Drying]
PR 4	Intensive glassware 90°C [+ Drying]
PR 5	Agar Intensive glassware 93°C [+ Drying]
PR 6	Thermal disinfection 93°C – 3 minutes [+ Drying]
PR 7	Thermal disinfection 93°C – 10 minutes [+ Drying]
PR 8	Thermal disinfection + intensive wash 93°C – 5 minutes [+ Drying]
PR 9	Thermal disinfection + intensive wash 93°C – 10 minutes [+ Drying]
PR 10	Custom 1 (default: Petroleum Grease, 93°C – 1 minute) [+ Drying]
PR 11	Custom 2 (default: Mineral Oil, 93°C – 10 minutes [+ Drying]
PR 12	Custom 3 (default: Petrol, 93°C – 5 minute) [+ Drying]
PR 13	Custom 4 (default: Diesel, 93°C – 1 minute) [+ Drying]
PR 14	Custom 5 (default: Intensive Petrol, 93°C – 10 minutes [+ Drying]
PR 15	Custom 6 (default: Appliance auto-disinfection 93°C – 1 minute) [+ Drying]
PR 16	Service - operation of peristaltic pumps, intake of detergent + rinsing

Tab. 1 – Brief table of the programs installed

DEFAULT DETERGENT DOSES AND WATER INTAKES	FOR SERVICE PROGRAM ONLY, AD HOC DETERGENT DOSES AND WATER INTAKES
P1 = 5ml/litre – alkaline detergent	P1 = 12ml/litre
P2 = 3ml/litre - acid	P2 = 12ml/litre
P3 = 10ml/litre - degreasing	P3 = 12ml/litre
P2 ^a = 4ml/litre (non-standard P2 doses used in special cases).	P2 ^a - Parameter not used in the Service program.
Water intake per phase, default = 9 litres	Water intake, phase 1 and phase 2 = 6 litres Water intake, phase 3 = 9 litres

Tab. 2 – Brief table of doses and water intakes per phase.

MAINS WATER - DEMINERALISED WATER in PROGRAMS

The enclosed table indicates when a specific phase takes in cold water and when demineralised water.



If demineralised water is deselected by means of the appliance setup parameters, the appliance will always take in cold water.

Whenever possible, connection to demineralised water is recommended, for better final rinsing results and to avoid the formation of limescale.

4.6.1 RESINS - SOFTENER, REGENERATION AND WASHING PHASE

The softener incorporated in the device removes calcium compounds from the intake water to reduce its hardness. The softener comprises a container of ion exchanger resins, which must be regenerated regularly. The appliance decides when a resin regeneration and washing cycle is needed depending on the hardness set and the quantity of water treated. Resin regeneration and washing cycles are performed at the start of a washing cycle, to avoid the risk that standing brine (salt and water) might cause corrosion.

4.6.2 DRYING PHASE – Only for GW4060, GW4190 series

The **drying** phase is included in the washing and thermal disinfection cycle, except in programs “1 Prewash”, “Auto-Disinfection” and “16 Service”. When a program which includes the drying phase is selected, the LED of the function illuminates.

ADJUSTING THE DRYING SETTINGS

The specific buttons for this function allow the effectiveness of the drying process to be regulated by increasing or decreasing the total time of this phase. 3 dashes light up on the display if the default drying settings are active, while 2 dashes or 1 dash indicate gradually shorter times.

STANDARD STRUCTURE OF THE DRYING PHASE, 3 STEPS

Drying total time: $t_d = t_1 + t_2 + t_3$	Step 1 [t ₁]	Step 2 [t ₂]	Step 3 [t ₃]
Characteristics:	The fan speed is reduced (1/4 - 1/2 of full speed) to aid the controlled venting of steam from the well.	Time and temperature equal to target values set by the program. Fan at full speed.	Cooling phase for load and heating elements.
Duration, with 3 dashes	$t_1 = 10 \text{ minutes}$	$t_2 = \text{target_time}$	$t_3 = 4 \text{ minutes}$
Duration, with 2 dashes	$t_1 = 10 \text{ minutes}$	$t_2 = 0.5 \times \text{target_time}$	$t_3 = 4 \text{ minutes}$
Duration, with 1 dash	$t_1 = 10 \text{ minutes}$	$t_2 = 0$	$t_3 = 4 \text{ minutes}$
Duration, with 0 dashes	$t_1 = 4 \text{ minutes}$	$t_2 = 0$	$t_3 = 4 \text{ minutes}$

NOTES

The “**target time**” stated in the table refers to the time specified in the annexed “PROGRAM TABLE” document. In a standard program, e.g. with drying having $t_{\text{target}}=25\text{min}$, the total time of the phase is calculated as follows: $t_d = t_1 + t_2 + t_3 = 10 + 25 + 4 = 39\text{minutes}$. If the final temperature is below 40°C during the last phase of the washing cycle, step 1 is completely omitted, since there is no need to expel steam: $t_1=0$.

Target temperature: the temperature value in the program table, which refers to the temperature reached by the monitoring probe TA in the drying duct. The temperature in the well during drying is about 30°C below the target value. For example, for a drying program with $T_{\text{target}}=110^\circ\text{C}$,

$$T_{\text{well}} \approx T_{\text{target}} - 30^\circ\text{C} = 110^\circ\text{C} - 30^\circ\text{C} = 80^\circ\text{C}.$$

SEPARATE DRYING PHASE

Users can also select a **SEPARATE DRYING CYCLE** by pressing the “**Drying Decrease**” button for 2 seconds: the message “Pdry” appears on the display. The drying time can still be adjusted as described above.




fig. 13 – Display screen: *separate drying* program.



fig. 14 – Buttons used to *Decrease* and *Increase* the *Drying Time*.




fig. 15 – Display during *Drying* setting adjustment.



N.B. - important for selecting the cycle and max. cycle temperature

The drying temperature stated in the programs table (e.g. 110°C) refers to the value measured at a specific point of the Dryer duct and does not correspond to the temperature reached by the load.

Temperatures in the chamber during drying never exceed 90°C.

4.7 A₀ THERMAL DISINFECTION PARAMETER

The A₀ parameter (introduced by the EN 15883 standard²) allows a numerical value to be assigned to the thermal disinfection carried out - shown by **Td** in the interface LEDs.

When calculating the parameter, only the time intervals during which the temperature is above 65°C are considered. For SMEG thermal disinfection programs, the calculation is simplified by only including the “extension” phase, when the temperature is kept constant at close to the target value set.

Programs that include thermal disinfection have therefore been designed to offer the following A₀ values:

Temperature [°C] and time [min]	A ₀
90°C - 1'	600
90°C - 5'	3000
93°C - 5'	6000
93°C - 10'	12000

The formula for calculating A₀ is given below.

$$A_0 = \tau \cdot 10^{\left(\frac{T-80}{10}\right)}$$

τ = Time in seconds for which the disinfection temperature must be maintained.

T = Disinfection temperature in °C.

If the temperature is 80°C, A₀ is equal to the temperature maintenance time.

²EUROPEAN STANDARD EN 15883 “Washer-Disinfectors”, with particular reference to section 3 *Terms and Definitions* and annex B, *A₀ concept*. of part 1 of the standard, 15883-1.

4.8 PREPARING THE LOAD FOR WASHING AND DISINFECTION

Loading instructions are provided in compliance with 5.4.4-k of IEC61010-2-040:2005.

The load for processing must be arranged appropriately on the most suitable supports for machine-washing.

Effective washing starts with preparation of the elements / glassware for processing.



Therefore, before placing the glassware and the other items of the load in the specific baskets, it is necessary to remove any large residues deriving from previous use, by soaking, treatment or rinsing.

Stainless steel items cannot be immersed in physiological solutions of sodium chloride, as prolonged contact causes corrosion and damage to the surface due to stress corrosion cracking.

Do not overload the washing baskets to ensure effective mechanical action of the water across the entire surface of the glassware.

To allow effective cleaning and prevent damage to the objects for processing, glassware must be properly arranged on the most suitable supports for machine-washing.

Contact Smeg for advice about the washing trolleys best suited to your needs instruments@smeg.it

Items which can be dismantled must be prepared and stowed as instructed by the manufacturer.

For effective cleaning, items with hinged joints must be opened to minimise the overlapped surfaces. The instrument holder devices used and the fixing devices must be designed so that no shadow zones are created during subsequent cleaning and disinfection.

Make sure that there are no labels on the glassware for washing which may detach during the process.



GLASSWARE WASHING SUPPORTS - TROLLEYS

The appliance comes without washing trolleys.

Refer to the washing trolley manuals for guidance on their correct use.



Avoid direct, repeated contact with dirty material.

Always take the greatest care and use suitable personal protection equipment, both before and after treatment.

Before processing glassware or any other object in the glassware washer, check on the manufacturer's instructions that the items are approved for automatic treatment in the Glassware Washer and also the recommended maximum washing temperature.

4.9 END OF CYCLE

Once the selected program has been successfully completed, the display shows **End**.

“End” only appears if the washing phases have been performed and been completed correctly; it does not appear if the program has been interrupted or if an error which may affect the cycles’ effectiveness has occurred.

A “beep” sounds to indicate the end of the cycle.

AUTOMATIC OPENING AT END OF CYCLE

If automatic opening at the end of the cycle is enabled, the door of the Glassware Washer opens automatically and only the 4 dashes appear on the display.

Automatic opening is not enabled if alarms occur during washing.

Drying LED: When a program which includes the drying phase is selected, at the end of the cycle, the relative LED is:

1. on and steady if the process has been completed correctly,
2. flashing if drying has been interrupted (e.g. due to a power blackout or forced interruption of the program by the user).





fig. 16 – End of cycle, display screen. The word End indicates that the washing process has been completed correctly. If automatic opening is enabled: the door is released automatically and the 4 dashes appear on the display.

GW - AUTOMATIC OPENING AT END OF CYCLE

At the end of the cycle, the door can open automatically to allow the remaining steam to disperse and help the processed material to cool and dry more quickly.

If the temperature inside the chamber is above 80°C, a delay of 10 minutes is automatically calculated before the door opens.

During this delay, the lock can be released by pressing the **Door opening** button.

Automatic door opening at the end of the cycle is a parameter which can be enabled or disabled: refer to the specific point: SETUP MODE - “FUNCTIONS”.

Wait about 10 minutes with the door open before removing the processed load, to allow the material in the washing chamber to cool, and if necessary to dry.

4.10 INTERRUPTING A PROGRAM IN EXECUTION



INTERRUPTING A PROGRAM

A cycle in progress can be interrupted by pressing the **Start/Stop** button.

During the “**Suspension**” of the program, the letter “**S**” appears on the display, followed by the code of the current program.




fig. 17 – Interrupting a Program: hold down the **Start/Stop** button

For GW - Laboratory Glassware Washer series products, SUSPENSION phase:

1. If the temperature in the chamber is below **40°C**: the door can be opened. If the door is closed again within 1 minute, the cycle can be restarted from the phase where it was interrupted. The cycle is always restarted by holding down the **Start/Stop** button for 2 seconds.
2. If the temperature in the chamber is above **40°C**: a reset cycle is required since the door cannot be opened.
3. If the cycle is interrupted and the user does not press any buttons within 1 minute: the appliance automatically switches to **Reset** mode



fig. 18 – Display screen, Suspension of program 5: **S 05** “.

4.11 PROGRAM SPECIAL OPTIONS

1. Temporary disabling of demi water

2. Delayed start (function introduced from firmware version 4.6.xx.yy).

The “Special Options” menu is accessed by pressing the “Start/Stop” and “Door Opening” buttons simultaneously. To scroll through the special functions: press the “Functions” button and release at once.

The special options only apply for 1 single cycle, executed after selection of the option. After execution of the cycle, the default values are restored.

Special options must be confirmed by pressing the “Start/Stop” button; after confirmation, the program display returns to the screen.

If the option is not confirmed, after 10 seconds without any action the system quits the option and the last parameter set is maintained.

N.B.: the current date and time must be set to allow the “delayed start” setting to function correctly.

4.11.1 DELAYED START - Setting the hours

Delayed start: initially, the display shows “dh:00”.

The selection keys can be pressed to set up to a max of 12 hours (“dh:12”).

The “clock” LED remains on while the selection is made and the value set flashes until confirmed.




fig. 19 – Delayed start setting screen

4.11.2 DELAYED START - Starting the program and display

To start a program with delayed start set: Hold down the “Start/Stop” (as when starting the cycle normally).

The display shows, in alternation (switching every 2 seconds):

the time remaining before the start of the program (format “hh:mm”) / **the program selected** (“Pr:XX”).

The “clock” LED remains steady when the time remaining is displayed.

At the end of the delay time: the program starts and from this point onward the contents of the display are the same as for a normal cycle.




fig. 20 – Once the cycle start command is given, the display alternates the program selected with the time remaining before it actually starts.




fig. 21 – Example of use of the delayed start to set a delay of 1 h

4.11.3 TEMPORARY DISABLING OF DEMI WATER



Press the “**Start/Stop**” and “**Door opening**” buttons simultaneously for a moment to display the status of the “demineralised water disabling” option.
The **Selection** buttons can be used to deselect the use of demineralised water.

If the option is active, the message “**no dE**”, meaning “**No demi**”, **temporary disabling of demineralised water**, appears on the screen.

If the option is disabled by pressing the two keys mentioned above, the message “**dE**” is shown on the screen.

This option can only be used if the use of demineralised water has been activated.



The **disabling option** can be activated before starting the selected program, and remains **active for only one cycle**, at the end of which it is deactivated automatically.

When the disabling option is active, demi water is replaced with cold water in all the phases where it would normally be used.




fig. 22 – Demi water disabling buttons. To access the function: hold down the **Functions + Door opening** buttons together. To select demi water disabling: **Selection** buttons




fig. 23 – Display screen: **TEMPORARY disabling of demi water**; the option can be enabled and disabled using the selection buttons

4.11.4 Deleting “Delayed start” - “No Demi” settings

To delete a setting:

1. return the parameter to zero (e.g. “dh:00” for delayed start, “ dE” for demi water).
2. by switching the device off and back on (the default value is restored: e.g. “dh:00”).
3. by opening the door of the device (for “delayed start” only) once the cycle start command has been given

4.12 RESET PROCEDURE

The **RESET** procedure performs water pump-out and intake cycles to return the device to conditions of safety in the event of a malfunction.

A **RESET** cycle can and should be started if an alarm occurs.

The Reset cycle can also be started independently, regardless of whether or not an alarm is present.

N.B.: the **RESET** procedure can only be started with the door **CLOSED**.



When an error occurs, generally the device manages the anomaly itself: in this case **the alarm code flashes on the display** until the problem has been dealt with: **no reset commands are accepted during this time**.

Any RESET forced by the user will not be accepted until the problem has been dealt with and the alarm code on the display is "steady".

If alarm "**AF:84**" (overheating during drying) occurs, the device manages its status independently and does not allow a **RESET** to be performed until the automatic process for dealing with the alarm has been completed.

STARTING A RESET

After an alarm or with the appliance at a standstill, press the **Functions** and **Start/Stop** buttons and hold them down for about 2 seconds, until a **beep** is heard.

"**P-**" appears on the display and the **RESET** procedure starts.

At the end of the procedure, "**E-**" flashes on the display, alternating with the program which was running and has been interrupted, if any.

Depending on the appliance's status, the **RESET** command may not be accepted: in this case, open and close the door and then repeat the procedure.



ATTENTION

If it is not possible to perform the **RESET** cycle, before calling the After-Sales Service switch the appliance off and back on and then try again.




fig. 24 – **RESET** cycle selection buttons: hold down the **Start/Stop + Functions** buttons together




fig. 25 – RESET: during the RESET the display shows “P-“. At the end of the RESET cycle, “E-“ appears on the display.

4.13 PRINTING THE CYCLE AND APPLIANCE PARAMETERS

If the Smeg “WD-PRINTE” accessory has been installed, it is possible to force printing of the data relating:

1. To the last cycle performed: by pressing the **1 On/Off** and **2 Functions** keys simultaneously.
2. To the setup parameters, by pressing the **3 Functions** and **4 LH Selection** keys simultaneously.

Print last cycle		+	
Print parameters		+	

Commands are only accepted if the door of the device is closed.




fig. 26 – Standard procedure for forcing printing.
Last cycle data: use “On/Off”+“Functions” keys; Parameters: use “Functions”+“LH Selection” keys.

5 SETUP MODE - “FUNCTIONS”

The appliance has setup procedures for modifying its operating parameters.

The setup mode can be accessed by pressing the “Functions” button and requires a password.

For reasons of security and responsibility, 4 levels of password are available to access different levels of functions and settings::

1. USER Level (1111)
2. SUPERUSER Level (provided by the installation engineer)
3. ENGINEER Level
4. SMEG Level

5.1 ACCESSING THE SETUP MODE / ENTERING THE PASSWORD



- Hold the **Functions** button (5) down for 5 seconds. 4 dashes appear on the display.
- Use the **Selection** buttons (4-8) to modify the individual characters of the password.
- Short press the **Functions** button (5) to move on to the next character
- Once the 4 characters of the password have been entered, confirm by holding down the **Start/Stop** button (3).
 - If the password entered is correct, the system accesses the first menu item
 - If the password is not correct, two **beeps** signal the error.
 - 5 seconds after the last button is pressed, the system automatically quits the Setup menu.
- Once the Setup mode has been accessed, the various **items** available can be scrolled using the **Functions** button (5).
- Note on GW4060, GW4190 series: a short press on the button **Drying +** (6) has the same meaning of the short press on button **Functions** (5) in the following cases:
 - When entering password, to switch from one character to the next.
 - Going into the setup mode, to scroll through the various items available.



fig. 27 – **Functions** button: hold down to access the Setup mode.

The **items** accessible from the **SETUP** menu are shown in the table below. To move on from one **item** to the next, press the “**Drying**” button.

Setup Item	Description
Set	Appliance parameter setup - detailed explanation provided below in “Set” Parameters
Cloc	Date and time - detailed explanation provided below in “Cloc” Parameters
Filt	Drying filter counter. (Only for GW4060, GW4190 series) Press the Start/Stop button to display the number of absolute filter hours left before the maintenance message appears. The initial default setting is 500. The value is in units. (Display only parameter, cannot be modified).
Tec	Appliance operating hour adds up towards maintenance. Press the Start/Stop button to display the number of cycles left before the maintenance message appears. The initial default setting is 900. The value is in units. (Display only parameter, cannot be modified).
Coun	Appliance total operation counter (display only parameter, cannot be modified)
PrEn	Setting which allows display and selection only of the washing programs actually of interest.

5.2 ACCESS AND PARAMETER MODIFICATION

The three **options** referred to above - **Set, Cloc and PrEn** – are used for accessing and modifying:

1. Appliance parameters, **Set**.
2. Date and Time parameters, **Cloc**.
3. Setting to allow display and use of a subset of programs only, **PrEn (PrEn: “Program Enable”)**.

The following basic procedure applies to every **item**.

Step	Button	Action
1		Start/Stop – hold down to access the option required . To access the Set parameters, hold down until the “ Set ” option appears on the display. To access the Cloc (date and time setting) parameters, hold down until the “ Cloc ” option appears on the display. Similarly, for the “PrEn” option
2		Functions - Selection of the parameter required (see detail of the parameters in the points which follow). Note on GW4060, GW4190 series: a short press on the button Drying + has the same meaning of the short press on button Functions .
3		Selection buttons - Modification of the selected parameter During modification, the parameter concerned flashes on the display.
4		Start/Stop – hold down to confirm the modification. Once the modification has been confirmed the parameter stops flashing and appears on the display with steady light.
5		Door opening – hold down to quit the screen for display/modification of the current parameter.




fig. 28 – Steps for display and modification of the **Set** and **Cloc** parameters:
Access the parameter: by pressing the **Start/Stop** button
Quit the parameter: by pressing the **Door opening** button.




fig. 29 – **Set** parameters. Hold down the **Start/Stop** button to access the parameters. Parameter -1 sets to the dose dispensed by peristaltic pump **P1**. Use the **Selection** buttons to modify the parameter. Use the **Functions** button to move to another parameter. Confirm changes using the **Start/Stop** button.




fig. 30 – **Cloc** parameters. Hold down the **Start/Stop** button to access the parameters, in this case identified as **C1**, **C2**, etc. Parameter **C1** sets the current year. Use the **Selection** buttons to modify the parameter. Use the “**Drying Decrease/Increase**” buttons to move to another parameter. Confirm changes using the **Start/Stop** button.

5.3 “SET” PARAMETERS

Set on the display.

To access this item press the **Start/Stop** button.

The **Set** item contains a number of setting parameters.

User level operators can only display the parameters and cannot modify any “Set” parameters.

Set PARAMETER	DESCRIPTION	POSSIBLE USER ACTION								
		With set parameter “-d=0”	With set parameter “-d=1”							
-1	<p>Detergent pump P1 dispensing rate. Setting in ml/litre, with maximum threshold of 10 ml/litre. The default value of this parameter is 5ml/litre. Maximum recommended quantity with Smeg detergents: 10ml/litre. Increase and decrease with the Selection buttons.</p>	Display and Modification	Display							
-2	<p>Detergent pump P2 dispensing rate. Setting in ml/litre, with maximum threshold of 20 ml/litre. The default value of this parameter is 3ml/litre. Maximum recommended quantity with Smeg detergents: 10ml/litre. Increase and decrease with the Selection buttons.</p>	Display and Modification	Display							
-3	<p>Detergent pump P3 dispensing rate. (not present on GW2145) Setting in ml/litre, with maximum threshold of 20 ml/litre. The default value of this parameter is 10ml/litre. Maximum recommended quantity with Smeg detergents: 15ml/litre. Increase and decrease with the Selection buttons.</p>	Display and Modification	Display							
-4	<p>Water hardness setting Set the hardness in °f, in steps of 5. Range: 5-60°f. The default setting is 40°f. Warning: Make sure that the value is actually the same as the hardness of the water used.</p>	Display and Modification	Display and Modification							
-5	<p>Automatic opening at end of cycle The default setting on GW appliances is “On”, but automatic door opening can be disabled by changing the parameter to “OF”. If the option is enabled, the device door is opened automatically at the end of the cycle.</p>	Display and Modification	Display and Modification							
-6	<p>Printer language English by default, the language can be chosen from the various options available</p>	Display and Modification	Display and Modification							
	<table border="1"> <thead> <tr> <th>German</th> <th>French</th> <th>Spanish</th> <th>Italian</th> <th>English</th> </tr> </thead> <tbody> <tr> <td>dE</td> <td>Fr</td> <td>ES</td> <td>It</td> <td>En</td> </tr> </tbody> </table>			German	French	Spanish	Italian	English	dE	Fr
German	French	Spanish	Italian	English						
dE	Fr	ES	It	En						
-7	<p>Demineralised water The use of demineralised water can be activated or deactivated depending on whether or not the relative connection has been made (Activated: dn, Deactivated: --)</p>	Display and Modification	Display and Modification							

Set PARAMETER	DESCRIPTION	POSSIBLE USER ACTION	
		With set parameter "-d=0"	With set parameter "-d=1"
-8	<p>Residual filter hours (for GW4060, GW4190) The parameter indicates, to the nearest ten, the number of absolute filter hours left before the maintenance message appears. The initial default setting is 500. (e.g. the display shows "-8:50"). The value can be reset by an authorised engineer.</p>	Display Only	Display Only
-9	<p>Residual cycles before maintenance Indicates, to the nearest ten, the number of cycles left before the maintenance message appears. The initial default setting is 900 (e.g. "-9:90"). The value can be reset by an authorised engineer.</p>	Display Only	Display Only
-A	<p>Power Supply The parameter indicates whether the connection is: - single-phase: parameter setting "1" - three-phase: parameter setting "3" N.B.: the value "3" refers to both 400V 3N~ / 50 Hz and 230V 3~ / 50Hz connections. Warning: this parameter is only used for the correct setting of times remaining.</p>	Display Only	Display Only
-b	<p>Memory Data Overwriting - Overwriting enabled: parameter setting "1" - Overwriting disabled: parameter setting "0" The parameter refers to whether or not the device's internal memory, which keeps a record of all the cycles performed and the alarms which have occurred, can be overwritten. WARNING If overwriting is disabled, alarm AF:91 is triggered once the memory is full. In this case the appliance will not perform any more cycles until the data from the memory are downloaded; this can only be done using the TRACELOG software.</p>	Display and Modification	Display and Modification
-C	<p>Powder dispenser Presence (Parameter Introduced from the Main Firmware version 1.11.0.33) Parameter for use by authorized technicians only. This parameter is available only on certain models, where the detergent dispenser DD is installed on the inner door. It can assume the following values: "-C: dd", it means presence and activation of the dispenser powder (DD), it is neglected the dosage set by parameter P1: "Set -1". "-C: P1" for the presence of the peristaltic pump P1. The dosage of the detergent is determined by the value assigned to the parameter "Set -1".</p>	Display Only	Display Only

Set PARAMETER	DESCRIPTION	POSSIBLE USER ACTION	
		With set parameter "-d=0"	With set parameter "-d=1"
-d	<p>Superuser dosage setting disabling parameter. (Parameter introduced from Main Firmware version 4.6.xx.yy). If "-d: 0" the superuser is able to modify the dosages If "-d: 1" the superuser is not able to modify the dosages but can display them. Default value: "-d: 0". This Parameter can only be modified by the authorised engineer.</p>	Display Only	Display Only
-E	<p>Parameter for setting presence or absence of the WD-LANE accessory. (Parameter introduced from Main Firmware version 4.6.xx.yy). Sets the serial port communication speed "-E: Pr" = 9600 bps (for Printer connection) "-E: Ln" = 11500 bps (for connection to WD-LANE - LAN accessory).</p>	Display and Modification	Display and Modification
-t	<p>Test cycle (Parameter Introduced from the Main Firmware version 1.11.0.33) The parameter can be used by authorized technicians only; it allows the display and the start of a Test program, on the display "Pr t". This test program performs a short washing cycle and activates all the loads of the device. The cycle is defined to be used for machine verification after components replacement. - Parameter value "--", you do not have access to the cycle "Pr t" - Parameter Value "on", the technician gets access to the Test Program "Pr t". In order not to run the Test cycle "Pr t": switch Off the device to reset its condition.</p>	Display Only	Display Only

In the following pictures some examples of "Set" parameters are featured.




fig. 31 – "Set - C" parameter. The user can fix the use of the powder dispenser (= "dd") or the peristaltic pump (= "P1"). When "dd" is chosen, a possible value saved in "Set - 1" will be ignored by the system.




fig. 32 – "Set - t" value can be "--" or "on". In this latter case a test cycle will be activated.

5.3.1 APPLIANCE DATA STORAGE, OVERWRITING, "Set -b"

The appliance stores the data relating to the **cycles performed** and all the **alarms which have occurred** in its own **internal memory**, on a PCB.

The data in the appliance’s memory can be read using the **Smeg TRACELOG** software.

Overwriting of the data in the appliance’s memory can be enabled or disabled using a setup parameter, **Set: -b –**; for details, refer to the **SETUP MODE - "FUNCTIONS"** section.

Parameter **Set -b** has two possible values:

- "0" meaning "overwriting disabled"
- "1" meaning "overwriting enabled" – 1 is the **default value** for the parameter; overwriting is permitted and no alarm is triggered when the appliance’s memory is full.

WARNING



If overwriting is disabled, parameter "-b=0", alarm AF:91 is triggered once the memory is full. In this case the appliance is cut out and will not perform any more cycles until the data from the memory are downloaded to free spaces; this can only be done using the TRACELOG software.




fig. 33 – SET appliance parameters; parameter -b is used to enable and disable overwriting of the internal data memory.

The following is an example of data stored in the internal memory, with a key to their interpretation. The data can be accessed using the TRACELOG software.

=====		
2012/10/17 10:39:46	START PROGRAM: 7	Oper: 255
	Cycle N.: 288	
2012/10/17 10:40:17	START PHASE N.: 0	
2012/10/17 10:41:12	Inflow Detergents	P1:0.40%
2012/10/17 10:42:26	TL1: 22.7	TCL: 22.5
2012/10/17 10:42:26	Start heating	Target temp. 93
2012/10/17 10:56:12	TL1: 94.1	TCL: 93.6
2012/10/17 10:57:36	TL1: 93.9	TCL: 93.4
2012/10/17 10:57:52	TL1: 94.1	TCL: 93.6
2012/10/17 10:59:43	TL1: 93.9	TCL: 93.4
2012/10/17 10:59:58	TL1: 94.1	TCL: 93.6
2012/10/17 11:01:53	TL1: 93.9	TCL: 93.4
2012/10/17 11:02:08	TL1: 94.1	TCL: 93.6
2012/10/17 11:04:12	TL1: 93.9	TCL: 93.4
2012/10/17 11:04:26	TL1: 94.1	TCL: 93.5
2012/10/17 11:05:53	Ao>=15983	
2012/10/17 11:05:53	Start heating	Target temp. 0
2012/10/17 11:05:53	TL1: 94.1	TCL: 93.6
2012/10/17 11:07:49	START PHASE N.: 1	
2012/10/17 11:08:41	Inflow Detergents	P2:0.20%
2012/10/17 11:09:34	TL1: 43.3	TCL: 42.3
2012/10/17 11:09:34	Start heating	Target temp. 0
2012/10/17 11:11:34	TL1: 59.7	TCL: 59.2
2012/10/17 11:11:34	Start heating	Target temp. 0
2012/10/17 11:11:34	TL1: 59.7	TCL: 59.2
2012/10/17 11:12:33	START PHASE N.: 2	
2012/10/17 11:13:19	Inflow Detergents	
2012/10/17 11:13:41	TL1: 36.9	TCL: 36.5
2012/10/17 11:13:41	Start heating	Target temp. 0
2012/10/17 11:15:41	TL1: 49.9	TCL: 49.5
2012/10/17 11:15:41	Start heating	Target temp. 0
2012/10/17 11:15:41	TL1: 49.9	TCL: 49.5
2012/10/17 11:16:41	START PHASE N.: 3	
2012/10/17 11:17:26	Inflow Detergents	
2012/10/17 11:17:49	TL1: 33.3	TCL: 32.9
2012/10/17 11:17:49	Start heating	Target temp. 75
2012/10/17 11:25:47	TL1: 76.2	TCL: 75.7
2012/10/17 11:26:30	Start heating	Target temp. 0
2012/10/17 11:26:30	TL1: 76.9	TCL: 76.5
2012/10/17 11:27:29	Start drying: TA_TARGET: 120	TA1: 73
2012/10/17 12:04:47	End drying: TA1: 89	
2012/10/17 12:09:07	END PROGRAM: 7	
=====		

The double line marks the start and end of a program

Cons. number of CYCLE

First phase, called "0"

Detergent dispensing

Phase Target temperature

A0 value reached

Second phase, called "1"

The two tank temperature probes are called:
TL1 – working probe
TCL – control probe

Drying phase, with Target temperature in Dryer duct

fig. 34 – Example of cycle data stored in the appliance's memory.

5.4 “CLOC” – DATE AND TIME PARAMETERS

The **Cloc** parameter (an abbreviation of “**CLOCK**”), is used to **modify the current date and time**. This parameter can be modified with just the user password.



The modification procedure is:

1. To select the **Cloc** parameter, hold pressed the **Start/Stop** button (3) when the indication “**Cloc**” is on the display (ref. also to the par . “**ACCESSING THE SETUP MODE**”).
2. **Selection** buttons (4-8): for modifying the value of the parameter.
3. **Functions** button (5) – quick press: to move from one parameter to the next. (on GW4060, GW4190 the operator can also use the “**Drying +**” key (6) with the same function).
4. **Start/Stop** button (3): hold pressed to confirm the modification
5. **Door opening** button (7): hold pressed to quit the parameter being modified and the Cloc function.

During modification, the parameter concerned flashes.

CLOC PARAMETER	DESCRIPTION
C1	Year, last two figures only (00-99), e.g.: for “2012” enter “12”
C2	Month (1-12)
C3	Day (1-31)
C4	Hour (1-24)
C5	Minutes (0-59)
C6	Seconds (0-59)




fig. 35 – **Cloc** parameters. Hold down the **Start/Stop** button to access the parameters, in this case identified as **C1**, **C2**, etc. Parameter **C1** sets the current year, **C2** the month, etc. Use the **Selection** buttons to modify the parameter. Use the “**Functions**” button to move to another parameter. Confirm changes using the **Start/Stop** button.

5.5 SETTING THE “PrEn” “PROGRAM ENABLING” FUNCTION

By means of the appropriate settings, the “**PrEn**” parameter allows the user to only display and start the programs actually of interest. This setting can be made from “superuser” level upward. (Function introduced from Main Firmware version 4.6.xx.yy).






Once the menu has been selected, the operator can scroll through the user programs (P01-P15).

The first 2 digits of the display show the program number and the second 2 the “on” or “oF” value (e.g. “06:on”) depending on whether or not the program currently selected is enabled.

By default: all programs are enabled, “on”. To disable a program, set “oF” next to its number (e.g. “06:oF”) using the **Modify** buttons.

At least 1 program is always enabled; the superuser can disable up to 15 of the 16 programs available.

The procedure is similar to that in use for the setup parameters.

Step	Button	Action
1		Start / stop – hold down to access the option . To access the PrEn parameters, hold down until the “ PrEn ” setup option appears on the display.
2		Functions – quick press - Selects the programs to be enabled/disabled.
3		Buttons for Selection - Modification of the selected parameter (switch from “on” to “oF” to disable display of the program on the screen).
4		Start / Stop – hold down to confirm the modification .
5		Door opening – hold down: to quit the screen used for display/modification.




fig. 36 - Example of use of the "PrEn" function to disable display of program "06" on the screen.

6 INSTALLATION INSTRUCTIONS

Once the appliance has been installed correctly, it must be prepared for use as follows. The authorised engineer commissions the device in accordance with the set procedures and records its serial number.

The basic operations are as follows:

- Fill up with regenerating salt.
- Add detergent and neutralising agent and any other chemicals used.
- At first installation, program “16 Service” must be run, preferably twice in succession, to operate the peristaltic pumps and prime the detergent inlet lines correctly.
- Run a complete program including a thermal disinfection phase with no load for processing in the tank.

The device is now ready for use.

6.1 USE OF THE WATER SOFTENER

The quantity of limescale in the water (which depends on its hardness) may leave a white veil on the instruments after being washed. These instruments may become dull over few washing cycles. To avoid this from happening, the appliance is equipped with an automatic water softener which removes the substances that cause hardness from the water with the aid of regenerating salt.




fig. 37 – The salt tank is accessible with the door open; filling the softener tank with regenerating salt.

If water of average hardness is used, fresh salt will have to be added about every 20 washes. The softener tank contains **about 1 kg of coarse salt**.

The salt tank is in the bottom of the glassware washer’s chamber.

Remove the bottom basket, unscrew the cap of the tank by turning it anticlockwise and **pour in the salt** using the funnel supplied with the glassware washer; take care not to spill salt in the chamber. Remove any residual of salt around the hole and in the tank before screwing the lid back on.



ATTENTION

- **The first time the glassware washer is used, a litre of water must be poured into the tank with the salt.** Whenever the tank is filled, always take care to ensure that the cap is firmly closed. The mixture of water and detergent must not get into the salt tank, as this would impair operation of the regeneration system. In this case the warranty cover will cease.
- **Only use coarse salt for domestic dishwashers.**
- **Do not use edible salt**, since it contains insoluble substances which might damage the water softening system over time.
- **Fill up with salt when necessary before starting a washing program**, so that any excess salt solution will immediately be removed by the water; if salty water is left inside the chamber, it may cause corrosion. Carry out a prewash to prevent this if necessary.

6.2 USING DETERGENT AND NEUTRALISING AGENT

The appliance is equipped with detergent dispensing devices.

1. **Pump P1**, dispensing of neutral or weak alkaline liquid detergent (optional on **GW** series).
2. **DD powder dispenser** (alternative to P1): powder detergent dispenser, on the inner door.
3. **Pump P2**: dispensing of acid neutralising agent.
4. **Pump P3**: optional, never installed by default. Auxiliary detergent/lubricant pump.

SYMBOL	MEANING
●	feature present
○	optional accessory, installable on the model in question.
-	feature not present and not installable on the model in question.

Dispensing components fitted on the various models (OPT = optional accessory).

GW Models	GW2145	GW1160 GW4060 GW4190	GW1160S GW4060S GW4190S	GW1160SC GW4060SC GW4190SC	GW1160C GW4060C GW4190C
DD inner door dispenser	●	●	-	-	●
P1 peristaltic detergent pump	○	○	●	●	○
P2 peristaltic neutralising agent pump	●	●	●	●	●
P3 extra peristaltic pump	-	○	○	○	○

Except for "PREWASH" programs only, dispenser **DD** (or pump **P1** if installed) is activated automatically during washing.

The neutralising agent for rinsing is automatically dispensed by pump **P2** during the phase after washing.

It is important to use good quality detergents to obtain good results. Keep detergent bottles closed in a dry place to prevent the formation of lumps which might lead to poorer washing results.

Once bottles have been opened they must not be kept for too long, since the detergents tend to lose their properties.



ATTENTION

If the level sensors are not fitted: check the levels of product in the cans / bottles to ensure that programs are not performed without detergent or neutralising agent.



During installation, or when a completely empty jerry can of product is replaced, run program no. 16 SERVICE to charge the system with liquid. The pipeline from the jerry can to the pump must be filled to ensure that the product is dispensed properly during subsequent washing cycles.

6.2.1 LIQUID DETERGENT INTAKE SYSTEM

Each peristaltic pump is supplied with a detergent intake system, in two possible versions.

a. STANDARD version (Intake system installed by default on the **GW** range):

1. one stainless steel intake pipe, to suck in detergent from the can;
2. a support, a conical rubber-fitting support, to position the intake pipe correctly;
3. the silicone hose for connection between stainless steel intake pipe and peristaltic pump.
4. Detergent intake filter, fitted directly on the stainless steel pipe.

b. Configuration with LEVEL SENSOR (optional on **GW** range):

1. Suction nozzle with integral level sensor and tapered rubber support for fitting onto the can.

2. Silicone hose for connection between intake pipe and peristaltic pump.
3. Detergent intake filter, fitted directly on the suction nozzle pipe.

WARNING

It is essential that the label on the intake pipe corresponds to the type of detergent being used.

Refer to the following colour code:

1. White/clear: P1, alkaline detergent
2. Red: P2 acid neutralising agent
3. Blue: P3, additive

Incorrect connections (e.g. if detergents P1 and P2 are exchanged) will reduce effectiveness of the process.




Connection errors leading to the mixing of different chemicals in the intake systems may cause irreparable damage to parts of the circuit.

Detergent intake system connection errors lead to the loss of Warranty cover for the parts concerned (peristaltic pumps, detergent suction and intake circuit, and level and flow sensors if installed).

The intake pipe is supplied complete with intake filter. The filter pushes into place on the pipe. Make sure that the filter is always fitted and correctly positioned to keep the detergent dispensing system in good working order.

Check regularly that the silicone hose is properly fixed to the stainless steel intake pipe and that there are no leaks. Use suitable fittings (e.g. plastic ties) to secure the silicone hose to the stainless steel pipe.



*fig. 38 – Detergent intake **STANDARD** version. Fitting the detergent intake pipe to the can. Fit the rubber cap over the top of the can for a perfect, secure connection. The pipe is supplied complete with intake filter.*

6.2.2 DETERGENT LEVEL SENSORS - OPTIONAL

The appliance can be fitted with level sensors integral in the detergent intake pipe, Smeg product "WD-LS xxx" (*Level Sensor*; the value of xxx may vary depending on the suction nozzle configuration required, for GW1160, GW4060: WD-LS3060).

Contact Smeg for advice about the accessory best suited to your needs and your device.



WARNING

The electrical connection to level sensors must be made using the terminal board provided inside the rear crosspiece (each peristaltic pump has a terminal for connection of the relative level sensor).

This procedure must be carried out by Smeg authorised technical staff.




fig. 39 – Configuration with LEVEL SENSOR: Detergent suction nozzle with integral level sensor. The suction nozzle is complete with filter: make sure that the filter is always correctly positioned to ensure no lumps are drawn in.

6.3 USING THE DETERGENT DISPENSER – GW SERIES ONLY

The “DD” powder detergent dispenser is on the inside of the door, integral in the “inner door”.

Except for the “Prewash” program, the user must place the right amount of detergent in the dispenser before each washing.

Press the catch to open the detergent dispenser lid.

Then add the detergent and close the lid carefully.




fig. 40 – DD powder detergent dispenser, accessible with door open.

During washing, the dispenser opens automatically.

When **Smeg Deterglass** detergent is used, the recommended dose, for the standard cycle, 9 water liters per phase, is **27-45 g** (3 – 5 g / liter), the quantity shown by the ridges on the inside rim of the dispenser.

To establish the right amount, use the **measuring beaker provided: the use of a measuring device will help to ensure the repeatability of the process.**



fig. 41 – Detergent measuring beaker provided.

Refer to the product technical and safety data sheets.

6.4 RECOMMENDED DETERGENTS

One of the fundamental factors in achieving a good washing process is the type of detergents used.

Smeg guarantees excellent washing results with the use of the products recommended below.

Contact your Smeg dealer for advice about the types of detergents suitable for your application, and their modes of use.



ATTENTION

Always comply with the INSTRUCTIONS provided by detergent PRODUCERS, especially with regard to the RECOMMENDED DOSES and the correct TEMPERATURES for their use.

Safety information about doses supplied in compliance with 5.4.4.s of IEC61010-1-040.

Refer to the product instructions and technical safety data sheets.

Technical data sheets are available on request from SMEG S.p.A. instruments@smeg.it

N.B.: the THIRD pump, P3, is an optional accessory.

Recommended detergents for GW series

The following is a list of the products with which Smeg guarantees the specified washing results.

Alkaline detergents:

- Smeg DETERLIQUID D and D2
- Smeg DETERLIQUID SP

Powder detergent

- Smeg DETERGLASS







Acid neutralising agents:

- Smeg ACIDGLASS P and P2

Special additives:

- Contact Smeg
-

GENERAL PRECAUTIONS FOR DETERGENTS

	<p>HANDLE DETERGENT CANS WITH CARE Warning: detergents may be TOXIC. Refer to the product safety data sheets. Once a product has run out, replace the empty can with a full one of the same product. If the product left in the old jerry can is poured into the new one, take care not to overfill the new cans to ensure that they will not overflow when the suction nozzles are inserted. Gloves must be worn during transfer of detergents, during fill-ups and during insertion of the suction nozzles. Information supplied in compliance with 5.4.3.m, 5.4.4.n and 5.4.4.q IEC61010-2-040:2005</p>
	<p>FIRST AID MEASURES IN CASE OF CONTACT WITH DETERGENTS Instructions provided in compliance with point 5.4.4.p IEC61010-2-040. Remove contaminated clothing and put it in a safe place. Contact with the skin or eyes: wash immediately and thoroughly with water. If possible, apply a sterile gauze dressing. Seek medical advice. Swallowing: rinse out the mouth with plenty of water. Seek for medical advice immediately.</p>
	<p>DETERGENT SAFETY DATA SHEETS. De detergent SAFETY DATA SHEETS must be kept:</p> <ol style="list-style-type: none"> 1. close to the place where detergents are stored; 2. close to the appliance. <p>In easily accessible positions. The latest versions of the safety data sheets should be requested regularly (e.g. once a year). Safety data sheets will be supplied by Smeg S.p.A. on request.</p>
 	<p>DISPOSAL Information supplied in compliance with point 5.4.4.L of IEC61010-2-040:2005. DISPOSAL of any product residues and containers (cans and bottles): refer to the "DISPOSAL ADVICE" section of product safety data sheets. If mislaid, data sheets may be requested from SMEG S.p.A. instruments@smeg.it The appliance's responsible authority must dispose of detergent residues and containers in accordance with the relevant national or local legislation.</p>
	<p>FLAMMABILITY Always refer to detergent technical data sheets to assess the products' flammability. Do not use flammable products in the appliance.</p>

7 ALARMS and WARNING MESSAGES

This section provides instructions for interpreting alarm messages and taking any appropriate measures. Information supplied in compliance with point 5.4.4.j of IEC61010-2-040:2005.

The appliance’s display has backlit LED symbols.

See also the explanations provided in the INDICATOR LEDS - DETAILED EXPLANATION section.

The symbols at the top of the display contain **WARNINGS** which the appliance provides to its users, as detailed below (For full information, also see the “DESCRIPTION OF CONTROLS AND PROGRAMS” section).










fig. 42 – Central display, detail of backlit symbols.

SYMBOL	MEANING
	<p>No detergent Lights up if the detergent level sensors are installed and one of them is at "minimum" level. The error is displayed at the end of the cycle and in combination with the LED the segment display shows the code of the detergent to which the warning refers: P1: code on display “A-:68” P2: code “A-:69” P3: code “A-:70”</p> <p>Activated at the end of the cycle and when the user attempts to start a new program. The washing cycle can be started with this message present (i.e. overriding the warning) by pressing the Start/Stop button again. If two or more similar alarms (e.g. cans P1 and P2 empty), the Start/Stop button will have to be pressed twice or more before the cycle can be started. (The device’s internal memory records the event but allows execution of the cycle).</p>
	<p>Maintenance: filter replacement (only on GW4060, GW4190 series) Active if the dryer is present and the absolute filter is installed. (The absolute filter is optional on the “GW” series). The LED indicates that the absolute filter needs to be replaced. Illumination of the symbol is triggered by the number of drying system operating hours, set at 500h when the filter is new. The filter must be replaced by Smeg authorised staff. Keeping the filter in operation after the end of its useful life will lead to poor drying results. If the filter is full, the drying air flow rate is reduced. The filter’s useful life may be less than the 500 h set, depending on the amount of dirt in the environment. If drying performance deteriorates, users should have the filter changed earlier. The appliance does not signal this type of wear, but only wear arising from the hours actually worked.</p>

SYMBOL	MEANING
	<p>Appliance maintenance The appliance counts the number of cycles performed and the LED comes on to alert the user that maintenance is required. These are inspection and maintenance operations, scheduled every 1000 cycles, which must be performed by Smeg authorised staff to keep the appliance safe and in good working order.</p> <p>These maintenance operations are not covered by the product Warranty, which does not include the replacement of components which deteriorate due to normal wear and tear.</p>
	<p>No Salt Salt must be added to the water softener reservoir in the chamber. Salt is necessary to enable the softener to reduce the hardness of the intake water. The intake water hardness is set on the appliance parameters on the basis of the information provided by the user. Regular checks (every 6 months/year) to ensure that the initial characteristics are retained are recommended. If the characteristics change, contact Smeg for modification of the parameters accordingly.</p>
	<p>Alarm The machine displays an anomaly, which can be generated by a fault or by the detection of abnormal conditions. The segment display shows "AF" plus a code number identifying the alarm triggered. The alarms table below details the components which may have caused the event for each alarm code, to simplify diagnostics and troubleshooting.</p>



The following is a list of the alarms which may occur on the device, they are of two types: "FATAL" and "NON FATAL".

- In the first case, "**FATAL ALARMS**", the message which appears is "AF" followed by the alarm number.
- The second category, "**NON FATAL ALARMS**", are warning messages and not malfunctions; the display shows "A-" followed by the alarm number. **Press the "Start/Stop" button to continue using the device.**



HOW TO DEAL WITH ALARMS - READ CAREFULLY

When an alarm occurs, the device's electronics will react to set itself in condition of safety..

When an error occurs, generally the device manages the anomaly itself: in this case **the alarm code flashes on the display** until the problem has been dealt with: **no reset commands are accepted during this time.**

Any RESET forced by the user will not be accepted until the problem has been dealt with and the alarm code on the display is "steady".

When an alarm message is noticed: Refer to the table below to find out what the code means and the relevant countermeasures. Apply the countermeasures recommended for the specific circumstances.

For convenience, the procedures that the user is normally required to carry out are summarised below:

1. First, follow the **DEFAULT PROCEDURE**.
2. If the alarm does not disappear: proceed as described for the **RESET PROCEDURE**.

1. DEFAULT PROCEDURE

If the event which triggered the alarm is due to a temporary anomaly, the situation can be resolved as follows:

Switch off using the "On/Off" button and then switch on again using the same button. The display shows a flashing "OFF" during automatic management of the anomaly. The word "OFF" becomes steady once the situation has been resolved.

2. RESET PROCEDURE for AF alarms:

1. Perform the **RESET** cycle.
2. If the alarm message disappears after the **RESET** cycle: normal use of the device can be restarted by switching it off and back on with the **On/Off** button.
3. If the alarm message persists after the **RESET** cycle: switch the device off and back on using the "On/Off" button and by **disconnecting it from the electricity supply** (Wait at least 10 seconds after switching off before switching back on again).
4. If the alarm occurs when the device is restarted: retry the **RESET**.

If the alarm persists even after the RESET procedure:

1. **Turn off the water supply taps.**
2. **Disconnect the appliance from the electricity supply.**
3. **Check that all the device's connections** (electricity and water) are correct and there have been no changes from the initial installation conditions.
4. **Contact the Smeg After-Sales Service.**




fig. 43 – RESET cycle selection buttons ("Start/Stop" + "Functions") buttons

ALARM ID	DESCRIPTION	USER ACTION
AF:01	Water not heated within time allowed.	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:02	The temperature difference between the two probes, "TL1" and "TC", is more than 2°C (only active on WD series)	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:04	Probe "TL1" (chamber temperature) reads a higher temperature than the effective one.	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:05	Working temperature probe "TL1" is generating an abnormal signal (probe "open").	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:07	Drying temperature probe TA1 is generating an abnormal signal (probe "open").	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:10	Control temperature probe TLC is generating an abnormal signal (probe "open"). (only active on WD series)	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:11	No cold water during filling.	<p>Check the water supply:</p> <ol style="list-style-type: none"> 1. Check that the intake tap is open. 2. Check the water supply pressure. 3. Check that the filler hoses are correctly positioned. <p>Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p>
AF:13	No demineralised water.	<p>Check the water and demi water supply:</p> <ol style="list-style-type: none"> 1. Check that the intake tap is open. 2. If using water from a tank (with PAD accessory), check that the tank is not empty or placed too low down. 3. Check the water supply pressure. 4. Check that the display settings actually correspond to the water connections (demi water present or not). <p>Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p>
AF:17	Cold water filling time not correct.	Checks and procedures as for AF:11
AF:19	Demi filling time too long. The appliance is taking too long to take in demineralised water.	Checks and procedures as for AF:13
AF:23	Not enough water. Chamber water level too low.	<p>Checks and procedures as for AF:11</p> <p>Also check that there are no leaks from the appliance (no water around it). If there are leaks: turn off all water supply taps at once and contact the After-Sales Service.</p>
AF:25	Water circuit malfunction. Anomaly related to the circulating pump "ML". Circulating pump pressure too low.	<p>There may be foam in the chamber. Check the type of detergent used.</p> <p>For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p>

ALARM ID	DESCRIPTION	USER ACTION
AF:26	Water being taken into chamber at wrong time. Cold water intake valve "EVF" may have failed.	<p>Check the water supply: Check that the water supply pressure is within the permitted range.</p> <p>Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p> <p>Also check that there are no leaks from the appliance (no water around it).</p> <p>If there are leaks: turn off all water supply taps at once and contact the After-Sales Service.</p>
AF:28	Demineralised water being taken into chamber at wrong time. Cold water intake valve "EVF" may have failed.	<p>Check the water supply: Check that the demineralised water supply pressure is within the permitted range.</p> <p>Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p> <p>Also check that there are no leaks from the appliance (no water around it).</p> <p>If there are leaks: turn off all water supply taps at once and contact the After-Sales Service.</p>
AF:29	Chamber does not empty. No drainage.	<p>Check the connection to the water drain, especially that the height of the drain connections complies with the specified requirements and that there are no restrictions in the drain hoses.</p>
AF:30	During the working cycle, the chamber water level exceeds the safety level. Water safety level.	<p>Check the water supply to the appliance:</p> <ol style="list-style-type: none"> 1. Intake pressure. 2. That the connections are correct, as specified in this manual. <p>For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p> <p>If there are leaks: turn off all water supply taps at once and contact the After-Sales Service.</p>
AF:32	Water standing in washing chamber with appliance in standby.	<p>For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p> <p>If there are leaks: turn off all water supply taps at once and contact the After-Sales Service.</p>
AF:33	Steam condenser water intake anomaly. No water in steam condenser.	<p>Indicates that there is no water in the steam condenser when there should be.</p> <p>Check the water supply:</p> <ol style="list-style-type: none"> 1. Intake water pressure 2. That the connections are correct, as specified in this manual. 3. Blockages or restrictions in hoses. <p>Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p>
AF:34	Steam condenser water drainage failure.	<p>Check that the drain connections are correct:</p> <ol style="list-style-type: none"> 1. Height and position of the drain connection 2. That the connections are correct, as specified in this manual. 3. Blockages or restrictions in hoses. <p>For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE.</p>

ALARM ID	DESCRIPTION	USER ACTION
AF:37	Drainage problems when using a mixture of water.	The use of a mixture of water is an option adopted in special cases to cool the water discharged into the drain. Check: 1. The intake water temperature; if the water is not cold enough (recommended $T < 25^{\circ}\text{C}$) there may be difficult in achieving the target temperature. 2. Problems with drain; check that the hoses and drain connections meet the specified requirements.
A-:41	P1 detergent intake anomaly	1. Check that the detergent intake nozzle is correctly positioned in the can. 2. Check for detergent leaks (pools of detergent around the device). This code provides a warning ; it is not an alarm which cuts out the appliance. Operation of the device can be continued by pressing the Start/Stop button.
A-:42	P2 detergent intake anomaly	Proceed as for AF:41
A-:43	P3 detergent intake anomaly	Proceed as for AF:41
AF:54	Door opening detected with cycle in progress. Door interlock microswitch malfunction.	1. Make sure the door is closed properly before starting a cycle. 2. Do not force the door open with a cycle in progress; always use the buttons provided on the appliance to stop a cycle and open the door. 3. Check there is nothing between the door and the appliance's chamber preventing the door from closing properly. 3. Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:56	Door locking device malfunction, automatic opening failure	Always make sure the door is closed properly before starting a cycle. Try to open the appliance again by pressing the Door opening button. If necessary, release the door lock by hand as described in this manual. Follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:58	Dryer heating failure. (Only on GW4060, GW4190).	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:67	Dryer motor "cooling" malfunction. Cooling is included at the end of the drying phase to bring the load processed and the heating elements to a safe temperature. (Only on GW4060, GW4190).	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .

ALARM ID	DESCRIPTION	USER ACTION
A-:68	P1 can empty	<p>Check that there is detergent in the jerry can or that the level sensor is operating correctly.</p> <p>This is a warning and not an alarm as such; the washing cycle can be started with this message present (i.e. overriding the warning) by pressing the Start/Stop button again (hold it down for 2 seconds).</p> <p>Activated at the end of the cycle and when the user attempts to start a new program. If two or more similar alarms (e.g. cans P1 and P2 empty), the Start/Stop button will have to be pressed twice or more before the cycle can be started. (The device's internal memory records the event but allows execution of the cycle).</p>
A-:69	P2 can empty	As for "A-:68"
A-:70	P3 can empty	As for "A-:68"
AF:73	Internal memory data storage error	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:74	Water leak from chamber. Alarm only enabled with the Aquastop accessory fitted.	Turn off the water supply taps. Contact the After-Sales Service.
AF:75	No salt in softener. This alarm does not appear on the Display but is recorded in the appliance memory; the relative LED illuminates on the display.	Add salt to the softener salt tank inside the chamber, accessible with the door open.
AF:77	Intake water temperature over 45°C; prewash temperature must be below 45°C. Alarm generally disabled on GW series products.	The alarm occurs when the temperature is over 45°C at the start of the cycle. Wait for the appliance to cool before starting a new cycle.
AF:78	<i>Restore fail.</i> Problem of motherboard.	For the user: follow the DEFAULT PROCEDURE described above. If the alarm does not disappear, follow the RESET PROCEDURE .
AF:79	Program not compatible.	<p>This alarm is triggered if a program is created using the TRACELOG software using procedures not compatible for the proper performance of the cycle. Reconsider the program created phase by phase; it may be useful to compare it with an original Smeg program provided as reference (see attached PROGRAMS TABLE doc.).</p> <p>Pay special attention to the version of the appliance being used; e.g. do not activate peristaltic pump P3 when writing the program if this pump is not installed.</p>
AF:84	Dryer temperature reading higher than actual temperature. (Only on GW4060, GW4190).	<p>Wait for the alarm to be resolved automatically.</p> <p>During automatic management of the AF:84 alarm the RESET cycle cannot be started.</p> <p>Do not disconnect the electricity supply to the device: the overheating problem is being dealt with automatically, with a heating element cooling cycle.</p> <p>The DEFAULT PROCEDURE described above cannot be implemented until the end of the automatic process.</p>
AF:91	Internal memory full.	<p>This alarm only occurs if the appliance parameter that prevents overwriting in the memory is enabled.</p> <p>With the default setting, this alarm will not occur.</p> <p>Memory space must be cleared to allow use of the appliance: this can be done by connecting to the appliance by means of the RS-232 serial port and using the Smeg TRACELOG software.</p>

ALARM ID	DESCRIPTION	USER ACTION
AF:92	Drying filter maintenance	<p>The absolute filter in the drying system has exceeded the set number of working hours. Contact the After-Sales Service to have the filter changed. The filter is a consumable and is not covered by the warranty.</p> <p>Press the Start/Stop button to override the message and continue using the device.</p> <p>The device's internal memory records the event.</p> <p>With a fouled filter, drying will be poor.</p>
AF:93	Appliance maintenance	<p>The appliance has exceeded the set number of operating hours since installation or since its last service: contact the After-Sales Service for routine maintenance.</p> <p>Press the Start/Stop button to override the message and continue using the device.</p> <p>The service must be carried out as soon as possible to keep the device in good working order.</p> <p>The device's internal memory records the event.</p>
AF:94	Temperature below actual temperature during extension phase. T not stable	<p>For the user: follow the DEFAULT PROCEDURE described above.</p> <p>If the alarm does not disappear, follow the RESET PROCEDURE.</p>
AF:96	Chamber water level incorrect.	<p>Check the intake and drain water connections and that the appliance has been installed as specified.</p> <p>If leaks are noticed around the device: turn off the water supply taps.</p> <p>For the user: follow the DEFAULT PROCEDURE described above.</p> <p>If the alarm does not disappear, follow the RESET PROCEDURE.</p>

8 CLEANING AND MAINTENANCE

DISCONNECTING THE ELECTRICITY AND WATER SUPPLIES

Before doing any work: disconnect the power supply using the appropriate switch on the panel or disconnecting the power cable, and turn off the water intake tap.



WORKING SPACE REQUIRED

A space of approx 1m² is required in front of the appliance for cleaning and maintenance to be carried out without difficulty.

POWER CABLE REPLACEMENT

If damaged power cables have to be replaced, **type FROR 450/750** cables with the same gauge as the cables mounted in the factory must be used (the cables fitted are marked with their types and characteristics). This replacement may only be made by authorised technical staff. Use genuine Smeg spare parts.



WARNING

Any work done on the appliance by unauthorised staff is not covered by the warranty and is performed under the user's responsibility.



During maintenance and cleaning: if appropriate, use **personal protection equipment**.

8.1 CLEANING THE DEVICE AND ITS PARTS

General cleaning

The external surfaces and inner door of the device must be cleaned at regular intervals (once a month is recommended) with a soft cloth soaked in water or a normal detergent for steel surfaces.

The door seals must be cleaned with a damp sponge.

After cleaning, a washing cycle should be carried out without a load for processing, to remove any detergent residues.




fig. 44 – The inside and outside of the device must be cleaned.

Cleaning the water intake filter

The water intake filter "A" on the outlet of the tap requires regular cleaning; the recommended interval is **once every 2 – 6 months, depending on the intake water** quality. Turn off the water supply tap, unscrew the end of the water intake hose, remove the filter "A" and clean it gently under running water. Fit the filter "A" back in place and reconnect the water intake hose with care.

Keep an eye on the free end of the hose to prevent water spills.




fig. 45 – water intake filter "A"

Cleaning the spray arms

The spray arms can easily be removed by unscrewing the knurled nut that fixes them to the pivot pin, so that the nozzles can be periodically cleaned to prevent fouling.

Wash the spray arms under running water and then fit them back in place. Make sure that they are free to turn unimpeded.

Recommended spray arm cleaning interval: **once a week**.

For trolleys with fixed spray nozzles: refer to the trolley manual for cleaning instructions.




fig. 46 – Appliance bottom spray arm: remove and clean regularly to maintain optimal washing results.


Cleaning the filter unit

The filter unit consists of a coarse filter, a conical fine filter, and a larger external filter. To ensure that the appliance operates effectively, it is extremely important to keep the filters clean. They should therefore be inspected frequently (e.g. if glassware with paper labels is washed, inspect after each cycle - **in normal conditions the filters should be cleaned once a week**) to remove deposits which may adversely affect operation.

Recommendations for correct maintenance

- The filters should be cleaned under running water using a stiff brush.
- It is essential to clean the filters carefully in accordance with the instructions given above: the appliance washer will be unable to operate correctly if the filters are fouled.
- Fit the filters back in place with care before starting a washing program.

8.1.1 GW2145 - filter elements



The central filter "D" should be inspected regularly and cleaned if necessary. To remove it, take hold of the handle, turn anti-clockwise and lift.

Push the central filter "D" up from below to extract it from the micro-filter and separate the two parts of the plastic filter by pressing the body of the filter in the zone shown by the arrows. Lift out the central filter.

8.1.2 GW1160 - GW4060 – GW4190 - filter elements




fig. 47 – Draining filters in the chamber; dismantling procedure and detailed view.

Coarse conical filter

To remove the coarse filter, press the "tabs" of the filter and pull upwards (detail A, see the figure above). Clean the filter and replace it.

Conical fine filter

This is located underneath the conical coarse filter (B) and should be checked and cleaned whenever the coarse filter is inspected. Clean with a brush and hot water for effective removal of deposits.

External circular filter

To remove this filter:

- Hold of the "tabs" of the coarse conical filter and turn it anticlockwise (C).
- Without pressing the tabs, lift the entire unit (D).

When cleaning this filter, it is best also to clean the others.

8.2 IF THE DEVICE IS TO BE OUT OF USE

If the appliance is to be left unused for a long time, follow these recommendations. The procedure is recommended for periods of 24 hours or more out of use.

- Perform the PREWASH program 1 with the appliance empty.
- If contaminated materials are processed: run program 15 (or a program with a thermal disinfection Td phase) with no load in the appliance.
- Disconnect the electricity supply.
- Leave the door slightly open to prevent unpleasant odours from forming inside the washing chamber.
- Turn off the water intake tap.

8.3 REUSING THE DEVICE AFTER A PERIOD OUT OF USE

If the appliance has been unused for a long time, before starting a cycle, follow these recommendations.

- Inspect the filters at the ends of the water intake hoses and check that no sludge or rust has deposited in the pipelines; if this has occurred, allow water to run from the water supply tap for a few minutes.
- Restore the electricity supply (if disconnected).
- Reconnect the water supply hose and turn on the tap.
- Run PROGRAM 15 (or another program with a thermal disinfection Td phase) with no load for processing.



A thermal disinfection cycle should be performed before using the device if it has been out of use for 24 hours or more.

8.4 TROUBLESHOOTING

Slight faults can sometimes be eliminated by the user with the aid of the following instructions.

1. If the programme fails to start, make sure that:

- The appliance is connected to the electrical mains.
- The appliance is receiving power.
- The water tap is turned on.
- The door of the appliance has been closed properly.

2. If water stagnates in the device, make sure that:

- There are no tight bends in the drain hose.
- The drain trap is not clogged.
- The appliance filters are not clogged.

3. If the items in the load for washing in general are not cleaned properly, make sure that:

- The right amount of detergent has been added.
- There is regenerating salt in the tank.
- The instruments have been positioned correctly.
- The program selected is suitable for the type and degree of dirt on the instruments.
- All the filters are clean and correctly fitted.
- The spray arm water outlet holes are not clogged.
- Nothing is preventing the spray arms from turning.
- Make sure that the salt tank cap is firmly closed.

4. If the items in the load fail to dry or are dull, make sure that:

- Absolute filter condition: operating hours can be viewed as an appliance parameter (the default operating hours are set for normal conditions of use; the drying filter's working life will be reduced by dirtier than average environments).
- There is neutralising agent in the container.
- The neutralising agent dispenser is set correctly.
- The detergent used is of good quality and has not lost its characteristics (e.g. owing to incorrect storage, pack left open, etc.).

- Make sure that the salt tank cap is firmly closed.
- 5. If the elements processed are stained make sure that:**
- The amount of neutralising agent dispensed is not excessive.
- 6. If there are visible traces of rust in the chamber:**
- The chamber is made of corrosion-proof steel, so any rust marks are due to external factors (pieces of rust from the water pipes, etc.). Specific products for removal of these stains are commercially available.
- Make sure that the detergent doses are correct. Some detergents can be more corrosive than others.
- Make sure that the salt tank cap is firmly closed.
- 7. If the optional printer does not work:**
- Check that the accessory is fitted with suitable thermal paper.
- Check that the device is connected properly (both electricity supply and data connection).
- On GW4190: check that the data output switch (next to the on-off switch) is set correctly for operation of the printer.

If the faults persist after the instructions given above have been followed, contact your nearest authorised service centre.

9 ROUTINE INSPECTIONS AND MAINTENANCE

9.1 DAILY

1. Check the detergent and neutralising agent levels in the cans and fill up if necessary.
2. Check spray arm rotation and inspect the water outlet holes visually to check that they are clean.

9.2 WEEKLY

1. Clean the sump filters, following the instructions provided.
2. Run PROGRAM 6, or another program with a thermal disinfection phase, without instruments for processing in the appliance, for precautionary cleaning/disinfection of the washing chamber.

9.3 EVERY SIX MONTHS

1. Inspect the water solenoid valve intake filters; if necessary, clean them by pouring hot water through them in the opposite direction to the operating flow.
2. Inspect the intake and outlet hoses of the detergent and neutralising agent pumps.

9.4 EVERY YEAR

Contact your nearest authorised Smeg service centre to have the appliance completely checked over at the end of the warranty period and annually for successive years, or, if this occurs before the end of the year, when the “**Appliance maintenance**” LED illuminates.



Maintenance operations are not covered by the product Warranty, which does not include the replacement of components which deteriorate due to normal wear and tear.

The operations which Smeg authorised staff will carry out are:

1. Inspection and if necessary replacement of worn **peristaltic pump** components (especially the internal hose)
2. Inspection of the **detergent intake** pipes and replacement if necessary
3. Inspection of the **door gasket** and replacement if necessary
4. Inspection of the **drying system filters** (relative and absolute) and replacement if necessary
5. Inspection and if necessary cleaning/replacement of the **filters** (water intake filters on the filler pipes, detergent filters on the suction systems)
6. Checking of the intake **water hardness setting** (the user must first have the intake water analysed to allow the correct settings to be made)
7. Checking of the **detergent dosage settings**.
8. Inspection of the **steam condenser** assembly
 - a. Inspection of the **nozzles**, check that water is flowing correctly
 - b. Inspection of the connection, water intake, water drainage and level pressure switch connection hoses
9. **Performance of a complete operating cycle**, including the drying phase, to check for any leaks or malfunctions.



WARNING

SMEG declines all responsibility in the event of appliance malfunctions, or injury or damage, arising from the failure to comply with the above recommendations.

10 TECHNICAL DATA

The documentation includes a description of the rated characteristics of the electrical connection, all the intake and outlet connections and the ambient conditions for which the appliance was designed (ref. point 5.4.2 of IEC61010-1:2001). For “non electric” intakes, the pressure and flow ranges are provided (ref. point 5.4.2 IEC61010-2-040:2005).

WATER SUPPLY REQUIRED			
Type of water required	1 - Mains (required)	2 - Demineralised (recommended)	3 - hot
Pressure	200 kPa – 500 kPa (2 bar – 5 bar)		
Type of threaded connection	3 / 4 ”		
Required flow [min - max]	2 - 12 litres / min		
Mains water hardness [max] <small>(for models equipped with a steam condenser, the use of purified water with hardness below 15°F is recommended).</small>	42 °f	-	42 °f
Fe2+ / Fe3+ [max]	0.5 ppm		
water temperature [max]	35°C	35°C	35°C
Demineralised water conductivity [max]	-	30µS/cm	-
Limits depending on microbial contamination	Minimum microbiological quality required: “potable water” (ref. limits in accordance with Italian Legislative Decree 31/2001)		
DIMENSIONS and WEIGHT			
Series	GW2145	GW1160, GW4060	GW4190
Standard External [Height x Width x Depth]	850 x 450 x 620	850 x 600 x 600	850 x 900 x 600
Standard External - aquastop accessory [Height x Width x Depth] - overall height is increased by +7mm	-	857 x 600 x 600	857 x 900 x 600
Standard External - PAD1 accessory [Height x Width x Depth] - overall depth is increased by +110mm / +150mm (depending on the model of PAD installed).	850 x 450 x 730 / 770	850 x 600 x 710 / 760	850 x 900 x 710 / 760
Standard External - peristaltic pump P3 accessory [Height x Width x Depth] – Only for GW1160, GW4060: overall depth is increased by +60mm	-	850 x 600 x 660	850 x 900 x 600
Net weight	56 kg	72 kg	105 kg
Weight of appliance + packaging	66 kg	90 kg	120 kg
MATERIALS			
Washing chamber	AISI 316L		
External surface	AISI 304		

ELECTRICITY SUPPLY			
MODELS - comments	SUFFIX TO THE MODEL NAME	ELECTRICITY SUPPLY	INTEGRAL PROTECTIVE DEVICES
GW2145	-	230V ~ / 50Hz / 14A / 3300W	FUSES 10.3x38 16 A ON POWER SUPPLY [L1, N]
GW1160, GW4060, GW4190 Max voltage/power Possible electrical versions: The names of the models have different suffixes to identify the power supply rating. E.g. If the model does not have a suffix, e.g. "GW4060" it will be three-phase, with neutral, 400V between phases. . Suffix "-3": "GW4060-3" indicates three-phase model, 230V between phases, no neutral. .Suffix "-1": "GW4060-1" indicates single-phase model. . The numerical suffix "6" identifies the 60 Hz version of 230V models. . Suffix "U": North American models, 208V convertible (* default connection)	standard	400V 3N~ / 50Hz / 12A / 7000W	FUSES 10.3x38 16 A ON POWER SUPPLY [L1, L2, L3]
	-6	400V 3N~ / 60Hz / 12A / 7000W	FUSES 10.3x38 16 A ON POWER SUPPLY [L1, L2, L3]
	-3	230V 3~ / 50Hz / 19A / 7000W	FUSES 10.3x38 20 A ON POWER SUPPLY [L1, L2, L3]
	-36	230V 3~ / 60Hz / 19A / 7000W	FUSES 10.3x38 20 A ON POWER SUPPLY [L1, L2, L3]
	-1	230V 1N~ / 50Hz / 12A / 2800W	FUSES 10.3x38 16 A ON POWER SUPPLY [L1, N]
	-16	230V 1N~ / 60Hz / 12A / 2800W	FUSES 10.3x38 16 A ON POWER SUPPLY [L1, N]
	-U	208V 3~ / 60Hz / 17A / 6000W *	FUSES 10.3x38 30 A ON POWER SUPPLY [L1, L2, L3]
		208V ~ / 60Hz / 29A / 6000W	FUSES 10.3x38 30 A ON POWER SUPPLY [L1, L2]
OTHER DATA			
DATA COMMUNICATION - PRINTER CONNECTION	RS-232 PORT		
MAX. NOISE LEVEL	50dB		
AMBIENT CONDITIONS			
USE	Indoor		
ALTITUDE	Up to 1000m		
AMBIENT TEMPERATURE	From 5°C to 40°C		
MAX RELATIVE HUMIDITY	80% for temperatures up to 31°C with linear decrease to 50% at the temperature of 40°C		
INSTALLATION CATEGORY (OR SURGE CATEGORY or <i>INSTALLATION CATEGORY</i>)	II		
ELECTRICAL INSULATION CLASS (ref. IEC 61140)	I		
DEGREE OF POLLUTION	2 (ref. 61010-1, point 3.6.6.2)		
STANDARDS AND CLASSIFICATION			
CONFORMITY - North American models	IEC61010-1(ed.2), IEC61010-2-040(ed.1), CAN/CSA C22.2 No. 61010-1-04 / UL 61010-1, 2nd Edition		
CONFORMITY - Other models	- European directive 2006/95/CE [Safety]: European Standard EN61010-1:2001, EN61010-2-040:2005 - European directive 2004/108/CE [Emc] : European Standard EN61326:2006 - European directive 2011/65/CE (RoHS 2) : European Standard EN 50581 :2012 - Other: European Standard EN 61770:2009		


11 PRODUCT DIMENSIONS - MEASUREMENTS IN mm

Preliminary note

The dimensions given below are nominal ones; when providing a housing space **always consider at least + 5mm (0.2 inch) per dimension** to create a minimum space necessary, for the extraction of the machine from the housing, for maintenance operations.

11.1 GW2145

The appliance can be requested with top for built-in installation for special requirements. In this case, the product is 20 mm shorter than the stated height.



CONFIGURAZIONE	A	B	C	D	E	F
STANDARD	620	600	30	450	850	> 10
STANDARD + PAD1	730 / 770	600	30	450	850	>10

11.2 GW1160, GW4060 SERIES

The appliance can be requested with top for built-in installation for special requirements. In this case, the product is 30 mm shorter than the stated height.





fig. 48 – Product overall dimension diagram

VERSION	A	B	C	D	E	F
STANDARD	600	600	30	600	850	> 10
STANDARD + AQUASTOP	600	600	30	600	857	> 10
STANDARD + PAD1	710 / 750	600	30	600	850	>10
STANDARD + P3	660	600	30	600	850	> 10

11.3 GW4190 SERIES

The appliance can be requested with top for built-in installation for special requirements. In this case, the product is 30 mm shorter than the stated height.



N.B.: on the GW4190 series, the optional peristaltic pump P3 does not cause an increase in depth.

VERSION	A	B	C	D	E	F
STANDARD	600	600	30	900	850	> 10
STANDARD + ACQUASTOP	600	600	30	900	857	> 10
STANDARD + PAD1	710 / 750	600	30	900	850	>10
STANDARD + P3	600	600	30	900	850	> 10

12 POSITIONING THE APPLIANCE

IMPORTANT

The appliance must be positioned with its back near a wall (with a gap of at least 10 mm) and must be installed by an authorised SMEG engineer.

The installer engineer is responsible for the correct operation of the appliance after installation and must also provide the user with all information required for its correct use.

During installation, the protective film must be removed from the steel outer surfaces.

The set of installation accessories (seals and hose ties) is inside the chamber.

The appliance can be installed with its sides flush against those of the adjacent cabinets, taking care not to obstruct the steam vent on the rear. The wall at the back should be in masonry or another waterproof material.

Also take care that the heat does not reach any electrical circuits or sockets behind the appliance.

The appliance has water intake and drain hoses which can be positioned leading to the right or left, depending on installation requirements.

If ordered in the appropriate version, the appliance can be installed built under a worktop; installation must be carried out by skilled staff.

LEVELLING

After positioning the appliance, screw the feet down or unscrew them to adjust its height and level with the aid of a spirit level so that it will be perfectly horizontal (max angle tolerance allowed: 0.5 °, corresponding to a maximum permissible difference in height between opposite corners of the appliance of about 5 mm).

Correct levelling will ensure that the appliance operates correctly.

WARNING

Any adjustments, maintenance, etc. must be performed with the device switched off and disconnected from power sources.


LIFTING AND TRANSPORT

Before shipment from the factory, the base of the appliance is fixed to a pallet used for both lifting and transport. When moving the appliance around, a forklift truck or pallet truck must be used.

CONNECTIONS ON THE BACK OF THE APPLIANCE

ID.	DESCRIPTION
A	Demi water intake hose
B	Mains water intake hose
C	Electricity supply cable
D	Steam condenser exhaust hose
E	Appliance drain hose
F	Load hot water pipe (for models with suffix "W", only es. GW4060WS1)
P1	Peristaltic pump P1 intake hose and P1 can detergent level sensor cable
P2	Peristaltic pump P2 intake hose and P2 can detergent level sensor cable
H	RS-232 port for PC or printer connection Cables with connectors for correct connection are provided with the following accessories: 1. "Smeg WD-PRINTE" printer 2. PC connection cable, for communication of data by means of Smeg TRACELOG software.

GW2145



GW1160, GW4060 Series




fig. 49 – Diagram of connections on rear of appliance.

GW4190 Series

On the GW4190 series, the detergent intake connections “P1” and “P2” are inside the side cabinet and are accessed by opening the hatch. The connections are marked for identification with P1 (detergent) and P2 (acid neutraliser).




fig. 50 – Rear of GW4190: Simplified drawing, connections on rear of the device.




fig. 51 – GW4190 series: the detergent intake hose connections are accessible in the side cabinet. P3 is always optional.

INSTALLATION ON STAND

If the device is installed on a stand (e.g. for GW1160, GW4060: Smeg product “B6040”) with having height (base height) “HB”, all the measurements of height above the supporting surface provided in this document must be increased by the said height “HB”.

$$HMi = Hi + HB$$

Where:

Hi = normal height stated in the document



HB = height of appliance stand

HMi = height recalculated for appliance installed on stand



fig. 52 – Smeg B6040 Stand.


13 ELECTRICAL SYSTEM REQUIREMENTS

	<p>WARNING <i>The electrical system to which the appliance is connected must comply with the relevant regulations in the place of installation.</i></p> <p><i>All electrical checks and system installation must be carried out to the proper standard by skilled staff with proven experience, qualified to work on electrical systems.</i></p> <p><i>The skilled staff are responsible for ensuring that the ground connection is in good working order.</i></p> <p><i>These devices are fitted with power surge protection devices: however, where possible a specific overload cutout suitably rated to protect the equipment should be installed for each device.</i></p> <p><i>For any spare parts: use only genuine parts, ordering them from your service centre.</i></p> <p><i>Please refer to the “TECHNICAL CHARACTERISTICS” section for all specifications regarding the device's voltage requirements and power drawdown.</i></p>
	<p>POWER CABLE REPLACEMENT <i>If damaged power cables have to be replaced, cables with the same characteristics and gauge as the cables mounted in the factory must be used (the cables fitted are marked with their types and characteristics). This replacement may only be made by authorised technical staff. Use genuine spare parts.</i></p>

13.1 GW2145 - POWER SUPPLY CABLE AND DISCONNECTION DEVICE

Characteristics of the power supply cable supplied with the appliance:

- H05V2V2-F 3 x 2.5 mm², (for 230V single-phase version), Schuko plug.

	<p>WARNING <i>The device is supplied complete with power cable with plug for connection to the electrical mains: the plug constitutes its disconnection device and must be easily accessible to the user.</i></p>
---	--

13.2 GW1160, GW4060, GW4190 - POWER SUPPLY CABLE AND DISCONNECTION DEVICE


Characteristics of the power supply cable supplied with the appliance:

- FROR 5 x 2.5 mm² (5G2.5), 450/750 V, IMQ mark (three-phase 440V version)
- FROR 4 x 2.5 mm² (4G2.5), 450/750 V, IMQ mark (three-phase 230V version, no neutral)
- FROR 3 x 2.5 mm² (3G2.5), 450/750 V, IMQ mark (single-phase 230V version)
- SJT 4 x 10 AWG, 300V (208V version, North American models).

The device is supplied without a plug, with a cable with insulated wire terminals.

The device's **electrical connection** must be made **by means of an industrial plug**.

The plug must be supplied and fitted by the user. The plug must be of suitable size for the electrical cable and the electrical rating of the device.

	<p>A CIRCUIT-BREAKER must be installed for each appliance.</p> <p>Circuit-breaker characteristics</p> <ol style="list-style-type: none"> <i>Omnipolar: must break all live conductors;</i> <i>Easily accessible to the user;</i> <i>Easily operated (no tool must be required);</i> <i>Located in close proximity to the appliance;</i> <i>Clearly marked as the appliance circuit-breaker.</i> <p><i>The appliance is fitted with overcurrent protection devices.</i></p> <p><i>However, a protective device specifically for the appliance (e.g. magnetothermal breaker or fuse on every phase, suitably rated for the electrical data stated above) should be installed in the room's electrical panel.</i></p>
---	---

NORTH AMERICAN MODELS - 208V

*The device's default connection is:
208V 3~ / 60Hz / 17A / 6000W*

***The appliance can be electrically converted into the single-phase version
208V ~ / 60Hz / 29A / 6000W.***

The electrical conversion must be carried out by authorised technical staff.

Reference document: device electrical wiring diagram.

Any components must always be replaced with genuine Smeg spare parts.

14 WATER CONNECTION FITTINGS

Key to abbreviations used for water connections.

CODE	FILLING/DISCHARGE	WATER TYPE
cw	INTAKE	cold water hose - cw (<i>cold water</i>)
dw	INTAKE	pressurised demineralised water hose - dw (<i>demi water</i>)
d	DRAIN	appliance drain hose - d (<i>drain</i>)
cd	DRAIN	steam condenser drain hose - cd (<i>condenser drain</i>)

14.1 WATER INTAKE

14.1.1 WATER INTAKE CONNECTION

The appliance has **two intake hoses, for mains and demineralised water.**

The hoses are designed to be connected to taps with 3/4" gas threaded attachment.

Use the filter provided, "A" in the illustration, when connecting the end of the intake hose.




fig. 53 – Fit the filters provided when connecting the intake hose.

MODELS	CONNECTIONS NUMBER	CONNECTION TYPE		
		cw	dw	hw
GW2145, GW1160, GW4060, GW4190	2	1	1	---
GW1160W, GW4060W, GW4190W	3	1	1	1

14.1.2 WATER SUPPLY TAP POSITION

Water supply taps must be positioned close to the appliance, in a position accessible to the user. With reference to the diagram below, the recommended measurements are:

$$l < 50 \text{ cm} - HC < 100 \text{ cm}$$

There are no particular height requirements for the position of the water supply taps, but bear in mind that the hoses supplied are about 2 metres long.




fig. 54 – Diagram. The water supply connections may be on the right or left of the appliance, always considering the maximum distance "l".

Note to prevent the risk of clogging or damage: if the water pipe is new or has been unused for a long time, before connecting the water supply make sure that the water is clear and free of impurities.

14.1.3 WATER SUPPLY SYSTEM REQUIREMENTS

Preliminary checks

1. **The water supply pressure must be within the required limits: min. 2 bar - max 5 bar.**
2. **The hot water temperature must not exceed 50°C** - Higher temperatures might reduce the efficiency of the integral water softener by damaging the resins it contains.
3. **The water supply taps must be accessible.**

WATER SUPPLY CHARACTERISTICS

Cw: MAINS WATER - essential:

The appliance must be connected to "potable" mains water having hardness 42°f max. and with total dissolved iron content, Fe^{2+} e Fe^{3+} , not exceeding 0.5 ppm.

For models equipped with a steam condenser, the use of purified water, preferably with hardness below 20°f, is recommended.

N.B.:

If the water supply contains more than 0.5 ppm of iron (Fe^{2+}/Fe^{3+}) and/or the hardness of the water supply exceeds 42°f (French degrees), the water must be pretreated by installing an iron removal and/or softening system upstream of the appliance.



WARNING - MAINS WATER HARDNESS

User's responsibility: the hardness of the mains water must be checked prior to installation.

The Smeg engineer will not be able to set the parameters of the appliance's integral water softener unless the hardness of the intake water is known.

Regular checks must be made (e.g. every 6 months/year) on the intake water to ensure that the initial values are maintained, to allow the appliance parameters to be corrected if the characteristics change.

Dw: DEMINERALISED WATER - optional connection but strongly recommended

If available, the use of demineralised water (**conductivity < 30 µS**) is recommended to ensure that washing is optimal from the chemical point of view, through better removal of the chemical residues in the water supply; however, dirt residues will still be effectively removed even if ordinary mains water is used.



WARNING - NO DEMI WATER

If demineralised water is not available, do not connect the relative hose to hot and/or cold water intakes. Leave the "demineralised water" hose disconnected.

The appliance's SETTINGS must be corrected accordingly by the Authorised Service Centre.

14.1.4 PAD ACCESSORY FOR UNPRESSURISED DEMI WATER

Connection to a non-pressurised demineralised water line (e.g. a tank operating by gravity) is only possible if a special demineralised water pump is installed. Contact Smeg for details about accessories and optionals.

WARNING

- WHEN INSTALLING THE "PAD" ACCESSORY, MAKE SURE THAT:

The supply pressure of the PAD is less than 1 bar and greater than 0.1 bar.

- "PAD" IN COMBINATION WITH AN UNPRESSURISED TANK

With reference to the illustration below, the tank must be at a height L such that:

$L > 100 \text{ cm}$



The height is measured between the surface on which the appliance stands and the bottom of the tank.

This height ensures that water is supplied to the PAD at a pressure of about 0.1 bar.

SMEG declines all responsibility for any damage or injury caused by incorrect installation of the pump by unauthorised technical staff.




fig. 55 – Unpressurised tank, height above appliance installation surface.

14.2 WATER DRAIN

Appliance drain hoses:

Rubber fitting for hose connector diameter 21mm (1/2").

WATER DRAIN CONNECTION

The appliance is equipped with one or more drain hoses, depending on the model.

The drain hoses are identified as:

- **d** - appliance drain hose – **d** (*drain*)
- **cd** - steam condenser drain hose – **cd** (*condenser drain*)

The diagram and summary table are provided below.




fig. 56 – Diagram. The water drain connections may be on the right or left of the appliance, always considering the maximum distance “l_max”.

MODELS	No. OF DRAIN HOSES	Distance "l_max" Maximum Distance between drain and side of appliance [cm]	DRAIN TYPE AND HEIGHT "H" ABOVE APPLIANCE INSTALLATION SURFACE [cm]		
			Type	Hmin	Hmax
GW2145 GW1160 xxx GW4060 xxx GW4090 xxx	1	50	d	65	80
GW1160 Cxx GW4060 Cxx GW4090 Cxx	2	50	d	65	80
			cd	65	80

Note about product "families"

Suffix "C" indicates "equipped with steam condenser."

Any additional suffixes ".." may vary depending on the appliance version.




fig. 57 – Drain connection. $65\text{ cm} < H < 80\text{ cm}$, $X < 20\text{ cm}$. Suffix "C" on the name of the model indicates "equipped with steam condenser", in this case n.ro 2 drain pipes are provided.



WARNING

The drain must comply with international regulations: our company accepts no liability for pollution caused by improper use of the appliance.



fig. 58 – Adapter with connector for 1/2" hose.

General guidelines for installing the drain

The use of a drain with syphon trap is recommended. During installation, the following precautions should be complied with:

- Since the drain water temperature is 95°C , **the end of the drain hose must be permanently connected to the hose connector**, using the **ties** provided.
- There must be no sharp bends in the drain hose which might cause restrictions.
- The difference in height between the end of the drain hose and the surface on which the appliance is installed **must comply with the stated specifications**.
- The end of the hose must never be submerged in water.
- The **inside diameter** of the drain pipe to which the hose is connected must be **at least 40 mm**.
- A drain pipe **50 mm in diameter** is recommended.
- No extensions must be added to the drain hose supplied with the appliance. Any extensions might cause problems of flow-back into the chamber.

15 DATA CONNECTION

The device has a data communication port which can be used for one of the following functions:

1. **Connection** to an optional **external “Smeg WD-PRINTE” printer**. The printer produces a brief report with the key cycle data.
2. **PC connection** for digital data storage and interaction with the appliance by means of the “Smeg TRACELOG” program.

Reference should be made to the specific instructions of the above accessories for any connection details.

GW2145, GW1160, GW4060 Series




fig. 59 – Position of the RS-232 serial port on the device

GW4190 Series

The following are located above the panel that carries the main switch and the printer (if installed):

- **RS232 port**: for data transmission – connection by means of Smeg TRACELOG software.
- **“Data output” switch**: used to activate the panel printer or data transmission on the RS232 port. These two devices can never operate simultaneously. The data output switch activates the printer when turned to “II” or the RS-232 port when turned to “I”.




fig. 60 – GW4190: Position of the RS-232 port and the data output switch which activates the printer or the RS-232 port.

The RS-232 cable used for the connection must have the following characteristics:

CONNECTION TYPE	CONNECTOR
Crossover (Null Modem)	Female/Female



fig. 61 – RS-232 female connector.

Not all PC models have serial ports.

Smeg offers a USB/SERIAL converter accessory; contact your dealer for technical details or offers, or write to instruments@smeg.it



fig. 62 – Smeg USB/serial converter.

Smeg S.p.A.

Instruments Division

Via Leonardo da Vinci, 4 – 42016 Guastalla (RE) Italy

Tel. +39 0522 8211 – Fax +39 0522 821 592

E-Mail: instruments@smeg.it – service.instruments@smeg.it

www.smeg-instruments.com