PIPETEADOR ELECTRÓNICO DIGIPETTE FLOWASSIST
ELECTRONIC PIPETTE CONTROLLER DIGIPETTE FLOWASSIST
PIPETTEUR ÉLECTRONIQUE DIGIPETTE FLOWASSIST

REF. - CODE - RÉF. MBD006

DIGIPETTE

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INDEX

1. INTENDED USE OF THE INSTRUMENT	10
2. SAFETY INSTRUCTIONS.....	10
3. OVERVIEW	11
4. CHARACTERISTICS.....	11
5. CONFIGURATION	12
5.1. BATTERY.....	12
5.2. CHARGING THE UNIT.....	12
5.3. STAND	12
5.4. PIPETTE PLACEMENT.....	12
5.5. SELECTION OF THE PIPETTE	12
6. OPERATION	13
6.1. FILLING AND DISPENSING	13
6.2. FILLING/DISPENSING SPEED	13
6.3. MODES OF DISPENSING.....	13
7. MAINTENANCE/CLEANING AND AUTOCLAVING	13
7.1. DISASSEMBLY OF THE CONICAL ADAPTER.....	13
7.2. ASSEMBLY OF THE CONICAL ADAPTER.....	14
8. RECOMMENDATIONS FOR USE	14
9. TROUBLESHOOTING	15

1. INTENDED USE OF THE INSTRUMENT

The DIGIPETTE FlowAssist electronic pipette controller is an electronic instrument for filling graduated and one-mark pipettes, in glass or plastic, with volumes from 0.1 to 100 ml. Battery operated, it guarantees precision, accuracy and reliability.

When used correctly, the liquid only meets the pipette.

2. SAFETY INSTRUCTIONS

The instrument can be used with hazardous materials, but the user is responsible for following safety rules and checking applicable regulations before use, as this manual does not cover all possible risks.

- Every user must read and understand this operating manual before using the instrument.
- Follow the general risk instructions and safety regulations.
- Comply with all specifications provided by the reagent manufacturers.
- Caution! Do not operate the instrument in partially explosive atmospheres and do not pipette highly flammable media.
- Always carry out work in a manner that does not endanger your safety or the safety of others. Avoid splashing. Use only suitable containers.
- Never use excessive force.
- Use only original accessories and spare parts. Do not make any technical modifications. Do not disassemble the instrument beyond what is indicated in the operating manual.
- Always check that the instrument is in good working order before use. Contact the manufacturer if necessary.
- Voltage/Electrical Warning! Only use the original power adapter to charge the lithium-ion battery.
- Do not use batteries from other manufacturers.
- Protect the power adapter from moisture and use it only with that instrument.
- Discontinue use of the power adapter if the connection cable is damaged.
- Only authorised service technicians are allowed to carry out repairs and internal maintenance of the instrument.
- During charging, ensure that the battery charger is set to the correct voltage according to the line voltage.

Note: The instrument cannot be used while charging.

- Do not allow the battery to discharge completely
- Make sure that the micropipette is properly secured to avoid accidental dropping during liquid handling.
- Be careful during suction to prevent liquid from entering the cone adapter assembly.

Note: Do not overfill the pipette and remove or replace the filter if liquid accidentally enters.

- Caution! Avoid liquids and solvents that emit vapours. These may cause dripping or their vapours to damage the internal components of the instrument.
- Store the instrument in the stand provided when not in use.
- During charging the device must be switched off.

3. OVERVIEW

4. CHARACTERISTICS

- Digipette FlowAssist is compatible with all standard pipettes between 1 and 100 ml.
- A hydrophobic 0.2-micron PTFE membrane filter integrated into the assembly prevents the accidental ingress of liquids into the instrument.
- The user can control the speed of suction and dispensing by adjusting the pressure applied to the plunger buttons.
- It uses a solenoid valve for precise liquid control.
- The instrument offers gravity and blow-out dispensing options.
- It can fill a 25 ml pipette in 3 seconds.
- It offers direct charging via the micro-USB port with Type-C charging cable.
- Allows up to 8 hours of continuous operation with 3.4 V, Li-ion 1400 mAh battery.
- It has fully autoclavable pipette adapters and filter holders.
- LED indicators for low battery (red) and charging (yellow). The yellow LED turns off when the battery is fully charged.

5. CONFIGURATION

5.1. BATTERY

- To insert the battery, press and slide the top cover while pushing it slightly.
- Carefully remove the cover and insert the battery.
- Replace the top cover on the housing, making sure that it fits into the slots on the display.

5.2. CHARGING THE UNIT

- Once the battery is installed, a welcome screen will appear.
Note: To use the instrument at its maximum performance, it is recommended that the battery be fully charged before first use.
- The yellow LED will light up when the device is charging.
- The red LED will light up when the battery is low.
- The recommended output voltage is 5V, 3A.

5.3. STAND

- The pipette controller is supplied with a desktop stand that can be used for storage or when not in use.

5.4. PIPETTE PLACEMENT

- To fit the pipette, hold the pipette controller by the cone and press the pipette into the fitting cone assembly.

Note: Never use force to fix the pipette, especially for glass or thin-walled pipettes.

Capillary phenomenon

During aspiration and dispensing with glass serological pipettes, the liquid may continue to move even after the plungers are released, due to capillary action. This phenomenon is due to surface tension and adhesion between the liquid and the inner walls of the glass. This does not usually occur with plastic serological pipettes, as their slightly larger opening reduces this effect.

5.5. SELECTION OF THE PIPETTE

- To optimise the aspiration and dispensing speed, click on the pipette selection and select the capacity of the pipette being used.

6. OPERATION

6.1. FILLING AND DISPENSING

To fill the pipette, press the aspiration button to draw the liquid into the pipette and the dispensing button to expel the liquid.

Note: Do not overfill the pipette.

Remove or replace the filter if liquid is accidentally aspirated.

6.2. FILLING/DISPENSING SPEED

The user can control the speed of aspiration and dispensing by adjusting the pressure applied to the plunger buttons.

Pipette capacity-Maximum speed:

- ✓ 1 ml-0.3 ml/s
- ✓ 2 ml-0.6 ml/s
- ✓ 5 ml-1.5 ml/s
- ✓ 10 ml-3 ml/s
- ✓ 25 ml-8 ml/s
- ✓ 50 ml-8ml/s
- ✓ 100 ml-8ml/s

6.3. MODES OF DISPENSING

The instrument offers gravity and blow-out dispensing options.

■ Gravity dispensing: Lightly press the dispense button. Gravity mode dispenses liquid by the action of gravity only.

■ Blow-out dispensing: Press the dispensing button deeper to activate the motor. Once the motor is activated, the dispensing speed can be adjusted by varying the pressure applied to the button.

7. MAINTENANCE/CLEANING AND AUTOCLAVING

If the instrument is used correctly, no maintenance or cleaning is required. However, if necessary, the conical adapter assembly can be autoclaved.

For cleaning, it is recommended to clean the outer surface with 70% ethanol.

Note: Do not use alkaline solutions as they may cause discolouration of the product material.

7.1. DISASSEMBLY OF THE CONICAL ADAPTER

- Hold the instrument at the top and turn the cone clockwise (Fig. 1).


(Fig. 1).

■ The bushing rubber will appear. Gently pull it out together with the filter to remove them (Fig. 2).

■ This bushing assembly can now be autoclaved at 121°C and 1 bar for 15-20 min.

Note: The instrument is equipped with a 0.2-micron membrane filter, which can be autoclaved 5-7 times. It is recommended to clean the filter only when the pipette is overfilled or when there is difficulty in aspirating the liquid.

(Fig. 2).

7.2. ASSEMBLY OF THE CONICAL ADAPTER

■ Push the filter from the narrow end into the internal grooves of the housing (Fig. 3).

■ Align the rubber bushing and insert it into the filter (Fig. 4).

■ Finally, complete the assembly by turning the bushing cover.

(Fig. 3).

(Fig. 4).

8. RECOMMENDATIONS FOR USE

	GOOD PRACTICE	PRACTICES TO AVOID
1	When attaching the pipette, hold it as close to the upper end as possible.	Do not use a damaged charging adapter or cable.
2	Always hold the instrument with the pipette inserted vertically, with the tip pointing downwards.	Do not use the instrument while it is charging or if the battery is too low.
3	Aspirate a little more liquid than the required volume and slowly press the lower pipetting button to dispense the liquid until the meniscus is adjusted exactly to the desired volume.	Do not overfill the pipette, as this may clog the filter and affect the filling process.
4	Always use the storage stand to store	The instrument must never be used or charged in potentially explosive atmospheres.

9. TROUBLESHOOTING

Problem	Cause	Solution
Pipette is loose	The rubber may be dirty or clogged.	Clean the rubber by autoclaving. Replace it if it is broken.
Liquid drops from the pipette	Pipette is not positioned correctly	Replace the pipette
	Fluid leakage at the conical adapter	Correct positioning of the conical adapter
	Fluid leaks from internal connections	Send the instrument to the supplier
Pump does not work	Battery is discharged	Charge the battery
	The connection is broken	Send the instrument to the supplier
	The button is not well assembled	Check the button, if the problem persists send the instrument to the supplier.
Charging LED does not flash when power is supplied	Incorrect power adapter	Use the original power adapter
	Battery is not connected correctly	Reinsert the battery correctly and charge it.
	Damaged battery/power adapter/charging cable	Replace battery/power adapter/charging cable
	Damaged LED	Send the instrument to the supplier
No display	Damaged screen or loose internal component	Send the instrument to the supplier
Liquid does not flow into the pipette	The filter is clogged due to chemical ingress.	Autoclave or change filter

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