DISCO DE HARTL COMPLETO
ROTABLE OPTICAL DISC HARTELEY

REF :QLG005

INDEX OF LANGUAGES- ÍNDICE DE IDIOMAS

| | |
|---------------|---|
| SPANISH | 2 |
| ENGLISH..... | 4 |

DISCO DE HARTL COMPLETO

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ROTABLE OPTICAL DISC HARTELEY

The instrument is useful to learn optical concepts including laws of reflection and laws of refraction that can be easily verified and demonstrated. It is composed by a steel disc white painted diameter 30cm graduated with steps of 1° until 90, 0-90 and 0-90, it has four marked areas. In the middle of the disc there are marks to help you to locate the prisms and they have blocks with powerful magnet attached. It has a stable cast iron base. Also there is a laser that allow you to adjust 1, 3 or 5 beam rays. Including 6 optic elements:

1. Acrylic prisms (refraction): semicircular, biconvex, biconcave, right angle prism, triangle equilateral and trapezoidal prism.
2. Mirrors (reflection): plane and semicircular.

The base of all this prisms is ground glass.

The main disc and lamp house with 12V, 24W bulb can be rotated and held in any position desired.

With routable disc with single, double, triple, quadruple slits with lens for divergent, parallel or convergent beam is provided.

Working instructions

Place the plane mirror on 90-270 line. It will be fixed with the help of fixing blocks having powerful magnet attached. Bring lamp house 24w bulb zero of optical disc and adjustable for parallel beam and single slit. Now check the ray to find that mirror is exactly perpendicular to steel disc.

- a) Now slit beam at 30° the reflected beam will be reflected will also be reflected by 30°.


- b) Now fix semicircular plate of glass with its plane facing PQ, fixing is done as above by powerful strong magnet system and 0-180 line passing normally through center.

Adjust light source and parallel beam. The incident ray AO falls on plane surface of the plate and suffers refraction on entering it. The refracted ray travels along the radius of the plate and does not further bend on leaving. OB represents the refracted ray.

By turning the main disc and lamp house slit system measure the values of angle of incidence and corresponding angles of red refraction the result will.

$$\frac{\sin \gamma}{\sin \gamma} \text{ is constant hence}$$

Hence laws of refraction are verified

- c) Place right-angled prism against single parallel beam. The beam will be totally internally reflected. Similarly it can be performed with double slit. The evident rays will be as shown. Functions can

be shown biconvex and biconcave with help of quadruple slits and path rays will be visible and demonstrated.

Refraction can be shown of trapezoidal acrylic by throwing proper light with the help of lamp house. Slits and divergent and convergent rays



Illustration by double slit parallel beam

Illuminated by quadruple parallel beam.