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1. INTENDED USE

This device can be widely used in pathology study and diagnosis, scientific research and experiment in hospitals, medical colleges, animals and plants R&D institution and food inspection department, for embedding the biological tissue with paraffin.

2. CAUTIONS

1. Before starting the machine, wax pan must have the right medical grade paraffin.
2. Before heating, please pay attention to the preset temperature, it should be within the range of melting point of the paraffin.
3. Before melting totally, please don't press "paraffin drain switch" (manual switch or foot switch), so as not to damage the paraffin spout.
4. Forbid water or other liquid flowing into machine case and the paraffin bath.
5. In the course of heating, don't open the cover of the bath.
6. While embedding, prevent the impurity from entering the paraffin bath, in case that the paraffin spout is filled.


2. Back to MENU interface

Image #2

3. Choose AUTOMATICALLY to access the interface shown in the following image; stay at this interface and don't press any button if you want the automatic operation; the device will be in stand-by state until reaching the OPEN time (automatic working time). During stand-by state, the machine does not work

Image #3

Note: Please don't press OK key at AUTOMATICALLY interface (above image); if you press OK key, the machine will jump to manual operation state and starts working directly (the elements will be heated according to the settings). If you want to back to automatic working state, please press RE-SET key to view the MENU interface and access AUTOMATICALLY interface again.

6. MAINTENANCE

1. Keep the equipment on a flat surface, away from a fire source.
2. The device must have good grounding and good ventilation.
3. After using, remove the scraps and clean the apparatus.


7. POSSIBLE PROBLEMS

Phenomenon	Reason
Incessant paraffin flow	Damaged electromagnetic valve
Temperature has not shown	7805 shorts out
Temperature of the paraffin bath has not risen	Temperature relay of the paraffin bath open circuit
Temperature of the working surface has not risen	Temperature relay of the working surface open circuit
Temperature shows unstable	The device is not well grounded
Temperature do not rise to the preset temperature	Damaged temperature sensor

8. TIPS FOR AUTOMATIC OPERATION

1. First, set all data, such as temperature of different elements, system clock time, etc., and then access the TIME interface to set the automatic power on and automatic power off time and date.


Image #1

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Max. temp. set point of wax pan	99°C, step of $\pm 1^\circ\text{C}$
Max. temp. set point of right and left boxes	99°C, step of $\pm 1^\circ\text{C}$
Max. temp. set point of wax spout	99°C, step of $\pm 1^\circ\text{C}$
Max. temp. set point of working platform	99°C, step of $\pm 1^\circ\text{C}$
Temperature range of cold spot	-5°C~-20°C, step of $\pm 1^\circ\text{C}$
Wax pan capacity	5L
Activation of the paraffin flow	-foot switch -metal proximity sensor (approximately 4mm)
Operating time programming	start and end time of operation can be set within a 24-hour period.
Nominal supply voltage	AC220V $\pm 10\%$
Nominal frequency	50 ± 1 Hz

4. STRUCTURE OF THE DEVICE

5. OPERATION

Turn on the power switch; the LCD display shows the following screen:

WELCOME EMBEDDING SYSTEM SYSTEM TIME XX:XX:XX

Press "OK" key in control panel to access the MENU.

MENU	Date XX:XX:XX
1.TIME	4.TEMPERATURE
2.MANUAL	5.AUTOMATICALLY
3.CLOCK	

Press the ↑, ↓, ← and → to move the cursor.

A. CLOCK

Firstly, set the system time same with local time.

Press the ↑, ↓, ← and → to move the cursor to "3.CLOCK". Press OK key to access SYSTEM CLOCK interface.

SYSTEM CLOCK		
DAY XX	MON XX	YEAR XX
SEC XX	MIU XX	HOUR XX
WEEK XX		

Press OK key to save settings and back to MENU.

Press RESET key to abort operation and back to MENU

B. TEMPERATURE

Press the ↑, ↓, ← and → to move the cursor to "4.TEMPERATURE".

Press OK key to access TEMPERATURE SETUP interface

TEMPERATURE SETUP (°C)			
WAX PAN	XX	PRES1	XX
WAX SPOUT	XX	PRES1	XX
TABLE	XX	COLD	XX

Set the required temperature for wax pan, wax spout, Pres1, Pres1, working platform and cold spot.

Press OK key to save settings and back to MENU.

Press RESET key to abort operation and back to MENU.

C. TIME

Press the ↑, ↓, ← and → to move the cursor to "1.TIME".

Press OK.

Set the required time and day of the week for auto power on or auto power off.

OPEN: mean the time to auto power on

EXIT: mean the time to auto power off

For the day of the week, "√" mean you choose this day of week and embedding center will auto power on or auto power off on chose day . Press ↑ to choose the day of the week, press ↓ to deselect the date of week, press → to next day of the week.

Press OK key to save settings and back to MENU.

Press RESET key to abort operation and back to MENU

TIME	DATE XX:XX:XX
OPEN XX:XX	EXIT XX:XX
MON <input checked="" type="checkbox"/>	TUE <input type="checkbox"/> WEN <input type="checkbox"/> SUN <input type="checkbox"/>
THU <input type="checkbox"/>	FRI <input type="checkbox"/> SAT <input type="checkbox"/>

D. AUTOMATICALLY

Press the ↑, ↓, ← and → to move the cursor to "5.AUTOMATICALLY".

Press OK.

The device will start running when the system time reach the OPEN time (auto power on) (set OPEN time at TIME interface).

WAITING	DATE XX:XX:XX
OPEN XX:XX :XX	
WEEK 1 2 3 4 5 6 7	

D. MANUAL

Press the ↑, ↓, ← and → to move the cursor to "2.MANUAL".

Press OK.

The device starts running immediately if into this MANUAL interface and stays at this interface (the device stops running if out of this MANUAL interface).

RUN (°C)	DATE XX:XX:XX
WAX PAN	XX↑ PRESS1 XX
WAX SPOUT	XX PRESS1 XX
TABLE	XX COLD XX

If "↑" appears means the temperature will keep rising