Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung

ECO 8 & ECO 25 Thermoreactors

F101A0127, F101A0125

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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1. Introduction

Thermo reaction can be carried out quickly and safely by placing the tubes into the thermoreactor heating block. The Start-Temp key offers the possibility to set 5 different work temperatures (70, 100, 120, 150 and 160 °C) whilst the Time key to set 4 different work times at the selected temperatures (30, 60, 120 or infinite). ECO 8 and ECO 25 offer the possibility to thermoreact respectively 8 and 25 tubes simultaneously (Ø 16 mm). ECO 8 model allows to hold also 1 test tube with Ø 22 mm.

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- ECO 8 or ECO 25
- Safety cover
- Power supply cord
- Instruction manual

2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

2.2 Start-up

Use the on-off switch on the front of the unit to turn the unit on (position “-”). When switched on the unit all the leds blink for few seconds. When the leds are off the unit is ready to be programmed at different work temperatures and times.

3. Operating controls

ON - OFF SWITCH	It turns the unit on and off. If the switch is in the “0” position the unit is off; “-” means it is on.
START-TEMP	When Start-Temperature key is pressed, a work cycle will start with the pre-set values and the corresponding leds will light: Temperature 160 °C and Time 30 min. The next pressing permits to select a different temperature: 160, 70, 100, 120 and 150 °C. It's possible to select different work temperature until the heating block has not reached the set temperature and the count-down is not working. After that the key doesn't work anymore. If the selected temperature is lower than the running temperature of the heating block, (when the Start-Temp key is pressed) the work cycle will not start and all temperature leds will light intermittently, signalling an error. To correct the error, select a temperature that is higher than the running temperature or wait until the heating block cools down to a lower temperature than the temperature to be programmed.
TIME	When the Time key is pressed it's possible to select the following different time: 120, 60, 30 minutes or ∞ (infinite). When a time is selected, the corresponding led will light. After switching on the instrument the Time key permits to set the time at the set temperature but not to start the reaction cycle. The run time will start when the heating block reaches the set temperature.
TEMP. LED	When a led lights it means that a work cycle is running at the selected temperature. If the leds are off a work cycle is not in progress.
TIME LED	When a pre-set time is selected the corresponding led will light. When the led blinks it means that the selected temperature has been reached and the count-down has started.

3.1 Heating block temperature verification

The temperature is controlled electronically which ensures stability without temperature oscillations around the set point. The Pt 100 probe provides high precision and does not require calibration because the internal software performs a self-calibration every time the instrument is turned on. Check the heating block temperature with VELP Scientifica thermometer.

3.2 Power interruptions

The unit will register a black out or a lack of power supply as a manual stop. The work cycle will have to be restarted.

3.3 Safety features

The unit is protected against overheating by a thermostat that turns off the power supply when the temperature exceeds the safety temperature. Temperature setting is continuously controlled by electronic equipment that stops thermoregulation if anomalous measured are sensed. In case, all the leds start blinking and an acoustical warning signal sounds.

3.4 End of work cycle

At the end of the count-down the work cycle will stop automatically. To stop in every moment the work cycle turn off the general switch. When the work session is completed, should the unit remain connected to the power supply, switch it off using the on-off switch on the frontal panel.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Code Model	F101A0127 ECO 8	F101A0125 ECO 25
Admitted power supply	115 V - 230 V / 50-60 Hz	115 V - 230 V / 50-60 Hz
Power	135 W	400 W
Dimensions (WxHxD)	135x95x230 mm	155x95x275 mm
Weight	2 Kg	3,6 Kg
Holes and diameter	8 (Ø16 mm) + 1 (Ø22 mm)	25 (Ø16 mm)
Construction material	Epoxy painted aluminum cast	
Thermoregulation	Pid microprocessor	
Selectable temperatures	70 - 100 - 120 - 150 - 160 °C	
Selectable times	30 - 60 - 120 min. - infinite	
Set temperature reached	Acoustic signal	
Count-down	Visual signal	
End of cycle	Acoustic signal	
Overtemperature	On the block at 180 °C	
Heating block temperature stability and homogeneity	± 0.5 °C	
Temperature precision	± 1.0 °C	
Time to reach 150 °C	12 min.	
Enviromental ambient temperature	5 - 40 °C (41 - 104°F)	
Storage temperature range	-10... + 60 °C (14 - 140°F)	
Max humidity	85%	
Pollution degree	2	

6. Accessories / Spare parts

CA0091666	Test tube for sample decomposition 22 mm external diameter, NS 19/26 cone with glass cap (ECO 8)
CA0091667	Condenser 200 mm type KS with 3 meters polyethylene tube (ECO 8)
10000002	Absorption attachment for condenser (ECO 8)
A00001039	Test tube extractor (ECO 25)
CM0091680	Set of 20 glass test tubes with internal Ø 14 mm and external Ø 16 mm (ECO 8 and ECO 25)
CA0091636	Holder for 12 round glass cells (ECO 8 and ECO 25)
10000282	Thermometer 0 - 250 °C
10000667	ECO 8 Retarded fuse 5x20 3A
10000670	ECO 25 Retarded fuse 5x20 8A
10000232	Foot

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)

2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2006/95/EC / Bassa tensione 2006/95/CE / Basse tension 2006/95/CE / Baja tensión 2006/95/CE / Niederspannung 2006/95/EG

- Electromagnetic compatibility directive 2004/108/EC / Compatibilità elettromagnetica 2004/108/CE / Compatibilité électromagnétique 2004/108/CE / Compatibilidad electromagnética 2004/108/CE / Elektromagnetische Verträglichkeit 2004/108/EG

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

Thank you for having chosen VELP!

Since 1983 VELP has offered to professionals in the sector a range of sophisticated and reliable equipment.

VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Our product lines:

Food&Feed Line

DK and DKL Digesters
JP Recirculating water pump for fumes aspiration
SMS Scrubber
UDK Distillation units
SER Solvent extractors
FIWE Raw fiber extractors
Dietary fiber extractors
Oxitest Oxidation test reactor
NDA 701 Elemental analyzer
Consumables

Environment Line

ECO Thermoreactors
BOD. Determination systems
Refrigerated thermostats and incubators
Flocculators
Overhead mixer
Mineralization unit for trace heavy metals determination
Turbidimeter
Radiation detector

Stirring Line

Heating magnetic stirrers
Vertex digital thermoregulator
Ultraflat magnetic stirrer
Magnetic stirrers
Overhead stirrers
Heating plates
Vortex mixers
Homogenizer

Other Lab Solution

Open circulating bath
Recirculating water vacuum pump
Peristaltic pump

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