
# ESTEREOMICROSCOPIOS STEREOMICROSCOPES

## UNA TECHNIC

### Serie 200

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APROBADO POR: \_\_\_\_\_

Thanks for choosing a professional stereomicroscope URA TECHNIC, its high mechanical and optical quality will allow you to have this equipment available for many years. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

In order to use this equipment, it is advisable to follow these points:

- ◆ This manual should be available for all the users of these equipments.  
Any doubt can be solved getting into touch with the URA TECHNIC TECHNICAL Service.
- ◆ As an accuracy instrument, it should be used carefully avoiding the sudden movements and knocks.
- ◆ Never dismantle the different pieces of the stereomicroscope to repair yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.


**¡ATTENTION! IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

## INDEX

1- DESCRIPTION .....	14
2- APLICATIONS OF THE INSTRUMENTS .....	15
3- INSTALLATION .....	15
4- CHARACTERISTICS .....	16
4.1- TECNICAL CHARACTERISTICS .....	16
4.2- OPTICAL CHARACTERISTICS .....	17
5- ASSEMBLY .....	18
6- OPERATING INSTRUCTIONS .....	18-19
7- MAITENANCE .....	19
8- CAUSES OF A DEFECTIVE IMAGE .....	20
9- ACCESSORIES .....	20
WARRANTY INFORMATION.....	21
SYMILARITY DECLARATION "CE".....	22
STANDARD CONFIGURATION.....	13

**1- DESCRIPCIÓN**

- 1- SUBJECTION CLAMPS.
- 2- MECHANICAL STAGE.
- 3- BASE.
- 4- POWER "ON/OFF".
- 5- INCIDENT ILLUMINATOR.
- 6- CONTROL OF APPROACH.
- 7- CONTROL OF HEIGHT ADJUSTMENT BINOCULAR.
- 8- RING OF ADJUSTMENT OF APPROACH TENSION.
- 9- EYEPIECES.
- 10- CONTROL OF ADJUSTMENT DIOPTRIC.
- 11- HEAD.
- 12- OBJETIVE
- 13- SWITCH OF INCIDENT OR TRANSMITED ILLUMINATION

## 2- APPLICATIONS OF THE INSTRUMENT

The stereomicroscope serve to observe images expanded of small objects, being the image created, isn't invest and conserving the original colour of the object observed.

Applications of stereomicroscope:

- ◆ Teaching: basic, secondary, university and professional formation.
- ◆ Industry.
- ◆ Agriculture y stockbreeding.

These stereomicroscopes are easy to service and excellent quality of image.

## 3-INSTALLATION

- ◆ Remove with care the stereomicroscope, lending special attention of catching it for the base or for the backup, never by the head neither by focusing control.
- ◆ Put the stereomicroscope stand on a steady table, maintaining it far away of the direct light of the sun, high temperatures, humidity, dust and vibrations.
- ◆ Remove the paper and the rubbers that protect it of the dust.
- ◆ Remove the protective cover of the objective (model 220).

These stereomicroscope have a Schuko plug. This plug must be connected to a plugging device connected to ground. This grounded device must be at hand to be able to disconnect the plug quickly in case of emergency.

Used always the stereomicroscope dust cover when it is out of use in order to avoid that the dust settles on the optical parts. When it is going to be out of use for a long period (weekend, holidays,...) protect it in its wooden box.

Used always halogen lamps of 12V-10W URA TECHNIC.

Neither the manufacturer nor the supplier will assume any responsibility for the damages caused to the stereomicroscope or for any other damage suffered by any person due to the neglect of the right process of electric connection. The AC must be 220-230V,50 Hz.

## 4- CHARACTERISTICS

### 4.1- TECNICAL CHARACTERISTICS

The series 200 of stereomicroscopes URA TECHNIC is composed for six models that with characteristics:

- ◆ Model 210: The simple model of the range, small stereomicroscope, carried out integrally in resistant plastic, operating with batteries 1.5V, lamp incident light 2.5V, 0.5A, objective 2X, right head.
- ◆ Model 214: A couple objectives interchangeable, right head. Two models, 214/1 with incident light and model 214/2 with incident and transmitted light, halogen lamp of 12 V 10 W.
- ◆ Model 216: A couple objectives interchangeable, inclined 45° head. Two models, 216/1 with incident light and model 216/2 with incident and transmitted light, halogen lamp of 12 V 10 W.
- ◆ Model 218: Two couple of objectives in turret, inclined 45° head. Two models, 218/1 with incident light and model 218/2 with incident and transmitted light, halogen lamp of 12 V 10 W.
- ◆ Model 219: Two couple of objectives in turret, inclined 45° head and rotate 360°. Two models, 219/1 with incident light and model 219/2 with incident and transmitted light, halogen lamp of 12 V 10 W.
- ◆ Model 220: Two couple of objectives in turret, inclined 45° head and rotate 360°, intensity knob, it has light knob in the base where we can select: incident light, transmitted light or incident and transmitted simultaneously, halogen lamp of 12 V 10 W
- ◆ Model 201: objectives couple interchangeable, inclined 45° head.

Others characteristics:

- ◆ BASE: metal, very firm, with transformer for illumination.
- ◆ POWER "ON/OFF".
- ◆ BINOCULAR HEAD: variable in height upon the backup.
- ◆ INTERPUPILARY GRADUATION: with scale between 51 and 75 mm.
- ◆ CORRECTION DIOPTRIC:  $\pm 5$  mm in ocular tube.
- ◆ Two EYEPIECES (10x) With system of fixation.
- ◆ FOCUSING SYSTEM: vertically variable by pinion and rack, with stop of security and regulation of tension.
- ◆ CONTROL OF APPROACH: both-way, with control of friction.
- ◆ STAGE of glass ground of 95 mm de diameter.
- ◆ Two BLINKERS.
- ◆ Couple or two couples objectives standard (2x or 2x and 4x), according to model.
- ◆ INTENSITY OF LIGHT KNOB, according to model.
- ◆ POWER OF LIGHT, according to model.

**NOTE:** According to model, some specifications, cannot be of application in its stereomicroscope.

## 4.2 OPTICAL CHARACTERISTICS:

The following tables show the different combinations that can be obtained with eyepieces and objectives optional.

- ◆ Table 1 models 214 and 216.
- ◆ Table 2 models 218, 219 and 220 the torret are: (2x y 4x).
- ◆ Table 3 model 201.

The working distance is to 57 mm for all models.

### MODELS 214 and 216

OBJECTIVES	EYEPIECES								WORKING DISTANCE
	5X		10X		15X		20X		
	Mag.	Field Ø	Mag.	Field Ø	Mag.	Field Ø	Mag.	Field Ø	
1X	5X	22mm	10X	20mm	15X	15mm	20X	10mm	57mm
2X	10X	12mm	20X	10mm	30X	7,5mm	40X	6,5mm	80mm
3X	15X	8mm	30X	6,7mm	45X	5mm	60X	4,3mm	61mm
4X	20X	6mm	40X	5mm	60X	3,7mm	80X	3,2mm	57mm
6X	30X	4mm	60X	3,3mm	90X	2,5mm	120X	1,8mm	66mm

Ø: Dimeter; Mag.: Total magnification

TABLE 1

### MODELS 218, 219 and 220

OBJECTIVES	EYEPIECES								WORKING DISTANCE
	5X		10X		15X		20X		
	Mag.	Field Ø	Mag.	Field Ø	Mag.	Field Ø	Mag.	Field Ø	
1X	5X	11mm	10X	20mm	15X	15mm	20X	10mm	57mm
2X	10X	11mm	20X	10mm	30X	7,5mm	40X	6,5mm	57mm
1X	5X	22mm	10X	20mm	15X	15mm	20X	10mm	57mm
3X	15X	8mm	30X	6,7mm	45X	5mm	60X	4,3mm	57mm
2X	10X	11mm	20X	10mm	30X	7,5mm	40X	6,5mm	57mm
4X	20X	5,6mm	40X	5mm	60X	3,7mm	80X	3,2mm	57mm

Ø: Dimeter; Mag.: Total magnification

TABLE 2

### MODEL 201

OBJECTIVES	EYEPIECES								WORKING DISTANCE
	5X		10X		15X		20X		
	Mag.	Field Ø	Mag.	Field Ø	Mag.	Field Ø	Mag.	Field Ø	
2X	10X	10mm	20X	10mm	30X	7,5mm	40X	6,5mm	113mm

Ø: Dimeter; Mag.: Total magnification

TABLE 3

## 5- ASSEMBLY

With the stereomicroscope has received a couple of eyepieces (normally 10x). Eyepieces are subjects for a screw. If you want to change the eyepieces by other that are available as option, before should loosen said screw and to substitute them. Return to press it once substituted the eyepieces.

To change the objective of the models 214 and 216, first one must free two small screw that has the head: one in the front and the other in the left side. Later the objective toward the right displaces himself until this is disjointed. The objective is placed desired and return to press the two screws.

To change the objective of the models 218 and 219 only one must rotate the turret.

The objective of the model 220 changes, freeing the four screws that are found in the head. The objective of the model 201 changes of the same form that the models 214 and 216 only that instead of having two screws has one in the left side of the head. But this model comes without mounting. To install it, read attentively the instructions that come subsequently.

### Installation model 201

- ◆ This model has a base in which are inserted the interrupter (On / Off) and the electrical socket of the illuminator.
- ◆ This base has a vertical rod, in him two are inserted ring: one ring that has in its splits internal a furrow by the one that goes lodged the cable of the illuminator. The other ring is the one that has a screw and serves to adjust the height of the body binocular.
- ◆ Upon this second ring is placed a nose that holds the horizontal arm. This nose serves for regular the height and to rotate the body binocular. The horizontal arm is inserted in this nose and is screwed strongly by means of 2 small screw with a key ALLEN nº2.
- ◆ The horizontal arm is of similar form to the vertical rod. In this arm is introduced the other ring with the stub cable.
- ◆ This arm in its extreme has an hole, this hole serves for join its prolongation to do it flexible. This prolongation has a screw which is introduced in the hole of the arm, and is curled with the nut with handle, without forgetting us of placing the to obtain greater washers subjection.
- ◆ In the extreme from the prolongation is placed the head, which is adjusted also by means of two screw ALLEN nº2.
- ◆ The iluminador is inserted in an hole that it has the head. This is a black piece in form of cylinder with different width.

## 6- OPERATING INSTRUTIONS

- ◆ Connect correctly the stereoscope to an electrical socket of 220 V, 50-60 Hz provided of takes of ground.
- ◆ Connect the type of lighting, incident or reflected, more convenient to its observation, optionally exists stage of contrast write /black for use with light reflected.
- ◆ Center the sample in the stage. Rotate the control of approach slowly until obtaining an image sharpness observing for the eyepiece left done not provide of ring of adjustment of dioptric, once focused look at for the eyepiece right, if the image

- is not clear rotate the ring of adjustment of dioptic until obtaining a good focus.
- ◆ With the tool that supplies can adjust the tension of the focusing controls, introducing in the small hole that are found in the axle of the focusing control the extreme of said tool, maintain the motionless command and rotate the tool to enlarge or to diminish the friction that opposes the focusing controls.
  - ◆ For large samples or of more thickness, can be adjusted the height of binocular. It free the screw of the back part and once placed to the necessary height return to set it strongly.
  - ◆ Observe for both eyes at the same time and graduate the distance interpupillary until the images of both eyes are founded correctly in an alone vision

## 7- MAINTENANCE

The maintenance of these stereomicroscopes is the same as any other optical instrument.

- ◆ They must be always clean and covered with the dust cover to avoid that the dust settles on them.
- ◆ The user must never dismantle the lens, if there would be any dirtiness on the outside surfaces of the lens, clean them with a soft cloth without fluff. It can clean it with toluene or xilol.
- ◆ To eliminate the dust settled on the lens, blow with a rubber pear or clean it with a small soft brush.
- ◆ For the cleaning of the mechanical pieces use non-corrosive lubricants, so that the paint does not become damaged. Be careful not to touch the optical pieces. Use always a neutral detergent.
- ◆ If the lamp burns out, change it with an URA halogen lamp 12V 10W, being careful not to touch the bulb with your hands.

**IMPORTANT:** Before changing the lamp bulb or the fuse, make sure that the microscope is disconnected. Make also sure of using halogen lamps 12V 10W, since the use of any other kind of lamps will produce a defective functioning of the instrument.

- ◆ For the change of incident lighting lamp loosen the screw, and replacing it for a new lamp. It be assured to adjust all correctly.
- ◆ For the change of transmitted lighting lamp, remove the screw, loosen the screw, replacing the lamp for a new lamp.
- ◆ If it is necessary to change the fuse you must follow this process: Loosen the protective lid that is on the base of the stereomicroscope (FUSE), change the fuse and place it in the same position, putting the protective lid again. the fuse are of 0.5A, 250 V.

## 8- CAUSES OF A DEFECTIVE IMAGE

In case of a defective image you must check that:

- ◆ The illumination is well made, the light intensity is not excessive enough or too weak.
- ◆ The eyepieces must fit well and the objectives must be well tightened.
- ◆ The cleaning of the whole optical system, following this process: make the eyepieces rotate and if the specks move, clean them. Make the whole of the head rotate. NEVER DISMANTLE THE HEAD. However you can clean carefully those accessible surfaces of the prisms blowing with a rubber pear. Rotate the objective, if the parasitic images turn around clean it with a small dry brush trying to remove the dust, observe the frontal surface with a magnifying glass or an inverted eyepiece.

## 9- ACCESSORIES

The stereomicroscope can be completed with several eyepieces, objectives and others accessories.

Code	Description
<b>Eyepieces</b>	
90200143	Eyepiece 5x
90200146	Eyepiece 15x
90200148	Eyepiece 20x
<b>Objectives</b>	
90200201	Objective 1x
90200203	Objective 3x
90200204	Objective 4x

Others accessories:

Code	Description
90200901	Stage of contrast Write /Black
90200991	Plastic dust cover
90200351	Halogen Lamp 12v 10w
90200902	Two blinkers.
90200903	Stage glass groud
Adjustment of tension tool	

## WARRANTY INFORMATION

AUXILAB, S.L. warrants this equipment for a period of twelve months from the date of purchasing to be free from any defects of manufacturing that affects its functioning, we will replace these pieces which break due to a manufacturing defect. The change of these pieces will be free during the warranty period and will not extend it.

This warranty does not extend to any products which have been subject to misuse, have been manipulated by somebody else other than our Technical Service, or have not been installed properly. This also affects to equipments that have not used the original accessories.

This limited warranty does not extend to any equipments which have been damaged as a result of catastrophic causes ( fire, flood,...), atmospheric, blows or falls.

The costs caused by this repair, the transport of the equipments to our Technical Services and its return will be paid by the user.

It is necessary the presentation of the Warranty Certificate properly filled (without strokes or crossing outs) in order to achieve any free of charge operation of our Technical Service during the warranty period. This warranty does not extend to the maintenance of the equipment or its pieces, or the damages produced by the everyday use of the apparatus.

To provide a better service, when you send an equipment to our Technical Service we would be grateful if you could provide us with the following information.

<b>Supplier:</b>			
<b>Equipment:</b>		<b>Model:</b>	
<b>Serie number:</b>			
<b>Date of purchase:</b>			
<b>Brief description of the damage:</b>			

AUXILAB, S.L. Polígono Morea Norte, 8-31191 BERIAIN (Navarra)  
 Tel(948) 310 513 Fax(948) 310 500 - 312 071  
 e-mail: asistencia@auxilab.com <http://www.auxilab.es>

## **SIMILARITY DECLARATION "CE"**

AUXILAB, S.L. as the owner of the trade mark URA TECHNIC

Declares that the equipments:

Model: PROFESSIONAL MICROSCOPES URA TECHNIC have been manufactured under our checking and they fulfill the following safety norm:

73/23/CEE Low tension norm

By means of the fulfillment of the following international norms:

EN 61010-1: 1993 + A2: 1995  
(IEC 1010-1: 1990 with the correction 1 of 1992 and the correction 2 of 1995)

Having available the similarity declaration approved by a registered organization (SEMKO AB N° 9803068/01 approved in Stockholm 31st March 1998).

Alfonso Ainciburu  
Quality Manager

ENGLISH

**STANDARD CONFIGURATION**

Main unit, European plug, binocular head.

According to the model will be composed by:

Model	214/1	214/2	216/1	216/2	218/1	218/2	219/1	219/2	220/2
Head	Right	Right	Incl.	Incl.	Incl.	Incl.	Incl.	Incl.	Incl.
Torreta	No	No	No	No	Yes	Yes	Yes	Yes	Yes
Objective	2x	2x	2x	2x	2x y 4x	2x y 4x	2x y 4x	2x y 4x	2x y 4x
Illumination*	1	2	1	2	1	2	1	2	2

Incl: Inclined

\*Illumination:

1: Incident (Upper).

2: Double, Incident y transmit (upper and lower).

All Models:

QUANTILY	DESCRIPTION
2	Eyepieces 10x
1	Stage of contrast Write /Black.
1	Adjustment of tension tool.
2	Halogen Lamp 12v 10w (Colocadas en el equipo)
1	Plastic dust cover.
1	Instruction manual, Warranty and CE certificate.
1	Microscopy manual
1	Packing: wooden box and card box.

Others

	Objectives:
	Eyepieces:
	Lamps:
	Others:

MADE BY: \_\_\_\_\_